

**DICOM CONFORMANCE STATEMENT**  
**FOR**  
**VANTAGE-ORIAN™ AND VANTAGE-GALAN™**  
**VERSION V4.5 AND V5.0**

**CANON MEDICAL SYSTEMS CORPORATION**

## **Trademarks**

VANTAGE-ORIAN and VANTAGE-GALAN are trademarks of Canon Medical Systems Corporation.

This document may include trademarks or registered trademarks of other companies.

### **IMPORTANT!**

- (1) No part of this document may be copied or reprinted, in whole or in part, without written permission.
- (2) The contents of this document are subject to change without prior notice and without our legal obligation.
- (3) Please refer to the Canon Medical Systems Corporation website for the most recent version of this conformance statement.

Global: <https://www.medical.canon/Interoperability/DICOM/EN>

Japan: <https://www.medical.canon/Interoperability/DICOM/JP>

## 1. CONFORMANCE STATEMENT OVERVIEW

Table 1-1 provides an overview of the network services supported by VANTAGE-ORIAN™ and VANTAGE-GALAN™

**Table 1-1  
NETWORK SERVICES**

| <b>SOP Classes</b>                            | <b>User of Service (SCU)</b> | <b>Provider of Service (SCP)</b> |
|---|------------------------------|----------------------------------|
| <b>Workflow Management</b>                    |                              |                                  |
| Modality Worklist Information Model – Find    | Yes                          | No                               |
| Modality Performed Procedure Step             | Yes*                         | No                               |
| <b>Transfer</b>                               |                              |                                  |
| MR Image Storage                              | Yes                          | Yes                              |
| Secondary Capture Image Storage               | Yes                          | Yes                              |
| Enhanced MR Image Storage                     | Yes                          | No                               |
| MR Spectroscopy Storage                       | Yes                          | No                               |
| Grayscale Softcopy Presentation State Storage | Yes                          | No                               |
| <b>Storage Commitment</b>                     |                              |                                  |
| Storage Commitment Push Model                 | Yes*                         | Yes                              |
| <b>Query/Retrieve</b>                         |                              |                                  |
| Patient Root Q/R Information Model – Find     | Yes*                         | Yes*                             |
| Patient Root Q/R Information Model – Move     | Yes*                         | Yes*                             |
| Study Root Q/R Information Model – Find       | Yes*                         | Yes*                             |
| Study Root Q/R Information Model – Move       | Yes*                         | Yes*                             |
| <b>Print Management</b>                       |                              |                                  |
| Basic Grayscale Print Management              | Yes                          | No                               |

\*:Option

Table 1-2 provides an overview of the Media Storage Application Profiles supported by VANTAGE-ORIAN™ and VANTAGE-GALAN™.

**Table 1-2  
MEDIA SERVICES**

| <b>Media Storage Application Profile</b> | <b>Write Files (FSC)</b> | <b>Update Files (FSU)</b> | <b>Read Files (FSR)</b> |
|--|--------------------------|---------------------------|-------------------------|
| <b>Compact Disk – Recordable</b>         |                          |                           |                         |
| MR Image CD-R                            | No                       | No                        | Yes                     |
| General Purpose CD-R                     | No                       | No                        | Yes                     |
| <b>DVD – Recordable</b>                  |                          |                           |                         |
| MR Image DVD-R                           | Yes                      | No                        | Yes                     |
| General Purpose DVD-R                    | Yes                      | No                        | Yes                     |
| <b>DVD – Random Access</b>               |                          |                           |                         |
| MR Image DVD-RAM                         | Yes                      | Yes                       | Yes                     |
| General Purpose DVD-RAM                  | Yes                      | No                        | Yes                     |

## 2. TABLE OF CONTENTS

|           |   |           |
|-----------|---|-----------|
| <b>1.</b> | <b>CONFORMANCE STATEMENT OVERVIEW .....</b>                     | <b>i</b>  |
| <b>2.</b> | <b>TABLE OF CONTENTS .....</b>                                  | <b>a</b>  |
| <b>3.</b> | <b>INTRODUCTION .....</b>                                       | <b>1</b>  |
| 3.1       | <b>REVISION HISTORY .....</b>                                   | <b>1</b>  |
| 3.2       | <b>AUDIENCE.....</b>  | <b>1</b>  |
| 3.3       | <b>REMARKS .....</b>  | <b>1</b>  |
| 3.4       | <b>DEFINITIONS, TERMS AND ABBREVIATIONS .....</b>               | <b>2</b>  |
| 3.5       | <b>REFERENCES.....</b>  | <b>2</b>  |
| <b>4.</b> | <b>NETWORKING .....</b>   | <b>3</b>  |
| 4.1       | <b>IMPLEMENTATION MODEL.....</b>                                | <b>3</b>  |
| 4.1.1     | Application Data Flow .....                                     | 3         |
| 4.1.2     | Functional Definition of AEs.....                               | 5         |
| 4.1.3     | Sequencing of Real-World Activities.....                        | 7         |
| 4.2       | <b>AE SPECIFICATIONS .....</b>                                  | <b>9</b>  |
| 4.2.1     | Verification SCU AE Specifications .....                        | 9         |
| 4.2.2     | Verification SCP AE Specifications.....                         | 12        |
| 4.2.3     | MWM SCU AE Specification.....                                   | 16        |
| 4.2.4     | MPPS SCU AE Specification .....                                 | 25        |
| 4.2.5     | Storage SCU AE Specification .....                              | 33        |
| 4.2.6     | Storage Commitment SCU AE Specification .....                   | 37        |
| 4.2.7     | Q/R SCP AE Specification .....                                  | 44        |
| 4.2.8     | Q/R SCU AE Specification.....                                   | 53        |
| 4.2.9     | Storage SCP AE Specification.....                               | 61        |
| 4.2.10    | Storage Commitment SCP AE Specification .....                   | 65        |
| 4.2.11    | Print SCU AE Specification.....                                 | 72        |
| 4.3       | <b>NETWORK INTERFACES.....</b>                                  | <b>81</b> |
| 4.3.1     | Physical Network Interface .....                                | 81        |
| 4.3.2     | Additional Protocols .....                                      | 81        |
| 4.4       | <b>CONFIGURATION.....</b>                                       | <b>82</b> |
| 4.4.1     | AE Title/Presentation Address Mapping .....                     | 82        |
| 4.4.2     | Parameters.....   | 83        |
| <b>5.</b> | <b>MEDIA INTERCHANGE .....</b>                                  | <b>87</b> |
| 5.1       | <b>IMPLEMENTATION MODEL.....</b>                                | <b>87</b> |
| 5.1.1     | Application Data Flow .....                                     | 87        |
| 5.1.2     | Functional Definition of AEs.....                               | 88        |
| 5.1.3     | Sequencing of Real-World Activities.....                        | 89        |
| 5.1.4     | File Meta Information for Implementation Class and Version..... | 89        |
| 5.2       | <b>AE SPECIFICATIONS .....</b>                                  | <b>90</b> |
| 5.2.1     | Offline-Media AE Specification .....                            | 90        |
| 5.3       | <b>MEDIA CONFIGURATION .....</b>                                | <b>91</b> |
| <b>6.</b> | <b>SUPPORT OF CHARACTER SETS .....</b>                          | <b>92</b> |
| <b>7.</b> | <b>SECURITY.....</b>  | <b>93</b> |
| 7.1       | <b>DE-IDENTIFICATION .....</b>                                  | <b>93</b> |

|            |  |            |
|------------|--|------------|
| <b>8.</b>  | <b>ANNEXES .....</b>   | <b>99</b>  |
| <b>8.1</b> | <b>IOD CONTENTS .....</b>                                      | <b>99</b>  |
| 8.1.1      | Created SOP Instances .....                                    | 99         |
| 8.1.2      | Usage of Attributes from received IOD's .....                  | 146        |
| 8.1.3      | Attribute Mapping .....  | 146        |
| 8.1.4      | Coerced/Modified Fields .....                                  | 146        |
| <b>8.2</b> | <b>DATA DICTIONARY OF PRIVATE ATTRIBUTES .....</b>             | <b>147</b> |
| <b>8.3</b> | <b>CONTROLLED TERMINOLOGY AND TEMPLATES.....</b>               | <b>149</b> |
| <b>8.4</b> | <b>GRAYSCALE IMAGE CONSISTENCY .....</b>                       | <b>149</b> |
| <b>8.5</b> | <b>STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES .....</b> | <b>149</b> |
| 8.5.1      | Standard Extended SOP Class - MR Image Storage .....           | 149        |
| <b>8.6</b> | <b>PRIVATE TRANSFER SYNTAXES .....</b>                         | <b>149</b> |

### 3. INTRODUCTION

#### 3.1 REVISION HISTORY

| REV | Date of Issue    | Author                | Description        |
|-----|------------------|-----------------------|--------------------|
| Nil | January 04, 2018 | Canon Medical Systems | Initial version    |
| A   | March 14, 2018   | Canon Medical Systems | Add target version |

#### 3.2 AUDIENCE

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

#### 3.3 REMARKS

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality.

This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended.

The scope of this Conformance Statement is to facilitate communication with Canon Medical Systems and other vendors' Medical equipment. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [DICOM]. However, by itself it is not guaranteed to ensure the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity between Canon Medical Systems and non- Canon Medical Systems equipment.
- Test procedures should be defined to validate the desired level of connectivity.
- The DICOM standard will evolve to meet the users' future requirements. Canon Medical Systems is actively involved in developing the standard further and therefore reserves the right to make changes to its products or to discontinue its delivery.

### 3.4 DEFINITIONS, TERMS AND ABBREVIATIONS

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard.

Abbreviations and terms are as follows:

|              |  |
|--------------|--|
| <b>AE</b>    | Application Entity                                       |
| <b>ASCE</b>  | Association Control Service Element                      |
| <b>CD-R</b>  | Compact Disk Recordable                                  |
| <b>DIMSE</b> | DICOM Message Service Element                            |
| <b>DVD</b>   | A trademark of the DVD Forum that is not an abbreviation |
| <b>FSC</b>   | File-Set Creator   |
| <b>FSU</b>   | File-Set Updater   |
| <b>FSR</b>   | File-Set Reader  |
| <b>GSPS</b>  | Grayscale Softcopy Presentation State                    |
| <b>IE</b>    | Information Entity                                       |
| <b>IOD</b>   | Information Object Definition                            |
| <b>MPPS</b>  | Modality Performed Procedure Step                        |
| <b>MPPSR</b> | Modality Performed Procedure Step Retrieve               |
| <b>MSPS</b>  | Modality Scheduled Procedure Step                        |
| <b>MWM</b>   | Modality Worklist Management                             |
| <b>R</b>     | Required Key Attribute                                   |
| <b>O</b>     | Optional Key Attribute                                   |
| <b>PDU</b>   | Protocol Data Unit                                       |
| <b>SCU</b>   | Service Class User (DICOM client)                        |
| <b>SCP</b>   | Service Class Provider (DICOM server)                    |
| <b>SOP</b>   | Service-Object Pair                                      |
| <b>U</b>     | Unique Key Attribute                                     |
| <b>UID</b>   | Unique Identifier  |

### 3.5 REFERENCES

[DICOM] Digital Imaging and Communications in Medicine (DICOM), NEMA PS 3.1-3.18, 2009



## 4. NETWORKING

### 4.1 IMPLEMENTATION MODEL

#### 4.1.1 Application Data Flow




Figure 4.1-1 APPLICATION DATA FLOW DIAGRAM

- The Verification SCU AE can send verification requests to the specified DICOM destination.
- The Verification SCP AE can respond to C-ECHO requests from the specified DICOM destination.
- The MWM SCU AE receives Worklist information from a remote AE. It is associated with the local real-world activity "Update Worklist". When the "Update Worklist" is performed the MWM SCU AE queries a remote AE for worklist items and provides the set of worklist items matching the query request. "Update Worklist" is performed as a result of an operator request.
- The MPPS SCU AE sends MPPS information to a remote AE. It is associated with the local real-world activity "Acquire Images". When the "Acquire Images" is performed the MPPS SCU AE creates and updates Modality Performed Procedure Step instances managed by a remote AE. Acquisition of images will result in automated creation of an MPPS Instance. Completion of the MPPS is performed as the result of an operator action. If the remote AE is configured as an MPPSR SCP, the MPPS SCU AE can receive MPPSR information.
- The Storage SCU AE sends images, spectroscopy data or GSPS to a remote AE. It is associated with the local real-world activity "Send objects". "Send objects" is performed upon user request for specific images selected.  
If the remote AE is configured as an archive device, the Storage SCU AE will send a storage commitment request to the Storage Commitment SCU AE. But Enhanced MR image and MR Spectroscopy is not sent.
- Receiving the storage commitment request from the Storage SCU AE, the Storage Commitment SCU AE will request Storage Commitment and if a commitment is successfully obtained will record this information in the local database.
- The Q/R SCP AE handles incoming query and retrieve requests issued by a remote AE. It is associated with the local real-world activity "Handle Q/R Requests". "Handle Q/R Requests" handles retrieval requests by issuing a command to the Storage SCU AE to send the requested Images to the destination specified by the remote AE. The Q/R SCP AE functions as an SCP for C-FIND and C-MOVE requests.
- The Q/R SCU AE queries a remote AE for lists of studies, series and images and retrieves selected studies, series or images. It is associated with the local real-world activity "Query and Retrieve Images".
- The Storage SCP AE receives incoming images. It is associated with the local real-world activity "Store Images to the Local File System". "Store Images to the Local File System" stores the received images to the local file system.
- If the Storage Commitment request is received then Storage Commitment SCP AE will immediately check if the referenced Composite SOP Instances are in the local database and return Report Notification.
- The Print SCU AE prints images on a remote AE (Printer). It is associated with the local real-world activity "Send Images & Print Management Information". "Send Images & Print Management Information" creates a print-job within the print queue containing one or more virtual film sheets composed from images selected by the user.

## **4.1.2 Functional Definition of AEs**

### **4.1.2.1 Functional Definition of Verification SCU AE**

The existence of a verification request with associated network destination will activate the Verification SCU AE. An Association request is sent to the destination AE and upon successful negotiation of a Presentation Context the verification is started. If the verification fails, the Verification SCU AE will show the error message.

### **4.1.2.2 Functional Definition of Verification SCP AE**

The Verification SCP AE waits for another application to connect at the presentation address configured for its Application Entity Title. The Verification SCP AE will accept Associations with Presentation Contexts for SOP Classes of the Verification Service Classes.

### **4.1.2.3 Functional Definition of MWM SCU AE**

The MWM SCU AE attempts to download a worklist from a remote node. If the MWM SCU AE establishes an Association to a remote AE, it will transfer patient's information and worklist items via the open Association. The results will be displayed in a separate list. The patient's information will be used for the patient registration.

### **4.1.2.4 Functional Definition of MPPS SCU AE**

The MPPS SCU AE performs the creation of an MPPS Instance automatically whenever images are acquired. Further updates on the MPPS data can be performed automatically or interactively. If the remote AE is configured as an MPPSR SCP, the MPPS SCU AE can receive MPPSR information.

### **4.1.2.5 Functional Definition of Storage SCU AE**

The existence of a send-job queue entry with associated network destination will activate the Storage SCU AE. An Association request is sent to the destination AE and upon successful negotiation of a Presentation Context the image, MRSpectroscopy data & GSPS transfer is started. If the image, MRSpectroscopy data & GSPS transfer fails, the Storage SCU AE will retry this send-job automatically. If the remote AE is configured as an archive device, the storage SCU AE will send a storage commitment request to the Storage Commitment SCU AE.

### **4.1.2.6 Functional Definition of Storage Commitment SCU AE**

Receiving the storage commitment request from the Storage SCU AE, the Storage Commitment SCU AE will request Storage Commitment and if a commitment is successfully obtained will record this information in the local database.

### **4.1.2.7 Functional Definition of Q/R SCP AE**

The Q/R SCP AE waits for another application to connect at the presentation address configured for its Application Entity Title. The Q/R SCP AE will accept Associations with Presentation Contexts for SOP Class of the Query/Retrieve Service Class. It will handle query and retrieve requests on these Presentation Contexts and respond with data objects with values corresponding to the contents of the local file system. When a retrieval request is received, the Q/R SCP AE issues a command to the Storage SCU AE to send the specified images to the destination.

### **4.1.2.8 Functional Definition of Q/R SCU AE**

The Q/R SCU AE is activated when the user selects a remote node to query and enters some key information, Patient's Name, Patient ID and/or Study Date. The user can select studies, series and images to be retrieved. The images will be received at the Storage SCP AE.

### **4.1.2.9 Functional Definition of Storage SCP AE**

The Storage SCP AE waits for another application to connect at the presentation address configured for its Application Entity Title. The Storage SCP AE will accept Associations with Presentation Contexts for SOP Classes of the Storage Service Classes. Any images received on such Presentation Contexts will be stored to the local file system.


#### **4.1.2.10 Functional Definition of Storage Commitment SCP AE**

If the Storage Commitment request is received then Storage Commitment SCP AE will immediately check if the referenced Composite SOP Instances are in the local database and return Report Notification.

#### **4.1.2.11 Functional Definition of Print SCU AE**

The existence of a print-job in the print queue will activate the Print SCU AE. An Association is established with the printer and the printer's status determined. If the printer is operating normally, the film sheets described within the print-job will be printed. If the printer is not operating normally, an error message will be displayed and this print-job can be canceled or restarted by the user operations.

### 4.1.3 Sequencing of Real-World Activities



**Figure 4.1-2  
SEQUENCING CONSTRAINTS**

Under typical scheduled workflow conditions the sequencing constraints illustrated in Figure 4.1-2 apply:

1. Query Worklist
2. Receive Worklist of Modality Scheduled Procedure Steps (MSPS)
3. Select Workitem (MSPS) from Worklist
4. Start Acquisition and Create MPPS
5. Acquire Images
6. Complete Acquisition and Finalize MPPS
7. Store Acquired Images, Spectroscopy data and any associated Grayscale Softcopy Presentation State (GSPS) instances
8. If the Image Manager is configured as an archive device the AE will request Storage Commitment for the images and associated GSPS instances
9. Q/R Requests by a Remote AE
10. Send Images
11. Query/Retrieve Images
12. Receive Images
13. Receive Commitment Request
14. Return Commitment Report
15. Print Images

Other workflow situations (e.g. unscheduled procedure steps) will have other sequencing constraints. Some activities may be omitted according to situations.

## 4.2 AE SPECIFICATIONS

### 4.2.1 Verification SCU AE Specifications

#### 4.2.1.1 SOP Class

The Verification SCU AE provides Standard Conformance to the following DICOM SOP classes:

**Table 4.2-1**

**SOP CLASSES FOR VERIFICATION SCU AE**

| SOP Class Name | SOP Class UID     | SCU | SCP |
|----------------|-------------------|-----|-----|
| Verification   | 1.2.840.10008.1.1 | Yes | No  |

#### 4.2.1.2 Association Policies

##### 4.2.1.2.1 General

The Verification SCU AE can form associations via user control. The Verification SCU AE can only request the starting of an association. It cannot accept association start requests from external application entities. The DICOM standard application context name for DICOM is always accepted.

**Table 4.2-2**

**DICOM APPLICATION CONTEXT FOR VERIFICATION SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

##### 4.2.1.2.2 Number of Associations

The Verification SCU AE initiates one Association at a time.

**Table 4.2-3**

**NUMBER OF ASSOCIATIONS INITIATED FOR VERIFICATION SCU AE**

|   |   |
|---|---|
| Maximum Number of Simultaneous Associations | 1 |
|---|---|

##### 4.2.1.2.3 Asynchronous Nature

The Verification SCU AE does not support asynchronous communication (multiple incomplete transactions on a single association). All association requests must be completed and confirmed before new actions can be performed.

**Table 4.2-4**

**ASYNCHRONOUS NATURE FOR VERIFICATION SCU AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

##### 4.2.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-5**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR VERIFICATION SCU AE**


|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.1.3 Association Initiation Policy

#### 4.2.1.3.1 Activity - Verify Connectivity

##### 4.2.1.3.1.1 Description and Sequencing of Activity

The Verification SCU AE initiates association through user control.



**Figure 4.2-1**  
**SEQUENCING OF ACTIVITY – VERIFY CONNECTIVITY**

The following sequencing restrictions, illustrated in figure above, apply when the Verification SCU AE:

1. The Verification SCU AE opens a new association with the specified destination AE.
2. The Verification SCU AE sends C-ECHO requests.
3. The Verification SCU AE closes the Association.



#### 4.2.1.3.1.2 Proposed Presentation Contexts

The Verification SCU AE is capable of proposing the Presentation Contexts shown in the following table:

**Table 4.2-6  
PROPOSED PRESENTATION CONTEXTS FOR VERIFICATION SCU AE**

| Presentation Context Table |                   |                           |                     |      |           |
|----------------------------|-------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                   | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID               | Name                      | UID                 |      |           |
| Verification               | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                            |                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

#### 4.2.1.3.1.3 SOP Specific Conformance for Verification SOP Class

The Verification SCU AE provides standard conformance to the Verification Service Class as an SCU.

The Verification SCU AE performs the following actions based on the status code values in the C-ECHO responses from the destination C-ECHO SCP:

**Table 4.2-7  
VERIFICATION SCU AE C-ECHO RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Status Code                       | Behavior  |
|----------------|-----------------|-----------------------------------|---|
| Success        | Success         | 0000                              | A message is sent to the User Interface.  |
| Error          | Failure         | Status codes other than the above | The association is terminated when the error occurs. A failure message is output to the Service Log. A message is sent to the user interface. |

#### 4.2.1.4 Association Acceptance Policy

The Verification SCU AE does not accept associations.

## 4.2.2 Verification SCP AE Specifications

### 4.2.2.1 SOP Class

The Verification SCP AE provides Standard Conformance to the following DICOM SOP classes:

**Table 4.2-8**  
**SOP CLASSES FOR VERIFICATION SCP AE**

| SOP Class Name | SOP Class UID     | SCU | SCP |
|----------------|-------------------|-----|-----|
| Verification   | 1.2.840.10008.1.1 | No  | Yes |

### 4.2.2.2 Association Policies

#### 4.2.2.2.1 General

The DICOM standard application context name for DICOM is always accepted.

**Table 4.2-9**  
**DICOM APPLICATION CONTEXT FOR VERIFICATION SCP AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.2.2.2 Number of Associations

The Verification SCP AE accepts one Association at a time.

**Table 4.2-10**  
**NUMBER OF ASSOCIATIONS ACCEPTED FOR VERIFICATION SCP AE**

|   |   |
|---|---|
| Maximum Number of Simultaneous Associations | 1 |
|---|---|

#### 4.2.2.2.3 Asynchronous Nature

The Verification SCP AE does not support asynchronous communication (multiple incomplete transactions on a single association).

**Table 4.2-11**  
**ASYNCHRONOUS NATURE FOR VERIFICATION SCP AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.2.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-12**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR VERIFICATION SCP AE**

|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.2.3 Association Initiation Policy


The Verification SCP AE does not initiate Associations.

## 4.2.2.4 Association Acceptance Policy

### 4.2.2.4.1 Activity - Verify Connectivity

#### 4.2.2.4.1.1 Description and Sequencing of Activity

The Verification SCP AE accepts Associations only if they have valid Presentation Contexts. If none of the requested Presentation Contexts are accepted then the Association Request itself is rejected. It can be configured to only accept Associations with certain hosts (using TCP/IP address) and/or Application Entity Titles.



**Figure 4.2-2**  
**SEQUENCING OF ACTIVITY – VERIFY CONNECTIVITY**

A possible sequence of interactions between the Verification SCP AE and a Verification SCU AE is illustrated in the Figure above:

- The Verification SCU AE opens a new association with the Verification SCP AE.
- The Verification SCU AE sends C-ECHO requests. The Verification SCP AE replies with a C-ECHO response (status success).
- The Verification SCU AE closes the Association.

The Verification SCP AE may reject association attempts as shown in the Table below. The Result, Source and Reason/Diag columns represent the values returned in the appropriate fields of an ASSOCIATE-RJ PDU. The contents of the Source column is abbreviated to save space and the meaning of the abbreviations are:

- DICOM UL service-user
- DICOM UL service-provider (ASCE related function)
- DICOM UL service-provider (Presentation related function)

**Table 4.2-13  
ASSOCIATION REJECTION REASONS**

| <b>Result</b>             | <b>Source</b> | <b>Reason/Diag</b>                         | <b>Explanation</b>   |
|---------------------------|---------------|--|--|
| 2 –<br>rejected-transient | c             | 2 – local-limit-exceeded                   | The (configurable) maximum number of simultaneous Associations has been reached. An Association request with the same parameters may succeed at a later time.  |
| 2 –<br>rejected-transient | c             | 1 – temporary-congestion                   | No Associations can be accepted at this time due to the real-time requirements of higher priority activities (e.g. during image acquisition no Associations will be accepted) or because insufficient resources are available (e.g. memory, processes, threads). An Association request with the same parameters may succeed at a later time.                        |
| 1 –<br>rejected-permanent | a             | 2 – application-context-name-not-supported | The Association request contained an unsupported Application Context Name. An association request with the same parameters will not succeed at a later time.   |
| 1 –<br>rejected-permanent | a             | 7 – called-AE-title-not-recognized         | The Association request contained an unrecognized Called AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association initiator is incorrectly configured and attempts to address the Association acceptor using the wrong AE Title. |
| 1 –<br>rejected-permanent | a             | 3 – calling-AE-title-not-recognized        | The Association request contained an unrecognized Calling AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association acceptor has not been configured to recognize the AE Title of the Association initiator.                      |
| 1 –<br>rejected-permanent | b             | 1 – no-reason-given                        | The Association request could not be parsed. An Association request with the same format will not succeed at a later time.   |

#### 4.2.2.4.1.2 Accepted Presentation Contexts

The Verification SCP AE will prefer to select the Explicit VR Little Endian Transfer Syntax if multiple transfer syntaxes are offered.

Any of the Presentation Contexts shown in the following table are acceptable to the Verification SCP AE.

**Table 4.2-14**

#### **ACCEPTED PRESENTATION CONTEXTS BY VERIFICATION SCP AE**

| <b>Presentation Context Table</b> |                   |                           |                     |             |                  |
|-----------------------------------|-------------------|---------------------------|---------------------|-------------|------------------|
| <b>Abstract Syntax</b>            |                   | <b>Transfer Syntax</b>    |                     | <b>Role</b> | <b>Ext. Neg.</b> |
| <b>Name</b>                       | <b>UID</b>        | <b>Name</b>               | <b>UID</b>          |             |                  |
| Verification                      | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCP         | None             |
|                                   |                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |             |                  |

#### 4.2.2.4.1.3 SOP Specific Conformance for Verification SOP Class

The Verification SCP AE provides standard conformance to the Verification SOP Class as an SCP.

## 4.2.3 MWM SCU AE Specification

### 4.2.3.1 SOP Classes

The MWM SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-15**  
**SOP CLASSES FOR MWM SCU AE**

| SOP Class Name                             | SOP Class UID          | SCU | SCP |
|--|------------------------|-----|-----|
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Yes | No  |

### 4.2.3.2 Association Policies

#### 4.2.3.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-16**  
**DICOM APPLICATION CONTEXT FOR MWM SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.3.2.2 Number of Associations

The MWM SCU AE initiates one Association at a time for a Worklist request.

**Table 4.2-17**  
**NUMBER OF ASSOCIATIONS INITIATED FOR MWM SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.3.2.3 Asynchronous Nature

The MWM SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-18**  
**ASYNCHRONOUS NATURE FOR MWM SCU AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.3.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-19**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR MWM SCU AE**

|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.3.3 Association Initiation Policy


#### 4.2.3.3.1 Activity – Update Worklist

##### 4.2.3.3.1.1 Description and Sequencing of Activities

The request for an “Update Worklist” is initiated by user interaction, i.e. pressing the buttons “Get Worklist” or automatically at the time of key-word change.

Upon initiation of the request, the MWM SCU AE will build an Identifier for the C-FIND request, will initiate an Association to send the request and will wait for Worklist responses. After retrieval of all responses, the MWM SCU AE will access the local database to add or update patient demographic data. The results will be displayed in a separate list.

The MWM SCU AE will initiate an Association in order to issue a C-FIND request according to the Modality Worklist Information Model.



**Figure 4.2-3**  
**SEQUENCING OF ACTIVITY – UPDATE WORKLIST**

A possible sequence of interactions between the MWM SCU AE and a Department Scheduler (e.g. a device such as a RIS or HIS which supports the Modality Worklist SOP Class as an SCP) is illustrated in the Figure above:

1. The MWM SCU AE opens an association with the Department Scheduler
2. The MWM SCU AE sends a C-FIND request to the Department Scheduler containing the Worklist Query attributes.
3. The Department Scheduler returns a C-FIND response containing the requested attributes of the first matching Worklist Item.
4. The Department Scheduler returns another C-FIND response containing the requested attributes of the second matching Worklist Item.
5. The Department Scheduler returns another C-FIND response with status Success indicating that no further matching Worklist Items exist. This example assumes that only 2 Worklist items match the Worklist Query.
6. The MWM SCU AE closes the association with the Department Scheduler.
7. The user selects a Worklist Item from the Worklist and prepares to acquire new images.

#### 4.2.3.3.1.2 Proposed Presentation Contexts

The MWM SCU AE will propose Presentation Contexts as shown in the following table:

**Table 4.2-20**

#### **PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY UPDATE WORKLIST**

| <b>Presentation Context Table</b>          |                 |                           |                     |             |                  |
|--|-----------------|---------------------------|---------------------|-------------|------------------|
| <b>Abstract Syntax</b>                     |                 | <b>Transfer Syntax</b>    |                     | <b>Role</b> | <b>Ext. Neg.</b> |
| <b>Name</b>                                | <b>UID</b>      | <b>Name</b>               | <b>UID</b>          |             |                  |
| Modality Worklist Information Model – FIND | 1.2.840.10008.5 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU         | None             |
|  | .1.4.31         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |             |                  |



#### 4.2.3.3.1.3 SOP Specific Conformance for Modality Worklist SOP Class

The MWM SCU AE provides standard conformance to the Modality Worklist SOP Class as an SCU.

The behavior of the MWM SCU when encountering status codes in a Modality Worklist C-FIND response is summarized in the Table below.

**Table 4.2-21  
MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning  | Status Code            | Behavior   |
|----------------|--|------------------------|--|
| Success        | Matching is complete   | 0000                   | The SCP has completed the matches. Worklist items are available for display or further processing.   |
| Refused        | Out of Resources   | A700                   | The Association is aborted using A-ABORT and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged.   |
| Failed         | Identifier does not match SOP Class  | A900                   | The Association is aborted using A-ABORT and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged.   |
| Failed         | Unable to Process  | Cxxx                   | The Association is aborted using A-ABORT and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged.   |
| Cancel         | Matching terminated due to Cancel request  | FE00                   | If the query was cancelled due to too many worklist items then the SCP has completed the matches. Worklist items are available for display or further processing. Otherwise, the Association is aborted using A-ABORT and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. |
| Pending        | Matches are continuing   | FF00                   | The worklist item contained in the Identifier is collected for later display or further processing.  |
| Pending        | Matches are continuing – Warning that one or more Optional Keys were not supported | FF01                   | The worklist item contained in the Identifier is collected for later display or further processing. The status meaning is logged only once for each C-FIND operation.  |
| *              | *  | Any other status code. | The Association is aborted using A-ABORT and the worklist is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged.   |

The behavior of the MWM SCU AE during communication failure is summarized in the Table below.

**Table 4.2-22**  
**MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior  |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and the worklist query is marked as failed. The reason is logged and reported to the user if an interactive query. |
| Association aborted by the SCP or network layers | The worklist query is marked as failed. The reason is logged and reported to the user if an interactive query.  |

Acquired images will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

The Table below provides a description of the MWM SCU AE Worklist Request Identifier and specifies the attributes that are copied into the images. Unexpected attributes returned in a C-FIND response are ignored.

**Table 4.2-23**  
**WORKLIST REQUEST IDENTIFIER**

| Module Name<br>Attribute Name           | Tag         | VR | M   | R | D | IOD |
|---|-------------|----|-----|---|---|-----|
| <b>SOP Common</b>                       |             |    |     |   |   |     |
| Specific Character Set                  | (0008,0005) | CS |     | x |   | x   |
| <b>Scheduled Procedure Step</b>         |             |    |     |   |   |     |
| Scheduled Procedure Step Sequence       | (0040,0100) | SQ |     |   |   |     |
| > Scheduled Station AE Title            | (0040,0001) | AE | S   | x |   |     |
| > Scheduled Station Name                | (0040,0010) | SH |     |   |   |     |
| > Scheduled Procedure Step Location     | (0040,0011) | SH |     |   |   |     |
| > Scheduled Procedure Step Start Date   | (0040,0002) | DA | R,S | x | x |     |
| > Scheduled Procedure Step Start Time   | (0040,0003) | TM | R   | x | x |     |
| > Scheduled Procedure Step End Date     | (0040,0004) | DA |     |   |   |     |
| > Scheduled Procedure Step End Time     | (0040,0005) | TM |     |   |   |     |
| > Scheduled Performing Physician's Name | (0040,0006) | PN |     | x |   |     |
| > Scheduled Procedure Step Description  | (0040,0007) | LO |     | x | x | x   |
| > Scheduled Protocol Code Sequence      | (0040,0008) | SQ |     |   | x | x   |
| >>Code Value                            | (0008,0100) | SH |     |   |   |     |
| >>Coding Scheme Designator              | (0008,0102) | SH |     |   |   |     |
| >>Coding Scheme Version                 | (0008,0103) | SH |     |   |   |     |
| >>Code Meaning                          | (0008,0104) | LO |     |   |   |     |
| > Scheduled Procedure Step ID           | (0040,0009) | SH |     | x | x | x   |
| > Scheduled Procedure Step Status       | (0040,0020) | CS |     |   |   |     |
| > Comments on Scheduled Procedure Step  | (0040,0400) | LT |     |   |   |     |
| > Modality                              | (0008,0060) | CS | S   | x | x | x   |
| > Requested Contrast Agent              | (0032,1070) | LO |     |   |   |     |
| > Pre-Medication                        | (0040,0012) | LO |     |   |   |     |

| <b>Requested Procedure</b>                 |             |    |     |   |   |   |
|--|-------------|----|-----|---|---|---|
| Requested Procedure ID                     | (0040,1001) | SH | S,* | x | x | x |
| Reason for the Requested Procedure         | (0040,1002) | LO |     |   |   |   |
| Requested Procedure Comments               | (0040,1400) | LT |     |   |   |   |
| Requested Procedure Code Sequence          | (0032,1064) | SQ |     |   |   |   |
| >Code Value                                | (0008,0100) | SH |     |   |   |   |
| >Coding Scheme Designator                  | (0008,0102) | SH |     |   |   |   |
| >Coding Scheme Version                     | (0008,0103) | SH |     |   |   |   |
| >Code Meaning                              | (0008,0104) | LO |     |   |   |   |
| Requested Procedure Description            | (0032,1060) | LO |     | x | x |   |
| Study Instance UID                         | (0020,000D) | UI |     | x |   | x |
| Referenced Study Sequence                  | (0008,1110) | SQ |     |   |   | x |
| >Referenced SOP Class UID                  | (0008,1150) | UI |     | x |   | x |
| >Referenced SOP Instance UID               | (0008,1155) | UI |     | x |   | x |
| Requested Procedure Priority               | (0040,1003) | SH |     |   |   |   |
| Patient Transport Arrangements             | (0040,1004) | LO |     |   |   |   |
| Requested Procedure Location               | (0040,1005) | LO |     |   |   |   |
| Confidentiality Code                       | (0040,1008) | LO |     |   |   |   |
| Reporting Priority                         | (0040,1009) | SH |     |   |   |   |
| Names of Intended Recipients of Results    | (0040,1010) | PN |     |   |   |   |
| <b>Imaging Service Request</b>             |             |    |     |   |   |   |
| Reason for the Imaging Service Request     | (0040,2001) | LO |     |   |   |   |
| Imaging Service Request Comments           | (0040,2400) | LT |     |   |   |   |
| Requesting Physician                       | (0032,1032) | PN |     |   | x | x |
| Referring Physician's Name                 | (0008,0090) | PN |     | x |   |   |
| Requesting Service                         | (0032,1033) | LO |     |   |   |   |
| Accession Number                           | (0008,0050) | SH | S,* | x | x | x |
| Issuer of Accession Number Sequence        | (0008,0051) | SQ |     |   |   |   |
| >Local Namespace Entity ID                 | (0040,0031) | UT |     |   |   |   |
| >Universal Entity ID                       | (0040,0032) | UT |     |   |   |   |
| >Universal Entity ID Type                  | (0040,0033) | CS |     |   |   |   |
| Issue Date of Imaging Service Request      | (0040,2004) | DA |     |   |   |   |
| Issue Time of Imaging Service Request      | (0040,2005) | TM |     |   |   |   |
| Order Entered By                           | (0040,2008) | PN |     |   |   |   |
| Order Enterer's Location                   | (0040,2009) | SH |     |   |   |   |
| Order Callback Phone Number                | (0040,2010) | SH |     |   |   |   |
| Placer Order Number                        | (0040,2016) | LO |     |   |   |   |
| Filler Order Number                        | (0040,2017) | LO |     |   |   |   |
| Confidentiality Constraint on Patient Data | (0040,3001) | LO |     |   |   |   |
| <b>Visit Relationship</b>                  |             |    |     |   |   |   |
| Referenced Patient Sequence                | (0008,1120) | SQ |     |   |   |   |
| >Referenced SOP Class UID                  | (0008,1150) | UI |     |   |   |   |
| >Referenced SOP Instance UID               | (0008,1155) | UI |     |   |   |   |
| <b>Visit Identification</b>                |             |    |     |   |   |   |
| Institution Name                           | (0008,0080) | LO |     |   |   | x |
| Institution Address                        | (0008,0081) | ST |     |   |   |   |
| Institution Code Sequence                  | (0008,0082) | SQ |     |   |   |   |
| >Code Value                                | (0008,0100) | SH |     |   |   |   |
| >Coding Scheme Designator                  | (0008,0102) | SH |     |   |   |   |
| >Coding Scheme Version                     | (0008,0103) | SH |     |   |   |   |
| >Code Meaning                              | (0008,0104) | LO |     |   |   |   |
| Admission ID                               | (0038,0010) | LO |     |   |   |   |
| Issuer of Admission ID                     | (0038,0011) | LO |     |   |   |   |

| <b>Visit Status</b>                      |             |    |   |   |   |   |
|--|-------------|----|---|---|---|---|
| Visit Status ID                          | (0038,0008) | CS |   |   |   |   |
| Current Patient Location                 | (0038,0300) | LO |   |   |   |   |
| Patient's Institution Residence          | (0038,0400) | LO |   |   |   |   |
| Visit Comments                           | (0038,4000) | LT |   |   |   |   |
| <b>Visit Admission</b>                   |             |    |   |   |   |   |
| Referring Physician's Address            | (0008,0092) | ST |   |   |   |   |
| Referring Physician's Telephone Number   | (0008,0094) | SH |   |   |   |   |
| Admitting Diagnosis Description          | (0008,1080) | LO |   |   |   |   |
| Admitting Diagnosis Code Sequence        | (0008,1084) | SQ |   |   |   |   |
| >Code Value                              | (0008,0100) | SH |   |   |   |   |
| >Coding Scheme Designator                | (0008,0102) | SH |   |   |   |   |
| >Coding Scheme Version                   | (0008,0103) | SH |   |   |   |   |
| >Code Meaning                            | (0008,0104) | LO |   |   |   |   |
| Route of Admissions                      | (0038,0016) | LO |   |   |   |   |
| Admitting Date                           | (0038,0020) | DA |   |   |   |   |
| Admitting Time                           | (0038,0021) | TM |   |   |   |   |
| <b>Patient Relationship</b>              |             |    |   |   |   |   |
| Referenced Patient Alias Sequence        | (0038,0004) | SQ |   |   |   |   |
| >Referenced SOP Class UID                | (0008,1150) | UI |   |   |   |   |
| >Referenced SOP Instance UID             | (0008,1155) | UI |   |   |   |   |
| <b>Patient Identification</b>            |             |    |   |   |   |   |
| Patient's Name                           | (0010,0010) | PN | * | x | x | x |
| Patient ID                               | (0010,0020) | LO | S | x | x | x |
| Issuer of Patient ID                     | (0010,0021) | LO |   |   |   |   |
| Issuer of Patient ID Qualifiers Sequence | (0010,0024) | SQ |   |   |   |   |
| >Universal Entity ID                     | (0040,0032) | UT |   |   |   |   |
| >Universal Entity ID Type                | (0040,0033) | CS |   |   |   |   |
| >Identifier Type Code                    | (0040,0035) | CS |   |   |   |   |
| Other Patient IDs                        | (0010,1000) | LO |   |   |   |   |
| Other Patient Names                      | (0010,1001) | PN |   |   |   |   |
| Patient's Birth Name                     | (0010,1005) | PN |   |   |   |   |
| Patient's Mother's Birth Name            | (0010,1060) | PN |   |   |   |   |
| Medical Record Locator                   | (0010,1090) | LO |   |   |   |   |
| <b>Patient Demographic</b>               |             |    |   |   |   |   |

|  |             |    |  |   |   |   |
|--|-------------|----|--|---|---|---|
| Patient's Age  | (0010,1010) | AS |  |   | x | x |
| Occupation   | (0010,2180) | SH |  |   |   |   |
| Confidentiality Constraint on Patient Data Description | (0040,3001) | LO |  |   |   |   |
| Patient's Birth Date                                   | (0010,0030) | DA |  | x | x | x |
| Patient's Birth Time                                   | (0010,0032) | TM |  |   |   |   |
| Patient's Sex  | (0010,0040) | CS |  | x | x | x |
| Patient's Insurance Plan Code Sequence                 | (0010,0050) | SQ |  |   |   |   |
| >Code Value  | (0008,0100) | SH |  |   |   |   |
| >Coding Scheme Designator                              | (0008,0102) | SH |  |   |   |   |
| >Coding Scheme Version                                 | (0008,0103) | SH |  |   |   |   |
| >Code Meaning  | (0008,0104) | LO |  |   |   |   |
| Patient's Size   | (0010,1020) | DS |  | x | x | x |
| Patient's Weight                                       | (0010,1030) | DS |  | x | x | x |
| Patient's Address                                      | (0010,1040) | LO |  |   |   |   |
| Military Rank  | (0010,1080) | LO |  |   |   |   |
| Branch of Service                                      | (0010,1081) | LO |  |   |   |   |
| Country Residence                                      | (0010,2150) | LO |  |   |   |   |
| Region of Residence                                    | (0010,2152) | LO |  |   |   |   |
| Patient's Telephone Numbers                            | (0010,2154) | SH |  |   |   |   |
| Ethnic Group   | (0010,2160) | SH |  |   |   |   |
| Patient's Religious Reference                          | (0010,21F0) | LO |  |   |   |   |
| Patient Comment  | (0010,4000) | LT |  | x | x | x |
| <b>Patient Medical</b>                                 |             |    |  |   |   |   |
| Medical Alerts   | (0010,2000) | LO |  | x | x |   |
| Allergies  | (0010,2110) | LO |  |   |   |   |
| Smoking Status   | (0010,21A0) | CS |  |   |   |   |
| Additional Patient History                             | (0010,21B0) | LT |  |   |   |   |
| Pregnancy Status                                       | (0010,21C0) | US |  | x | x | x |
| Last Menstrual Date                                    | (0010,21D0) | DA |  |   |   |   |
| Special Needs  | (0038,0050) | LO |  |   |   |   |
| Patient State  | (0038,0500) | LO |  |   |   |   |
| <b>Other Extended Attributes</b>                       |             |    |  |   |   |   |
| Study Description                                      | (0008,1030) | LO |  |   | x | x |
| Institutional Department Name                          | (0008,1040) | LO |  |   | x | x |
| Performing Physician's Name                            | (0008,1050) | PN |  |   | x | x |
| Name of Physician(s) Reading Study                     | (0008,1060) | PN |  |   | x | x |
| Operator's Name  | (0008,1070) | PN |  |   | x | x |

The above table should be read as follows:

- Module Name: The name of the associated module for supported worklist attributes.
- Attribute Name: Attributes supported to build the MWM SCU AE Worklist Request Identifier.
- Tag: DICOM tag for this attribute.
- VR: DICOM VR for this attribute.
- M: Matching keys for (automatic) Worklist Update. An "S" will indicate that the MWM SCU AE will supply an attribute value for Single Value Matching, a "R" will indicate Range Matching, an "\*" will indicate Wildcard Matching.
- R: Return keys. An "x" will indicate that the MWM SCU AE will supply this attribute as Return Key with zero length for Universal Matching. This setting can be configured using the service tool. The system's default setting is described in the above table.
- D: Displayed keys. An "x" indicates that this worklist attribute is displayed to the user during a patient registration dialog. For example, Patient Name will be displayed when registering the patient prior to an examination. This setting can be configured using the service tool. The system's default setting is described in the above table.
- IOD: An "x" indicates that this Worklist attribute is included into all Object Instances created during performance of the related Procedure Step. This setting can be configured using the service tool. The system's default setting is described in the above table.

#### **4.2.3.4 Association Acceptance Policy**

The MWM SCU AE does not accept Associations.

## 4.2.4 MPPS SCU AE Specification

### 4.2.4.1 SOP Classes

The MPPS SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-24**  
**SOP CLASSES FOR MPPS SCU AE**

| SOP Class Name                    | SOP Class UID           | SCU | SCP |
|-----------------------------------|-------------------------|-----|-----|
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Yes | No  |

### 4.2.4.2 Association Policies

#### 4.2.4.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-25**  
**DICOM APPLICATION CONTEXT FOR MPPS SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.4.2.2 Number of Associations

The MPPS SCU AE initiates one Association at a time.

**Table 4.2-26**  
**NUMBER OF ASSOCIATIONS INITIATED FOR MPPS SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.4.2.3 Asynchronous Nature

The MPPS SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-27**  
**ASYNCHRONOUS NATURE FOR MPPS SCU AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.4.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-28**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR MPPS SCU AE**

|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.4.3 Association Initiation Policy


#### 4.2.4.3.1 Activity – Acquire Images

##### 4.2.4.3.1.1 Description and Sequencing of Activities

The MPPS SCU AE performs the creation of a MPPS Instance automatically whenever images are acquired. Further updates on the MPPS data can be performed automatically or interactively.

The MPPS SCU AE will initiate an Association to issue an:

- N-CREATE request according to the CREATE Modality Performed Procedure Step SOP Instance operation, or an:
- N-SET request to update the contents and state of the MPPS according to the SET Modality Performed Procedure Step Information operation.



**Figure 4.2-4**  
**SEQUENCING OF ACTIVITY – ACQUIRE IMAGES**



A possible sequence of interactions between the MPPS SCU AE and a Department Scheduler (e.g. a device such as a RIS or HIS which supports the MPPS SOP Class and MPPSR SCP Class as an SCP) is illustrated in the Figure above:

1. The MPPS SCU AE opens an association with the Department Scheduler
2. The MPPS SCU AE sends an N-CREATE request to the Department Scheduler to create an MPPS instance with status of "IN PROGRESS" and create all necessary attributes. The Department Scheduler acknowledges the MPPS creation with an N-CREATE response (status success).
3. The MPPS SCU AE closes the association with the Department Scheduler.
4. All images are acquired and stored in the local database. (Figure 4.2-4)
5. The MPPS SCU AE opens an association with the Department Scheduler.
6. The MPPS SCU AE sends an N-SET request to the Department Scheduler to update the MPPS instance with status of "COMPLETED" and set all necessary attributes. The Department Scheduler acknowledges the MPPS update with an N-SET response (status success).
7. The MPPS SCU AE closes the association with the Department Scheduler.

#### 4.2.4.3.1.2 Proposed Presentation Contexts

The MPPS SCU AE will propose Presentation Contexts as shown in the following Table:

**Table 4.2-29  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY ACQUIRE IMAGES**

| Presentation Context Table        |                 |                           |                     |      |           |
|-----------------------------------|-----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                   |                 | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                              | UID             | Name                      | UID                 |      |           |
| Modality Performed Procedure Step | 1.2.840.10008.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                                   | .1.2.3.3        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

#### 4.2.4.3.1.3 SOP Specific Conformance for MPPS SOP Class

The MPPS SCU AE provides standard conformance to the Modality Performed Procedure Step SOP Class as an SCU.

The behavior of the MPPS SCU AE when encountering status codes in an MPPS N-CREATE or N-SET response is summarized in the Table below.

**Table 4.2-30  
MPPS N-CREATE / N-SET RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning   | Status Code            | Behavior   |
|----------------|---|------------------------|--|
| Success        | Success   | 0000                   | The SCP has completed the operation successfully.  |
| Failure        | Processing Failure – Performed Procedure Step Object may no longer be updated | 0110                   | The Association is aborted using A-ABORT and the MPPS is marked as failed. The status meaning is logged and reported to the user. Additional information in the Response will be logged (i.e. Error Comment and Error ID).       |
| Warning        | Attribute Value Out of Range  | 0116                   | The MPPS operation is considered successful but the status meaning is logged. Additional information in the response identifying the attributes out of range is logged. (I.e., Elements in the Modification List/Attribute List) |
| *              | *   | Any other status code. | The Association is aborted using A-ABORT and the MPPS is marked as failed. The status meaning is logged and reported to the user.  |

The behavior of the MPPS SCU AE during communication failure is summarized in the Table below:

**Table 4.2-31  
MPPS COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior  |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and MPPS is marked as failed. The reason is logged and reported to the user. |
| Association aborted by the SCP or network layers | The MPPS is marked as failed. The reason is logged and reported to the user.  |

The Table below provides a description of the MPPS N-CREATE and N-SET request identifiers sent by the MPPS SCU AE. Empty cells in the N-CREATE and N-SET columns indicate that the attribute is not sent. An "x" indicates that an appropriate value will be sent. A "Zero length" attribute will be sent with zero length.

**Table 4.2-32**  
**MPPS N-CREATE / N-SET REQUEST IDENTIFIER**

| Attribute Name                               | Tag         | VR | N-CREATE   | N-SET  |
|--|-------------|----|--|--|
| Specific Character Set                       | (0008,0005) | CS | Created, if an extended or replacement character set is used. Refer to 6.SUPPORT OF CHARACTER SETS | Created, if an extended or replacement character set is used. Refer to 6.SUPPORT OF CHARACTER SETS |
| <b>Performed Procedure Step Relationship</b> |             |    |  |  |
| Scheduled Step Attributes Sequence           | (0040,0270) | SQ | Always Set   |  |
| > Study Instance UID                         | (0020,000D) | UI | From Modality Worklist or autosect   |  |
| >Issuer of Accession Number Sequence         | (0008,0051) | SQ | From Modality Worklist   |  |
| >>Local Namespace Entity ID                  | (0040,0031) | UT | From Modality Worklist   |  |
| >>Universal Entity ID                        | (0040,0032) | UT | From Modality Worklist   |  |
| >>Universal Entity ID Type                   | (0040,0033) | CS | From Modality Worklist   |  |
| > Referenced Study Sequence                  | (0008,1110) | SQ | From Modality Worklist   |  |
| >> Referenced SOP Class UID                  | (0008,1150) | UI | From Modality Worklist   |  |
| >> Referenced SOP Instance UID               | (0008,1155) | UI | From Modality Worklist   |  |
| > Accession Number                           | (0008,0050) | SH | From Modality Worklist or user input   |  |
| > Requested Procedure ID                     | (0040,1001) | SH | From Modality Worklist   |  |
| > Requested Procedure Description            | (0032,1060) | LO | From Modality Worklist   |  |
| > Scheduled Procedure Step ID                | (0040,0009) | SH | From Modality Worklist   |  |
| > Scheduled Procedure Step Description       | (0040,0007) | LO | From Modality Worklist   |  |
| > Scheduled Protocol Code Sequence           | (0040,0008) | SQ | From Modality Worklist   |  |
| >>Code Value                                 | (0008,0100) | SH | From Modality Worklist   |  |
| >>Coding Scheme Designator                   | (0008,0102) | SH | From Modality Worklist   |  |
| >>Coding Scheme Version                      | (0008,0103) | SH | From Modality Worklist   |  |
| >>Code Meaning                               | (0008,0104) | LO | From Modality Worklist   |  |

|   |             |    |                                      |                               |
|---|-------------|----|--------------------------------------|-------------------------------|
| > Placer Order Number / Imaging Service Request | (0040,2016) | LO | Zero length                          |                               |
| > Filler Order Number / Imaging Service Request | (0040,2017) | LO | Zero length                          |                               |
| Patient's Name                                  | (0010,0010) | PN | From Modality Worklist or user input |                               |
| Patient ID                                      | (0010,0020) | LO | From Modality Worklist or user input |                               |
| Issuer of Patient ID                            | (0010,0021) | LO | From Modality Worklist               |                               |
| Issuer of Patient ID Qualifiers Sequence        | (0010,0024) | SQ | From Modality Worklist               |                               |
| >Universal Entity ID                            | (0040,0032) | UT | From Modality Worklist               |                               |
| >Universal Entity ID Type                       | (0040,0033) | CS | From Modality Worklist               |                               |
| >Identifier Type Code                           | (0040,0035) | CS | From Modality Worklist               |                               |
| Patient's Birth Date                            | (0010,0030) | DA | From Modality Worklist or user input |                               |
| Patient's Sex                                   | (0010,0040) | CS | From Modality Worklist or user input |                               |
| Referenced Patient Sequence                     | (0008,1120) | SQ | Zero length                          |                               |
| > Referenced SOP Class UID                      | (0008,1150) | UI |                                      |                               |
| > Referenced SOP Instance UID                   | (0008,1155) | UI |                                      |                               |
| <b>Performed Procedure Step Information</b>     |             |    |                                      |                               |
| Performed Procedure Step ID                     | (0040,0253) | SH | Automatically created.               |                               |
| Performed Station AE Title                      | (0040,0241) | AE | MPPS AE Title                        |                               |
| Performed Station Name                          | (0040,0242) | SH | From configuration                   |                               |
| Performed Location                              | (0040,0243) | SH | From configuration                   |                               |
| Performed Procedure Step Start Date             | (0040,0244) | DA | Actual start date                    |                               |
| Performed Procedure Step Start Time             | (0040,0245) | TM | Actual start time                    |                               |
| Performed Procedure Step Status                 | (0040,0252) | CS | "IN PROGRESS"                        | "COMPLETED" or "DISCONTINUED" |
| Performed Procedure Step Description            | (0040,0254) | LO | Zero or more items                   | Zero or more items            |
| Performed Procedure Type Description            | (0040,0255) | LO | Zero length                          | Zero or more items            |
| Procedure Code Sequence                         | (0008,1032) | SQ | Zero or more items                   | Zero or more items            |
| >Code Value                                     | (0008,0100) | SH | x                                    | x                             |
| >Coding Scheme Designator                       | (0008,0102) | SH | x                                    | x                             |
| >Coding Scheme Version                          | (0008,0103) | SH | x                                    | x                             |
| >Code Meaning                                   | (0008,0104) | LO | x                                    | x                             |
| Performed Procedure Step End Date               | (0040,0250) | DA | Zero length                          | Actual end date               |

|   |             |    |                                      |                               |
|---|-------------|----|--------------------------------------|-------------------------------|
| Performed Procedure Step End Time                             | (0040,0251) | TM | Zero length                          | Actual end time               |
| Comments on the Performed Procedure Step                      | (0040,0280) | ST | Zero length                          | Zero length                   |
| Performed Procedure Step Discontinuation Reason Code Sequence | (0040,0281) | SQ | Zero length                          | “COMPLETED” or “DISCONTINUED” |
| >Code Value   | (0008,0100) | SH |                                      | x                             |
| >Coding Scheme Designator                                     | (0008,0102) | SH |                                      | x                             |
| >Coding Scheme Version  | (0008,0103) | SH |                                      | x                             |
| >Code Meaning   | (0008,0104) | LO |                                      | x                             |
| <b>Image Acquisition Results</b>                              |             |    |                                      |                               |
| Modality  | (0008,0060) | CS | “MR”                                 |                               |
| Study ID  | (0020,0010) | SH | Automatically created.               |                               |
| Performed Protocol Code Sequence                              | (0040,0260) | SQ | Zero or more items                   | Zero or more items            |
| >Code Value   | (0008,0100) | SH | From Modality Worklist               |                               |
| >Coding Scheme Designator                                     | (0008,0102) | SH | From configuration                   |                               |
| >Coding Scheme Version  | (0008,0103) | SH | From configuration                   |                               |
| >Code Meaning   | (0008,0104) | LO | From Modality Worklist               |                               |
| Performed Series Sequence                                     | (0040,0340) | SQ | One item                             | One or more items             |
| > Performing Physician’s Name                                 | (0008,1050) | PN | From Modality Worklist or user input | x                             |
| > Protocol Name   | (0018,1030) | LO | User input                           | x                             |
| > Operator’s Name   | (0008,1070) | PN | From Modality Worklist or user input | x                             |
| > Series Instance UID   | (0020,000E) | UI | Automatically created                | x                             |
| > Series Description  | (0008,103E) | LO | User input                           | x                             |
| > Retrieve AE Title   | (0008,0054) | AE | Zero length                          | x                             |
| > Referenced Image Sequence                                   | (0008,1140) | SQ | Zero length                          | One or more items             |
| >> Referenced SOP Class UID                                   | (0008,1150) | UI |                                      | x                             |
| >> Referenced SOP Instance UID                                | (0008,1155) | UI |                                      | x                             |
| > Referenced Non-Image Composite SOP Instance Sequence        | (0040,0220) | SQ | Zero length                          | x                             |
| >> Referenced SOP Class UID                                   | (0008,1150) | UI |                                      | x                             |
| >> Referenced SOP Instance UID                                | (0008,1155) | UI |                                      | x                             |
| <b>Billing and Material Code</b>                              |             |    |                                      |                               |
| Billing Procedure Step Sequence                               | (0040,0320) | SQ | Zero Length                          | Zero Length                   |
| Film Consumption Sequence                                     | (0040,0321) | SQ | Zero or more items                   | Zero or more items            |
| >Number of Films  | (2100,0170) | IS |                                      | x                             |
| >Medium Type  | (2000,0030) | CS |                                      | x                             |

|                                       |             |    |             |             |
|---------------------------------------|-------------|----|-------------|-------------|
| >Film Size ID                         | (2010,0050) | CS |             | x           |
| Billing Supplies and Devices Sequence | (0040,0324) | SQ | Zero Length | Zero Length |

#### 4.2.4.4 Association Acceptance Policy

The MPPS SCU AE does not accept Associations.

## 4.2.5 Storage SCU AE Specification

### 4.2.5.1 SOP Classes

The Storage SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-33**  
**SOP CLASSES FOR STORAGE SCU AE**

| SOP Class Name                                | SOP Class UID                | SCU | SCP |
|---|------------------------------|-----|-----|
| MR Image Storage                              | 1.2.840.10008.5.1.4.1.1.4    | Yes | No  |
| Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7    | Yes | No  |
| Enhanced MR Image Storage                     | 1.2.840.10008.5.1.4.1.1.4.1  | Yes | No  |
| MR Spectroscopy Storage                       | 1.2.840.10008.5.1.4.1.1.4.2  | Yes | No  |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 | Yes | No  |

### 4.2.5.2 Association Policies

#### 4.2.5.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-34**  
**DICOM APPLICATION CONTEXT FOR STORAGE SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.5.2.2 Number of Associations

The Storage SCU AE initiates one Association at a time.

**Table 4.2-35**  
**NUMBER OF ASSOCIATIONS INITIATED FOR STORAGE SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.5.2.3 Asynchronous Nature

The Storage SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-36**  
**ASYNCHRONOUS NATURE FOR STORAGE SCU AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.5.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-37**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR STORAGE SCU AE**


|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.5.3 Association Initiation Policy

#### 4.2.5.3.1 Activity – Send objects

##### 4.2.5.3.1.1 Description and Sequencing of Activities

The Storage SCU AE attempts to initiate a new Association in order to issue a Storage request (C-STORE). If the job contains multiple objects then multiple C-STORE requests will be issued over the same Association. If the object transfer fails, the Storage SCU AE will retry this send-job automatically.



**Figure 4.2-5**  
**SEQUENCING OF ACTIVITY – SEND OBJECTS**

A possible sequence of interactions between the Storage SCU AE and an Image Server (e.g. a storage or archive device supporting the Storage SOP Classes as an SCP) is illustrated in the Figure above:

1. The Storage SCU AE opens an Association with the Image Server
2. Acquired images, spectroscopy data or generated GSPS are transmitted to the Image Server using a Storage request (C-STORE) and the Image Server replies with a C-STORE response (status success).
3. The Storage SCU AE closes the Association with the Image Server.



#### 4.2.5.3.1.2 Proposed Presentation Contexts

The Storage SCU AE is capable of proposing the Presentation Contexts shown in the following table:

**Table 4.2-38  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND OBJECTS**

| Presentation Context Table                    |                                  |                           |                     |      |           |
|---|----------------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                               |                                  | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name  | UID                              | Name                      | UID                 |      |           |
| MR Image Storage                              | 1.2.840.10008.5<br>.1.4.1.1.4    | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   |                                  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| Secondary Capture Image Storage               | 1.2.840.10008.5<br>.1.4.1.1.7    | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   |                                  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| Enhanced MR Image Storage                     | 1.2.840.10008.5<br>.1.4.1.1.4.1  | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   |                                  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| MR Spectroscopy Storage                       | 1.2.840.10008.5<br>.1.4.1.1.4.2  | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   |                                  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5<br>.1.4.1.1.11.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   |                                  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

#### 4.2.5.3.1.3 SOP Specific Conformance for Storage SOP Classes

The Storage SCU AE provides standard conformance to the Storage Service Class as an SCU.

The behavior of Storage SCU AE when encountering status codes in a C-STORE response is summarized in the Table below:

**Table 4.2-39  
STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning                   | Status Code            | Behavior  |
|----------------|-----------------------------------|------------------------|---|
| Success        | Success                           | 0000                   | The SCP has successfully stored the SOP Instance. If all SOP Instances in a send job have status success then the job is marked as complete.  |
| Refused        | Out of Resources                  | A7xx                   | The Association is aborted using A-ABORT and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application. This is a transient failure. |
| Error          | Data Set does not match SOP Class | A9xx                   | The Association is aborted using A-ABORT and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.                              |
| Error          | Cannot Understand                 | Cxxx                   | The Association is aborted using A-ABORT and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.                              |
| *              | *                                 | Any other status code. | The Association is aborted using A-ABORT and the send job is marked as failed. The status code is logged and the job failure is reported to the user via the job control application.                                 |

The behavior of Storage SCU AE during communication failure is summarized in the Table below:

**Table 4.2-40**  
**STORAGE COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>  |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |
| Association aborted by the SCP or network layers | The send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application.  |

If the object transfer fails, the Storage SCU AE will retry this send-job automatically. The interval of time to resend a failed job and the number of times of retries are also configurable.

The contents of Storage SOP Instances created by the Storage SCU AE conform to the DICOM Image IOD definitions and are described in section 8.1.

#### **4.2.5.4 Association Acceptance Policy**

The Storage SCU AE does not accept Associations.

## 4.2.6 Storage Commitment SCU AE Specification

### 4.2.6.1 SOP Classes

The Storage Commitment SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-41**

**SOP CLASSES FOR THE STORAGE COMMITMENT SCU AE**

| SOP Class Name                | SOP Class UID        | SCU | SCP |
|-------------------------------|----------------------|-----|-----|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Yes | No  |

### 4.2.6.2 Association Policies

#### 4.2.6.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-42**

**DICOM APPLICATION CONTEXT FOR STORAGE COMMITMENT SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.6.2.2 Number of Associations

The Storage Commitment SCU AE initiates one Association at a time.

**Table 4.2-43**

**NUMBER OF ASSOCIATIONS INITIATED FOR STORAGE COMMITMENT SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

The Storage Commitment SCU AE accepts Associations to receive N-EVENT-REPORT notifications for the Storage Commitment Push Model SOP Class.

**Table 4.2-44**

**NUMBER OF ASSOCIATIONS ACCEPTED FOR STORAGE COMMITMENT SCU AE**

|   |    |
|---|----|
| Maximum number of simultaneous Associations | 10 |
|---|----|

#### 4.2.6.2.3 Asynchronous Nature

The Storage Commitment SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-45**

**ASYNCHRONOUS NATURE FOR STORAGE COMMITMENT SCU AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.6.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-46**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR STORAGE COMMITMENT SCU AE**


|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.6.3 Association Initiation Policy

#### 4.2.6.3.1 Activity – Commit Sent Objects

##### 4.2.6.3.1.1 Description and Sequencing of Activities

If the remote AE is configured as an archive device the Storage Commitment SCU AE will, after all objects have been sent, transmit a single Storage Commitment request (N-ACTION). Upon receiving the N-ACTION response the Storage Commitment SCU AE will release the Association. The notification of Storage commitment (N-EVENT-REPORT) will be received over a separate Association.



**Figure 4.2-6**  
**SEQUENCING OF ACTIVITY – COMMIT SENT OBJECTS**

A possible sequence of interactions between the Storage Commitment SCU AE and an Image Server (e.g. a storage or archive device supporting the Storage Commitment SOP Classes as an SCP) is illustrated in the Figure above:

1. The Storage Commitment SCU AE opens an Association with the Image Server.
2. A Storage Commitment request (N-ACTION) is transmitted to the Image Server to obtain Storage Commitment of previously transmitted objects. The Image Server replies with an N-ACTION response indicating the request has been received and is being processed.
3. The Storage Commitment AE closes the Association with the Image Server.

NOTE: The N-EVENT-REPORT will be sent over a separate Association initiated by the Image Server (see Section 4.2.6.4.1).

##### 4.2.6.3.1.2 Proposed Presentation Contexts

The Storage Commitment SCU AE is capable of proposing the Presentation Contexts shown in the following table:

**Table 4.2-47**  
**PROPOSED PRESENTATION CONTEXTS BY STORAGE COMMITMENT SCU AE**

| Presentation Context Table    |                      |                           |                     |      |           |
|-------------------------------|----------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax               |                      | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                          | UID                  | Name                      | UID                 |      |           |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                               |                      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

A Presentation Context for the Storage Commitment Push Model will only be proposed if the remote AE is configured as an archive device.

### 4.2.6.3.1.3 SOP Specific Conformance for Storage Commitment SOP Class

#### 4.2.6.3.1.3.1 Storage Commitment Operations (N-ACTION)

The Storage Commitment SCU AE provides standard conformance to the Storage Commitment Service Class as an SCU.

The Storage Commitment SCU AE will request storage commitment for instances of the Storage SOP Classes if the remote AE is configured as an archive device and a presentation context for the Storage Commitment Push Model has been accepted.

The behavior of Storage SCU Commitment AE when encountering status codes in an N-ACTION response is summarized in the Table below:

**Table 4.2-48  
STORAGE COMMITMENT N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b> | <b>Status Code</b>     | <b>Behavior</b>   |
|-----------------------|------------------------|------------------------|---|
| Success               | Success                | 0000                   | The request for storage commitment is considered successfully sent. A timer is started which will expire if no N-EVENT-REPORT for the Transaction UID is received within a configurable timeout period. |
| *                     | *                      | Any other status code. | The Association is aborted using A-ABORT and the request for storage commitment is marked as failed.  |

The behavior of Storage Commitment AE during communication failure is summarized in the Table below:

**Table 4.2-49  
STORAGE COMMITMENT COMMUNICATION FAILURE BEHAVIOR**


| <b>Exception</b>                                 | <b>Behavior</b>  |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |
| Association aborted by the SCP or network layers | The send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application.  |

## 4.2.6.4 Association Acceptance Policy

### 4.2.6.4.1 Activity – Receive Storage Commitment Response

#### 4.2.6.4.1.1 Description and Sequencing of Activities

The Storage Commitment SCU AE will accept Associations in order to receive responses to a Storage Commitment Request.



**Figure 4.2-7**  
**SEQUENCING OF ACTIVITY - RECEIVE STORAGE COMMITMENT RESPONSE**

A possible sequence of interactions between the Storage Commitment SCU AE and an Image Server (e.g. a storage or archive device supporting Storage Commitment SOP Classes as an SCP) is illustrated in the Figure above:

1. The Image Server opens an Association with the Storage Commitment SCU AE.
2. The Image Server sends an N-EVENT-REPORT request notifying the Storage Commitment SCU AE of the status of a previous Storage Commitment Request. The Storage SCU AE replies with an N-EVENT-REPORT response confirming receipt.
3. The Image Server closes the Association with the Storage Commitment SCU AE.

The Storage Commitment SCU AE may reject Association attempts as shown in the Table below. The Result, Source and Reason/Diag columns represent the values returned in the appropriate fields of an ASSOCIATE-RJ PDU (see PS 3.8, Section 9.3.4). The contents of the Source column is abbreviated to save space and the meaning of the abbreviations are:

- DICOM UL service-user
- DICOM UL service-provider (ASCE related function)
- DICOM UL service-provider (Presentation related function)

**Table 4.2-50**  
**ASSOCIATION REJECTION REASONS**

| <b>Result</b>             | <b>Source</b> | <b>Reason/Diag</b>                                | <b>Explanation</b>   |
|---------------------------|---------------|---|--|
| 2 –<br>rejected-transient | c             | 2 –<br>local-limit-exceeded                       | The (configurable) maximum number of simultaneous Associations has been reached. An Association request with the same parameters may succeed at a later time.  |
| 2 –<br>rejected-transient | c             | 1 –<br>temporary-congestion                       | No Associations can be accepted at this time due to the real-time requirements of higher priority activities (e.g. during image acquisition no Associations will be accepted) or because insufficient resources are available (e.g. memory, processes, threads). An Association request with the same parameters may succeed at a later time.                        |
| 1 –<br>rejected-permanent | a             | 2 –<br>application-context-<br>name-not-supported | The Association request contained an unsupported Application Context Name. An Association request with the same parameters will not succeed at a later time.   |
| 1 –<br>rejected-permanent | a             | 7 – called-AE-title-not-<br>recognized            | The Association request contained an unrecognized Called AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association initiator is incorrectly configured and attempts to address the Association acceptor using the wrong AE Title. |
| 1 –<br>rejected-permanent | a             | 3 – calling-AE-title-not-<br>recognized           | The Association request contained an unrecognized Calling AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association acceptor has not been configured to recognize the AE Title of the Association initiator.                      |
| 1 –<br>rejected-permanent | b             | 1 – no-reason-given                               | The Association request could not be parsed. An Association request with the same format will not succeed at a later time.   |

#### 4.2.6.4.1.2 Accepted Presentation Contexts

The Storage Commitment SCU AE will prefer to select the Explicit VR Little Endian Transfer Syntax if multiple transfer syntaxes are offered.

Any of the Presentation Contexts shown in the following table are acceptable to the Storage Commitment SCU AE.

**Table 4.2-51  
ACCEPTABLE PRESENTATION CONTEXTS BY  
ACTIVITY RECEIVE STORAGE COMMITMENT RESPONSE**

| Presentation Context Table    |                      |  |  |      |           |
|-------------------------------|----------------------|--|--|------|-----------|
| Abstract Syntax               |                      | Transfer Syntax  |  | Role | Ext. Neg. |
| Name                          | UID                  | Name   | UID                                      |      |           |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |



### 4.2.6.4.1.3 SOP Specific Conformance for Storage Commitment SOP Class

#### 4.2.6.4.1.3.1 Storage Commitment Notifications (N-EVENT-REPORT)

The Storage Commitment SCU AE provides standard conformance to the Storage Commitment Service Class as an SCU.

The behavior of Storage Commitment SCU AE when receiving Event Types within the N-EVENT-REPORT is summarized in the Table below.

**Table 4.2-52  
STORAGE COMMITMENT N-EVENT-REPORT BEHAVIOUR**

| Event Type Name                                      | Event Type ID | Behavior  |
|--|---------------|---|
| Storage Commitment Request Successful                | 1             | The Storage Commitment SCU AE permits the operator(s) to delete the Referenced SOP Instances under Referenced SOP Sequence (0018,1199). |
| Storage Commitment Request Complete – Failures Exist | 2             | The Storage Commitment SCU AE requests the Storage SCU AE to send the Referenced SOP Instances under Failed SOP Sequence (0018,1198).   |

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in the Table below.

**Table 4.2-53  
STORAGE COMMITMENT N-EVENT-REPORT RESPONSE STATUS REASONS**

| Service Status | Further Meaning        | Status Code | Reasons   |
|----------------|------------------------|-------------|---|
| Success        | Success                | 0000        | The storage commitment result has been successfully received.   |
| Failure        | Unrecognized Operation | 0211        | The Transaction UID in the N-EVENT-REPORT request is not recognized (was never issued within an N-ACTION request).  |
| Failure        | Resource Limitation    | 0213        | The Transaction UID in the N-EVENT-REPORT request has expired (no N-EVENT-REPORT was received within a configurable time limit).  |
| Failure        | No Such Event Type     | 0113        | An invalid Event Type ID was supplied in the N-EVENT-REPORT request.  |
| Failure        | Processing Failure     | 0110        | An internal error occurred during processing of the N-EVENT-REPORT. A short description of the error will be returned in Error Comment (0000,0902).   |
| Failure        | Invalid Argument Value | 0115        | One or more SOP Instance UIDs with the Referenced SOP Sequence (0008,1199) or Failed SOP Sequence (0008,1198) was not included in the Storage Commitment Request associated with this Transaction UID. The unrecognized SOP Instance UIDs will be returned within the Event Information of the N-EVENT-REPORT response. |

## 4.2.7 Q/R SCP AE Specification

### 4.2.7.1 SOP Classes

The Q/R SCP AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-54**  
**SOP CLASSES FOR Q/R SCP AE**

| SOP Class Name                            | SOP Class UID               | SCU | SCP |
|---|-----------------------------|-----|-----|
| Patient Root Q/R Information Model – Find | 1.2.840.10008.5.1.4.1.2.1.1 | No  | Yes |
| Patient Root Q/R Information Model – Move | 1.2.840.10008.5.1.4.1.2.1.2 | No  | Yes |
| Study Root Q/R Information Model – Find   | 1.2.840.10008.5.1.4.1.2.2.1 | No  | Yes |
| Study Root Q/R Information Model – Move   | 1.2.840.10008.5.1.4.1.2.2.2 | No  | Yes |

### 4.2.7.2 Association Policies

#### 4.2.7.2.1 General

The DICOM standard application context name for DICOM 3.0 is always accepted:

**Table 4.2-55**  
**DICOM APPLICATION CONTEXT FOR Q/R SCP AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.7.2.2 Number of Associations

The Q/R SCP AE accepts one Association at a time.

**Table 4.2-56**  
**NUMBER OF ASSOCIATIONS ACCEPTED FOR Q/R SCP AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.7.2.3 Asynchronous Nature

The Q/R SCP AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-57**  
**ASYNCHRONOUS NATURE FOR Q/R SCP AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.7.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-58**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR Q/R SCP AE**

|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.7.3 Association Initiation Policy

The Q/R SCP AE does not initiate Associations.

### 4.2.7.4 Association Acceptance Policy


#### 4.2.7.4.1 Activity – Handle Q/R Requests

##### 4.2.7.4.1.1 Description and Sequencing of Activities

The Q/R SCP AE accepts Associations only if they have valid Presentation Contexts. If none of the requested Presentation Contexts are accepted then the Association Request itself is rejected. It can be configured to only accept Associations with certain hosts (using TCP/IP address) and/or Application Entity Titles.

If the Q/R SCP AE receives a query (C-FIND) request then the response(s) will be sent over the same Association used to send the C-FIND-Request.

If the Q/R SCP AE receives a retrieval (C-MOVE) request then the responses will be sent over the same Association used to send the C-MOVE-Request. The Q/R SCP AE will notify the Storage SCU AE to send the requested SOP Instances to the C-MOVE Destination AE. The Storage SCU AE notifies the Q/R SCP AE of the success or failure of each attempt to send a Composite SOP Instance to the peer C-MOVE Destination AE. The Q/R SCP AE then sends a C-MOVE Response indicating this status after each attempt. Once the Storage SCU AE has finished attempting to transfer all the requested SOP Instances, the Q/R SCP AE sends a final C-MOVE Response indicating the overall status of the attempted retrieval.



**Figure 4.2-8**  
**SEQUENCING OF ACTIVITY – HANDLE Q/R REQUESTS**

The following sequencing constraints illustrated in the Figure above:

1. The Q/R SCU AE opens an Association with the Q/R SCP AE.
2. The Q/R SCU AE sends a C-FIND-RQ Message
3. The Q/R SCP AE returns a C-FIND-RSP Message to the Q/R SCU AE with matching information. A C-FIND-RSP is sent for each entity matching the identifier specified in the C-FIND-RQ. A final C-FIND-RSP is sent indicating that the matching is complete.
4. The Q/R SCU AE closes the Association.
5. The Q/R SCU AE opens an Association with the Q/R SCP AE.
6. The Q/R SCU AE sends a C-MOVE-RQ Message. The Q/R SCP AE notifies the Storage SCU AE to send the Composite SOP Instances to the peer C-MOVE Destination AE as indicated in the C-MOVE-RQ.
7. The Storage SCU AE opens an Association with the C-MOVE Destination AE.
8. The Storage SCU AE sends images to the C-MOVE Destination AE. The Storage SCU AE indicates to the Q/R SCP AE whether the transfer succeeded or failed.
9. The Q/R SCP AE then returns a C-MOVE-RSP indicating this success or failure.
10. The Storage SCU AE closes the Association.
11. The Q/R SCP AE sends a final C-MOVE-RSP indicating the overall success or failure of the retrieval.
12. The Q/R SCU AE closes the Association.

The Q/R SCP AE may reject Association attempts as shown in the table below. The Result, Source and Reason/Diag columns represent the values returned in the corresponding fields of an ASSOCIATE-RJ PDU (see PS 3.8, Section 9.3.4). The following abbreviations are used in the Source column:

- DICOM UL service-user
- DICOM UL service-provider (ASCE related function)
- DICOM UL service-provider (Presentation related function)

**Table 4.2-59**  
**ASSOCIATION REJECTION REASONS**

| <b>Result</b>             | <b>Source</b> | <b>Reason/Diag</b>                         | <b>Explanation</b>   |
|---------------------------|---------------|--|--|
| 2 –<br>rejected-transient | c             | 2 – local-limit-exceeded                   | The (configurable) maximum number of simultaneous Associations has been reached. An Association request with the same parameters may succeed at a later time.  |
| 2 –<br>rejected-transient | c             | 1 – temporary-congestion                   | No Associations can be accepted at this time due to the real-time requirements of higher priority activities (e.g. during image acquisition no Associations will be accepted) or because insufficient resources are available (e.g. memory, processes, threads). An Association request with the same parameters may succeed at a later time.                        |
| 1 –<br>rejected-permanent | a             | 2 – application-context-name-not-supported | The Association request contained an unsupported Application Context Name. An association request with the same parameters will not succeed at a later time.   |
| 1 –<br>rejected-permanent | a             | 7 – called-AE-title-not-recognized         | The Association request contained an unrecognized Called AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association initiator is incorrectly configured and attempts to address the Association acceptor using the wrong AE Title. |
| 1 –<br>rejected-permanent | a             | 3 – calling-AE-title-not-recognized        | The Association request contained an unrecognized Calling AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association acceptor has not been configured to recognize the AE Title of the Association initiator.                      |
| 1 –<br>rejected-permanent | b             | 1 – no-reason-given                        | The Association request could not be parsed. An Association request with the same format will not succeed at a later time.   |

#### 4.2.7.4.1.2 Accepted Presentation Contexts

The Q/R SCP AE will prefer to select the Explicit VR Little Endian Transfer Syntax if multiple transfer syntaxes are offered.

Any of the Presentation Contexts shown in the following table are acceptable to the Q/R SCP AE.

**Table 4.2-60**  
**ACCEPTED PRESENTATION CONTEXTS BY Q/R SCP AE**

| <b>Presentation Context Table</b>               |                                 |                           |                     |             |                  |
|---|---------------------------------|---------------------------|---------------------|-------------|------------------|
| <b>Abstract Syntax</b>                          |                                 | <b>Transfer Syntax</b>    |                     | <b>Role</b> | <b>Ext. Neg.</b> |
| <b>Name</b>                                     | <b>UID</b>                      | <b>Name</b>               | <b>UID</b>          |             |                  |
| Patient Root<br>Q/R Information<br>Model - FIND | 1.2.840.10008.5<br>.1.4.1.2.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU         | None             |
|   |                                 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |             |                  |
| Patient Root<br>Q/R Information<br>Model - MOVE | 1.2.840.10008.5<br>.1.4.1.2.1.2 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU         | None             |
|   |                                 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |             |                  |
| Study Root<br>Q/R Information<br>Model - FIND   | 1.2.840.10008.5<br>.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU         | None             |
|   |                                 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |             |                  |
| Study Root<br>Q/R Information<br>Model - MOVE   | 1.2.840.10008.5<br>.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU         | None             |
|   |                                 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |             |                  |

#### 4.2.7.4.1.3 SOP Specific Conformance for Q/R Find SOP Classes

The Q/R SCP AE provides standard conformance to the Query/Retrieve Find SOP Class as an SCP. It supports hierarchical queries and not relational queries.

**Table 4.2-61**  
**PATIENT ROOT C-FIND SCP SUPPORTED ELEMENTS**

| Level Name<br>Attribute Name      | Tag         | VR | Types of<br>Matching |
|-----------------------------------|-------------|----|----------------------|
| <b>Patient Level</b>              |             |    |                      |
| Patient's Name                    | (0010,0010) | PN | S,*,U                |
| Patient ID                        | (0010,0020) | LO | S,*,U                |
| <b>Study Level</b>                |             |    |                      |
| Study Date                        | (0008,0020) | DA | S,R,U                |
| Study Time                        | (0008,0030) | TM | R,U                  |
| Accession Number                  | (0008,0050) | SH | S,*,U                |
| Study Instance UID                | (0020,000D) | UI | S,*,U                |
| Study ID                          | (0020,0010) | SH | S,*,U                |
| <b>Series Level</b>               |             |    |                      |
| Modality                          | (0008,0060) | CS | S,*,U                |
| Series Number                     | (0020,0011) | IS | S,*,U                |
| Series Instance UID               | (0020,000E) | UI | S,*,U                |
| Number of Series RelatedInstances | (0020,1209) | IS | U                    |
| <b>Instance Level</b>             |             |    |                      |
| SOP Instance UID                  | (0008,0018) | UI | S,*,U                |
| Instance Number                   | (0020,0013) | IS | S,*,U                |

**Table 4.2-62  
STUDY ROOT C-FIND SCP SUPPORTED ELEMENTS**

| <b>Level Name<br/>Attribute Name</b> | <b>Tag</b>  | <b>VR</b> | <b>Types of<br/>Matching</b> |
|--------------------------------------|-------------|-----------|------------------------------|
| <b>Study Level</b>                   |             |           |                              |
| Study Date                           | (0008,0020) | DA        | S,R,U                        |
| Study Time                           | (0008,0030) | TM        | R,U                          |
| Accession Number                     | (0008,0050) | SH        | S,*,U                        |
| Patient's Name                       | (0010,0010) | PN        | S,*,U                        |
| Patient ID                           | (0010,0020) | LO        | S,*,U                        |
| Study Instance UID                   | (0020,000D) | UI        | S,*,U                        |
| Study ID                             | (0020,0010) | SH        | S,*,U                        |
| <b>Series Level</b>                  |             |           |                              |
| Modality                             | (0008,0060) | CS        | S,*,U                        |
| Series Number                        | (0020,0011) | IS        | S,*,U                        |
| Series Instance UID                  | (0020,000E) | UI        | S,*,U                        |
| <b>Instance Level</b>                |             |           |                              |
| SOP Instance UID                     | (0008,0018) | UI        | S,*,U                        |
| Instance Number                      | (0020,0013) | IS        | S,*,U                        |

The tables should be read as follows:

Attribute Name: Attributes supported for returned C-FIND Responses.

Tag: Appropriate DICOM tag for this attribute.

VR: Appropriate DICOM VR for this attribute.

Types of Matching: The types of Matching supported by the C-FIND SCP.

A "S" indicates the identifier attribute can specify Single Value Matching, a "R" will indicate Range Matching, an "\*" will denote wildcard matching, and a "U" will indicate universal matching.



The Q/R SCP AE returns C-FIND response status as specified below.

**Table 4.2-63**  
**Q/R SCP AE C-FIND RESPONSE STATUS RETURN REASONS**

| <b>Service Status</b> | <b>Further Meaning</b>   | <b>Status Code</b> | <b>Reasons</b>   |
|-----------------------|--|--------------------|--|
| Success               | Success  | 0000               | Matching is complete. No final identifier is supplied.   |
| Refused               | Out of Resources   | A700               | System reached the limit in disk space or memory usage.<br>Error message is output to the Service Log.   |
| Failed                | Identifier does not match SOP Class                                      | A900               | The C-FIND query identifier contains invalid Elements or values, or is missing mandatory Elements or values for the specified SOP Class.<br>Error message is output to the Service Log.  |
|                       | Unable to process  | C000<br>C001       | The C-FIND query identifier is valid for the specified SOP Class but cannot be used to query the database. For example, this can occur if received data contains unsupported character sets. (See section 6 'SUPPORT OF CHARACTER SETS'.)                              |
| Cancel                | Matching terminated due to Cancel Request                                | FE00               | The C-FIND SCU sent a Cancel Request. This has been acknowledged and the search for matches has been halted.   |
| Pending               | Matches are continuing and current match is supplied.                    | FF00               | Indicates that the search for further matches is continuing. This is returned when each successful match is returned and when further matches are forthcoming. This status code is returned if all Optional keys in the query identifier are actually supported.       |
|                       | Matches are continuing but one or more Optional Keys were not supported. | FF01               | Indicates that the search for further matches is continuing. This is returned when each successful match is returned and when further matches are forthcoming. This status code is returned if there are Optional keys in the query identifier that are not supported. |

#### 4.2.7.4.1.4 SOP Specific Conformance for Q/R Move SOP Classes

The Q/R SCP AE provides standard conformance to the Query/Retrieve Move SOP Classes as an SCP.

The Q/R SCP AE will convey to the Storage SCU AE that an Association with a DICOM Application Entity named by the external C-MOVE SCU (through a MOVE Destination AE Title) should be established. It will also convey to the Storage SCU AE to perform C-STORE operations on specific images requested by the external C-MOVE SCU.

The Q/R SCP AE returns C-MOVE response status as specified below.

**Table 4.2-64**  
**Q/R SCP AE C-MOVE RESPONSE STATUS RETURN REASONS**

| Service Status | Further Meaning  | Status Code | Reasons  |
|----------------|--|-------------|--|
| Success        | Sub-operations complete – No Failures                    | 0000        | All the Composite SOP Instances have been successfully sent to the C-MOVE Destination AE.  |
| Refused        | Out of Resources – Unable to calculate number of matches | A701        | Number of matches cannot be determined due to system failure. Returned if the server's database is not functioning so the search for matches to the C-MOVE Request cannot be found.<br>Error message is output to the Service Log. |
|                | Out of Resources – Unable to perform sub-operations      | A702        | C-STORE sub-operations cannot be performed due to failure to access Composite SOP Instances in archive, or failure of a C-STORE Request.<br>Error message is output to the Service Log.  |
|                | Move destination unknown                                 | A801        | The Destination Application Entity named in the C-MOVE Request is unknown to Q/R SCP AE.<br>Error message is output to the Service Log.  |
| Failed         | Identifier does not match SOP Class                      | A900        | The C-MOVE identifier contains invalid Elements or values, or is missing mandatory Elements or values for the specified SOP Class or retrieval level.<br>Error message is output to the Service Log.                               |
| Cancel         | Matching terminated due to Cancel Request                | FE00        | The C-MOVE SCU sent a Cancel Request. This has been acknowledged and the export of Composite SOP Instances to the C-MOVE Destination AE has been halted.   |

## 4.2.8 Q/R SCU AE Specification

### 4.2.8.1 SOP Classes

The Q/R SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-65**  
**SOP CLASSES FOR Q/R SCU AE**

| SOP Class Name                          | SOP Class UID               | SCU | SCP |
|---|-----------------------------|-----|-----|
| Study Root Q/R Information Model – Find | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | No  |
| Study Root Q/R Information Model – Move | 1.2.840.10008.5.1.4.1.2.2.2 | Yes | No  |

### 4.2.8.2 Association Policies

#### 4.2.8.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-66**  
**DICOM APPLICATION CONTEXT FOR Q/R SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.8.2.2 Number of Associations

The Q/R SCU AE initiates one Association at a time.

**Table 4.2-67**  
**NUMBER OF ASSOCIATIONS INITIATED FOR Q/R SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.8.2.3 Asynchronous Nature

The Q/R SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-68**  
**ASYNCHRONOUS NATURE FOR Q/R SCU AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.8.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-69**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR Q/R SCU AE**


|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.8.3 Association Initiation Policy

#### 4.2.8.3.1 Activity – Query and Retrieve Images

##### 4.2.8.3.1.1 Description and Sequencing of Activities

The Q/R SCU AE is activated when the user selects a remote node to query and enters some key information, Patient's Name, Patient ID and/or Study Date. The user can select studies, series and images to be retrieved. The images will be received at the Storage SCP AE.



**Figure 4.2-9**  
**SEQUENCING OF ACTIVITY – QUERY AND RETRIEVE IMAGES**

The following sequencing constraints illustrated in the Figure above:

1. The Q/R SCU AE opens an Association with the Q/R SCP AE.
2. The Q/R SCU AE sends a C-FIND-RQ Message
3. The Q/R SCP AE returns a C-FIND-RSP Message to the Q/R SCU AE with matching information. A C-FIND-RSP is sent for each entity matching the identifier specified in the C-FIND-RQ. A final C-FIND-RSP is sent indicating that the matching is complete.
4. The Q/R SCU AE closes the Association.
5. The Q/R SCU AE opens an Association with the Q/R SCP AE.
6. The Q/R SCU AE sends a C-MOVE-RQ Message. The Q/R SCP AE notifies the Storage SCU AE to send the Composite SOP Instances to the peer C-MOVE Destination AE as indicated in the C-MOVE-RQ.
7. The Storage SCU AE opens an Association with the C-MOVE Destination AE.
8. The Storage SCU AE sends images to the C-MOVE Destination AE. The Storage SCU AE indicates to the Q/R SCP AE whether the transfer succeeded or failed.
9. The Q/R SCP AE then returns a C-MOVE-RSP indicating this success or failure.
10. The Storage SCU AE closes the Association.
11. The Q/R SCU AE closes the Association.

#### 4.2.8.3.1.2 Proposed Presentation Contexts

The Q/R SCU AE will propose Presentation Contexts as shown in the following Table:

**Table 4.2-70**  
**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY**  
**QUERY AND RETRIEVE IMAGES**

| Presentation Context Table              |                 |                           |                     |      |           |
|---|-----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                         |                 | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                                    | UID             | Name                      | UID                 |      |           |
| Study Root Q/R Information Model – Find | 1.2.840.10008.5 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   | .1.4.1.2.2.1    | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| Study Root Q/R Information Model – Move | 1.2.840.10008.5 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   | .1.4.1.2.2.2    | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

#### 4.2.8.3.1.3 SOP Specific Conformance for Q/R Find SOP Classes

The Q/R SCU AE provides standard conformance to the Query/Retrieve Find SOP Classes as an SCU.

The behavior of the Q/R SCU AE when encountering status codes in a Q/R C-FIND response is summarized in the Table below:

**Table 4.2-71  
Q/R SCU AE C-FIND RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>   | <b>Status Code</b>     | <b>Behavior</b>   |
|-----------------------|--|------------------------|---|
| Success               | Matching is complete   | 0000                   | The SCP has completed the matches. Study, Series or Image information items are available for display or further processing.  |
| Refused               | Out of Resources   | A700                   | The Association is aborted using A-ABORT and the Study, Series, or Image information query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Failed                | Identifier does not match SOP Class  | A900                   | The Association is aborted using A-ABORT and the Study, Series, or Image information query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Failed                | Unable to Process  | Cxxx                   | The Association is aborted using A-ABORT and the Study, Series, or Image information query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Cancel                | Matching terminated due to Cancel request  | FE00                   | The Association is aborted using A-ABORT and the Study, Series, or Image information query is marked as failed. The status meaning is logged and reported to the user if an interactive query.  |
| Pending               | Matches are continuing   | FF00                   | The Study, Series, or Image information items contained in the Identifier is collected for later display or further processing.   |
| Pending               | Matches are continuing – Warning that one or more Optional Keys were not supported | FF01                   | The Study, Series, or Image information items contained in the Identifier is collected for later display or further processing.   |
| *                     | *  | Any other status code. | The Association is aborted using A-ABORT and the Study, Series, or Image information is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged.       |

The behavior of the Q/R SCU AE during communication failure is summarized in the Table below.

**Table 4.2-72**  
**Q/R FIND COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>   |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and the study, series or image query is marked as failed. The reason is logged and reported to the user if an interactive query. |
| Association aborted by the SCP or network layers | The study, series or image query is marked as failed. The reason is logged and reported to the user if an interactive query.  |

All queries are initiated at the highest level of the information model (the STUDY level), and then for each response received, recursively repeated at the next lower levels (the SERIES and then IMAGE levels), in order to completely elucidate the “tree” of instances available on the remote AE.

The Table below provides a description of the Q/R SCU AE C-FIND Request Identifier.

**Table 4.2-73**  
**STUDY ROOT REQUEST IDENTIFIER FOR C-FIND**

| <b>Name</b>                        | <b>Tag</b>  | <b>Types of Matching</b> |
|------------------------------------|-------------|--------------------------|
| Specific Character Set             | (0008,0005) | N/A                      |
| <b>Study Level</b>                 |             |                          |
| Study Date                         | (0008,0020) | S,*,U,R                  |
| Study Time                         | (0008,0030) | U                        |
| Accession Number                   | (0008,0050) | U                        |
| Modality in Study                  | (0008,0061) | U                        |
| Study Description                  | (0008,1030) | U                        |
| Referring Physician's Name         | (0008,0090) | U                        |
| Study Description                  | (0008,1030) | U                        |
| Procedure Code Sequence            | (0008,1032) | U                        |
| >Code Value                        | (0008,0100) | U                        |
| >Coding Scheme Designator          | (0008,0102) | U                        |
| >Coding Scheme Version             | (0008,0103) | U                        |
| >Code Meaning                      | (0008,0104) | U                        |
| Name of Physician(s) Reading Study | (0008,1060) | U                        |
| Referenced Study Sequence          | (0008,1110) | U                        |
| >Referenced SOP Class UID          | (0008,1150) | U                        |
| > Referenced SOP Instance UID      | (0008,1155) | U                        |
| Referenced Patient Sequence        | (0008,1120) | U                        |
| >Referenced SOP Class UID          | (0008,1150) | U                        |
| > Referenced SOP Instance UID      | (0008,1155) | U                        |
| Patient's Name                     | (0010,0010) | U                        |
| Patient's ID                       | (0010,0020) | S,*,U                    |
| Patient's Birth Date               | (0010,0030) | U                        |
| Patient's Birth Time               | (0010,0032) | U                        |
| Patient's Sex                      | (0010,0040) | U                        |
| Other Patient IDs                  | (0010,1000) | U                        |
| Other Patient Names                | (0010,1001) | U                        |
| Patient's Age                      | (0010,1010) | U                        |
| Patient's Size                     | (0010,1020) | U                        |
| Patient's Weight                   | (0010,1030) | U                        |
| Ethnic Group                       | (0010,2160) | U                        |
| Occupation                         | (0010,2180) | U                        |
| Additional Patient History         | (0010,21B0) | U                        |
| Patient Comments                   | (0010,4000) | U                        |
| Study Instance UID                 | (0020,000D) | UNIQUE                   |
| Study ID                           | (0020,0010) | U                        |
| Other Study Numbers (RET)          | (0020,1070) | U                        |
| Number of Study Related Series     | (0020,1206) | U                        |
| Number of Study Related Instances  | (0020,1208) | U                        |
| Interpretation Author (RET)        | (4008,010C) | U                        |



| <b>Series Level</b>               |              |        |
|-----------------------------------|--------------|--------|
| Series Date                       | (0008,0021)  | U      |
| Series Time                       | (0008,0031)  | U      |
| Modality                          | (0008,0060)  | S,*,U  |
| Station Name                      | (0008,1010)  | S,*,U  |
| Institutional Department Name     | (0008,1040)  | U      |
| Performing Physician's Name       | (0008,1050)  | U      |
| Operator's Name                   | (0008,1070)  | U      |
| Body Part Examined                | (0018,0015)  | U      |
| Series Instance UID               | (0020,000E)  | UNIQUE |
| Series Number                     | (0020,0011)  | S,*,U  |
| Number of Series RelatedInstances | (0020,1209)  | U      |
| <b>Instance Level</b>             |              |        |
| SOP Class UID                     | (0008,0016)  | U      |
| SOP Instance UID                  | (0008,0018)  | UNIQUE |
| Content Date                      | (0008,0023)  | U      |
| Content Time                      | (0008,0033)  | U      |
| Referenced SOP Class UID          | >(0008,1150) | U      |
| Referenced SOP Instance UID       | >(0008,1155) | U      |
| Instance Number                   | (0020,0013)  | U      |
| Overlay Number (RET)              | (0020,0022)  | U      |
| Curve Number (RET)                | (0020,0024)  | U      |
| LUT Number (RET)                  | (0020,0026)  | U      |

#### Types of Matching:

The types of Matching supported by the Q/R SCU AE. An "S" indicates the identifier attribute uses Single Value Matching, an "R" indicates Range Matching, an "\*" indicates wildcard matching, and a 'U' indicates Universal Matching. "UNIQUE" indicates that this is the Unique Key for that query level, in which case Universal Matching or Single Value Matching is used depending on the query level.

The Q/R SCU AE supports Study and Series level as a Query Level.

#### 4.2.8.3.1.4 SOP Specific Conformance for Q/R Move SOP Classes

The Q/R SCU AE provides standard conformance to the Query/Retrieve Move SOP Classes as an SCU.

The behavior of the Q/R SCU AE when encountering status codes in a Q/R C-MOVE response is summarized in the Table below:

**Table 4.2-74  
Q/R SCU AE C-MOVE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning  | Status Code | Behavior  |
|----------------|--|-------------|---|
| Success        | Sub-operations complete – No Failures                    | 0000        | The Storage SCP AE has successfully received the SOP Instance. If all SOP Instances in a move job have status success then the job is marked as complete.   |
| Refused        | Out of Resources – Unable to calculate number of matches | A701        | The Association is aborted using A-ABORT and the move job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application. This is a transient failure. |
|                | Out of Resources – Unable to perform sub-operations      | A702        | The Association is aborted using A-ABORT and the move job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.                              |
|                | Move destination unknown                                 | A801        | The Association is aborted using A-ABORT and the move job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.                              |
| Failed         | Identifier does not match SOP Class                      | A900        | The Association is aborted using A-ABORT and the move job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.                              |
| Warning        | Sub-operations complete but one or more failures.        | B000        | The Association is aborted using A-ABORT and the move job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.                              |

The behavior of the Q/R SCU AE during communication failure is summarized in the Table below.

**Table 4.2-75  
Q/R MOVE COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior  |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and the retrieve is marked as failed. The reason is logged and reported to the user if an interactive query. |
| Association aborted by the SCP or network layers | The retrieve is marked as failed. The reason is logged and reported to the user if an interactive query.  |

The system requests Image Level Move only.

#### 4.2.8.4 Association Acceptance Policy

The Q/R SCU AE does not accept Associations.

## 4.2.9 Storage SCP AE Specification

### 4.2.9.1 SOP Classes

The Storage SCP AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-76**  
**SOP CLASSES FOR STORAGE SCP AE**

| SOP Class Name                  | SOP Class UID             | SCU | SCP |
|---------------------------------|---------------------------|-----|-----|
| MR Image Storage                | 1.2.840.10008.5.1.4.1.1.4 | No  | Yes |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | No  | Yes |

### 4.2.9.2 Association Policies

#### 4.2.9.2.1 General

The DICOM standard application context name for DICOM 3.0 is always accepted:

**Table 4.2-77**  
**DICOM APPLICATION CONTEXT FOR STORAGE SCP AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.9.2.2 Number of Associations

The Storage SCP AE accepts Associations to receive C-STORE requests.

**Table 4.2-78**  
**NUMBER OF ASSOCIATIONS ACCEPTED FOR STORAGE SCP AE**

|   |    |
|---|----|
| Maximum number of simultaneous Associations | 10 |
|---|----|

#### 4.2.9.2.3 Asynchronous Nature

The Storage SCP AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-79**  
**ASYNCHRONOUS NATURE FOR STORAGE SCP AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.9.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-80**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR STORAGE SCP AE**

|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.9.3 Association Initiation Policy

The Storage SCP AE does not initiate Associations.


### 4.2.9.4 Association Acceptance Policy

#### 4.2.9.4.1 Activity – Store Images to the local file system

##### 4.2.9.4.1.1 Description and Sequencing of Activities

The Storage SCP AE accepts Associations only if they have valid Presentation Contexts. If none of the requested Presentation Contexts are accepted then the Association Request itself is rejected. It can be

configured to only accept Associations with certain hosts (using TCP/IP address) and/or Application Entity Titles.



**Figure 4.2-10**  
**SEQUENCING OF ACTIVITY – STORE IMAGES TO THE LOCAL FILE SYSTEM**

A possible sequence of interactions between the Storage SCP AE and a Storage SCU AE is illustrated in the Figure above:

1. The Storage SCU AE opens an Association with the Storage SCP AE.
2. The Storage SCU AE sends images to the Storage SCP AE using a Storage request (C-STORE) and the Storage SCP AE replies with a C-STORE response (status success).
3. The Storage SCU AE closes the Association with the Storage SCP AE.

The Storage SCP AE may reject association attempts as shown in the Table below. The Result, Source and Reason/Diag columns represent the values returned in the appropriate fields of an ASSOCIATE-RJ PDU. The contents of the Source column is abbreviated to save space and the meaning of the abbreviations are:

- DICOM UL service-user
- DICOM UL service-provider (ASCE related function)
- DICOM UL service-provider (Presentation related function)

**Table 4.2-81  
ASSOCIATION REJECTION REASONS**

| <b>Result</b>             | <b>Source</b> | <b>Reason/Diag</b>                         | <b>Explanation</b>   |
|---------------------------|---------------|--|--|
| 2 –<br>rejected-transient | c             | 2 – local-limit-exceeded                   | The (configurable) maximum number of simultaneous Associations has been reached. An Association request with the same parameters may succeed at a later time.  |
| 2 –<br>rejected-transient | c             | 1 – temporary-congestion                   | No Associations can be accepted at this time due to the real-time requirements of higher priority activities (e.g. during image acquisition no Associations will be accepted) or because insufficient resources are available (e.g. memory, processes, threads). An Association request with the same parameters may succeed at a later time.                        |
| 1 –<br>rejected-permanent | a             | 2 – application-context-name-not-supported | The Association request contained an unsupported Application Context Name. An association request with the same parameters will not succeed at a later time.   |
| 1 –<br>rejected-permanent | a             | 7 – called-AE-title-not-recognized         | The Association request contained an unrecognized Called AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association initiator is incorrectly configured and attempts to address the Association acceptor using the wrong AE Title. |
| 1 –<br>rejected-permanent | a             | 3 –<br>calling-AE-title-not-recognized     | The Association request contained an unrecognized Calling AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association acceptor has not been configured to recognize the AE Title of the Association initiator.                      |
| 1 –<br>rejected-permanent | b             | 1 – no-reason-given                        | The Association request could not be parsed. An Association request with the same format will not succeed at a later time.   |

#### 4.2.9.4.1.2 Accepted Presentation Contexts

The Storage SCP AE will prefer to select the Explicit VR Little Endian Transfer Syntax if multiple transfer syntaxes are offered.

Any of the Presentation Contexts shown in the following table are acceptable to the Storage SCP AE.

**Table 4.2-82**  
**ACCEPTED PRESENTATION CONTEXTS BY STORAGE SCP AE**

| Presentation Context Table      |                 |                           |                     |      |           |
|---------------------------------|-----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                 |                 | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                            | UID             | Name                      | UID                 |      |           |
| MR Image Storage                | 1.2.840.10008.5 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCP  | None      |
|                                 | .1.4.1.1.4      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
| Secondary Capture Image Storage | 1.2.840.10008.5 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCP  | None      |
|                                 | .1.4.1.1.7      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

#### 4.2.9.4.1.3 SOP Specific Conformance for Storage SOP Classes

The associated Activity with the Storage service is the storage of medical image data received over the network on a designated hard disk. The Storage SCP AE will return a failure status if it is unable to store the images on to the hard disk.

The Storage SCP AE is Level 0 conformant as a Storage SCP.

**Table 4.2-83**  
**STORAGE SCP AE C-STORE RESPONSE STATUS RETURN REASONS**

| Service Status | Further Meaning                   | Status Code | Reasons  |
|----------------|-----------------------------------|-------------|--|
| Success        | Success                           | 0000        | The Composite SOP Instance was successfully received, verified, and stored in the system database. |
| Refused        | Out of Resources                  | A7xx        | Indicates that there were not enough local resources.  |
| Error          | Data Set does not match SOP Class | A9xx        | Indicates that the Data Set does not encode a valid instance of the SOP Class specified.           |
|                | Processing Failed                 | B006        | Indicates that some elements discarded.  |
|                | Processing failed                 | B007        | Indicates that the Data Set does not match a supported SOP Class.                                  |
|                | Cannot understand                 | C0xx        | Indicates that the Storage SCP AE cannot parse the Data Set into Elements.                         |

## 4.2.10 Storage Commitment SCP AE Specification

### 4.2.10.1 SOP Classes

The Storage Commitment SCP AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-84**

**SOP CLASSES FOR THE STORAGE COMMITMENT SCP AE**

| SOP Class Name                | SOP Class UID        | SCU | SCP |
|-------------------------------|----------------------|-----|-----|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | No  | Yes |

### 4.2.10.2 Association Policies

#### 4.2.10.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-85**

**DICOM APPLICATION CONTEXT FOR STORAGE COMMITMENT SCP AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.10.2.2 Number of Associations

The Storage Commitment SCP AE accepts one Association at a time

**Table 4.2-86**

**NUMBER OF ASSOCIATIONS INITIATED FOR STORAGE COMMITMENT SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

The Storage Commitment SCP AE initiates Associations to send N-EVENT-REPORT notifications for the Storage Commitment Push Model SOP Class.

**Table 4.2-87**

**NUMBER OF ASSOCIATIONS ACCEPTED FOR STORAGE COMMITMENT SCU AE**

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.10.2.3 Asynchronous Nature

The Storage Commitment SCP AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-88**

**ASYNCHRONOUS NATURE FOR STORAGE COMMITMENT SCP AE**

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.10.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-89**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR STORAGE COMMITMENT SCP AE**


|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

### 4.2.10.3 Association Initiation Policy

#### 4.2.10.3.1 Activity – Commit Sent Objects

##### 4.2.10.3.1.1 Description and Sequencing of Activities

The Storage Commitment SCP AE will initiate Associations in order to send responses to a Storage Commitment Response(N-EVENT-REPORT).



**Figure 4.2-11**  
**SEQUENCING OF ACTIVITY - SEND STORAGE COMMITMENT RESPONSE**

A possible sequence of interactions between the Storage Commitment SCP AE and Storage Commitment SCU AE is illustrated in the Figure above:

1. The Storage Commitment SCP AE opens an Association with the Storage Commitment SCU AE.
2. The Storage Commitment SCP AE sends an N-EVENT-REPORT request notifying the Storage Commitment SCU AE of the status of a previous Storage Commitment Request. The Storage SCU AE replies with an N-EVENT-REPORT response confirming receipt.
3. The Storage Commitment SCP AE closes the Association with the Storage Commitment SCU AE.

##### 4.2.10.3.1.2 Proposed Presentation Contexts

The Storage Commitment SCP AE is capable of proposing the Presentation Contexts shown in the following table:

**Table 4.2-90**  
**PROPOSED PRESENTATION CONTEXTS BY STORAGE COMMITMENT SCP AE**

| Presentation Context Table    |                      |  |  |      |           |
|-------------------------------|----------------------|--|--|------|-----------|
| Abstract Syntax               |                      | Transfer Syntax  |  | Role | Ext. Neg. |
| Name                          | UID                  | Name   | UID                                      |      |           |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None      |

A Presentation Context for the Storage Commitment Push Model will only be proposed if the remote AE is configured as an archive device.



### 4.2.10.3.1.3 SOP Specific Conformance for Storage Commitment SOP Class

#### 4.2.10.3.1.3.1 Storage Commitment Notifications (N-EVENT-REPORT)

The Storage Commitment SCP AE provides standard conformance to the Storage Commitment Service Class as an SCP.

The behavior of Storage SCP Commitment AE when encountering status codes in an N-EVENT-REPORT response is summarized in the Table below:

**Table 4.2-91  
STORAGE COMMITMENT N-EVENT-REPORT RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning      | Status Code            | Behavior  |
|----------------|----------------------|------------------------|---|
| Success        | Success              | 0000                   | The SCU has successfully received the Storage Commitment PushModel N-EVENT-REPORT Request. Success indication message is output to the Service Logs. No message is posted to the User Interface.    |
| Warning        | Attribute List Error | 0107                   | Transmission of Storage Commitment Push Model N-EVENT-REPORTRequest is considered successful. Warning indication message is output to the Service Logs. No message is posted to the User Interface. |
| *              | *                    | Any other status code. | This is treated as a permanent Failure. Error indication message is output to the Service Logs. No message is posted to the User Interface.   |

: The behavior of Storage Commitment AE during communication failure is summarized in the Table below:

**Table 4.2-92  
STORAGE COMMITMENT COMMUNICATION FAILURE BEHAVIOR**


| Exception  | Behavior   |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |
| Association aborted by the SCU or network layers | The send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application.  |

## 4.2.10.4 Association Acceptance Policy

### 4.2.10.4.1 Activity – Receive Storage Commitment Response

#### 4.2.10.4.1.1 Description and Sequencing of Activities

The Storage Commitment SCP AE will accept Associations in order to receive responses to a Storage Commitment Request(N-ACTION).



**Figure 4.2-12**  
**SEQUENCING OF ACTIVITY – COMMIT RECEIVE OBJECTS**

A possible sequence of interactions between the Storage Commitment SCU AE and the Storage Commitment SCP AE is illustrated in the Figure above:

1. The Storage Commitment SCU AE opens an Association with the Storage Commitment SCP AE.
2. A Storage Commitment request (N-ACTION) is transmitted to the Storage Commitment SCP AE to obtain Storage Commitment of previously transmitted objects. The Storage Commitment SCP AE replies with an N-ACTION response indicating the request has been received and is being processed.
3. The Storage Commitment SCU AE closes the Association with the Storage Commitment SCP AE.

The Storage Commitment SCP AE may reject Association attempts as shown in the Table below. The Result, Source and Reason/Diag columns represent the values returned in the appropriate fields of an ASSOCIATE-RJ PDU (see PS 3.8, Section 9.3.4). The contents of the Source column is abbreviated to save space and the meaning of the abbreviations are:

- DICOM UL service-user
- DICOM UL service-provider (ASCE related function)
- DICOM UL service-provider (Presentation related function)

**Table 4.2-93**  
**ASSOCIATION REJECTION REASONS**

| <b>Result</b>             | <b>Source</b> | <b>Reason/Diag</b>                                | <b>Explanation</b>   |
|---------------------------|---------------|---|--|
| 2 –<br>rejected-transient | c             | 2 –<br>local-limit-exceeded                       | The (configurable) maximum number of simultaneous Associations has been reached. An Association request with the same parameters may succeed at a later time.  |
| 2 –<br>rejected-transient | c             | 1 –<br>temporary-congestion                       | No Associations can be accepted at this time due to the real-time requirements of higher priority activities (e.g. during image acquisition no Associations will be accepted) or because insufficient resources are available (e.g. memory, processes, threads). An Association request with the same parameters may succeed at a later time.                        |
| 1 –<br>rejected-permanent | a             | 2 –<br>application-context-<br>name-not-supported | The Association request contained an unsupported Application Context Name. An Association request with the same parameters will not succeed at a later time.   |
| 1 –<br>rejected-permanent | a             | 7 – called-AE-title-not-<br>recognized            | The Association request contained an unrecognized Called AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association initiator is incorrectly configured and attempts to address the Association acceptor using the wrong AE Title. |
| 1 –<br>rejected-permanent | a             | 3 – calling-AE-title-not-<br>recognized           | The Association request contained an unrecognized Calling AE Title. An Association request with the same parameters will not succeed at a later time unless configuration changes are made. This rejection reason normally occurs when the Association acceptor has not been configured to recognize the AE Title of the Association initiator.                      |
| 1 –<br>rejected-permanent | b             | 1 – no-reason-given                               | The Association request could not be parsed. An Association request with the same format will not succeed at a later time.   |

#### 4.2.10.4.1.2 Accepted Presentation Contexts

The Storage Commitment SCP AE will prefer to select the Explicit VR Little Endian Transfer Syntax if multiple transfer syntaxes are offered.

Any of the Presentation Contexts shown in the following table are acceptable to the Storage Commitment SCP AE.

**Table 4.2-94  
ACCEPTABLE PRESENTATION CONTEXTS BY  
ACTIVITY RECEIVE STORAGE COMMITMENT REQUEST**

| Presentation Context Table    |                      |  |  |      |           |
|-------------------------------|----------------------|--|--|------|-----------|
| Abstract Syntax               |                      | Transfer Syntax  |  | Role | Ext. Neg. |
| Name                          | UID                  | Name   | UID                                      |      |           |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None      |

#### 4.2.10.4.1.3 SOP Specific Conformance for Storage Commitment SOP Class

##### 4.2.10.4.1.3.1 Storage Commitment Operation (N-ACTION)

The Storage Commitment SCP AE provides standard conformance to the Storage Commitment Service Class as an SCP

The behavior of Storage Commitment SCP AE when receiving Service Status within the N-ACTION is summarized in the Table below.

**Table 4.2-95  
STORAGE COMMITMENT N-ACTION RESPONSE STATUS RETURN BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>  | <b>Status Code</b> | <b>Behavior</b>   |
|-----------------------|-------------------------|--------------------|---|
| Success               | Success                 | 0000               | The SCP has successfully received the Storage Commitment Push Model N-ACTION Request and can process the commitment request for the indicated SOP Instances.                |
| Error                 | Processing Failure      | 0110               | Indicates that the Storage Commitment Push Model N-ACTION Request cannot be parsed or fully processed due to a database or system failure.                                  |
| Error                 | Missing Attribute       | 0120               | Indicates that the Storage Commitment Push Model N-ACTION Request cannot be processed because a required attribute is missing from the N-ACTION Request Data Set.           |
| Error                 | Missing Attribute Value | 0121               | Indicates that the Storage Commitment Push Model N-ACTION Request cannot be processed because a Type 1 attribute in the N-ACTION Request Data Set does not specify a value. |

The behavior of Storage Commitment AE during communication failure is summarized in the Table below:

**Table 4.2-96  
STORAGE COMMITMENT COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>  |
|--|--|
| Timeout  | The Association is aborted using A-ABORT and the send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |
| Association aborted by the SCP or network layers | The send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application.  |

## 4.2.11 Print SCU AE Specification

### 4.2.11.1 SOP Classes

The Print SCU AE provides Standard Conformance to the following Meta SOP Classes:

**Table 4.2-97**

#### META SOP CLASSES FOR PRINT SCU AE

| SOP Class Name                        | SOP Class UID         | SCU | SCP |
|---------------------------------------|-----------------------|-----|-----|
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9 | Yes | No  |

The above Meta SOP Classes are defined by the following set of supported SOP Classes:

**Table 4.2-98**

#### SOP CLASSES FOR PRINT SCU AE

| SOP Class Name                      | SOP Class UID          | SCU | SCP |
|-------------------------------------|------------------------|-----|-----|
| Basic Film Session SOP Class        | 1.2.840.10008.5.1.1.1  | Yes | No  |
| Basic film Box SOP Class            | 1.2.840.10008.5.1.1.2  | Yes | No  |
| Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4  | Yes | No  |
| Printer SOP Class                   | 1.2.840.10008.5.1.1.16 | Yes | No  |

### 4.2.11.2 Association Policies

#### 4.2.11.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-99**

#### DICOM APPLICATION CONTEXT FOR PRINT SCU AE

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.11.2.2 Number of Associations

The Print SCU AE initiates one Association at a time.

**Table 4.2-100**

#### NUMBER OF ASSOCIATIONS INITIATED FOR PRINT SCU AE

|   |   |
|---|---|
| Maximum number of simultaneous Associations | 1 |
|---|---|

#### 4.2.11.2.3 Asynchronous Nature

The Print SCU AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 4.2-101**

#### ASYNCHRONOUS NATURE FOR PRINT SCU AE

|   |   |
|---|---|
| Maximum Number of Outstanding Asynchronous Transactions | 1 |
|---|---|

#### 4.2.11.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-102**

**DICOM IMPLEMENTATION CLASS AND VERSION FOR PRINT SCU AE**


|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name | TM_MR_DCM_V3.0             |

#### 4.2.11.3 Association Initiation Policy

##### 4.2.11.3.1 Activity – Send Images & Print Management Information

###### 4.2.11.3.1.1 Description and Sequencing of Activities

A user composes images onto film sheets and requests them to be sent to a specific hardcopy device. The user can select the desired film format and number of copies.



**Figure 4.2-13**  
**SEQUENCING OF ACTIVITY – SEND IMAGES & PRINT MANAGEMENT INFORMATION**

A typical sequence of DIMSE messages sent over an association between the Print SCU AE and a Printer is illustrated in the Figure above:

1. The Print SCU AE opens an Association with the Printer.
2. N-GET on the Printer SOP Class is used to obtain current printer status information.
3. The Print SCU AE closes the Association with the Printer.
4. The Print SCU AE opens an Association with the Printer.
5. N-CREATE on the Film Session SOP Class creates a Film Session.
6. N-CREATE on the Film Box SOP Class creates a Film Box linked to the Film Session.
7. N-SET on the Image Box SOP Class transfers the contents of the film sheet to the printer.
8. N-ACTION on the Film Box SOP Class instructs the Printer to print the Film Box.
9. The Printer prints the requested number of film sheets.
10. The Printer asynchronously reports its status via N-EVENT-REPORT notification (Printer SOP Class). The printer can send this message at any time. The Print SCU AE does not require the N-EVENT-REPORT to be sent. The Print SCU AE is capable of receiving an N-EVENT-REPORT notification at any time during an association.
11. The Print SCU AE closes the Association with the Printer.

#### 4.2.11.3.1.2 Proposed Presentation Contexts

The Print SCU AE is capable of proposing the Presentation Contexts shown in the Table below:

**Table 4.2-103**  
**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY**  
**SEND IMAGES & PRINT MANAGEMENT INFORMATION**

| Presentation Context Table                  |                 |                           |                     |      |           |
|---|-----------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                             |                 | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name  | UID             | Name                      | UID                 |      |           |
| Basic Grayscale<br>Print Management<br>Meta | 1.2.840.10008.5 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|   | .1.1.9          | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |



#### 4.2.11.3.1.3 Common SOP Specific Conformance for all Print SOP Classes

The general behavior of the Print SCU AE during communication failure is summarized in the Table below. This behavior is common for all SOP Classes supported by the Print SCU AE.

**Table 4.2-104**  
**PRINT COMMUNICATION FAILURE BEHAVIOR**

| Exception  | Behavior  |
|--|---|
| Timeout  | The Association is aborted using A-ABORT and the print-job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |
| Association aborted by the SCP or network layers | The print-job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application.  |

#### 4.2.11.3.1.4 SOP Specific Conformance for Printer SOP Class

The Print SCU AE supports the following DIMSE operations and notifications for the Printer SOP Class:

- N-GET

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.11.3.1.4.1 Printer SOP Class Operations (N-GET)

The Print SCU AE uses the Printer SOP Class N-GET operation to obtain information about the current printer status. The attributes obtained via N-GET are listed in the Table below:

**Table 4.2-105**  
**PRINTER SOP CLASS N-GET REQUEST ATTRIBUTES**

| Attribute Name            | Tag         | VR | Value               | Source  |
|---------------------------|-------------|----|---------------------|---------|
| Printer Status            | (2110,0010) | CS | Provided by Printer | Printer |
| Printer Status Info       | (2110,0020) | CS | Provided by Printer | Printer |
| Printer Name              | (2110,0030) | LO | Provided by Printer | Printer |
| Manufacturer              | (0008,0070) | LO | Provided by Printer | Printer |
| Manufacturer's Model Name | (0008,1090) | LO | Provided by Printer | Printer |
| Device Serial Number      | (0018,1000) | LO | Provided by Printer | Printer |
| Software Version          | (0018,1020) | LO | Provided by Printer | Printer |
| Date of Last Calibration  | (0018,1200) | DA | Provided by Printer | Printer |
| Time of Last Calibration  | (0018,1201) | TM | Provided by Printer | Printer |

The Printer Status information is evaluated as follows:

- If Printer status (2110,0010) is NORMAL, the print-job continues to be printed.
- If Printer status (2110,0010) is FAILURE, the print-job is marked as failed. The contents of Printer Status Info (2110,0020) is logged and reported to the user via the job control application.
- If Printer status (2110,0010) is WARNING, the print-job continues to be printed. The contents of Printer Status Info (2110,0020) is logged and reported to the user via the job control application.

The behavior of The Print SCU AE when encountering status codes in an N-GET response is summarized in the Table below:

**Table 4.2-106  
PRINTER SOP CLASS N-GET RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Status Code            | Behavior   |
|----------------|-----------------|------------------------|--|
| Success        | Success         | 0000                   | The request to get printer status information was success.   |
| *              | *               | Any other status code. | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |

#### 4.2.11.3.1.5 SOP Specific Conformance for the Film Session SOP Class

The Print SCU AE supports the following DIMSE operations for the Film Session SOP Class:

- N-CREATE

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.11.3.1.5.1 Film Session SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below:

**Table 4.2-107  
FILM SESSION SOP CLASS N-CREATE REQUEST ATTRIBUTES**

| Attribute Name     | Tag         | VR | Value                                | Source |
|--------------------|-------------|----|--------------------------------------|--------|
| Number of Copies   | (2000,0010) | IS | 1 .. 99                              | User   |
| Print Priority     | (2000,0020) | CS | "MED"                                | Auto   |
| Medium Type        | (2000,0030) | CS | "BLUE FILM", "CLEAR FILM" or "PAPER" | User   |
| Film Destination   | (2000,0040) | CS | "MAGAZINE" or "PROCESSOR"            | User   |
| Film Session Label | (2000,0050) | LO | "TOSHIBA_MRI"                        | Auto   |

The behavior of The Print SCU AE when encountering status codes in a N-CREATE response is summarized in the Table below:

**Table 4.2-108  
FILM SESSION SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning              | Status Code            | Behavior   |
|----------------|------------------------------|------------------------|--|
| Success        | Success                      | 0000                   | The SCP has completed the operation successfully.  |
| Warning        | Attribute Value Out of Range | 0116                   | The N-CREATE operation is considered successful but the status meaning is logged. Additional information in the response identifying the attribute out of range is logged.(I.e., Elements in the modification list/Attribute list) |
| Warning        | Attribute List Error         | 0107                   | The N-CREATE operation is considered successful but the status meaning is logged. Additional information in the response identifying the attribute out of range is logged.(I.e., Elements in the Attribute identifier list)        |
| *              | *                            | Any other status code. | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user.   |

#### 4.2.11.3.1.6 SOP Specific Conformance for the Film Box SOP Class

The Print SCU AE supports the following DIMSE operations for the Film Box SOP Class:

- N-CREATE
- N-ACTION

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.11.3.1.6.1 Film Box SOP Class Operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below:

**Table 4.2-109**  
**FILM BOX SOP CLASS N-CREATE REQUEST ATTRIBUTES**

| Attribute Name                   | Tag         | VR | Value                                      | Source |
|----------------------------------|-------------|----|--|--------|
| Image Display Format             | (2010,0010) | ST | STANDARD\C,R                               | User   |
| Referenced Film Session Sequence | (2010,0500) | SQ |  | Auto   |
| >Referenced SOP Class UID        | (0008,1150) | UI | 1.2.840.10008.5.1.1.1                      | Auto   |
| >Referenced SOP Instance UID     | (0008,1155) | UI | From created Film Session SOP Instance     | Auto   |
| Film Orientation                 | (2010,0040) | CS | "PORTRAIT" or "LANDSCAPE"                  | User   |
| Film Size ID                     | (2010,0050) | CS | 14INX17IN, 14INX14IN, 11INX14IN, 8INX10IN  | User   |
| Magnification Type               | (2010,0060) | CS | "REPLICATE", "BILINEAR", "CUBIC" or "NONE" | User   |
| Smoothing Type                   | (2010,0080) | CS |  | User   |
| Border Density                   | (2010,0100) | CS | "BLACK" or "WHITE"                         | User   |
| Empty Image Density              | (2010,0110) | CS | "BLACK" or "WHITE"                         | User   |
| Min Density                      | (2010,0120) | US | 0 .. 9999                                  | User   |
| Max Density                      | (2010,0130) | US | 0 .. 9999                                  | User   |
| Trim                             | (2010,0140) | CS | "YES" or "NO"                              | User   |
| Configuration Information        | (2010,0150) | ST |  | Auto   |

The behavior of the Print SCU AE when encountering status codes in a N-CREATE response is summarized in the Table below:

**Table 4.2-110  
FILM BOX SOP CLASS N-CREATE RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>  | <b>Status Code</b>     | <b>Behavior</b>  |
|-----------------------|---|------------------------|--|
| Success               | Success   | 0000                   | The SCP has completed the operation successfully.  |
| Warning               | Requested Min Density or Max Density outside of printer's operating range | B605                   | The N-CREATE operation is considered successful but the status meaning is logged.  |
| *                     | *   | Any other status code. | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |

#### 4.2.11.3.1.6.2 Film Box SOP Class Operations (N-ACTION)

An N-ACTION Request is issued to instruct the Print SCP to print the contents of the Film Box. The Action Reply argument in an N-ACTION response is not evaluated.

The behavior of The Print SCU AE when encountering status codes in a N-ACTION response is summarized in the Table below:

**Table 4.2-111**  
**FILM BOX SOP CLASS N-ACTION RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning   | Status Code            | Behavior   |
|----------------|---|------------------------|--|
| Success        | Success   | 0000                   | The SCP has completed the operation successfully. The film has been accepted for printing.   |
| Warning        | Film Box SOP Instance hierarchy does not contain Image Box SOP Instances (empty page)   | B603                   | The N-ACTION operation is considered successful if it is configured that the status would be considered successful.                    |
| Warning        | Image size is larger than Image Box size. The image has been demagnified.   | B604                   | The N-ACTION operation is considered successful if it is configured that the status would be considered successful.                    |
| Warning        | Image size is larger than Image Box size. The image has been cropped to fit.  | B609                   | The N-ACTION operation is considered successful if it is configured that the status would be considered successful.                    |
| Warning        | Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit. | B60A                   | The N-ACTION operation is considered successful if it is configured that the status would be considered successful.                    |
| Failure        | Unable to create Print Job SOP Instance; print queue is full.   | C602                   | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure        | Image size is larger than Image Box size.   | C603                   | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure        | Combined Print Image Size is larger than Image Box size.  | C613                   | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| *              | *   | Any other status code. | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |

#### 4.2.11.3.1.7 SOP Specific Conformance for the Grayscale Image Box SOP Class

The Print SCU AE supports the following DIMSE operations for the Grayscale Image Box SOP Class:

- N-SET

Details of the supported attributes and status handling behavior are described in the following subsections.

##### 4.2.11.3.1.7.1 Grayscale Image Box SOP Class Operations (N-SET)

The attributes supplied in an N-SET Request are listed in the Table below:

**Table 4.2-112**  
**GRAYSCALE IMAGE BOX SOP CLASS N-SET REQUEST ATTRIBUTES**

| Attribute Name                 | Tag         | VR | Value                                      | Source |
|--------------------------------|-------------|----|--|--------|
| Image Position                 | (2020,0010) | US | 1 .. 36                                    | Auto   |
| Polarity                       | (2020,0020) | CS | "NORMAL" or "REVERSE"                      | User   |
| Magnification Type             | (2010,0060) | CS | "REPLICATE", "BILINEAR", "CUBIC" or "NONE" | User   |
| Smoothing Type                 | (2010,0080) | CS |  | User   |
| Basic Grayscale Image Sequence | (2020,0110) | SQ |  | Auto   |
| >Samples Per Pixel             | (0028,0002) | US | 1  | Auto   |
| >Photometric Interpretation    | (0028,0004) | CS | "MONOCHROME2"                              | Auto   |
| >Rows                          | (0028,0010) | US |  | Auto   |
| >Columns                       | (0028,0011) | US |  | Auto   |
| >Pixel Aspect Ratio            | (0028,0034) | IS | 1\1  | Auto   |
| >Bits Allocated                | (0028,0100) | US | 8  | Auto   |
| >Bits Stored                   | (0028,0101) | US | 8  | Auto   |
| >High Bit                      | (0028,0102) | US | 7  | Auto   |
| >Pixel Representation          | (0028,0103) | US | 1  | Auto   |
| >Pixel Data                    | (7FE0,0010) | OB |  | Auto   |

The behavior of the Print SCU AE when encountering status codes in a N-SET response is summarized in the Table below:

**Table 4.2-113  
GRAYSCALE IMAGE BOX SOP CLASS N-SET RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b>  | <b>Status Code</b>     | <b>Behavior</b>  |
|-----------------------|---|------------------------|--|
| Success               | Success   | 0000                   | The SCP has completed the operation successfully. Image successfully stored in Image Box.  |
| Warning               | Image size is larger than Image Box size. The image has been demagnified.   | B604                   | The N-SET operation is considered successful but the status meaning is logged.   |
| Warning               | Requested Min Density or Max Density outside of printer's operating range.  | B605                   | The N-SET operation is considered successful but the status meaning is logged.   |
| Warning               | Image size is larger than Image Box size. The image has been cropped to fit.  | B609                   | The N-SET operation is considered successful but the status meaning is logged.   |
| Warning               | Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit. | B60A                   | The N-SET operation is considered successful but the status meaning is logged.   |
| Failure               | Image size is larger than Image Box size.   | C603                   | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure               | Insufficient memory in printer to store the image.  | C605                   | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure               | Combined Print Image Size is larger than Image Box size.  | C613                   | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| *                     | *   | Any other status code. | The Association is aborted using A-ABORT and the print-job is marked as failed. The status meaning is logged and reported to the user. |

#### 4.2.11.4 Association Acceptance Policy

The Print SCU AE does not accept Associations.

### 4.3 NETWORK INTERFACES

#### 4.3.1 Physical Network Interface

This Product supports a single network interface. One of the following physical network interfaces will be available depending on installed hardware options:

**Table 4.3-1  
SUPPORTED PHYSICAL NETWORK INTERFACES**

|                    |
|--------------------|
| Ethernet 1000baseT |
| Ethernet 100baseT  |

#### 4.3.2 Additional Protocols

None.

## 4.4 CONFIGURATION

### 4.4.1 AE Title/Presentation Address Mapping

#### 4.4.1.1 Local AE Titles

All local applications use the AE Titles and TCP/IP Ports configured by the Field Service Engineer via the Service/Installation Tool.

**Table 4.4-1  
AE TITLE CONFIGURATION TABLE**

| Application Entity     | Default AE Title | Default TCP/IP Port |
|------------------------|------------------|---------------------|
| Verification SCU       | TM_MR_DCM_V1.2   | Not Applicable      |
| MWM SCU                |                  |                     |
| MPPS SCU               |                  |                     |
| Q/R SCU                |                  |                     |
| Print SCU              |                  |                     |
| Storage SCU            |                  |                     |
| Verification SCP       |                  | 104                 |
| Storage Commitment SCU |                  |                     |
| Storage SCP            |                  |                     |
| Q/R SCP                |                  |                     |
| Storage Commitment SCP |                  |                     |

#### 4.4.1.2 Remote AE Title/Presentation Address Mapping

The AE Titles, host names and port numbers of remote applications are configured using the Utility Tool. The system is case-insensitive for the AE Titles and configuration by the same AE Title is prohibited. Characters excluding shaded characters in the table below can be used as AE titles.

**Table 4.4-2  
AE TITLE CHARACTER REPERTOIRE**

|      | 0  | 1 | 2 | 3 | 4  | 5 | 6 | 7 | 8 | 9 | A | B | C | D | E | F |
|------|----|---|---|---|----|---|---|---|---|---|---|---|---|---|---|---|
| 0x20 | SP | ! | " | # | \$ | % | & | ' | ( | ) | * | + | , | - | . | / |
| 0x30 | 0  | 1 | 2 | 3 | 4  | 5 | 6 | 7 | 8 | 9 | : | ; | < | = | > | ? |
| 0x40 | @  | A | B | C | D  | E | F | G | H | I | J | K | L | M | N | O |
| 0x50 | P  | Q | R | S | T  | U | V | W | X | Y | Z | [ | \ | ] | ^ | _ |
| 0x60 | `  | a | b | c | d  | e | f | g | h | i | j | k | l | m | n | o |
| 0x70 | p  | q | r | s | t  | u | v | w | x | y | z | { |   | } | ~ |   |



## 4.4.2 Parameters

A large number of parameters related to acquisition and general operation can be configured using the Utility Tool. The Table below only shows those configuration parameters relevant to DICOM communication. See the Product's Service Manual for details on general configuration capabilities.

**Table 4.4-3  
CONFIGURATION PARAMETERS TABLE**

| Parameter   | Configurable<br>(Yes/No)<br>[Range]       | Default Value             |
|---|---|---------------------------|
| <b>General Parameters</b>   |   |                           |
| Max PDU Receive Size  | Yes                                       | 64 Kbytes                 |
| Max PDU Send Size   | [1-999999]                                |                           |
| Time-out waiting for a acceptance or rejection response to an Association Request (Application Level Timeout) | Yes<br>[1-999999]                         | 30 sec                    |
| Time-out waiting for a response to an Association release request (Application Level Timeout)                 | Yes<br>[1-999999]                         | 15 sec                    |
| Time-out waiting for completion of a TCP/IP connect request (Low-level timeout)                               | Yes<br>[1-999999]                         | 15 sec                    |
| Time-out awaiting a Response to a DIMSE Request (Low-Level Timeout)   | Yes<br>[1-999999]                         | 15 sec                    |
| Time-out for waiting for data between TCP/IP-packets (Low Level Timeout)                                      | Yes<br>[1-999999]                         | 15 sec                    |
| Supported Transfer Syntaxes (separately configurable for each service class and remote AE)                    | Yes<br>[ILE,ELE,EBE]                      | Implicit VR Little Endian |
| <b>Modality Worklist SCU Parameters</b>   |   |                           |
| Modality Worklist SCU time-out waiting for the final response to a C-FIND-RQ                                  | Yes<br>[1-999999]                         | 30sec                     |
| Maximum number of simultaneously initiated Associations by the MWM SCU AE                                     | No  | 1                         |
| <b>MPPS SCU Parameters</b>  |   |                           |
| MPPS SCU time-out waiting for a response to a N-CREATE-RQ   | Yes<br>[1-999999]                         | 30sec                     |
| MPPS SCU time-out waiting for a response to a N-SET-RQ  | Yes<br>[1-999999]                         | 30sec                     |
| MPPS SCU time-out waiting for a response to a N-GET-RQ  | Yes<br>[1-999999]                         | 30sec                     |
| Maximum number of simultaneously initiated Associations by the MPPS SCU AE                                    | No  | 1                         |
| Supported Transfer Syntaxes (separately configurable for each remote AE)                                      | Yes                                       | Implicit VR Little Endian |
| Behavior when receiving the Warning "Attribute Value Out of Range" as service status.                         | Yes<br>[Considered as Success or Failure] | Considered as Failure     |

| Parameter  | Configurable<br>(Yes/No)<br>[Range]       | Default Value                                   |
|--|---|---|
| <b>Storage SCU Parameters</b>  |   |   |
| Storage SCU time-out waiting for a response to a C-STORE-RQ  | Yes<br>[1-999999]                         | 30sec   |
| Number of times a failed send job may be retried   | Yes<br>[1-99999]                          | 10 times  |
| Delay between retrying failed send jobs  | Yes<br>[1-99999]                          | 60sec   |
| Maximum number of simultaneously initiated Associations by the Storage SCU AE  | No  | 1   |
| Behavior when receiving the Warning "Coercion of Data Elements" as service status.   | Yes<br>[Considered as Success or Failure] | Considered as Failure                           |
| Behavior when receiving the Warning "Data Set does not match SOP Class" as service status.   | Yes<br>[Considered as Success or Failure] | Considered as Failure                           |
| Behavior when receiving the Warning "Elements Discarded" as service status.  | Yes<br>[Considered as Success or Failure] | Considered as Failure                           |
| <b>Storage Commitment SCU Parameters</b>   |   |   |
| Storage Commitment SCU time-out waiting for a response to a N-ACTION-RQ  | Yes<br>[1-999999]                         | 30 sec  |
| Maximum number of simultaneously initiated Associations by the Storage Commitment SCU AE   | No  | 1   |
| Timeout waiting for a Storage Commitment Notification (maximum duration of applicability for a Storage Commitment Transaction UID).          | Yes<br>[1sec-10hrs]                       | 3 min   |
| Maximum number of simultaneously accepted Associations by the Storage Commitment SCU AE  | No  | 10  |
| Delay association release after sending a Storage Commitment Request (wait for a Storage Commitment Notification over the same association). | No  | 0   |
| Behavior when receiving N-EVENT Report - the Storage Commitment Request Successful.  | Yes                                       | Permit the operator (s) to delete the Instances |
| <b>Storage SCP parameters</b>  |   |   |
| Maximum number of simultaneously accepted Associations by the Storage SCP AE   | No  | 10  |


| Parameter   | Configurable<br>(Yes/No)<br>[Range]       | Default Value             |
|---|---|---------------------------|
| <b>Print SCU Parameters</b>   |   |                           |
| Print SCU time-out waiting for a response to a N-GET-RQ   | No  | 30sec                     |
| Print SCU time-out waiting for a response to a N-CREATE-RQ  | No  | 30sec                     |
| Print SCU time-out waiting for a response to a N-SET-RQ   | No  | 30sec                     |
| Print SCU time-out waiting for a response to a N-ACTION-RQ  | No  | 30sec                     |
| Maximum number of simultaneously initiated Associations by the Print SCU AE   | No  | 1                         |
| Supported Transfer Syntaxes (separately configurable for each remote printer)   | Yes                                       | Implicit VR Little Endian |
| Behavior when receiving the Warning "Attribute Value Out of Range" as service status of the Film Session N-CREATE.  | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Attribute List Error" as service status of the Film Session N-CREATE.  | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Requested Min Density or Max Density outside of printer's operating range" as service status of the Film Box N-CREATE.   | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Film Box SOP Instance hierarchy does not contain Image Box SOP Instances (empty page)" as service status of the Film Box N-ACTION.   | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Image size is larger than Image Box size. The image has been demagnified." as service status of the Film Box N-ACTION.   | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Image size is larger than Image Box size. The image has been cropped to fit." as service status of the Film Box N-ACTION.  | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit." as service status of the Film Box N-ACTION. | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Image size is larger than Image Box size. The image has been demagnified." as service status of the Grayscale Image Box N-SET.   | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Requested Min Density or Max Density outside of printer's operating range." as service status of the Grayscale Image Box N-SET.  | Yes<br>[Considered as Success or Failure] | Considered as Failure     |
| Behavior when receiving the Warning "Image size is larger than Image Box size. The image has been cropped to fit." as service status of the Grayscale Image Box N-SET.  | Yes<br>[Considered as Success or Failure] | Considered as Failure     |

| Parameter   | Configurable<br>(Yes/No)<br>[Range]             | Default Value            |
|---|---|--------------------------|
| Behavior when receiving the Warning "Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit." as service status of the Grayscale Image Box N-SET. | Yes<br>[Considered<br>as Success or<br>Failure] | Considered as<br>Failure |
| Behavior when receiving the Warning "Image size is larger than Image Box size. The image has been demagnified." as service status of the Color Image Box N-SET.   | Yes<br>[Considered<br>as Success or<br>Failure] | Considered as<br>Failure |
| Behavior when receiving the Warning "Requested Min Density or Max Density outside of printer's operating range." as service status of the Color Image Box N-SET.  | Yes<br>[Considered<br>as Success or<br>Failure] | Considered as<br>Failure |
| Behavior when receiving the Warning "Image size is larger than Image Box size. The image has been cropped to fit." as service status of the Color Image Box N-SET.  | Yes<br>[Considered<br>as Success or<br>Failure] | Considered as<br>Failure |
| Behavior when receiving the Warning "Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit." as service status of the Color Image Box N-SET.     | Yes<br>[Considered<br>as Success or<br>Failure] | Considered as<br>Failure |

## 5. MEDIA INTERCHANGE

### 5.1 IMPLEMENTATION MODEL

#### 5.1.1 Application Data Flow



**Figure 5.1-1**  
**APPLICATION DATA FLOW DIAGRAM FOR MEDIA STORAGE**

- The Offline-Media AE exports image files to an optical medium. It is associated with the local real-world activity “Export Image Files” performed upon user request.

## 5.1.2 Functional Definition of AEs

### 5.1.2.1 Functional Definition of Offline-Media AE

The Offline-Media AE is performed upon user request for selected studies /series/images to/from an offline DICOM CD-R medium. It therefore performs the following tasks:

Export:

- Builds DICOM Information Objects.
- Creates a DICOMDIR file that represents the contents of the DICOM Information Objects to be recorded.
- Records DICOM Information Objects and the DICOMDIR file to the DVD-R or the DVD-RAM medium.

Import:

- Reads the DICOMDIR file that represents the contents of the data as recorded.
- Displays the ordered list of studies/series/images, identifying information.
- Loads the selected studies/series/images from the optical medium and displays them on the screen.

Addition:

- Reads a File-set of the DVD-RAM medium and writes it to the local storage device.
- Adds the studies/series/images to the File-Set, and then writes it to the medium.
- Modifies the DICOMDIR file.

### 5.1.3 Sequencing of Real-World Activities

#### 5.1.3.1 Activity - Export Image Files

Operator requests to create new File-set(s) onto a new DVD-R or a new DVD-RAM. The requests are placed in a queue and are executed in the background.

The operations for “Export Image Files” are described below:

- Step-1: Select the instance(s), series or studies on the local storage device to be created to the optical medium.
- Step-2: Select the image archiving.
- Step-3: Select the DVD device as a destination.
- Step-4: Request to copy to the DVD-R.

#### 5.1.3.2 Activity - Import Image Files

Operator requests to retrieve File-set(s) on the optical medium. The requests are placed in a queue and are executed in the background.

The operations for “Import Image Files” are described below:

- Step-1: Select the MR image(s) and/or SC image(s), series or studies on the medium to be retrieved to the local storage device.
- Step-2: Select the data retrieval.
- Step-3: Select the local storage device as a destination.

#### 5.1.3.3 Activity – Add Image Files

Operator requests to add new objects to an already existing File-set on the DVD-RAM. The requests are placed in a queue and are executed in the background.

The operations for “Add Image Files” are described below:

- Step-1: Select the instance(s), series or studies on the local storage device to be added to the DVD-RAM medium.
- Step-2: Select the image archiving.
- Step-3: Select the DVD device as a destination.
- Step-4: Request to copy to the DVD-RAM.

### 5.1.4 File Meta Information for Implementation Class and Version

The implementation information written to the File Meta Header in each file is:

**Table 5.1-1  
DICOM IMPLEMENTATION CLASS AND VERSION FOR MEDIA STORAGE**

|                               |                            |
|-------------------------------|----------------------------|
| File Meta Information Version | 1                          |
| Implementation Class UID      | 1.2.392.200036.9116.4.2.10 |
| Implementation Version Name   | TM_MR_DCM_V3.0             |

## 5.2 AE SPECIFICATIONS

### 5.2.1 Offline-Media AE Specification

The Offline-Media AE provides standard conformance to the DICOM Interchange Option of the Media Storage Service Class. The Application Profiles and roles are listed below:

**Table 5.2-1**

**APPLICATION PROFILES, ACTIVITIES AND ROLES FOR OFFLINE-MEDIA**

| Application Profiles Supported | Real World Activity | Role | SC Option   |
|--------------------------------|---------------------|------|-------------|
| STD-CTMR-CD, STD-GEN-CD        | Export Image Files  | FSC  | Interchange |

#### 5.2.1.1 File Meta Information for the Application Entity

The Offline-Media AE does not set the Source Application Entity Title.

#### 5.2.1.2 Real-World Activities

##### 5.2.1.2.1.1 Activity – Export Image Files

The Offline-Media AE acts as an FSC using the interchange option when requested to export SOP Instances from the local database to a CD-R medium.

##### 5.2.1.2.1.2 Media Storage Application Profiles

The Offline-Media AE supports the STD-CTMR-CD and the STD-GEN-CD Application Profile.

##### 5.2.1.2.1.2.1 Options

The Offline-Media AE supports the SOP Classes and Transfer Syntaxes listed in the Table below for the STD-CTMR-CD Application Profile as an FSC.

**Table 5.2-2**

**IODS, SOP CLASSES AND TRANSFER SYNTAXES FOR STD-CTMR-CD PROFILE**

| Information Object Definition   | SOP Class UID               | Transfer Syntax                           | Transfer Syntax UID |
|---------------------------------|-----------------------------|---|---------------------|
| Media Storage Directory Storage | 1.2.840.10008.1.3.10        | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1 |
| MR Image Storage                | 1.2.840.10008.5.1.4.1.1.4   | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7   | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1 |
| Enhanced MR Image Storage       | 1.2.840.10008.5.1.4.1.1.4.1 | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1 |
| MR Spectroscopy                 | 1.2.840.10008.5.1.4.1.1.4.2 | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1 |



The Offline-Media AE supports the SOP Classes and Transfer Syntaxes listed in the Table below for the STD-GEN-CD Application Profile as an FSC.

**Table 5.2-3**  
**IODS, SOP CLASSES AND TRANSFER SYNTAXES FOR STD-GEN-CD PROFILE (FSC)**

| <b>Information Object Definition</b> | <b>SOP Class UID</b>        | <b>Transfer Syntax</b>                    | <b>Transfer Syntax UID</b> |
|--------------------------------------|-----------------------------|---|----------------------------|
| Media Storage Directory Storage      | 1.2.840.10008.1.3.10        | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1        |
| MR Image Storage                     | 1.2.840.10008.5.1.4.1.1.4   | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1        |
| Secondary Capture Image Storage      | 1.2.840.10008.5.1.4.1.1.7   | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1        |
| Enhanced MR Image Storage            | 1.2.840.10008.5.1.4.1.1.4.1 | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1        |
| MR Spectroscopy                      | 1.2.840.10008.5.1.4.1.1.4.2 | Explicit VR Little Endian<br>Uncompressed | 1.2.840.10008.1.2.1        |

### **5.3 MEDIA CONFIGURATION**

Not applicable to the Offline-Media AE.



## 7. SECURITY

This product does not support any specific security measures.

It is assumed that the product is used within a secured environment. It is assumed that a secured environment includes at a minimum:

- Firewall or router protections to ensure that only approved external hosts have network access to the product.
- Firewall or router protections to ensure that the product only has network access to approved external hosts and services.
- Any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g. such as a Virtual Private Network (VPN))

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

### 7.1 DE-IDENTIFICATION

This product supports the following requirements for de-identification. De-identification occurs, for example, when an anonymize patient operation is performed. This de-identification only covers standard DICOM attributes that contain protected patient information. The encrypted data capability is not supported, so once an image has been de-identified no facility is provided to recover the lost information.

De-identification supports the following object types:

- MR Image Objects
- Secondary capture images
- Enhanced MR Image Objects
- MR Spectroscopy Objects

De-identification does not support the following object types:

- Grayscale Softcopy Presentation State Objects

Table 7.1-1 provides the list of attributes and the expected action when de-identifying images.

**Table 7.1-1  
DICOM Attributes De-Identified(Patient/Study Module)**

| Attribute Name       | Tag         | Action |
|----------------------|-------------|--------|
| Patient's Name       | (0010,0010) | D      |
| Patient ID           | (0010,0020) | D      |
| Issuer of Patient ID | (0010,0021) | X      |
| Patient's Birth Date | (0010,0030) | D      |
| Patient's Birth Time | (0010,0032) | X      |
| Other Patient IDs    | (0010,1000) | X      |
| Other Patient Names  | (0010,1001) | X      |

|  |             |   |
|--|-------------|---|
| Patient Comments                                       | (0010,4000) | Z |
| Instance Creator UID                                   | (0008,0014) | U |
| Institution Code Sequence                              | (0008,0082) | X |
| Referring Physician's Address                          | (0008,0092) | X |
| Referring Physician's Telephone Numbers                | (0008,0094) | X |
| Admitting Diagnoses Description                        | (0008,1080) | X |
| Admitting Diagnoses Code Sequence                      | (0008,1084) | X |
| Patient's Insurance Plan Code Sequence                 | (0010,0050) | X |
| Patient's Birth Name                                   | (0010,1005) | X |
| Patient Address  | (0010,1040) | X |
| Patient's Mother's Birth Name                          | (0010,1060) | X |
| Military Rank  | (0010,1080) | X |
| Branch of Service                                      | (0010,1081) | X |
| Medical Record Locator                                 | (0010,1090) | X |
| Medical Alerts   | (0010,2000) | X |
| Allergies  | (0010,2110) | X |
| Country of Residence                                   | (0010,2150) | X |
| Region of Residence                                    | (0010,2152) | X |
| Patient's Telephone Number                             | (0010,2154) | X |
| Smoking Status   | (0010,21A0) | X |
| Pregnancy Status                                       | (0010,21C0) | X |
| Last Menstrual Date                                    | (0010,21D0) | X |
| Patient's Religious Preference                         | (0010,21F0) | X |
| Admission ID   | (0038,0010) | X |
| Issuer of Admission ID                                 | (0038,0011) | X |
| Admitting Date   | (0038,0020) | X |
| Admitting Time   | (0038,0021) | X |
| Special Needs  | (0038,0050) | X |
| Current Patient Location                               | (0038,0300) | X |
| Patient's Institution Residence                        | (0038,0400) | X |
| Patient State  | (0038,0500) | X |
| Visit Comments   | (0038,4000) | X |
| Confidentiality Constraint on Patient Data Description | (0040,3001) | X |
| Study Date   | (0008,0020) | Z |
| Study Time   | (0008,0030) | Z |
| Accession Number                                       | (0008,0050) | Z |
| Referring Physician's Name                             | (0008,0090) | Z |
| Referenced Study Sequence                              | (0008,1110) | X |
| Study ID   | (0020,0010) | D |
| Requesting Physician                                   | (0032,1032) | X |

|  |             |   |
|--|-------------|---|
| Requesting Service                             | (0032,1033) | X |
| Requested Procedure Description                | (0032,1060) | X |
| Requested Contrast Agent                       | (0032,1070) | X |
| Scheduled Station AE Title                     | (0040,0001) | X |
| Scheduled Procedure Step Start Date            | (0040,0002) | X |
| Scheduled Procedure Step Start Time            | (0040,0003) | X |
| Scheduled Procedure Step End Date              | (0040,0004) | X |
| Scheduled Procedure Step End Time              | (0040,0005) | X |
| Scheduled Performing Physician Name            | (0040,0006) | X |
| Scheduled Station Name                         | (0040,0010) | X |
| Scheduled Procedure Step Location              | (0040,0011) | X |
| Pre-Medication                                 | (0040,0012) | X |
| Patient Transport Arrangements                 | (0040,1004) | X |
| Requested Procedure Location                   | (0040,1005) | X |
| Names of Intended Recipient of Results         | (0040,1010) | X |
| Requested Procedure Comments                   | (0040,1400) | X |
| Reason for Imaging Service Request             | (0040,2001) | X |
| Order Entered By                               | (0040,2008) | X |
| Order Enterer Location                         | (0040,2009) | X |
| Order Callback Phone Number                    | (0040,2010) | X |
| Placer Order Number of Imaging Service Request | (0040,2016) | Z |
| Filler Order Number of Imaging Service Request | (0040,2017) | Z |
| Imaging Service Request Comments               | (0040,2400) | X |

**Table 7.1-2  
DICOM Attributes De-Identified(MR/SC Image IOD)**

| <b>Attribute Name</b>                        | <b>Tag</b>  | <b>Action</b> |
|--|-------------|---------------|
| Series Date                                  | (0008,0021) | D             |
| Performing Physician's Name                  | (0008,1050) | X             |
| Operator's Name                              | (0008,1070) | X             |
| Referenced Performed Procedure Step Sequence | (0008,1111) | X             |
| Performed Procedure Step Date                | (0040,0244) | X             |
| Performed Procedure Step Time                | (0040,0245) | X             |
| Performed Procedure Step ID                  | (0040,0253) | X             |
| Performed Procedure Step Description         | (0040,0254) | X             |
| Request Attributes Sequence                  | (0040,0275) | X             |
| Institution Name                             | (0008,0080) | D             |
| Institution Address                          | (0008,0081) | X             |
| Station Name                                 | (0008,1010) | X             |
| Institutional Department Name                | (0008,1040) | X             |
| Acquisition Date                             | (0008,0022) | X             |

**Table 7.1-3  
DICOM Attributes De-Identified(Enhanced MR Image IOD)**

| <b>Attribute Name</b>                        | <b>Tag</b>  | <b>Action</b> |
|--|-------------|---------------|
| Series Date                                  | (0008,0021) | D             |
| Performing Physician's Name                  | (0008,1050) | X             |
| Operator's Name                              | (0008,1070) | X             |
| Referenced Performed Procedure Step Sequence | (0008,1111) | X             |
| Performed Procedure Step Date                | (0040,0244) | X             |
| Performed Procedure Step Time                | (0040,0245) | X             |
| Performed Procedure Step ID                  | (0040,0253) | X             |
| Performed Procedure Step Description         | (0040,0254) | X             |
| Request Attributes Sequence                  | (0040,0275) | X             |
| Referenced Performed Procedure Step Sequence | (0008,1111) | X             |
| Institution Name                             | (0008,0080) | D             |
| Institution Address                          | (0008,0081) | X             |
| Station Name                                 | (0008,1010) | X             |
| Institutional Department Name                | (0008,1040) | X             |
| Institution Name                             | (0008,0080) | D             |
| Institution Address                          | (0008,0081) | X             |
| Station Name                                 | (0008,1010) | X             |
| Acquisition Context Sequence                 | (0040,0555) | X             |
| Icon Image Sequence                          | (0088,0200) | X             |
| Frame Comments                               | (0020,9158) | X             |
| Acquisition Date                             | (0008,0022) | X             |

**Table 7.1-4  
DICOM Attributes De-Identified(Enhanced MR Image IOD)**

| <b>Attribute Name</b>                        | <b>Tag</b>  | <b>Action</b> |
|--|-------------|---------------|
| Series Date                                  | (0008,0021) | D             |
| Performing Physician's Name                  | (0008,1050) | X             |
| Operator's Name                              | (0008,1070) | X             |
| Referenced Performed Procedure Step Sequence | (0008,1111) | X             |
| Performed Procedure Step Date                | (0040,0244) | X             |
| Performed Procedure Step Time                | (0040,0245) | X             |
| Performed Procedure Step ID                  | (0040,0253) | X             |
| Performed Procedure Step Description         | (0040,0254) | X             |
| Request Attributes Sequence                  | (0040,0275) | X             |
| Referenced Performed Procedure Step Sequence | (0008,1111) | X             |
| Institution Name                             | (0008,0080) | D             |
| Institution Address                          | (0008,0081) | X             |
| Station Name                                 | (0008,1010) | X             |
| Institutional Department Name                | (0008,1040) | X             |
| Institution Name                             | (0008,0080) | D             |
| Institution Address                          | (0008,0081) | X             |
| Station Name                                 | (0008,1010) | X             |
| Acquisition Context Sequence                 | (0040,0555) | X             |
| Frame Comments                               | (0020,9158) | X             |
| Acquisition Date                             | (0008,0022) | X             |

The above table should be read as follows:

Attribute Name: Attributes supported to De-identification

Tag: DICOM tag for this attribute.

Action: Action for De-identification

D: replace with a non-zero length value that may be a dummy value and consistent with the VR

Z: replace with a zero length value

X: remove

U: replace with a non-zero length UID that is internally consistent within a set of Instances



## 8. ANNEXES

### 8.1 IOD CONTENTS

#### 8.1.1 Created SOP Instances

Table 8.1-1 specifies the attributes of a MR Image transmitted by the Storage SCU AE.

Table 8.1-2 specifies the attributes of a Secondary Capture Image transmitted by the Storage SCU AE.

Table 8.1-3 specifies the attributes of an Enhanced MR Image transmitted by the Storage SCU AE.

Table 8.1-4 specifies the attributes of a MR Spectroscopy transmitted by the Storage SCU AE.

Table 8.1-5 specifies the attributes of a Grayscale Softcopy Presentation State transmitted by the Storage SCU AE.

The following tables use a number of abbreviations. The abbreviations used in the “Presence of ...” column are:

|             |  |
|-------------|--|
| VNAP        | Value Not Always Present (attribute sent zero length if no value is present) |
| ANAP        | Attribute Not Always Present   |
| ALWAYS      | Always Present   |
| EMPTY       | Attribute is sent without a value  |
| Not Set     | Attribute is Not Present   |
| CONDITIONAL | the module is used under specified condition                                 |

The abbreviations used in the “Source” column:

|        |   |
|--------|---|
| MWL    | the attribute value source is from Modality Worklist                              |
| USER   | the attribute value source is from User input                                     |
| AUTO   | the attribute value is generated automatically                                    |
| MPPS   | the attribute value is the same as that use for Modality Performed Procedure Step |
| CONFIG | the attribute value source is a configurable parameter                            |

## 8.1.1.1 MR Image IOD

**Table 8.1-1  
IOD OF CREATED MR IMAGE SOP INSTANCES**

| <b>IE</b>          | <b>Module</b>           | <b>Reference</b> | <b>Presence of Module</b>                                      |
|--------------------|-------------------------|------------------|--|
| Patient            | Patient                 | Table 8.1-6      | ALWAYS   |
|                    | Patient Identification  | Table 8.1-13     | ALWAYS   |
| Study              | General Study           | Table 8.1-7      | ALWAYS   |
|                    | Patient Study           | Table 8.1-8      | ALWAYS   |
|                    | Patient Demographic     | Table 8.1-14     | ALWAYS   |
|                    | Patient Medical         | Table 8.1-15     | ALWAYS   |
|                    | Visit Admission         | Table 8.1-16     | ALWAYS   |
|                    | Requested Procedure     | Table 8.1-17     | ALWAYS   |
|                    | Imaging Service Request | Table 8.1-18     | ALWAYS   |
| Series             | General Series          | Table 8.1-9      | ALWAYS   |
| Frame of Reference | Frame of Reference      | Table 8.1-10     | ALWAYS   |
| Equipment          | General Equipment       | Table 8.1-11     | ALWAYS   |
| Image              | General Image           | Table 8.1-12     | ALWAYS   |
|                    | Image Plane             | Table 8.1-21     | ALWAYS   |
|                    | Image Pixel             | Table 8.1-22     | ALWAYS   |
|                    | Contrast/Bolus          | Table 8.1-23     | CONDITIONAL<br>(Only if contrast media was used in this image) |
|                    | VOI LUT                 | Table 8.1-24     | ALWAYS   |
|                    | SOP Common              | Table 8.1-25     | ALWAYS   |
|                    | MR Image                | Table 8.1-26     | ALWAYS   |
|                    | Other Application       | Table 8.1-27     | ALWAYS   |

## 8.1.1.2 SC Image IOD

**Table 8.1-2  
IOD OF CREATED SC IMAGE SOP INSTANCES**

| <b>IE</b> | <b>Module</b>           | <b>Reference</b> | <b>Presence of Module</b> |
|-----------|-------------------------|------------------|---------------------------|
| Patient   | Patient                 | Table 8.1-6      | ALWAYS                    |
|           | Patient Identification  | Table 8.1-13     | ALWAYS                    |
| Study     | General Study           | Table 8.1-7      | ALWAYS                    |
|           | Patient Study           | Table 8.1-8      | ALWAYS                    |
|           | Patient Demographic     | Table 8.1-14     | ALWAYS                    |
|           | Patient Medical         | Table 8.1-15     | ALWAYS                    |
|           | Visit Admission         | Table 8.1-16     | ALWAYS                    |
|           | Requested Procedure     | Table 8.1-17     | ALWAYS                    |
|           | Imaging Service Request | Table 8.1-18     | ALWAYS                    |
| Series    | General Series          | Table 8.1-9      | ALWAYS                    |
| Equipment | General Equipment       | Table 8.1-11     | ALWAYS                    |
|           | SC Equipment            | Table 8.1-28     | ALWAYS                    |
| Image     | General Image           | Table 8.1-12     | ALWAYS                    |
|           | Image Pixel             | Table 8.1-29     | ALWAYS                    |
|           | SC Image                | Table 8.1-30     | ALWAYS                    |
|           | VOI LUT                 | Table 8.1-31     | ALWAYS                    |
|           | SOP Common              | Table 8.1-32     | ALWAYS                    |

## 8.1.1.3 Enhanced MR Image IOD

**Table 8.1-3**  
**IOD OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>IE</b>          | <b>Module</b>                 | <b>Reference</b> | <b>Presence of Module</b>  |
|--------------------|-------------------------------|------------------|--|
| Patient            | Patient                       | Table 8.1-6      | ALWAYS   |
|                    | Patient Identification        | Table 8.1-13     | ALWAYS   |
| Study              | General Study                 | Table 8.1-7      | ALWAYS   |
|                    | Patient Study                 | Table 8.1-8      | ALWAYS   |
|                    | Patient Demographic           | Table 8.1-14     | ALWAYS   |
|                    | Patient Medical               | Table 8.1-15     | ALWAYS   |
|                    | Visit Admission               | Table 8.1-16     | ALWAYS   |
|                    | Requested Procedure           | Table 8.1-17     | ALWAYS   |
|                    | Imaging Service Request       | Table 8.1-18     | ALWAYS   |
| Series             | General Series                | Table 8.1-9      | ALWAYS   |
|                    | MR Series                     | Table 8.1-33     | ALWAYS   |
| Frame of Reference | Frame of Reference            | Table 8.1-10     | ALWAYS   |
|                    | Synchronization               | Table 8.1-34     | CONDITIONAL<br>(Only if time synchronization was applied.)                                     |
| Equipment          | General Equipment             | Table 8.1-11     | ALWAYS   |
|                    | Enhanced General Equipment    | Table 8.1-35     | ALWAYS   |
| Image              | Image Pixel                   | Table 8.1-22     | ALWAYS   |
|                    | Enhanced Contrast/Bolus       | Table 8.1-36     | CONDITIONAL<br>(Only if contrast media were applied.)  |
|                    | Multi-frame Functional Groups | Table 8.1-37     | ALWAYS   |
|                    | Multi-frame Dimension         | Table 8.1-66     | ALWAYS   |
|                    | Cardiac Synchronization       | Table 8.1-67     | CONDITIONAL<br>(Only if cardiac synchronization was applied.)                                  |
|                    | Respiratory Synchronization   | Table 8.1-68     | CONDITIONAL<br>(Only if respiratory synchronization was applied.)                              |
|                    | Bulk Motion Synchronization   | Table 8.1-69     | CONDITIONAL<br>(Only if bulk motion synchronization was applied.)                              |
|                    | Acquisition Context           | Table 8.1-70     | ALWAYS   |
|                    | Enhanced MR Image             | Table 8.1-71     | ALWAYS   |
|                    | MR Pulse Sequence             | Table 8.1-72     | CONDITIONAL<br>(Only if Image Type (0008,0008) Value 1 is ORIGINAL. May be present otherwise.) |
|                    | SOP Common                    | Table 8.1-73     | ALWAYS   |

## 8.1.1.4 MR Spectroscopy IOD

**Table 8.1-4**  
**IOD OF CREATED MR SPECTROSCOPY SOP INSTANCES**

| IE                 | Module                         | Reference    | Presence of Module   |
|--------------------|--------------------------------|--------------|--|
| Patient            | Patient                        | Table 8.1-6  | ALWAYS   |
|                    | Patient Identification         | Table 8.1-13 | ALWAYS   |
| Study              | General Study                  | Table 8.1-7  | ALWAYS   |
|                    | Patient Study                  | Table 8.1-8  | ALWAYS   |
|                    | Patient Demographic            | Table 8.1-14 | ALWAYS   |
|                    | Patient Medical                | Table 8.1-15 | ALWAYS   |
|                    | Visit Admission                | Table 8.1-16 | ALWAYS   |
|                    | Requested Procedure            | Table 8.1-17 | ALWAYS   |
|                    | Imaging Service Request        | Table 8.1-18 | ALWAYS   |
| Series             | General Series                 | Table 8.1-9  | ALWAYS   |
|                    | MR Series                      | Table 8.1-33 | ALWAYS   |
| Frame of Reference | Frame of Reference             | Table 8.1-10 | ALWAYS   |
|                    | Synchronization                | Table 8.1-34 | CONDITIONAL<br>(Only if time synchronization was applied.)                                     |
| Equipment          | General Equipment              | Table 8.1-11 |  |
|                    | Enhanced General Equipment     | Table 8.1-35 | ALWAYS   |
| Image              | Enhanced Contrast/Bolus        | Table 8.1-36 | CONDITIONAL<br>(Only if contrast media were applied.)  |
|                    | Multi-frame Functional Groups  | Table 8.1-37 | ALWAYS   |
|                    | Multi-frame Dimension          | Table 8.1-66 | ALWAYS   |
|                    | Cardiac Synchronization        | Table 8.1-67 | CONDITIONAL<br>(Only if cardiac synchronization was applied.)                                  |
|                    | Respiratory Synchronization    | Table 8.1-68 | CONDITIONAL<br>(Only if respiratory synchronization was applied.)                              |
|                    | Bulk Motion Synchronization    | Table 8.1-69 | CONDITIONAL<br>(Only if bulk motion synchronization was applied.)                              |
|                    | Acquisition Context            | Table 8.1-70 | ALWAYS   |
|                    | MR Spectroscopy                | Table 8.1-74 | ALWAYS   |
|                    | MR Spectroscopy Pulse Sequence | Table 8.1-75 | CONDITIONAL<br>(Only if Image Type (0008,0008) Value 1 is ORIGINAL. May be present otherwise.) |
|                    | MR Spectroscopy Data           | Table 8.1-76 | ALWAYS   |
|                    | SOP Common                     | Table 8.1-77 | ALWAYS   |
|                    | Private Application            | Table 8.1-78 | CONDITIONAL<br>(Only if private data are present)  |

### 8.1.1.5 Grayscale Softcopy Presentation State IOD

Table 8.1-5

#### IOD OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES

| IE                 | Module                    | Reference    | Presence of Module  |
|--------------------|---------------------------|--------------|---|
| Patient            | Patient                   | Table 8.1-6  | ALWAYS  |
|                    | Patient Identification    | Table 8.1-13 | ALWAYS  |
| Study              | General Study             | Table 8.1-7  | ALWAYS  |
|                    | Patient Study             | Table 8.1-8  | ALWAYS  |
|                    | Patient Demographic       | Table 8.1-14 | ALWAYS  |
|                    | Patient Medical           | Table 8.1-15 | ALWAYS  |
|                    | Visit Admission           | Table 8.1-16 | ALWAYS  |
|                    | Requested Procedure       | Table 8.1-17 | ALWAYS  |
|                    | Imaging Service Request   | Table 8.1-18 | ALWAYS  |
| Series             | General Series            | Table 8.1-9  | ALWAYS  |
|                    | Presentation Series       | Table 8.1-79 | ALWAYS  |
| Equipment          | General Equipment         | Table 8.1-11 | ALWAYS  |
| Presentation State | Presentation State        | Table 8.1-80 | ALWAYS  |
|                    | Displayed Area            | Table 8.1-81 | ALWAYS  |
|                    | Spatial Transformation    | Table 8.1-82 | CONDITIONAL<br>(Only if Graphic Annotations are to be applied to referenced image(s)) |
|                    | Modality LUT              | Table 8.1-83 | ALWAYS  |
|                    | Softcopy VOI LUT          | Table 8.1-84 | CONDITIONAL<br>(Only if a VOI LUT is to be applied to referenced image(s))            |
|                    | Softcopy Presentation LUT | Table 8.1-85 | ALWAYS  |
|                    | SOP Common                | Table 8.1-86 | ALWAYS  |

## 8.1.1.6 Common Modules

**Table 8.1-6  
PATIENT MODULE OF CREATED SOP INSTANCES**

| Attribute Name                           | Tag         | VR | Value   | Presence of Value | Source       |
|--|-------------|----|---|-------------------|--------------|
| Patient's Name                           | (0010,0010) | PN | From Modality Worklist or user input. Values supplied via Modality Worklist will be entered as received. Maximum 64 characters. | VNAP              | MWL/<br>USER |
| Patient ID                               | (0010,0020) | LO | From Modality Worklist or user input. Maximum 64 characters.  | VNAP              | MWL/<br>USER |
| Issuer of Patient ID                     | (0010,0021) | LO | From Modality Worklist. Maximum 64 characters.  | VNAP              | MWL          |
| Issuer of Patient ID Qualifiers Sequence | (0010,0024) | SQ | From Modality Worklist  | ANAP              | MWL          |
| >Universal Entity ID                     | (0040,0032) | UT | From Modality Worklist  | ANAP              | MWL          |
| >Universal Entity ID Type                | (0040,0033) | CS | From Modality Worklist  | ANAP              | MWL          |
| >Identifier Type Code                    | (0040,0035) | CS | From Modality Worklist  | ANAP              | MWL          |
| Patient's Birth Date                     | (0010,0030) | DA | From Modality Worklist or user input  | VNAP              | MWL/<br>USER |
| Patient's Birth Time                     | (0010,0032) | TM | From Modality Worklist  | VNAP              | MWL          |
| Patient's Sex                            | (0010,0040) | CS | From Modality Worklist or user input  | VNAP              | MWL/<br>USER |
| Other Patient IDs                        | (0010,1000) | LO | From Modality Worklist. Maximum 64 characters.  | ANAP              | MWL          |
| Other Patient Names                      | (0010,1001) | PN | From Modality Worklist. Values supplied via Modality Worklist will be entered as received. Maximum 64 characters.               | ANAP              | MWL          |
| Patient Comments                         | (0010,4000) | LT | From Modality Worklist or user Input. Maximum 1024 characters.  | ANAP              | MWL/<br>USER |

**Table 8.1-7  
GENERAL STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name                      | Tag         | VR | Value  | Presence of Value | Source   |
|-------------------------------------|-------------|----|--|-------------------|----------|
| Study Date                          | (0008,0020) | DA | <yyyymmdd>   | ALWAYS            | AUTO     |
| Study Time                          | (0008,0030) | TM | <hhmmss.frac>  | ALWAYS            | AUTO     |
| Accession Number                    | (0008,0050) | SH | From Modality Worklist or user input. Maximum 16 characters. | VNAP              | MWL/USER |
| Issuer of Accession Number Sequence | (0008,0051) | SQ | From Modality Worklist                                       | ANAP              | MWL      |
| >Local Namespace Entity ID          | (0040,0031) | UT | From Modality Worklist                                       | ANAP              | MWL      |
| >Universal Entity ID                | (0040,0032) | UT | From Modality Worklist                                       | ANAP              | MWL      |
| >Universal Entity ID Type           | (0040,0033) | CS | From Modality Worklist                                       | ANAP              | MWL      |
| Referring Physician's Name          | (0008,0090) | PN | From Modality Worklist.                                      | VNAP              | MWL      |
| Study Description                   | (0008,1030) | LO | From Modality Worklist or user input. Maximum 64 characters. | ANAP              | MWL/USER |
| Physician(s) of Record              | (0008,1048) | PN |  | ANAP              | MWL/AUTO |
| Name of Physician(s) Reading Study  | (0008,1060) | PN |  | ANAP              | MWL/AUTO |
| Referenced Study Sequence           | (0008,1110) | SQ | From Modality Worklist                                       | ANAP              | MWL      |
| > Referenced SOP Class UID          | (0008,1150) | UI |  | ANAP              | MWL      |
| > Referenced SOP Instance UID       | (0008,1155) | UI |  | ANAP              | MWL      |
| Study Instance UID                  | (0020,000D) | UI | From Modality Worklist or generated by device                | ALWAYS            | MWL/AUTO |
| Study ID                            | (0020,0010) | SH | Requested Procedure ID from Worklist or generated by device  | VNAP              | MWL/AUTO |

**Table 8.1-8  
PATIENT STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name   | Tag         | VR | Value  | Presence of Value | Source    |
|------------------|-------------|----|--|-------------------|-----------|
| Patient's Age    | (0010,1010) | AS | Calculated from DoB input on base of actual Date | VNAP              | AUTO/USER |
| Patient's Size   | (0010,1020) | DS | From Modality Worklist or user input             | ALWAYS            | MWL/USER  |
| Patient's Weight | (0010,1030) | DS | From Modality Worklist or user input             | ALWAYS            | MWL/USER  |



**Table 8.1-9  
GENERAL SERIES MODULE OF CREATED SOP INSTANCES**

| Attribute Name                               | Tag         | VR | Value  | Presence of Value | Source       |
|--|-------------|----|--|-------------------|--------------|
| Series Date                                  | (0008,0021) | DA |  | ALWAYS            | AUTO         |
| Series Time                                  | (0008,0031) | TM |  | ALWAYS            | AUTO         |
| Modality                                     | (0008,0060) | CS | "MR"   | ALWAYS            | AUTO         |
| Series Description                           | (0008,103E) | LO | Scan Comment field in sequence queue. Maximum 64 characters.           | ANAP              | USER         |
| Performing Physician's Name                  | (0008,1050) | PN | Radiologist field in Study list. From Modality Worklist or user input. | ANAP              | MWL/<br>USER |
| Operator's Name                              | (0008,1070) | PN | Operator field in Study list. From Modality Worklist or user input.    | ANAP              | MWL/<br>USER |
| Referenced Performed Procedure Step Sequence | (0008,1111) | SQ |  | ANAP              | AUTO         |
| > Referenced SOP Class UID                   | (0008,1150) | UI |  | ANAP              | AUTO         |
| > Referenced SOP Instance UID                | (0008,1155) | UI |  | ANAP              | AUTO         |
| Body Part Examined                           | (0018,0015) | CS | User input. Maximum 16 characters.                                     | ALWAYS            | USER         |
| Protocol Name                                | (0018,1030) | LO | Scan ID field in sequence queue. Maximum 64 characters.                | ANAP              | USER         |
| Patient Position                             | (0018,5100) | CS | Generated by device  | ALWAYS            | AUTO         |
| Series Instance UID                          | (0020,000E) | UI | Generated by device  | ALWAYS            | AUTO         |
| Series Number                                | (0020,0011) | IS | Generated by device  | ALWAYS            | AUTO         |
| Laterality                                   | (0020,0060) | CS |  | ANAP              | AUTO         |
| Smallest Pixel Value in Series               | (0028,0108) | SS |  | ANAP              | AUTO         |
| Largest Pixel Value in Series                | (0028,0109) | SS |  | ANAP              | AUTO         |
| Performed Procedure Step Date                | (0040,0244) | DA |  | ANAP              | AUTO         |
| Performed Procedure Step Time                | (0040,0245) | TM |  | ANAP              | AUTO         |
| Performed Procedure Step ID                  | (0040,0253) | SH | From Modality Worklist or generated by device                          | ANAP              | MWL/<br>AUTO |
| Performed Procedure Step Description         | (0040,0254) | LO | From Modality Worklist   | ANAP              | MWL          |
| Performed Protocol Code Sequence             | (0040,0260) | SQ |  | ANAP              | AUTO         |
| >Code Value                                  | (0008,0100) | SH |  | ANAP              | AUTO         |
| >Coding Scheme Designator                    | (0008,0102) | SH |  | ANAP              | AUTO         |
| >Coding Scheme Version                       | (0008,0103) | SH |  | ANAP              | AUTO         |
| >Code Meaning                                | (0008,0104) | LO |  | ANAP              | AUTO         |
| Request Attributes Sequence                  | (0040,0275) | SQ | From Modality Worklist   | ANAP              | MWL          |

|                                       |             |    |                        |      |     |
|---------------------------------------|-------------|----|------------------------|------|-----|
| >Scheduled Procedure Step Description | (0040,0007) | LO | From Modality Worklist | ANAP | MWL |
| >Scheduled Protocol Code Sequence     | (0040,0008) | SQ | From Modality Worklist | ANAP | MWL |
| >>Code Value                          | (0008,0100) | SH |                        | ANAP | MWL |
| >>Coding Scheme Designator            | (0008,0102) | SH |                        | ANAP | MWL |
| >>Coding Scheme Version               | (0008,0103) | SH |                        | ANAP | MWL |
| >>Code Meaning                        | (0008,0104) | LO |                        | ANAP | MWL |
| >Scheduled Procedure Step ID          | (0040,0009) | SH | From Modality Worklist | ANAP | MWL |
| >Requested Procedure ID               | (0040,1001) | SH | From Modality Worklist | ANAP | MWL |
| >Reason for the Requested Procedure   | (0040,1002) | LO | From Modality Worklist | ANAP | MWL |

**Table 8.1-10  
FRAME OF REFERENCE MODULE OF CREATED SOP INSTANCES**

| Attribute Name               | Tag         | VR | Value | Presence of Value | Source |
|------------------------------|-------------|----|-------|-------------------|--------|
| Frame of Reference UID       | (0020,0052) | UI |       | ALWAYS            | AUTO   |
| Position Reference Indicator | (0020,1040) | LO |       | VNAP              | AUTO   |

**Table 8.1-11  
GENERAL EQUIPMENT MODULE OF CREATED SOP INSTANCES**

| Attribute Name                | Tag         | VR | Value  | Presence of Value | Source       |
|-------------------------------|-------------|----|--|-------------------|--------------|
| Manufacturer                  | (0008,0070) | LO | "TOSHIBA"  | ALWAYS            | AUTO         |
| Institution Name              | (0008,0080) | LO | From Configuration   | ALWAYS            | CONFIG       |
| Institution Address           | (0008,0081) | ST | From Configuration   | ANAP              | CONFIG       |
| Station Name                  | (0008,1010) | SH | From Configuration   | ALWAYS            | CONFIG       |
| Institutional Department Name | (0008,1040) | LO | From Modality Worklist or user input. Maximum 64 characters. | ANAP              | MWL/<br>USER |
| Manufacturer's Model Name     | (0008,1090) | LO |  | ALWAYS            | AUTO         |
| Device Serial Number          | (0018,1000) | LO |  | ANAP              | AUTO         |
| Software Version              | (0018,1020) | LO |  | ALWAYS            | AUTO         |

**Table 8.1-12  
GENERAL IMAGE MODULE OF CREATED SOP INSTANCES**

| Attribute Name               | Tag         | VR | Value  | Presence of Value | Source        |
|------------------------------|-------------|----|--|-------------------|---------------|
| Image Type                   | (0008,0008) | CS | If SOP Class UID(0008,0016) is "1.2.840.10008.5.1.4.1.1.4", see Table 8.1-19.<br>Else if SOP Class UID(0008,0016) is "1.2.840.10008.5.1.4.1.1.7", see Table 8.1-20 | ALWAYS            | AUTO          |
| Acquisition Date             | (0008,0022) | DA | <yyyymmdd>   | ALWAYS            | AUTO          |
| Content Date                 | (0008,0023) | DA | <yyyymmdd>   | ANAP              | AUTO          |
| Acquisition Time             | (0008,0032) | TM | <hhmmss.frac>  | ALWAYS            | AUTO          |
| Content Time                 | (0008,0033) | TM | <hhmmss.frac>  | ANAP              | AUTO          |
| Referenced Image Sequence    | (0008,1140) | SQ |  | ANAP              | AUTO          |
| >Referenced SOP Class UID    | (0008,1150) | UI |  | ANAP              | AUTO          |
| >Referenced SOP Instance UID | (0008,1155) | UI |  | ANAP              | AUTO          |
| Frame Type                   | (0008,9007) | CS |  | ANAP              | AUTO/<br>USER |
| Acquisition Number           | (0020,0012) | IS | Generated by device  | ANAP              | AUTO          |
| Instance Number              | (0020,0013) | IS | Generated by device  | ALWAYS            | AUTO          |
| Patient Orientation          | (0020,0020) | CS |  | ANAP              | AUTO          |
| Image in Acquisition         | (0020,1002) | IS |  | ANAP              | AUTO          |
| Image Comments               | (0020,4000) | LT | Generated by application or from user input. Maximum 44 characters.  | ANAP              | AUTO/<br>USER |

**Table 8.1-13  
PATIENT IDENTIFICATION MODULE ATTRIBUTES**

| Attribute Name                | Tag         | VR | Value | Presence of Value | Source       |
|-------------------------------|-------------|----|-------|-------------------|--------------|
| Patient's Birth Name          | (0010,1005) | PN |       | ANAP              | MWL/<br>AUTO |
| Patient's Mother's Birth Name | (0010,1060) | PN |       | ANAP              | MWL/<br>AUTO |
| Medical Record Locator        | (0010,1090) | LO |       | ANAP              | MWL/<br>AUTO |

**Table 8.1-14  
PATIENT DEMOGRAPHIC MODULE ATTRIBUTES**

| <b>Attribute Name</b>          | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|--------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Patient's Address              | (0010,1040) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Military Rank                  | (0010,1080) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Branch of Service              | (0010,1081) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Country of Residence           | (0010,2150) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Region of Residence            | (0010,2152) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Patient's Telephone Numbers    | (0010,2154) | SH        |              | ANAP                     | MWL/<br>AUTO  |
| Ethnic Group                   | (0010,2160) | SH        |              | ANAP                     | MWL/<br>AUTO  |
| Occupation                     | (0010,2180) | SH        |              | ANAP                     | MWL/<br>AUTO  |
| Patient's Religious Preference | (0010,21F0) | LO        |              | ANAP                     | MWL/<br>AUTO  |

**Table 8.1-15  
PATIENT MEDICAL MODULE ATTRIBUTES**

| <b>Attribute Name</b>      | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|----------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Medical Alerts             | (0010,2000) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Allergies                  | (0010,2110) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Smoking Status             | (0010,21A0) | CS        |              | ANAP                     | MWL/<br>AUTO  |
| Additional Patient History | (0010,21B0) | LT        |              | ANAP                     | MWL/<br>AUTO  |
| Last Menstrual Date        | (0010,21D0) | DA        |              | ANAP                     | MWL/<br>AUTO  |
| Special Needs              | (0038,0050) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Patient State              | (0038,0500) | LO        |              | ANAP                     | MWL/<br>AUTO  |

**Table 8.1-16  
VISIT ADMISSION MODULE ATTRIBUTES**

| <b>Attribute Name</b>             | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|-----------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Admitting<br>DiagnosesDescription | (0008,1080) | LO        |              | ANAP                     | MWL/<br>AUTO  |

**Table 8.1-17  
REQUESTED PROCEDURE MODULE ATTRIBUTES**

| <b>Attribute Name</b>              | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Requested Procedure<br>Description | (0032,1060) | LO        |              | ANAP                     | MWL           |

**Table 8.1-18  
IMAGING SERVICE REQUEST MODULE ATTRIBUTES**

| <b>Attribute Name</b>        | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Requesting Physician         | (0032,1032) | PN        |              | ANAP                     | MWL/<br>AUTO  |
| Requesting Service           | (0032,1033) | LO        |              | ANAP                     | MWL/<br>AUTO  |
| Scheduled StudyStart<br>Date | (0032,1000) | DA        |              | ANAP                     | MWL/<br>AUTO  |
| Scheduled StudyStart<br>Time | (0032,1001) | TM        |              | ANAP                     | MWL/<br>AUTO  |
| Study Comments               | (0032,4000) | LT        |              | ANAP                     | MWL/<br>AUTO  |

**Table 8.1-19  
IMAGE TYPE LIST OF MR IMAGE IOD**

| <b>Image Type</b>                  | <b>Explanation</b>                  |
|------------------------------------|-------------------------------------|
| ORIGINAL\PRIMARY\GDC               | Reconstruction image                |
| ORIGINAL\PRIMARY\OTHER             | Other original image                |
| DERIVED\SECONDARY\DWI_ALIGN        | DWI aligned image                   |
| DERIVED\SECONDARY\ISODWI           | Isotropic image                     |
| DERIVED\SECONDARY\ADC              | ADC image                           |
| DERIVED\SECONDARY\DCI              | Post processed image of DCI         |
| DERIVED\SECONDARY\FA               | Post processed image of Tensor      |
| DERIVED\SECONDARY\RA               |                                     |
| DERIVED\SECONDARY\VR               |                                     |
| DERIVED\SECONDARY\LAMBDA           |                                     |
| DERIVED\SECONDARY\TRACE            |                                     |
| DERIVED\SECONDARY\BOLD_MEAN        | Post processed image of BOLD        |
| DERIVED\SECONDARY\BOLD_CORR        |                                     |
| DERIVED\SECONDARY\BOLD_T           |                                     |
| DERIVED\SECONDARY\BOLD_PERCENT     |                                     |
| DERIVED\SECONDARY\CALC             | Post processed image of Calculation |
| DERIVED\SECONDARY\STITCH AUTO      | Post processed image of Stitching   |
| DERIVED\SECONDARY\MPR              | Post processed image of MPR         |
| DERIVED\SECONDARY\REFINE           | Post processed image of Filter      |
| DERIVED\SECONDARY\REFINE2          |                                     |
| DERIVED\SECONDARY\LSI              |                                     |
| DERIVED\SECONDARY\BEST             |                                     |
| DERIVED\SECONDARY\SINC             |                                     |
| DERIVED\SECONDARY\GOP              |                                     |
| DERIVED\SECONDARY\PROJECTION IMAGE | Auto MIP image                      |
| DERIVED\SECONDARY\PREVIEWMIP3      |                                     |
| DERIVED\SECONDARY\OTHER            | Other derived image                 |

**Table 8.1-20  
IMAGE TYPE LIST OF SC IMAGE IOD**

| <b>Image Type</b>                  | <b>Explanation</b>                  |
|------------------------------------|-------------------------------------|
| DERIVED\SECONDARY\SAVED IMAGE      | Captured image of any post process  |
| DERIVED\SECONDARY\APPDATA_PWI      | Job Data for DCI                    |
| DERIVED\SECONDARY\APPDATA_DTI      | Job Data for Tensor                 |
| DERIVED\SECONDARY\APPDATA_FMRI     | Job Data for BOLD                   |
| DERIVED\SECONDARY\MRS SINGLE VOXEL | Job Data for MRS                    |
| DERIVED\SECONDARY\MRS MULTI VOXEL  |                                     |
| DERIVED\SECONDARY\APPDATA_CARDIAC  | Job Data for Cardiac                |
| DERIVED\SECONDARY\CALC             | Post processed image of Calculation |
| DERIVED\SECONDARY\STITCH MANUAL    | Post processed image of Stitching   |
| DERIVED\SECONDARY\3D               | Captured image of 3D                |
| DERIVED\SECONDARY\MPR_SAVED_IMAGE  | Captured image of MPR               |
| DERIVED\SECONDARY\FUSION           | Captured image of Fusion            |
| DERIVED\SECONDARY\COMPOSITE        | Captured image of Composite         |

### 8.1.1.7 MR Image Modules

**Table 8.1-21  
IMAGE PLANE MODULE OF CREATED MR IMAGE SOP INSTANCES**

| Attribute Name              | Tag         | VR | Value | Presence of Value | Source |
|-----------------------------|-------------|----|-------|-------------------|--------|
| Slice Thickness             | (0018,0050) | DS |       | ALWAYS            | AUTO   |
| Image Position (Patient)    | (0020,0032) | DS |       | ALWAYS            | AUTO   |
| Image Orientation (Patient) | (0020,0037) | DS |       | ALWAYS            | AUTO   |
| Slice Location              | (0020,1041) | DS |       | ANAP              | AUTO   |
| Pixel Spacing               | (0028,0030) | DS |       | ALWAYS            | AUTO   |

**Table 8.1-22  
IMAGE PIXEL MODULE OF CREATED MR IMAGE SOP INSTANCES**

| Attribute Name             | Tag         | VR             | Value         | Presence of Value | Source |
|----------------------------|-------------|----------------|---------------|-------------------|--------|
| Samples per Pixel          | (0028,0002) | US             | 1             | ALWAYS            | AUTO   |
| Photometric Interpretation | (0028,0004) | CS             | "MONOCHROME2" | ALWAYS            | AUTO   |
| Rows                       | (0028,0010) | US             |               | ALWAYS            | AUTO   |
| Columns                    | (0028,0011) | US             |               | ALWAYS            | AUTO   |
| Bits Allocated             | (0028,0100) | US             | 16            | ALWAYS            | AUTO   |
| Bits Stored                | (0028,0101) | US             | 16            | ALWAYS            | AUTO   |
| High Bit                   | (0028,0102) | US             | 15            | ALWAYS            | AUTO   |
| Pixel Representation       | (0028,0103) | US             | 1             | ALWAYS            | AUTO   |
| Smallest Image Pixel Value | (0028,0106) | SS             |               | ANAP              | AUTO   |
| Largest Image Pixel Value  | (0028,0107) | SS             |               | ANAP              | AUTO   |
| Pixel Data                 | (7FE0,0010) | OB<br>or<br>OW |               | ALWAYS            | AUTO   |



**Table 8.1-23  
CONTRAST/BOLUS MODULE OF CREATED MR IMAGE SOP INSTANCES**

| Attribute Name       | Tag         | VR | Value | Presence of Value | Source |
|----------------------|-------------|----|-------|-------------------|--------|
| Contrast/Bolus Agent | (0018,0010) | LO |       | ALWAYS            | USER   |

**Table 8.1-24  
VOI/LUT MODULE OF CREATED MR IMAGE SOP INSTANCES**

| Attribute Name | Tag         | VR | Value | Presence of Value | Source       |
|----------------|-------------|----|-------|-------------------|--------------|
| Window Center  | (0028,1050) | DS |       | ALWAYS            | USER or AUTO |
| Window Width   | (0028,1051) | DS |       | ALWAYS            | USER or AUTO |

**Table 8.1-25  
SOP COMMON MODULE OF CREATED MR IMAGE SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value                       | Presence of Value | Source |
|------------------------|-------------|----|-----------------------------|-------------------|--------|
| Specific Character set | (0008,0005) | CS | Refer to Section 6          | ANAP              | CONFIG |
| Instance Creation Date | (0008,0012) | DA |                             | ANAP              | CONFIG |
| Instance Creation Time | (0008,0013) | TM |                             | ANAP              | CONFIG |
| Instance Creator UID   | (0008,0014) | UI |                             | ALWAYS            | AUTO   |
| SOP Class UID          | (0008,0016) | UI | "1.2.840.10008.5.1.4.1.1.4" | ALWAYS            | AUTO   |
| SOP Instance UID       | (0008,0018) | UI | Generated by device         | ALWAYS            | AUTO   |

**Table 8.1-26**  
**MR IMAGE MODULE OF CREATED MR IMAGE SOP INSTANCES**

| Attribute Name                | Tag         | VR | Value  | Presence of Value | Source |
|-------------------------------|-------------|----|--|-------------------|--------|
| Image Type                    | (0008,0008) | CS | See Table 8.1-19   | ALWAYS            | AUTO   |
| Scanning Sequence             | (0018,0020) | CS | Enumerated Values:<br>"SE":Spin Echo<br>"IR":Inversion Recovery<br>"GR":Gradient Recalled<br>"EP":Echo Planar<br>Multi-valued, but not all combinations are valid (e.g., "SE≠GR", etc.). | ALWAYS            | AUTO   |
| Sequence Variant              | (0018,0021) | CS | "NONE"   | ALWAYS            | AUTO   |
| Scan Options                  | (0018,0022) | CS | Gating Options:<br>"RG","CG","PPG","RCG","RPG"<br>Other Options:<br>"PER","FC","PFF","PFP","SP","FS"   | VNAP              | AUTO   |
| MR Acquisition Type           | (0018,0023) | CS |  | VNAP              | AUTO   |
| Sequence Name                 | (0018,0024) | SH |  | ANAP              | AUTO   |
| Repetition Time               | (0018,0080) | DS |  | ALWAYS            | AUTO   |
| Echo Time                     | (0018,0081) | DS |  | ALWAYS            | AUTO   |
| Inversion Time                | (0018,0082) | DS |  | ANAP              | AUTO   |
| Number of Averages            | (0018,0083) | DS |  | ALWAYS            | AUTO   |
| Imaging Frequency             | (0018,0084) | DS |  | ANAP              | AUTO   |
| Imaged Nucleus                | (0018,0085) | SH |  | ANAP              | AUTO   |
| Echo Number(s)                | (0018,0086) | IS |  | ALWAYS            | AUTO   |
| Magnetic Field Strength       | (0018,0087) | DS |  | ALWAYS            | AUTO   |
| Spacing Between Slices        | (0018,0088) | DS |  | ANAP              | AUTO   |
| Number of Phase Encoding Step | (0018,0089) | IS |  | ANAP              | AUTO   |
| Echo Train Length             | (0018,0091) | IS |  | VNAP              | AUTO   |
| Percent Phase Field of View   | (0018,0094) | DS |  | ANAP              | AUTO   |
| Pixel Bandwidth               | (0018,0095) | DS |  | ANAP              | AUTO   |
| Trigger Time                  | (0018,1060) | DS | CONDITIONAL<br>(Only if Scan Option (0018,0022) include heart gating (RG,CG,PPG,RCG,RPG).)   | ANAP              | AUTO   |
| Nominal Interval              | (0018,1062) | IS |  | ANAP              | AUTO   |
| Beat Rejection Flag           | (0018,1081) | CS |  | ANAP              | AUTO   |
| Low R-R Value                 | (0018,1082) | IS |  | ANAP              | AUTO   |
| High R-R Value                | (0018,1083) | IS |  | ANAP              | AUTO   |
| Heart Rate                    | (0018,1088) | IS |  | ANAP              | AUTO   |
| Cardiac Number of Images      | (0018,1090) | IS |  | ANAP              | AUTO   |
| Trigger Window                | (0018,1094) | IS |  | ANAP              | AUTO   |
| Reconstruction Diameter       | (0018,1100) | DS |  | ANAP              | AUTO   |

|                              |             |    |  |      |      |
|------------------------------|-------------|----|--|------|------|
| Receive Coil Name            | (0018,1250) | SH |  | ANAP | AUTO |
| Transmit Coil Name           | (0018,1251) | SH |  | ANAP | AUTO |
| Acquisition Matrix           | (0018,1310) | US |  | ANAP | AUTO |
| Phase Encoding Direction     | (0018,1312) | CS |  | ANAP | AUTO |
| Flip Angle                   | (0018,1314) | DS |  | ANAP | AUTO |
| SAR                          | (0018,1316) | DS |  | ANAP | AUTO |
| dB/dt                        | (0018,1318) | DS |  | ANAP | AUTO |
| Temporal Position Identifier | (0020,0100) | IS |  | ANAP | AUTO |
| Number of Temporal Position  | (0020,0105) | IS |  | ANAP | AUTO |
| Temporal Resolution          | (0020,0110) | DS |  | ANAP | AUTO |

Table 8.1-27

## OTHER APPLICATION MODULE OF CREATED MR IMAGE SOP INSTANCES

| Attribute Name                                   | Tag         | VR | Value                | Presence of Value | Source |
|--|-------------|----|----------------------|-------------------|--------|
| Acquisition Duration                             | (0018,9073) | FD |                      | ALWAYS            | AUTO   |
| Parallel Acquisition                             | (0018,9077) | CS |                      | ANAP              | AUTO   |
| Parallel Acquisition Technique                   | (0018,9078) | CS |                      | ANAP              | AUTO   |
| Private Creator Code                             | (0029,00xx) | LO | "TOSHIBA_MEC_MR3"    | VNAP              | AUTO   |
| Other Private Data                               | (0029,xx01) | SQ |                      | VNAP              | AUTO   |
| Private Creator Code                             | (700D,00xx) | LO | "TOSHIBA_MEC_MR3"    | ANAP              | AUTO   |
| Scale Factor                                     | (700D,xx00) | DS |                      | ANAP              | AUTO   |
| FOV  | (700D,xx05) | DS |                      | ANAP              | AUTO   |
| Receiver Gain Correction Check Flag              | (700D,xx0C) | CS |                      | ANAP              | AUTO   |
| Identification Flag of 3D GDC                    | (700D,xx20) | SH |                      | ANAP              | AUTO   |
| Private Creator Code                             | (700D,00yy) | LO | "TOSHIBA_MEC_MR3^10" | ANAP              | AUTO   |
| 2nd Flip Angle [degree]                          | (700D,yy10) | DS |                      | ANAP              | AUTO   |
| Acquisition Inner Matrix                         | (700D,yy11) | US |                      | ANAP              | AUTO   |
| MP2RAGE Flag                                     | (700D,yy12) | US |                      | ANAP              | AUTO   |
| Inversion efficiency of inversion recovery pulse | (700D,yy13) | FL |                      | ANAP              | AUTO   |
| Number of dummy shot                             | (700D,yy14) | SL |                      | ANAP              | AUTO   |
| FFE total repetition time[s]                     | (700D,yy15) | FL |                      | ANAP              | AUTO   |
| PAS Name   | (700D,yy16) | LO |                      | ANAP              | AUTO   |
| Intended Processing                              | (700D,yy17) | LT |                      | ANAP              | AUTO   |
| Scanned Orientation IDs                          | (700D,yy18) | SS |                      | ANAP              | AUTO   |
| PAS Reproduct Information                        | (700D,yy19) | OB |                      | ALWAYS            | AUTO   |

## 8.1.1.8 SC Image Modules

**Table 8.1-28**  
**SC EQUIPMENT MODULE OF CREATED SC IMAGE SOP INSTANCES**

| Attribute Name                              | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| Modality                                    | (0008,0060) | CS | "MR"  | ALWAYS            | AUTO   |
| Conversion Type                             | (0008,0064) | CS | "WSD" | ALWAYS            | AUTO   |
| Secondary Capture Device ID                 | (0018,1010) | LO |       | ANAP              | AUTO   |
| Secondary Capture Device Manufacture        | (0018,1016) | LO |       | ANAP              | AUTO   |
| Secondary Capture Manufacturer's Model Name | (0018,1018) | LO |       | ANAP              | AUTO   |
| Secondary Capture Device Software Version   | (0018,1019) | LO |       | ANAP              | AUTO   |

**Table 8.1-29**  
**IMAGE PIXEL MODULE OF CREATED SC IMAGE SOP INSTANCES**

| Attribute Name             | Tag         | VR       | Value                  | Presence of Value | Source |
|----------------------------|-------------|----------|------------------------|-------------------|--------|
| Samples per Pixel          | (0028,0002) | US       | 1 or 3                 | ALWAYS            | AUTO   |
| Photometric Interpretation | (0028,0004) | CS       | "MONOCHROME2" or "RGB" | ALWAYS            | AUTO   |
| Planar Configuration       | (0028,0006) | US       | 0                      | ANAP              | AUTO   |
| Rows                       | (0028,0010) | US       |                        | ALWAYS            | AUTO   |
| Columns                    | (0028,0011) | US       |                        | ALWAYS            | AUTO   |
| Bits Allocated             | (0028,0100) | US       | 16 or 8                | ALWAYS            | AUTO   |
| Bits Stored                | (0028,0101) | US       | 16 or 8                | ALWAYS            | AUTO   |
| High Bit                   | (0028,0102) | US       | 15 or 7                | ALWAYS            | AUTO   |
| Pixel Representation       | (0028,0103) | US       | 1                      | ALWAYS            | AUTO   |
| Pixel Data                 | (7FE0,0010) | OB or OW |                        | ALWAYS            | AUTO   |

**Table 8.1-30**  
**SC IMAGE MODULE OF CREATED SC IMAGE SOP INSTANCES**

| Attribute Name            | Tag         | VR | Value | Presence of Value | Source |
|---------------------------|-------------|----|-------|-------------------|--------|
| Date of Secondary Capture | (0018,1012) | DA |       | ANAP              | AUTO   |
| Date of Secondary Capture | (0018,1014) | TM |       | ANAP              | AUTO   |

**Table 8.1-31  
VOI/LUT MODULE OF CREATED SC IMAGE SOP INSTANCES**

| <b>Attribute Name</b> | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|-----------------------|-------------|-----------|--------------|--------------------------|---------------|
| Window Center         | (0028,1050) | DS        |              | ANAP                     | USER or AUTO  |
| Window Width          | (0028,1051) | DS        |              | ANAP                     | USER or AUTO  |

**Table 8.1-32  
SOP COMMON MODULE OF CREATED SC IMAGE SOP INSTANCES**

| <b>Attribute Name</b>  | <b>Tag</b>  | <b>VR</b> | <b>Value</b>                | <b>Presence of Value</b> | <b>Source</b> |
|------------------------|-------------|-----------|-----------------------------|--------------------------|---------------|
| Specific Character set | (0008,0005) | CS        | Refer to Section 6          | ANAP                     | CONFIG        |
| SOP Class UID          | (0008,0016) | UI        | "1.2.840.10008.5.1.4.1.1.7" | ALWAYS                   | AUTO          |
| SOP Instance UID       | (0008,0018) | UI        | Generated by device         | ALWAYS                   | AUTO          |

### 8.1.1.9 Enhanced MR Image Modules

**Table 8.1-33**  
**MR SERIES MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                               | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| Modality                                     | (0008,0060) | CS | "MR"  | ALWAYS            | AUTO   |
| Referenced Performed Procedure Step Sequence | (0008,1111) | SQ |       | ANAP              | AUTO   |
| >Referenced SOP Class UID                    | (0008,1150) | UI |       | ANAP              | AUTO   |
| >Referenced SOP Instance UID                 | (0008,1155) | UI |       | ANAP              | AUTO   |

**Table 8.1-34**  
**SYNCHRONIZATION MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                         | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| Trigger Source or Type                 | (0018,1061) | LO |       | Not Set           |        |
| Synchronization Trigger                | (0018,106A) | CS |       | ALWAYS            | AUTO   |
| Synchronization Channel                | (0018,106C) | US |       | Not Set           |        |
| Acquisition Time Synchronized          | (0018,1800) | CS |       | ALWAYS            | AUTO   |
| Time Source                            | (0018,1801) | SH |       | Not Set           |        |
| Time Distribution Protocol             | (0018,1802) | CS |       | Not Set           |        |
| NTP Source Address                     | (0018,1803) | LO |       | Not Set           |        |
| Synchronization Frame of Reference UID | (0020,0200) | UI |       | ALWAYS            | AUTO   |

**Table 8.1-35**  
**ENHANCED GENERAL EQUIPMENT MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name            | Tag         | VR | Value     | Presence of Value | Source |
|---------------------------|-------------|----|-----------|-------------------|--------|
| Manufacturer              | (0008,0070) | LO | "TOSHIBA" | ALWAYS            | AUTO   |
| Institution Name          | (0008,0080) | LO |           | ANAP              | CONFIG |
| Institution Address       | (0008,0081) | ST |           | ANAP              | CONFIG |
| Station Name              | (0008,1010) | SH |           | ANAP              | CONFIG |
| Manufacturer's Model Name | (0008,1090) | LO |           | ALWAYS            | AUTO   |
| Device Serial Number      | (0018,1000) | LO |           | ALWAYS            | AUTO   |
| Software Version          | (0018,1020) | LO |           | ALWAYS            | AUTO   |

Table 8.1-36

**ENHANCED CONTRAST/BOLUS MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                                | Tag         | VR | Value    | Presence of Value | Source |
|---|-------------|----|----------|-------------------|--------|
| Contrast/Bolus Agent Sequence                 | (0018,0012) | SQ |          | ANAP              | AUTO   |
| >Code Value                                   | (0008,0100) | SH |          | ANAP              | AUTO   |
| >Coding Scheme Designator                     | (0008,0102) | SH |          | ANAP              | AUTO   |
| >Coding Scheme Version                        | (0008,0103) | SH |          | Not Set           |        |
| >Code Meaning                                 | (0008,0104) | LO |          | ANAP              | AUTO   |
| >Contrast/Bolus Administration Route Sequence | (0018,0014) | SQ |          | ANAP              | AUTO   |
| >>Code Value                                  | (0008,0100) | SH |          | ANAP              | AUTO   |
| >>Coding Scheme Designator                    | (0008,0102) | SH |          | ANAP              | AUTO   |
| >>Coding Scheme Version                       | (0008,0103) | SH |          | Not Set           |        |
| >>Code Meaning                                | (0008,0104) | LO |          | ANAP              | AUTO   |
| >Contrast/Bolus Volume                        | (0018,1041) | DS | Length=0 | EMPTY             | AUTO   |
| >Contrast/Bolus Ingredient Concentration      | (0018,1049) | DS | Length=0 | EMPTY             | AUTO   |
| >Contrast/Bolus Agent Number                  | (0018,9337) | US |          | ANAP              | AUTO   |
| >Contrast/Bolus Ingredient Code Sequence      | (0018,9338) | SQ | Length=0 | EMPTY             | AUTO   |
| >Contrast Administration Profile Sequence     | (0018,9340) | SQ |          | Not Set           |        |
| >>Contrast/Bolus Volume                       | (0018,1041) | DS |          | Not Set           |        |
| >>Contrast/Bolus Start Time                   | (0018,1042) | TM |          | Not Set           |        |
| >>Contrast/Bolus Stop Time                    | (0018,1043) | TM |          | Not Set           |        |
| >>Contrast Flow Rate                          | (0018,1046) | DS |          | Not Set           |        |
| >>Contrast Flow Duration                      | (0018,1047) | DS |          | Not Set           |        |
| >Contrast/Bolus Ingredient Opaque             | (0018,9425) | CS |          | Not Set           |        |

**Table 8.1-37**  
**MULTI-FRAME FUNCTIONAL GROUP MODULE OF CREATED ENHANCED MR IMAGE SOP**  
**INSTANCES/MR SPECTROSCOPY SOP INSTANCES**

| Attribute Name                       | Tag         | VR | Value            | Presence of Value | Source |
|--------------------------------------|-------------|----|------------------|-------------------|--------|
| Shared Functional Groups Sequence    | (5200,9229) | SQ | See Table 8.1-38 | ALWAYS            | AUTO   |
| Per-frame Functional Groups Sequence | (5200,9230) | SQ | See Table 8.1-38 | ALWAYS            | AUTO   |
| Content Date                         | (0008,0023) | DA |                  | ALWAYS            | AUTO   |
| Content Time                         | (0008,0033) | TM |                  | ALWAYS            | AUTO   |
| Instance Number                      | (0020,0013) | IS |                  | ALWAYS            | AUTO   |
| Number of Frames                     | (0028,0008) | IS |                  | ALWAYS            | AUTO   |

**Table 8.1-38**  
**MULTI-FRAME FUNCTIONAL GROUP MACROS**

| Functional Group Macro             | Section      | Presence of Module For Enhanced MR Image | Presence of Module For MR Spectroscopy |
|------------------------------------|--------------|--|--|
| Pixel Measures                     | Table 8.1-39 | M  | M                                      |
| Frame Content                      | Table 8.1-40 | M  | M                                      |
| Plan Position                      | Table 8.1-41 | M  | M                                      |
| Plane Orientation                  | Table 8.1-42 | M  | M                                      |
| Referenced Image                   | Table 8.1-43 | C  | C                                      |
| Derivation Image                   | Table 8.1-44 | C  | C                                      |
| Cardiac Synchronization            | Table 8.1-45 | C  | C                                      |
| Frame Anatomy                      | Table 8.1-46 | M  | M                                      |
| Pixel value Transformation         | Table 8.1-47 | M  | M                                      |
| Frame VOI LUT                      | Table 8.1-48 | U  | Not Used                               |
| Real World Value Mapping           | Table 8.1-49 | Not Used                                 | Not Used                               |
| Contrast/Bolus Usage               | Table 8.1-50 | C  | C                                      |
| Respiratory Synchronization        |              | Not Used                                 | Not Used                               |
| MR Image Frame Type                | Table 8.1-51 | M  | Not Used                               |
| MR Spectroscopy Frame Type Macro   | Table 8.1-65 | Not Used                                 | M                                      |
| MR Timing and Related Parameters   | Table 8.1-52 | C  | C                                      |
| MR FOV/Geometry                    | Table 8.1-53 | C  | Not Used                               |
| MR Spectroscopy FOV/Geometry Macro | Table 8.1-64 | Not Used                                 | C                                      |
| MR Echo                            | Table 8.1-54 | C  | C                                      |
| MR Modifier                        | Table 8.1-55 | C  | C                                      |
| MR Imaging Modifier                | Table 8.1-56 | C  | Not Used                               |
| MR Receive Coil                    | Table 8.1-57 | C  | C                                      |
| MR Transmit Coil                   | Table 8.1-58 | C  | C                                      |
| MR Diffusion                       | Table 8.1-59 | C  | Not Used                               |
| MR Averages                        | Table 8.1-60 | C  | C                                      |



|                       |              |          |          |
|-----------------------|--------------|----------|----------|
| MR Spatial Saturation | Table 8.1-61 | C        | C        |
| MR Metabolite Map     | Table 8.1-62 | Not Used | Not Used |
| MR Velocity Encoding  | Table 8.1-63 | C        | Not Used |

**Table 8.1-39  
PIXEL MEASURES MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name          | Tag         | VR | Value | Presence of Value | Source |
|-------------------------|-------------|----|-------|-------------------|--------|
| Pixel Measures Sequence | (0028,9110) | SQ |       | ALWAYS            | AUTO   |
| >Slice Thickness        | (0018,0050) | DS |       | ALWAYS            | AUTO   |
| >Pixel Spacing          | (0028,0030) | DS |       | ALWAYS            | AUTO   |

**Table 8.1-40  
FRAME CONTENT MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name              | Tag         | VR | Value | Presence of Value | Source |
|-----------------------------|-------------|----|-------|-------------------|--------|
| Frame Content Sequence      | (0020,9111) | SQ |       | ALWAYS            | AUTO   |
| >Frame Acquisition Datetime | (0018,9074) | DT |       | VNAP              | AUTO   |
| >Frame Reference Datetime   | (0018,9151) | DT |       | VNAP              | AUTO   |
| >Respiratory Cycle Position | (0018,9214) | CS |       | Not Set           |        |
| >Frame Acquisition Duration | (0018,9220) | FD |       | VNAP              | AUTO   |
| >Cardiac Cycle Position     | (0018,9236) | CS |       | Not Set           |        |
| >Stack ID                   | (0020,9056) | SH |       | ANAP              | AUTO   |
| >In-Stack Position Number   | (0020,9057) | UL |       | ANAP              | AUTO   |
| >Temporal Position Index    | (0020,9128) | UL |       | ANAP              | AUTO   |
| >Frame Acquisition Number   | (0020,9156) | US |       | Not Set           |        |
| >Dimension Index Values     | (0020,9157) | UL |       | ANAP              | AUTO   |
| >Frame Comments             | (0020,9158) | LT |       | ANAP              | AUTO   |

**Table 8.1-41  
PIXEL POSITION MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name             | Tag         | VR | Value | Presence of Value | Source |
|----------------------------|-------------|----|-------|-------------------|--------|
| Plane Position Sequence    | (0020,9113) | SQ |       | ALWAYS            | AUTO   |
| > Image Position (Patient) | (0020,0032) | DS |       | ALWAYS            | AUTO   |

**Table 8.1-42  
PLANE ORIENTATION MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>         | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|-------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Plane Orientation Sequence    | (0020,9116) | SQ        |              | ALWAYS                   | AUTO          |
| > Image Orientation (Patient) | (0020,0037) | DS        |              | ALWAYS                   | AUTO          |

**Table 8.1-43  
REFERENCED IMAGE MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>               | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|-------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Referenced Image Sequence           | (0008,1140) | SQ        |              | ANAP                     | AUTO          |
| >Referenced SOP Class UID           | (0008,1150) | UI        |              | ANAP                     | AUTO          |
| >Referenced SOP Instance UID        | (0008,1155) | UI        |              | ANAP                     | AUTO          |
| >Referenced Frame Number            | (0008,1160) | IS        |              | ANAP                     | AUTO          |
| >Purpose of Reference Code Sequence | (0040,A170) | SQ        |              | ANAP                     | AUTO          |
| >>Code Value                        | (0008,0100) | SH        |              | Not Set                  |               |
| >>Coding Scheme Designator          | (0008,0102) | SH        |              | Not Set                  |               |
| >>Coding Scheme Version             | (0008,0103) | SH        |              | Not Set                  |               |
| >>Code Meaning                      | (0008,0104) | LO        |              | Not Set                  |               |

**Table 8.1-44  
DERIVATION IMAGE MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                       | Tag         | VR | Value | Presence of Value | Source |
|--------------------------------------|-------------|----|-------|-------------------|--------|
| Derivation Image Sequence            | (0008,9124) | SQ |       | EMPTY             | AUTO   |
| >Derivation Description              | (0008,2111) | ST |       | ANAP              | AUTO   |
| >Source Image Sequence               | (0008,2112) | SQ |       | ANAP              | AUTO   |
| >>Referenced SOP Class UID           | (0008,1150) | UI |       | ANAP              | AUTO   |
| >>Referenced SOP Instance UID        | (0008,1155) | UI |       | ANAP              | AUTO   |
| >>Referenced Frame Number            | (0008,1160) | IS |       | ANAP              | AUTO   |
| >>Purpose of Reference Code Sequence | (0040,A170) | SQ |       | ANAP              | AUTO   |
| >>>Code Value                        | (0008,0100) | SH |       | ANAP              |        |
| >>>Coding Scheme Designator          | (0008,0102) | SH |       | ANAP              |        |
| >>>Coding Scheme Version             | (0008,0103) | SH |       | ANAP              |        |
| >>>Code Meaning                      | (0008,0104) | LO |       | ANAP              |        |
| >Derivation Code Sequence            | (0008,9215) | SQ |       | ANAP              | AUTO   |
| >>Code Value                         | (0008,0100) | SH |       | ANAP              |        |
| >>Coding Scheme Designator           | (0008,0102) | SH |       | ANAP              |        |
| >>Coding Scheme Version              | (0008,0103) | SH |       | ANAP              |        |
| >>Code Meaning                       | (0008,0104) | LO |       | ANAP              |        |

**Table 8.1-45  
CARDIAC SYNCHRONIZATION MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                   | Tag         | VR | Value | Presence of Value | Source |
|----------------------------------|-------------|----|-------|-------------------|--------|
| Cardiac Synchronization Sequence | (0018,9118) | SQ |       | VNAP              | AUTO   |
| > Trigger Delay Time             | (0020,9153) | FD |       | VNAP              | AUTO   |

**Table 8.1-46  
FRAME ANATOMY MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name            | Tag         | VR | Value | Presence of Value | Source |
|---------------------------|-------------|----|-------|-------------------|--------|
| Frame Anatomy Sequence    | (0020,9071) | SQ |       | EMPTY             | AUTO   |
| >Anatomic Region Sequence | (0008,2218) | SQ |       | EMPTY             | AUTO   |
| >Frame Laterality         | (0020,9072) | CS |       | ANAP              | AUTO   |

Table 8.1-47

**PIXEL VALUE TRANSFORMATION MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                      | Tag         | VR | Value | Presence of Value | Source |
|-------------------------------------|-------------|----|-------|-------------------|--------|
| Pixel Value Transformation Sequence | (0028,9145) | SQ |       | ALWAYS            | AUTO   |
| >Rescale Intercept                  | (0028,1052) | DS | 0     | ALWAYS            | AUTO   |
| >Rescale Slope                      | (0028,1053) | DS | 1     | ALWAYS            | AUTO   |
| >Rescale Type                       | (0028,1054) | LO | "US"  | ALWAYS            | AUTO   |

Table 8.1-48

**FRAME VOI LUT MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value | Presence of Value | Source |
|------------------------|-------------|----|-------|-------------------|--------|
| Frame VOI LUT Sequence | (0028,9132) | SQ |       | VNAP              | AUTO   |
| > Window Center        | (0028,1050) | DS |       | VNAP              | AUTO   |
| > Window Width         | (0028,1051) | DS |       | VNAP              | AUTO   |

Table 8.1-49

**REAL WORLD VALUE MAPPING MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                       | Tag         | VR | Value | Presence of Value | Source |
|--------------------------------------|-------------|----|-------|-------------------|--------|
| Real World Value Mapping Sequence    | (0040,9096) | SQ |       | Not Set           |        |
| >LUT Explanation                     | (0028,3003) | LO |       | Not Set           |        |
| >Measurement Units Code Sequence     | (0040,08EA) | SQ |       | Not Set           |        |
| >LUT Label                           | (0040,9210) | SH |       | Not Set           |        |
| >Real World Value Last Value Mapped  | (0040,9211) | US |       | Not Set           |        |
| >Real World Value LUT Data           | (0040,9212) | FD |       | Not Set           |        |
| >Real World Value First Value Mapped | (0040,9216) | US |       | Not Set           |        |
| >Real World Value Intercept          | (0040,9224) | FD |       | Not Set           |        |
| >Real World Value Slope              | (0040,9225) | FD |       | Not Set           |        |

**Table 8.1-50  
CONTRAST/BOLUS USAGE MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>              | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Contrast/Bolus Usage Sequence      | (0018,9341) | SQ        |              | ANAP                     | AUTO          |
| >Contrast/Bolus Agent Number       | (0018,9337) | US        | "1"          | ANAP                     | AUTO          |
| >Contrast/Bolus Agent Administered | (0018,9342) | CS        |              | ANAP                     | AUTO          |
| >Contrast/Bolus Agent Detected     | (0018,9343) | CS        |              | VNAP                     | AUTO          |
| >Contrast/Bolus Agent Phase        | (0018,9344) | CS        |              | Not Set                  |               |

**Table 8.1-51  
MR IMAGE FRAME TYPE MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>               | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|-------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| MR Image Frame Type Sequence        | (0018,9226) | SQ        |              | ALWAYS                   | AUTO          |
| >Frame Type                         | (0008,9007) | CS        |              | ALWAYS                   | AUTO          |
| >Pixel Presentation                 | (0008,9205) | CS        | "MONOCHROME" | ALWAYS                   | AUTO          |
| >Volumetric Properties              | (0008,9206) | CS        |              | ALWAYS                   | AUTO          |
| >Volume Based Calculation Technique | (0008,9207) | CS        | "NONE"       | ALWAYS                   | AUTO          |
| >Complex Image Component            | (0008,9208) | CS        |              | ALWAYS                   | AUTO          |
| >Acquisition Contrast               | (0008,9209) | CS        |              | ALWAYS                   | AUTO          |

**Table 8.1-52**  
**MR TIMING AND RELATED PARAMETERS MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                            | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| MR Timing and Related Parameters Sequence | (0018,9112) | SQ |       | VNAP              | AUTO   |
| >Repetition Time                          | (0018,0080) | DS |       | VNAP              | AUTO   |
| >Echo Train Length                        | (0018,0091) | IS |       | VNAP              | AUTO   |
| >Flip Angle                               | (0018,1314) | DS |       | VNAP              | AUTO   |
| >Operating Mode Sequence                  | (0018,9176) | SQ |       | VNAP              | AUTO   |
| >>Operating Mode Type                     | (0018,9177) | CS |       | ALWAYS            | AUTO   |
| >>Operating Mode                          | (0018,9178) | CS |       | ALWAYS            | AUTO   |
| >Gradient Output Type                     | (0018,9180) | CS |       | VNAP              | AUTO   |
| >Gradient Output                          | (0018,9182) | FD |       | VNAP              | AUTO   |
| >Specific Absorption Rate Sequence        | (0018,9239) | SQ |       | VNAP              | AUTO   |
| >>Specific Absorption Rate Definition     | (0018,9179) | CS |       | ALWAYS            | AUTO   |
| >>Specific Absorption Rate Value          | (0018,9181) | FD |       | ALWAYS            | AUTO   |
| >RF Echo Train Length                     | (0018,9240) | US |       | VNAP              | AUTO   |
| >Gradient Echo Train Length               | (0018,9241) | US |       | VNAP              | AUTO   |

**Table 8.1-53**  
**MR FOV/GEOMETRY MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                                    | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| MR FOV/Geometry Sequence                          | (0018,9125) | SQ |       | VNAP              | AUTO   |
| >Percent Sampling                                 | (0018,0093) | DS |       | VNAP              | AUTO   |
| >Percent Phase Field of View                      | (0018,0094) | DS |       | VNAP              | AUTO   |
| >In-plane Phase Encoding Direction                | (0018,1312) | CS |       | VNAP              | AUTO   |
| >MR Acquisition Frequency Encoding Steps          | (0018,9058) | US |       | VNAP              | AUTO   |
| >MR Acquisition Phase Encoding Steps in-plane     | (0018,9231) | US |       | VNAP              | AUTO   |
| >MR Acquisition Phase Encoding Steps out-of-plane | (0018,9232) | US |       | VNAP              | AUTO   |

**Table 8.1-54**  
**MR ECHO MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name       | Tag         | VR | Value | Presence of Value | Source |
|----------------------|-------------|----|-------|-------------------|--------|
| MR Echo Sequence     | (0018,9114) | SQ |       | VNAP              | AUTO   |
| >Effective Echo Time | (0018,9082) | FD |       | VNAP              | AUTO   |

**Table 8.1-55**  
**MR MODIFIER MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                             | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| MR Modifier Sequence                       | (0018,9115) | SQ |       | VNAP              | AUTO   |
| >Inversion Recovery                        | (0018,9009) | CS |       | VNAP              | AUTO   |
| >Flow Compensation                         | (0018,9010) | CS |       | VNAP              | AUTO   |
| >Spoiling                                  | (0018,9016) | CS |       | VNAP              | AUTO   |
| >T2 Preparation                            | (0018,9021) | CS |       | VNAP              | AUTO   |
| >Spectrally Selected Excitation            | (0018,9026) | CS |       | VNAP              | AUTO   |
| >Spatial Pre-saturation                    | (0018,9027) | CS |       | VNAP              | AUTO   |
| >Partial Fourier Direction                 | (0018,9036) | CS |       | VNAP              | AUTO   |
| >Parallel Reduction Factor In-plane        | (0018,9069) | FD |       | VNAP              | AUTO   |
| >Parallel Acquisition                      | (0018,9077) | CS |       | VNAP              | AUTO   |
| >Parallel Acquisition Technique            | (0018,9078) | CS |       | VNAP              | AUTO   |
| >Inversion Times                           | (0018,9079) | FD |       | VNAP              | AUTO   |
| >Partial Fourier                           | (0018,9081) | CS |       | VNAP              | AUTO   |
| >Parallel Reduction Factor out-of plane    | (0018,9155) | FD |       | VNAP              | AUTO   |
| >Parallel Reduction Factor Second In-plane | (0018,9168) | FD |       | Not Set           |        |
| >Flow Compensation Direction               | (0018,9183) | CS |       | VNAP              | AUTO   |



**Table 8.1-56**  
**MR IMAGING MODIFIER MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                | Tag         | VR | Value  | Presence of Value | Source |
|-------------------------------|-------------|----|--------|-------------------|--------|
| MR Imaging Modifier Sequence  | (0018,9006) | SQ |        | VNAP              | AUTO   |
| >Tag Angle First Axis         | (0018,9019) | FD |        | ANAP              | AUTO   |
| >Magnetization Transfer       | (0018,9020) | CS | "NONE" | VNAP              | AUTO   |
| >Blood Signal Nulling         | (0018,9022) | CS | "NO"   | VNAP              | AUTO   |
| >Tagging                      | (0018,9028) | CS | "NONE" | VNAP              | AUTO   |
| >Tag Spacing First Dimension  | (0018,9030) | FD |        | ANAP              | AUTO   |
| >Tag Thickness                | (0018,9035) | FD |        | ANAP              | AUTO   |
| > Pixel Bandwidth             | (0018,0095) | DS |        | VNAP              | AUTO   |
| >Transmitter Frequency        | (0018,9098) | FD |        | VNAP              | AUTO   |
| >Tagging Delay                | (0018,9184) | FD |        | ANAP              | AUTO   |
| >Tag Spacing Second Dimension | (0018,9218) | FD |        | ANAP              | AUTO   |
| >Tag Angle Second Axis        | (0018,9219) | SS |        | ANAP              | AUTO   |

**Table 8.1-57**  
**MR RECEIVE COIL MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                  | Tag         | VR | Value     | Presence of Value | Source |
|---------------------------------|-------------|----|-----------|-------------------|--------|
| MR Receive Coil Sequence        | (0018,9042) | SQ |           | VNAP              | AUTO   |
| >Receive Coil Name              | (0018,1250) | SH |           | VNAP              | AUTO   |
| >Receive Coil Manufacturer Name | (0018,9041) | LO | "TOSHIBA" | ALWAYS            | AUTO   |
| >Receive Coil Type              | (0018,9043) | CS |           | VNAP              | AUTO   |
| >Quadrature Receive Coil        | (0018,9044) | CS |           | VNAP              | AUTO   |

**Table 8.1-58**  
**MR TRANSMIT COIL MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                   | Tag         | VR | Value     | Presence of Value | Source |
|----------------------------------|-------------|----|-----------|-------------------|--------|
| MR Transmit Coil Sequence        | (0018,9049) | SQ |           | VNAP              | AUTO   |
| >Transmit Coil Name              | (0018,1251) | SH |           | VNAP              | AUTO   |
| >Transmit Coil Manufacturer Name | (0018,9050) | LO | "TOSHIBA" | ALWAYS            | AUTO   |
| >Transmit Coil Type              | (0018,9051) | CS |           | VNAP              | AUTO   |

**Table 8.1-59**  
**MR DIFFUSION MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                         | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| MR Diffusion Sequence                  | (0018,9117) | SQ |       | VNAP              | AUTO   |
| >Diffusion Directionality              | (0018,9075) | CS |       | VNAP              | AUTO   |
| >Diffusion Gradient Direction Sequence | (0018,9076) | SQ |       | ANAP              | AUTO   |
| >>Diffusion Gradient Orientation       | (0018,9089) | FD |       | ANAP              | AUTO   |
| >Diffusion b-value                     | (0018,9087) | FD |       | VNAP              | AUTO   |
| >Diffusion Anisotropy Type             | (0018,9147) | CS |       | ANAP              | AUTO   |

**Table 8.1-60**  
**MR AVERAGES MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name       | Tag         | VR | Value | Presence of Value | Source |
|----------------------|-------------|----|-------|-------------------|--------|
| MR Averages Sequence | (0018,9119) | SQ |       | VNAP              | AUTO   |
| >Number of Averages  | (0018,0083) | DS |       | VNAP              | AUTO   |

**Table 8.1-61**  
**MR SPATIAL SATURATION MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                 | Tag         | VR | Value | Presence of Value | Source |
|--------------------------------|-------------|----|-------|-------------------|--------|
| MR Spatial Saturation Sequence | (0018,9107) | SQ |       | EMPTY             | AUTO   |
| >Slab Thickness                | (0018,9104) | FD |       | ANAP              | AUTO   |
| >Slab Orientation              | (0018,9105) | FD |       | ANAP              | AUTO   |
| >Mid Slab Position             | (0018,9106) | FD |       | ANAP              | AUTO   |

**Table 8.1-62  
MR METABOLITE MAP MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>                             | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|---|-------------|-----------|--------------|--------------------------|---------------|
| MR Metabolite Map Sequence                        | (0018,9152) | SQ        |              | Not Set                  |               |
| >Metabolite Map Description                       | (0018,9080) | ST        |              | Not Set                  |               |
| >Metabolite Map Code Sequence                     | (0018,9083) | SQ        |              | Not Set                  |               |
| >>Code Value                                      | (0008,0100) | SH        |              | Not Set                  |               |
| >>Coding Scheme Designator                        | (0008,0102) | SH        |              | Not Set                  |               |
| >>Coding Scheme Version                           | (0008,0103) | SH        |              | Not Set                  |               |
| >>Code Meaning                                    | (0008,0104) | LO        |              | Not Set                  |               |
| >Chemical Shift Sequence                          | (0018,9084) | SQ        |              | Not Set                  |               |
| >>Chemical Shift Minimum Integration Limit in ppm | (0018,9295) | FD        |              | Not Set                  |               |
| >>Chemical Shift Maximum Integration Limit in ppm | (0018,9296) | FD        |              | Not Set                  |               |

**Table 8.1-63  
MR VELOCITY ENCODING MACRO OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>            | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|----------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| MR Velocity Encoding Sequence    | (0018,9197) | SQ        |              | VNAP                     | AUTO          |
| >Velocity Encoding Direction     | (0018,9090) | FD        |              | VNAP                     | AUTO          |
| >Velocity Encoding Minimum Value | (0018,9091) | FD        |              | VNAP                     | AUTO          |
| >Velocity Encoding Maximum Value | (0018,9217) | FD        |              | VNAP                     | AUTO          |

Table 8.1-64

**MR SPECTROSCOPY FRAME TYPE MACRO OF CREATED MR SPECTROSCOPY SOP INSTANCES**

| Attribute Name                      | Tag         | VR | Value  | Presence of Value | Source |
|-------------------------------------|-------------|----|--------|-------------------|--------|
| MR Spectroscopy Frame TypeSequence  | (0018,9227) | SQ |        | ALWAYS            | AUTO   |
| >Frame Type                         | (0008,9007) | CS |        | ALWAYS            | AUTO   |
| >Volumetric Properties              | (0008,9206) | CS |        | ALWAYS            | AUTO   |
| >Volume Based Calculation Technique | (0008,9207) | CS | "NONE" | ALWAYS            | AUTO   |
| >Complex Image Component            | (0008,9208) | CS |        | ALWAYS            | AUTO   |
| >Acquisition Contrast               | (0008,9209) | CS |        | ALWAYS            | AUTO   |

Table 8.1-65

**MR SPECTROSCOPY FOV/GEOMETRY MACRO OF CREATED MR SPECTROSCOPY SOP INSTANCES**

| Attribute Name                                    | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| MR Spectroscopy FOV/Geometry Sequence             | (0018,9103) | SQ |       | VNAP              | AUTO   |
| >Spectroscopy Acquisition DataColumns             | (0018,9127) | UL |       | VNAP              | AUTO   |
| >Spectroscopy Acquisition PhaseRows               | (0018,0095) | DS |       | VNAP              | AUTO   |
| >Spectroscopy Acquisition PhaseColumns            | (0018,9234) | UL |       | VNAP              | AUTO   |
| >Spectroscopy AcquisitionOut-of-plane Phase Steps | (0018,9159) | UL |       | VNAP              | AUTO   |
| >Percent Sampling                                 | (0018,0093) | DS |       | VNAP              | AUTO   |
| >Percent Phase Field of View                      | (0018,0094) | DS |       | VNAP              | AUTO   |

**Table 8.1-66  
MULTI-FRAME DIMENSION MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                    | Tag         | VR | Value | Presence of Value | Source |
|-----------------------------------|-------------|----|-------|-------------------|--------|
| Dimension Organization Sequence   | (0020,9221) | SQ |       | VNAP              | AUTO   |
| >Dimension Organization UID       | (0020,9164) | UI |       | ALWAYS            | AUTO   |
| Dimension Index Sequence          | (0020,9222) | SQ |       | VNAP              | AUTO   |
| >Dimension Organization UID       | (0020,9164) | UI |       | ANAP              | AUTO   |
| >Dimension Index Pointer          | (0020,9165) | AT |       | ALWAYS            | AUTO   |
| >Functional Group Pointer         | (0020,9167) | AT |       | ANAP              | AUTO   |
| >Dimension Index Private Creator  | (0020,9213) | LO |       | Not Set           |        |
| >Functional Group Private Creator | (0020,9238) | LO |       | Not Set           |        |

**Table 8.1-67  
CARDIAC SYNCHRONIZATION MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                    | Tag         | VR | Value | Presence of Value | Source |
|-----------------------------------|-------------|----|-------|-------------------|--------|
| Low R-R Value                     | (0018,1081) | IS |       | VNAP              | AUTO   |
| High R-R Value                    | (0018,1082) | IS |       | VNAP              | AUTO   |
| Intervals Acquired                | (0018,1083) | IS |       | VNAP              | AUTO   |
| Intervals Rejected                | (0018,1084) | IS |       | VNAP              | AUTO   |
| Cardiac Framing Type              | (0018,1064) | LO |       | VNAP              | AUTO   |
| Cardiac Synchronization Technique | (0018,9037) | CS |       | VNAP              | AUTO   |
| Cardiac Beat Rejection Technique  | (0018,9169) | CS |       | VNAP              | AUTO   |
| Cardiac RR Interval Specified     | (0018,9070) | FD |       | VNAP              | AUTO   |
| Cardiac Signal Source             | (0018,9085) | CS |       | VNAP              | AUTO   |

**Table 8.1-68  
RESPIRATORY SYNCHRONIZATION MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                            | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| Respiratory Motion Compensation Technique | (0018,9170) | CS |       | VNAP              | AUTO   |
| Respiratory Signal Source                 | (0018,9171) | CS |       | ANAP              | AUTO   |
| Respiratory Trigger Type                  | (0020,9250) | CS |       | ANAP              | AUTO   |
| Respiratory Trigger Delay Threshold       | (0020,9256) | FD |       | VNAP              | AUTO   |

Table 8.1-69

**BULK MOTION SYNCHRONIZATION MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>              | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Bulk Motion Compensation Technique | (0018,9172) | CS        |              | VNAP                     | AUTO          |
| Bulk Motion Signal Source          | (0018,9173) | CS        |              | VNAP                     | AUTO          |

**Table 8.1-70**  
**ACQUISITION CONTEXT MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>            | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|----------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Acquisition Context Sequence     | (0040,0555) | SQ        |              | EMPTY                    | AUTO          |
| >Measurement Units Code Sequence | (0040,08EA) | SQ        |              | Not Set                  |               |
| >>Code Value                     | (0008,0100) | SH        |              | Not Set                  |               |
| >>Coding Scheme Designator       | (0008,0102) | SH        |              | Not Set                  |               |
| >>Coding Scheme Version          | (0008,0103) | SH        |              | Not Set                  |               |
| >>Code Meaning                   | (0008,0104) | LO        |              | Not Set                  |               |
| >Value Type                      | (0040,A040) | CS        |              | Not Set                  |               |
| >Concept Name Code Sequence      | (0040,A043) | SQ        |              | Not Set                  |               |
| >>Code Value                     | (0008,0100) | SH        |              | Not Set                  |               |
| >>Coding Scheme Designator       | (0008,0102) | SH        |              | Not Set                  |               |
| >>Coding Scheme Version          | (0008,0103) | SH        |              | Not Set                  |               |
| >>Code Meaning                   | (0008,0104) | LO        |              | Not Set                  |               |
| >Date                            | (0040,A121) | DA        |              | Not Set                  |               |
| >Time                            | (0040,A122) | TM        |              | Not Set                  |               |
| >Person Name                     | (0040,A123) | PN        |              | Not Set                  |               |
| >Referenced Frame Numbers        | (0040,A136) | US        |              | Not Set                  |               |
| >Text Value                      | (0040,A160) | UT        |              | Not Set                  |               |
| >Concept Code Sequence           | (0040,A168) | SQ        |              | Not Set                  |               |
| >>Code Value                     | (0008,0100) | SH        |              | Not Set                  |               |
| >>Coding Scheme Designator       | (0008,0102) | SH        |              | Not Set                  |               |
| >>Coding Scheme Version          | (0008,0103) | SH        |              | Not Set                  |               |
| >>Code Meaning                   | (0008,0104) | LO        |              | Not Set                  |               |
| >Numeric Value                   | (0040,A30A) | DS        |              | Not Set                  |               |
| Acquisition Context Description  | (0040,0556) | ST        |              | Not Set                  |               |

**Table 8.1-71  
ENHANCED MR IMAGE MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| <b>Attribute Name</b>                            | <b>Tag</b>  | <b>VR</b> | <b>Value</b>  | <b>Presence of Value</b> | <b>Source</b> |
|--|-------------|-----------|---------------|--------------------------|---------------|
| Image Type                                       | (0008,0008) | CS        |               | ALWAYS                   | AUTO          |
| Acquisition Datetime                             | (0008,002A) | DT        |               | VNAP                     | AUTO          |
| Referenced Waveform Sequence                     | (0008,113A) | SQ        |               | Not Set                  |               |
| Referenced Image Evidence Sequence               | (0008,9092) | SQ        |               | Not Set                  |               |
| Referenced Raw Data Sequence                     | (0008,9121) | SQ        |               | Not Set                  |               |
| Source Image Evidence Sequence                   | (0008,9154) | SQ        |               | Not Set                  |               |
| Pixel Presentation                               | (0008,9205) | CS        |               | ALWAYS                   | AUTO          |
| Volumetric Properties                            | (0008,9206) | CS        |               | ALWAYS                   | AUTO          |
| Volume Based Calculation Technique               | (0008,9207) | CS        |               | ALWAYS                   | AUTO          |
| Complex Image Component                          | (0008,9208) | CS        |               | ALWAYS                   | AUTO          |
| Acquisition Contrast                             | (0008,9209) | CS        |               | ALWAYS                   | AUTO          |
| Referenced Grayscale Presentation State Sequence | (0008,9237) | SQ        |               | Not Set                  |               |
| Spacing between Slices                           | (0018,0088) | DS        |               | ANAP                     | AUTO          |
| Content Qualification                            | (0018,9004) | CS        |               | VNAP                     | AUTO          |
| k-space Filtering                                | (0018,9064) | CS        |               | VNAP                     | AUTO          |
| Acquisition Duration                             | (0018,9073) | FD        |               | VNAP                     | AUTO          |
| Magnetic Field Strength                          | (0018,0087) | DS        |               | VNAP                     | AUTO          |
| Resonant Nucleus                                 | (0018,9100) | CS        |               | VNAP                     | AUTO          |
| Applicable Safety Standard Agency                | (0018,9174) | CS        |               | VNAP                     | AUTO          |
| Applicable Safety Standard Description           | (0018,9175) | LO        |               | ANAP                     | AUTO          |
| Acquisition Number                               | (0020,0012) | IS        |               | VNAP                     | AUTO          |
| Image Comments                                   | (0020,4000) | LT        |               | ANAP                     | AUTO          |
| Samples per Pixel                                | (0028,0002) | US        | "1"           | ALWAYS                   | AUTO          |
| Photometric Interpretation                       | (0028,0004) | CS        | "MONOCHROME2" | ALWAYS                   | AUTO          |
| Bits Allocated                                   | (0028,0100) | US        | "16"          | ALWAYS                   | AUTO          |
| Bits Stored                                      | (0028,0101) | US        | "16"          | ALWAYS                   | AUTO          |
| High Bit   | (0028,0102) | US        | "15"          | ALWAYS                   | AUTO          |
| Lossy Image Compression                          | (0028,2110) | CS        | "00"          | ALWAYS                   | AUTO          |
| Lossy Image Compression Ratio                    | (0028,2112) | DS        |               | ANAP                     | AUTO          |



|                     |             |    |  |      |      |
|---------------------|-------------|----|--|------|------|
| Icon Image Sequence | (0088,0200) | SQ |  | ANAP | AUTO |
|---------------------|-------------|----|--|------|------|

Table 8.1-72

**MR PULSE SEQUENCE MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name                         | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| Pulse Sequence Name                    | (0018,9005) | SH |       | VNAP              | AUTO   |
| Echo Pulse Sequence                    | (0018,9008) | CS |       | ALWAYS            | AUTO   |
| Multiple Spin Echo                     | (0018,9011) | CS |       | VNAP              | AUTO   |
| Multi-planar Excitation                | (0018,9012) | CS |       | ALWAYS            | AUTO   |
| Phase Contrast                         | (0018,9014) | CS |       | ALWAYS            | AUTO   |
| Time of Flight Contrast                | (0018,9015) | CS |       | ALWAYS            | AUTO   |
| Steady State Pulse Sequence            | (0018,9017) | CS |       | ALWAYS            | AUTO   |
| Echo Planar Pulse Sequence             | (0018,9018) | CS |       | ALWAYS            | AUTO   |
| MR Acquisition Type                    | (0018,0023) | CS |       | ALWAYS            | AUTO   |
| Saturation Recovery                    | (0018,9024) | CS |       | ALWAYS            | AUTO   |
| Spectrally Selected Suppression        | (0018,9025) | CS |       | ALWAYS            | AUTO   |
| Oversampling Phase                     | (0018,9029) | CS |       | ALWAYS            | AUTO   |
| Geometry of k-Space Traversal          | (0018,9032) | CS |       | ALWAYS            | AUTO   |
| Segmented k-Space Traversal            | (0018,9033) | CS |       | ALWAYS            | AUTO   |
| Rectilinear Phase Encode Reordering    | (0018,9034) | CS |       | ANAP              | AUTO   |
| Velocity Encoding Acquisition Sequence | (0018,9092) | SQ |       | ANAP              | AUTO   |
| >Velocity Encoding Direction           | (0018,9090) | FD |       | ANAP              | AUTO   |
| Number of k-Space Trajectories         | (0018,9093) | US |       | ALWAYS            | AUTO   |
| Coverage of k-Space                    | (0018,9094) | CS |       | ANAP              | AUTO   |

Table 8.1-73

**SOP COMMON MODULE OF CREATED ENHANCED MR IMAGE SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value                         | Presence of Value | Source |
|------------------------|-------------|----|-------------------------------|-------------------|--------|
| Specific Character set | (0008,0005) | CS | Refer to Section 6            | ANAP              | CONFIG |
| SOP Class UID          | (0008,0016) | UI | "1.2.840.10008.5.1.4.1.1.4.1" | ALWAYS            | AUTO   |
| SOP Instance UID       | (0008,0018) | UI | Generated by device           | ALWAYS            | AUTO   |

## 8.1.1.10 MR Spectroscopy Modules

Table 8.1-74

## MR SPECTROSCOPY MODULE OF CREATED MR SPECTROSCOPY SOP INSTANCES

| Attribute Name                       | Tag          | VR | Value | Presence of Value | Source |
|--------------------------------------|--------------|----|-------|-------------------|--------|
| Image Type                           | (0008,0008)  | CS |       | VNAP              | AUTO   |
| Acquisition DateTime                 | (0008,002A)  | DT |       | VNAP              | AUTO   |
| Volumetric Properties                | (0008,9206)  | CS |       | VNAP              | AUTO   |
| Volume Based Calculation Technique   | (0008,9207)  | CS |       | VNAP              | AUTO   |
| Complex Image Component              | (0008,9208)  | CS |       | VNAP              | AUTO   |
| Acquisition Contrast                 | (0008,9209)  | CS |       | VNAP              | AUTO   |
| Magnetic Field Strength              | (0018,0087)  | DS |       | VNAP              | AUTO   |
| Content Qualification                | (0018,9004)  | CS |       | VNAP              | AUTO   |
| Spectral Width                       | (0018,9052)  | FD |       | VNAP              | AUTO   |
| Chemical Shift Reference             | (0018,9053)  | FD |       | VNAP              | AUTO   |
| Volume Localization Technique        | (0018,9054)  | CS |       | VNAP              | AUTO   |
| De-Coupling                          | (0018,9059)  | CS |       | VNAP              | AUTO   |
| De-Coupled Nucleus                   | (0018,9060)  | CS |       | VNAP              | AUTO   |
| De-Coupling Frequency                | (0018,9061)  | FD |       | VNAP              | AUTO   |
| De-Coupling Method                   | (0018,9062)  | CS |       | VNAP              | AUTO   |
| De-Coupling Chemical Shift Reference | (0018,9063)  | FD |       | VNAP              | AUTO   |
| k-space Filtering                    | (0018,9064)  | CS |       | VNAP              | AUTO   |
| Time Domain Filtering                | (0018,9065)  | CS |       | VNAP              | AUTO   |
| Number of Zero Fills                 | (0018,9066)  | US |       | VNAP              | AUTO   |
| Baseline Correction                  | (0018,9067)  | CS |       | VNAP              | AUTO   |
| Acquisition Duration                 | (0018,9073)  | FD |       | VNAP              | AUTO   |
| Applicable Safety Standard Agency    | (0018,9174)  | CS |       |                   |        |
| Transmitter Frequency                | (0018, 9098) | FD |       | VNAP              | AUTO   |
| Resonant Nucleus                     | (0018,9100)  | CS |       | VNAP              | AUTO   |
| Frequency Correction                 | (0018, 9101) | CS |       | VNAP              | AUTO   |
| Volume Localization Sequence         | (0018,9126)  | SQ |       | VNAP              | AUTO   |
| >Slab Thickness                      | (0018,9104)  | FD |       | VNAP              | AUTO   |
| >Slab Orientation                    | (0018,9105)  | FD |       | VNAP              | AUTO   |
| >Mid Slab Position                   | (0018,9106)  | FD |       | VNAP              | AUTO   |
| First Order Phase Correction         | (0018,9198)  | CS |       | VNAP              | AUTO   |
| Water Referenced Phase Correction    | (0018,9199)  | CS |       | VNAP              | AUTO   |

**Table 8.1-75**  
**MR SPECTROSCOPY PULSE SEQUENCE MODULE OF CREATED MR SPECTROSCOPY SOP**  
**INSTANCES**

| <b>Attribute Name</b>               | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|-------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Pulse Sequence Name                 | (0018,9005) | SH        |              | ALWAYS                   | AUTO          |
| Echo Pulse Sequence                 | (0018,9008) | CS        |              | ALWAYS                   | AUTO          |
| Multiple Spin Echo                  | (0018,9011) | CS        |              | ALWAYS                   | AUTO          |
| Multi-planar Excitation             | (0018,9012) | CS        |              | ALWAYS                   | AUTO          |
| Steady State Pulse Sequence         | (0018,9017) | CS        |              | ALWAYS                   | AUTO          |
| Echo Planar Pulse Sequence          | (0018,9018) | CS        |              | ANAP                     | AUTO          |
| Spectrally Selected Suppression     | (0018,9025) | CS        |              | ANAP                     | AUTO          |
| Geometry of k-Space Traversal       | (0018,9032) | CS        |              | ANAP                     | AUTO          |
| Segmented k-Space Traversal         | (0018,9033) | CS        |              | ANAP                     | AUTO          |
| Rectilinear Phase Encode Reordering | (0018,9034) | CS        |              | ANAP                     | AUTO          |
| Number of k-Space Trajectories      | (0018,9093) | US        |              | ANAP                     | AUTO          |
| Coverage of k-Space                 | (0018,9094) | CS        |              | ANAP                     | AUTO          |
| MR Spectroscopy Acquisition Type    | (0018,9200) | CS        |              | ALWAYS                   | AUTO          |

**Table 8.1-76**  
**MR SPECTROSCOPY DATA MODULE OF CREATED MR SPECTROSCOPY SOP INSTANCES**

| <b>Attribute Name</b>              | <b>Tag</b>  | <b>VR</b> | <b>Value</b> | <b>Presence of Value</b> | <b>Source</b> |
|------------------------------------|-------------|-----------|--------------|--------------------------|---------------|
| Rows                               | (0028,0010) | US        |              | ALWAYS                   | AUTO          |
| Columns                            | (0028,0011) | US        |              | ALWAYS                   | AUTO          |
| Data Point Rows                    | (0028,9001) | UL        |              | ALWAYS                   | AUTO          |
| Data Point Columns                 | (0028,9002) | UL        |              | ALWAYS                   | AUTO          |
| Signal Domain Columns              | (0028,9003) | CS        |              | ALWAYS                   | AUTO          |
| Data Representation                | (0028,9108) | CS        |              | ALWAYS                   | AUTO          |
| Signal Domain Rows                 | (0028,9235) | CS        |              | ANAP                     | AUTO          |
| First Order Phase Correction Angle | (5600,0010) | OF        |              | ANAP                     | AUTO          |
| Spectroscopy Data                  | (5600,0020) | OF        |              | ALWAYS                   | AUTO          |

**Table 8.1-77**  
**SOP COMMON MODULE OF CREATED MR SPECTROSCOPY SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value                         | Presence of Value | Source |
|------------------------|-------------|----|-------------------------------|-------------------|--------|
| Specific Character set | (0008,0005) | CS | Refer to Section 6            | ANAP              | CONFIG |
| SOP Class UID          | (0008,0016) | UI | "1.2.840.10008.5.1.4.1.1.4.2" | ALWAYS            | AUTO   |
| SOP Instance UID       | (0008,0018) | UI | Generated by device           | ALWAYS            | AUTO   |

**Table 8.1-78**  
**PRIVATE APPLICATION MODULE OF CREATED MR SPECTROSCOPY SOP INSTANCES**

| Attribute Name                                  | Tag         | VR | Value             | Presence of Value | Source |
|---|-------------|----|-------------------|-------------------|--------|
| Private Creator Code                            | (700D,00xx) | LO | "TOSHIBA_MEC_MR3" | ANAP              | AUTO   |
| Number of PAC channel                           | (700D,xx80) | US |                   | ANAP              | AUTO   |
| Reference mode                                  | (700D,xx81) | US |                   | ANAP              | AUTO   |
| Gain value group for MRS                        | (700D,xx82) | SQ |                   | ANAP              | AUTO   |
| > Gain value of each channel for MRS            | (700D,xx83) | FL |                   | ANAP              | AUTO   |
| > Phase value of each channel for MRS           | (700D,xx84) | FL |                   | ANAP              | AUTO   |
| > Reference gain value of each channel for MRS  | (700D,xx85) | FL |                   | ANAP              | AUTO   |
| > Reference phase value of each channel for MRS | (700D,xx86) | FL |                   | ANAP              | AUTO   |
| Receiver gain of reference scan                 | (700D,xx87) | FL |                   | ANAP              | AUTO   |
| Flag of Water Sat pulse                         | (700D,xx88) | UL |                   | ANAP              | AUTO   |
| Selected contrast TE                            | (700D,xx89) | FL |                   | ANAP              | AUTO   |
| Raw Data Set Sequence                           | (700D,xx8A) | SQ |                   | ANAP              | AUTO   |
| >Multi Slice number                             | (700D,xx8B) | DS |                   | ANAP              | AUTO   |
| >Multi Coverage number                          | (700D,xx8C) | DS |                   | ANAP              | AUTO   |
| >Raw Data Sequence                              | (700D,xx8D) | SQ |                   | ANAP              | AUTO   |
| >>Coil Channel No.                              | (700D,xx8E) | DS |                   | ANAP              | AUTO   |
| >>MRS Raw Data                                  | (700D,xx8F) | OB |                   | ANAP              | AUTO   |
| >>MRS Ref Raw Data                              | (700D,xx90) | OF |                   | ANAP              | AUTO   |
| >>MRS Raw Data Type                             | (700D,xx92) | US |                   | ANAP              | AUTO   |
| Receiver gain of prescan                        | (700D,xx91) | FL | in dB             | ANAP              | AUTO   |

### 8.1.1.11 Grayscale Softcopy Presentation State Modules

Table 8.1-79

#### PRESENTATION SERIES MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES

| Attribute Name | Tag         | VR | Value | Presence of Value | Source |
|----------------|-------------|----|-------|-------------------|--------|
| Modality       | (0008,0060) | CS | "PR"  | ALWAYS            | AUTO   |

Table 8.1-80

#### PRESENTATION STATE MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES

| Attribute Name                | Tag         | VR | Value               | Presence of Value | Source |
|-------------------------------|-------------|----|---------------------|-------------------|--------|
| Presentation Series Sequence  | (0008,1115) | SQ |                     | ALWAYS            | AUTO   |
| >Referenced Image Sequence    | (0008,1140) | SQ |                     | ALWAYS            | AUTO   |
| >>Referenced SOP Class UID    | (0008,1150) | UI |                     | ALWAYS            | AUTO   |
| >>Referenced SOP Instance UID | (0008,1155) | UI |                     | ALWAYS            | AUTO   |
| >Series Instance UID          | (0020,000E) | UI |                     | ALWAYS            | AUTO   |
| Instance Number               | (0020,0013) | IS | Generated by device | ALWAYS            | AUTO   |
| Presentation Label            | (0070,0080) | CS | Generated by device | ALWAYS            | AUTO   |
| Presentation Description      | (0070,0081) | LO |                     | VNAP              | USER   |
| Presentation Creation Date    | (0070,0082) | DA | <yyyymmdd>          | ALWAYS            | AUTO   |
| Presentation Creation Time    | (0070,0083) | TM | <hhmmss.frac>       | ALWAYS            | AUTO   |
| Presentation Creator's Name   | (0070,0084) | PN |                     | VNAP              | USER   |

**Table 8.1-81  
DISPLAYED AREA MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP  
INSTANCES**

| Attribute Name                           | Tag         | VR | Value                   | Presence of Value | Source |
|--|-------------|----|-------------------------|-------------------|--------|
| Displayed Area Selection Sequence        | (0070,005A) | SQ |                         | ALWAYS            | AUTO   |
| >Referenced Image Sequence               | (0008,1140) | SQ |                         | ANAP              | AUTO   |
| >>Referenced SOP Class UID               | (0008,1150) | UI |                         | ANAP              | AUTO   |
| >>Referenced SOP Instance UID            | (0008,1155) | UI |                         | ANAP              | AUTO   |
| >Displayed Area Top Left Hand Corner     | (0070,0052) | IS |                         | ALWAYS            | AUTO   |
| >Displayed Area Bottom Right Hand Corner | (0070,0053) | IS |                         | ALWAYS            | AUTO   |
| >Presentation Size Model                 | (0070,0100) | CS | MAGNIFY or SCALE TO FIT | ALWAYS            | AUTO   |
| >Presentation PixelSpacing               | (0070,0101) | DS |                         | ANAP              | AUTO   |
| >Presentation Pixel Aspect Ratio         | (0070,0102) | IS |                         | ANAP              | AUTO   |
| >Presentation Pixel Magnification Ratio  | (0070,0103) | FL |                         | ANAP              | AUTO   |

**Table 8.1-82  
SPATIAL TRANSFORMATION MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION  
STATE SOP INSTANCES**

| Attribute Name        | Tag         | VR | Value | Presence of Value | Source |
|-----------------------|-------------|----|-------|-------------------|--------|
| Image Horizontal Flip | (0070,0041) | US |       | ALWAYS            | AUTO   |
| Image Rotation        | (0070,0042) | US |       | ALWAYS            | AUTO   |

**Table 8.1-83  
MODALITY LUT MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP  
INSTANCES**

| Attribute Name    | Tag         | VR | Value | Presence of Value | Source |
|-------------------|-------------|----|-------|-------------------|--------|
| Rescale Intercept | (0028,1052) | DS |       | ANAP              | AUTO   |
| Rescale Slope     | (0028,1053) | DS |       | ANAP              | AUTO   |
| Rescale Type      | (0028,1054) | LO |       | ANAP              | AUTO   |

**Table 8.1-84**  
**SOFTCOPY VOI LUT MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES**

| Attribute Name                     | Tag         | VR | Value | Presence of Value | Source |
|------------------------------------|-------------|----|-------|-------------------|--------|
| Softcopy VOI LUT Sequence          | (0028,3110) | SQ |       | ALWAYS            | AUTO   |
| >Referenced Image Sequence         | (0008,1140) | SQ |       | ALWAYS            | AUTO   |
| >>Referenced SOP Class UID         | (0008,1150) | UI |       | ALWAYS            | AUTO   |
| >>Referenced SOP Instance UID      | (0008,1155) | UI |       | ALWAYS            | AUTO   |
| >Window Center                     | (0028,1050) | DS |       | ALWAYS            | AUTO   |
| >Window Width                      | (0028,1051) | DS |       | ALWAYS            | AUTO   |
| >Window Center & Width Explanation | (0028,1055) | LO |       | ANAP              | AUTO   |

**Table 8.1-85**  
**SOFTCOPY PRESENTATION LUT MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value | Presence of Value | Source |
|------------------------|-------------|----|-------|-------------------|--------|
| Presentation LUT Shape | (2050,0020) | CS |       | ALWAYS            | AUTO   |

**Table 8.1-86**  
**SOP COMMON MODULE OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value                          | Presence of Value | Source |
|------------------------|-------------|----|--------------------------------|-------------------|--------|
| Specific Character Set | (0008,0005) | CS | Refer to Section 6             | ANAP              | AUTO   |
| Instance Creation Date | (0008,0012) | DA | <yyyymmdd>                     | ANAP              | AUTO   |
| Instance Creation Time | (0008,0013) | TM | <hhmmss.frac>                  | ANAP              | AUTO   |
| SOP Class UID          | (0008,0016) | UI | "1.2.840.10008.5.1.4.1.1.11.1" | ALWAYS            | AUTO   |
| SOP Instance UID       | (0008,0018) | UI | Generated by device            | ALWAYS            | AUTO   |

## 8.1.2 Usage of Attributes from received IOD's

No SOP Class specific fields are required.

## 8.1.3 Attribute Mapping

The relationships between attributes received via Modality Worklist, stored in acquired images and communicated via MPPS are summarized in

**Table 8.1-87**  
**ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST, IMAGE AND MPPS**

| Modality Worklist                    | Image IOD                             | MPPS IOD                              |
|--------------------------------------|---------------------------------------|---------------------------------------|
| --                                   | --                                    | Scheduled Step Attribute Sequence     |
| Study Instance UID                   | Study Instance UID                    | >Study Instance UID                   |
| Referenced Study Sequence            | Referenced Study Sequence             | >Referenced Study Sequence            |
| Accession Number                     | Accession Number                      | >Accession Number                     |
| --                                   | Request Attributes Sequence           | --                                    |
| Requested Procedure ID               | >Requested Procedure ID               | >Requested Procedure ID               |
| Scheduled Procedure Step ID          | >Scheduled Procedure Step ID          | >Scheduled Procedure Step ID          |
| Scheduled Procedure Step Description | >Scheduled Procedure Step Description | >Scheduled Procedure Step Description |
| --                                   | Study ID                              | Study ID                              |
| Requested Procedure Description      | Study Description                     | Requested Procedure Description       |
| Scheduled Procedure Step ID          | Performed Procedure Step ID           | Performed Procedure Step ID           |
| Scheduled Procedure Step Description | Performed Procedure Step Description  | Performed Procedure Step Description  |
| Requested Procedure Code Sequence    | Requested Procedure Code Sequence     | Requested Procedure Code Sequence     |

This table shows only typical data sets.

Other data sets are also set as default settings.

All map settings, including the default setting data sets, can be customized.

## 8.1.4 Coerced/Modified Fields

Not applicable to this product



## 8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

Table 8.2-1  
DATA DICTIONARY OF PRIVATE ATTRIBUTES

| Attribute Name                                   | Tag         | VR | VM  | Value                 |
|--|-------------|----|-----|-----------------------|
| Private Creator Code                             | (0029,00xx) | LO | 1   | "TOSHIBA_MEC_MR3"     |
| Other Private Data                               | (0029,xx01) | SQ | 1   |                       |
| Private Creator Code                             | (700D,00xx) | LO | 1   | "TOSHIBA_MEC_MR3"     |
| Scale Factor                                     | (700D,xx00) | DS | 1   | See TABLE 8.1-27      |
| FOV  | (700D,xx05) | DS | 1   | See Table 8.1-27      |
| Receiver Gain Correction Check Flag              | (700D,xx0C) | CS | 1   | See Table 8.1-27      |
| Identification Flag of 3D GDC                    | (700D,xx20) | SH | 1   | See Table 8.1-27      |
| Private Creator Code                             | (700D,00yy) | LO | 1   | "TOSHIBA_MEC_MR3^10". |
| 2nd Flip Angle [degree]                          | (700D,yy10) | DS | 1   | See Table 8.1-27      |
| Acquisition Inner Matrix                         | (700D,yy11) | US | 2   | See Table 8.1-27      |
| MP2RAGE Flag                                     | (700D,yy12) | US | 1   | See Table 8.1-27      |
| Inversion efficiency of inversion recovery pulse | (700D,yy13) | FL | 1   | See Table 8.1-27      |
| Number of dummy shot                             | (700D,yy14) | SL | 1   | See Table 8.1-27      |
| FFE total repetition time[s]                     | (700D,yy15) | FL | 1   | See Table 8.1-27      |
| PAS Name   | (700D,yy16) | LO | 3   | See Table 8.1-27      |
| Intended Processing                              | (700D,yy17) | LT | 1   | See Table 8.1-27      |
| Scanned Orientation IDs                          | (700D,yy18) | SS | 1-n | See Table 8.1-27      |
| PAS Reproduct Information                        | (700D,yy19) | OB | 1   | See Table 8.1-27      |
| Private Creator Code                             | (700D,00xx) | LO | 1   | "TOSHIBA_MEC_MR3"     |
| Number of PAC channel                            | (700D,xx80) | US | 1   | See Table 8.1-78      |
| Reference mode                                   | (700D,xx81) | US | 1   | See Table 8.1-78      |
| Gain value group for MRS                         | (700D,xx82) | SQ | 1-n | See Table 8.1-78      |
| Gain value of each channel for MRS               | (700D,xx83) | FL | 1   | See Table 8.1-78      |
| Phase value of each channel for MRS              | (700D,xx84) | FL | 1   | See Table 8.1-78      |
| Reference gain value of each channel for MRS     | (700D,xx85) | FL | 1   | See Table 8.1-78      |
| Reference phase value of each channel for MRS    | (700D,xx86) | FL | 1   | See Table 8.1-78      |
| Receiver gain of reference scan                  | (700D,xx87) | FL | 1   | See Table 8.1-78      |
| Flag of Water Sat pulse                          | (700D,xx88) | UL | 1   | See Table 8.1-78      |
| Selected contrast TE                             | (700D,xx89) | FL | 1   | See Table 8.1-78      |
| Raw Data Set Sequence                            | (700D,xx8A) | SQ | 1-n | See Table 8.1-78      |
| Multi Slice number                               | (700D,xx8B) | DS | 1   | See Table 8.1-78      |
| Multi Coverage number                            | (700D,xx8C) | DS | 1   | See Table 8.1-78      |
| Raw Data Sequence                                | (700D,xx8D) | SQ | 1-n | See Table 8.1-78      |
| Coil Channel No.                                 | (700D,xx8E) | DS | 1   | See Table 8.1-78      |
| MRS Raw Data                                     | (700D,xx8F) | OB | 1   | See Table 8.1-78      |

|                          |             |    |   |                  |
|--------------------------|-------------|----|---|------------------|
| MRS Ref Raw Data         | (700D,xx90) | OF | 1 | See Table 8.1-78 |
| MRS Raw Data Type        | (700D,xx92) | US | 1 | See Table 8.1-78 |
| Receiver gain of prescan | (700D,xx91) | FL | 1 | See Table 8.1-78 |

### **8.3 CONTROLLED TERMINOLOGY AND TEMPLATES**

Not applicable to this product

### **8.4 GRAYSCALE IMAGE CONSISTENCY**

Not applicable to this product

### **8.5 STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES**

#### **8.5.1 Standard Extended SOP Class - MR Image Storage**

The Storage SCU AE and the Offline-Media AE are making the following extensions to DICOM SOP Classes:

SOP Class : MR Image Storage  
Attribute : Acquisition Duration(0018,9073)  
Diffusion b-value (0018,9087)  
Diffusion Gradient Orientation (0018,9089)  
Stack ID (0020,9056)  
In-Stack Position Number (0020,9057)  
Temporal Position Index (0020,9128)

SOP Class : SC Image Storage  
Attribute : Acquisition Duration(0018,9073)  
Stack ID (0020,9056)  
In-Stack Position Number (0020,9057)  
Temporal Position Index (0020,9128)

### **8.6 PRIVATE TRANSFER SYNTAXES**

Not applicable to this product