



**DICOM CONFORMANCE STATEMENT  
FOR  
DIAGNOSTIC ULTRASOUND SYSTEM**

***Aplio air***

**MODEL CUS-AAR00 V1.0**

## Trademarks

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.

This document may include trademarks or registered trademarks of other companies.

### **IMPORTANT!**

- (1) No part of this document may be copied or reprinted, in whole or in part, without written permission.
- (2) The contents of this document are subject to change without prior notice and without our legal obligation.
- (3) Please refer to the CANON INC. website for the most recent version of this conformance statement.

Global: <https://www.medical.canon/Interoperability/DICOM/EN>

Japan: <https://www.medical.canon/Interoperability/DICOM/JP>

## 1. CONFORMANCE STATEMENT OVERVIEW

Table 1-1 provides an overview of the network services supported by *Aplio air*.

**Table 1-1  
NETWORK SERVICES**

| <b>SOP Classes</b>                         | <b>User of Service<br/>(SCU)</b> | <b>Provider of Service<br/>(SCP)</b> |
|--------------------------------------------|----------------------------------|--------------------------------------|
| <b>Transfer</b>                            |                                  |                                      |
| Ultrasound Image Storage                   | Yes                              | No                                   |
| Ultrasound Multi-frame Image Storage       | Yes                              | No                                   |
| <b>Workflow Management</b>                 |                                  |                                      |
| Modality Worklist Information Model – FIND | Yes                              | No                                   |

## 2. TABLE OF CONTENTS

|           |                                                         |           |
|-----------|---------------------------------------------------------|-----------|
| <b>1.</b> | <b>CONFORMANCE STATEMENT OVERVIEW .....</b>             | <b>i</b>  |
| <b>2.</b> | <b>TABLE OF CONTENTS .....</b>                          | <b>a</b>  |
| <b>3.</b> | <b>INTRODUCTION .....</b>                               | <b>1</b>  |
| 3.1       | REVISION HISTORY .....                                  | 1         |
| 3.2       | AUDIENCE.....                                           | 1         |
| 3.3       | REMARKS .....                                           | 1         |
| 3.4       | DEFINITIONS, TERMS AND ABBREVIATIONS .....              | 2         |
| 3.5       | REFERENCES.....                                         | 3         |
| <b>4.</b> | <b>NETWORKING.....</b>                                  | <b>4</b>  |
| 4.1       | IMPLEMENTATION MODEL.....                               | 4         |
| 4.1.1     | Application Data Flow .....                             | 4         |
| 4.1.2     | Functional Definition of AEs.....                       | 6         |
| 4.1.3     | Sequencing of Real-World Activities.....                | 7         |
| 4.2       | AE SPECIFICATIONS .....                                 | 8         |
| 4.2.1     | Verification SCU AE Specification .....                 | 8         |
| 4.2.2     | Storage SCU AE Specification.....                       | 11        |
| 4.2.3     | MWM SCU AE Specification.....                           | 15        |
| 4.3       | NETWORK INTERFACES.....                                 | 21        |
| 4.3.1     | Physical Network Interface .....                        | 21        |
| 4.3.2     | Additional Protocols .....                              | 21        |
| 4.3.3     | IPv4 and IPv6 Support.....                              | 21        |
| 4.4       | CONFIGURATION .....                                     | 22        |
| 4.4.1     | AE Title/Presentation Address Mapping .....             | 22        |
| 4.4.2     | Parameters.....                                         | 23        |
| <b>5.</b> | <b>MEDIA INTERCHANGE.....</b>                           | <b>24</b> |
| <b>6.</b> | <b>SUPPORT OF CHARACTER SETS.....</b>                   | <b>25</b> |
| <b>7.</b> | <b>SECURITY.....</b>                                    | <b>26</b> |
| <b>8.</b> | <b>ANNEXES .....</b>                                    | <b>27</b> |
| 8.1       | IOD CONTENTS .....                                      | 27        |
| 8.1.1     | Created SOP Instances .....                             | 27        |
| 8.1.2     | Usage of Attributes from received IOD's.....            | 34        |
| 8.1.3     | Attribute Mapping .....                                 | 34        |
| 8.1.4     | Coerced/Modified Fields .....                           | 36        |
| 8.2       | DATA DICTIONARY OF PRIVATE ATTRIBUTES .....             | 37        |
| 8.3       | CODED TERMINOLOGY AND TEMPLATES.....                    | 37        |
| 8.4       | GRAYSCALE IMAGE CONSISTENCY.....                        | 37        |
| 8.5       | STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES ..... | 37        |
| 8.6       | PRIVATE TRANSFER SYNTAXES.....                          | 37        |
| 8.7       | STANDARD EXTENDED AND PRIVATE TEMPLATES.....            | 37        |
| 8.8       | DICOM Security Profile Details .....                    | 37        |

### 3. INTRODUCTION

#### 3.1 REVISION HISTORY

**Table 3.1-1**  
**REVISION HISTORY**

| REV. | Date of Issue | Author                | Description             |
|------|---------------|-----------------------|-------------------------|
|      | April 5, 2024 | Canon Medical Systems | Initial Version         |
| *A   | Feb 13, 2025  | Canon Medical Systems | Fix VOI LUT description |
| *B   | Jun. 2026     | CANON INC.            | Change of company name  |

#### 3.2 AUDIENCE

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

#### 3.3 REMARKS

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality.

This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended.

The scope of this Conformance Statement is to facilitate communication with CANON INC. and other vendors' Medical equipment. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [DICOM]. However, by itself it is not guaranteed to ensure the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity between CANON INC. and non-CANON INC. equipment.
- Test procedures should be defined to validate the desired level of connectivity.
- The DICOM standard will evolve to meet the users' future requirements. CANON INC. is actively involved in developing the standard further and therefore reserves the right to make changes to its products or to discontinue its delivery.

### 3.4 DEFINITIONS, TERMS AND ABBREVIATIONS

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard.

Abbreviations and terms are as follows:

|                  |                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------|
| <b>AE</b>        | Application Entity                                                                                            |
| <b>ACSE</b>      | Association Control Service Element                                                                           |
| <b>BCP195</b>    | Recommendations for Secure Use of Transport Layer Security (TLS) and Datagram Transport Layer Security (DTLS) |
| <b>CD-R</b>      | Compact Disk Recordable                                                                                       |
| <b>CM</b>        | Code Meaning (0008,0104)                                                                                      |
| <b>CSD</b>       | Coding Scheme Designator (0008,0102)                                                                          |
| <b>CV</b>        | Code Value (0008,0100)                                                                                        |
| <b>DHCP</b>      | Dynamic Host Configuration Protocol                                                                           |
| <b>DIMSE</b>     | DICOM Message Service Element                                                                                 |
| <b>DNS</b>       | Domain Name System                                                                                            |
| <b>DVD</b>       | A trademark of the DVD forum that is not an abbreviation                                                      |
| <b>DVD±R</b>     | DVD Plus/Dash Recordable                                                                                      |
| <b>EV</b>        | Enumerated Value                                                                                              |
| <b>FSC</b>       | File-Set Creator                                                                                              |
| <b>FSR</b>       | File-Set Reader                                                                                               |
| <b>FSU</b>       | File-Set Updater                                                                                              |
| <b>IE</b>        | Information Entity                                                                                            |
| <b>IEEE</b>      | Institute of Electrical and Electronics Engineers                                                             |
| <b>IOD</b>       | Information Object Definition                                                                                 |
| <b>IPv4</b>      | Internet Protocol version 4                                                                                   |
| <b>IPv6</b>      | Internet Protocol version 6                                                                                   |
| <b>ISO</b>       | International Standard Organization                                                                           |
| <b>JPEG</b>      | Joint Photographic Experts Group                                                                              |
| <b>MPPS</b>      | Modality Performed Procedure Step                                                                             |
| <b>MSPS</b>      | Modality Scheduled Procedure Step                                                                             |
| <b>MWM</b>       | Modality Worklist Management                                                                                  |
| <b>NTP</b>       | Network Time Protocol                                                                                         |
| <b>PDU</b>       | Protocol Data Unit                                                                                            |
| <b>Q/R</b>       | Query/Retrieve                                                                                                |
| <b>RLE</b>       | Run Length Encoding                                                                                           |
| <b>SC Option</b> | Service Class Option                                                                                          |
| <b>SCU</b>       | Service Class User (DICOM client)                                                                             |
| <b>SCP</b>       | Service Class Provider (DICOM server)                                                                         |
| <b>SOP</b>       | Service-Object Pair                                                                                           |
| <b>SR</b>        | Structured Reporting                                                                                          |
| <b>TCP/IP</b>    | Transmission Control Protocol/Internet Protocol                                                               |

|            |                          |
|------------|--------------------------|
| <b>TLS</b> | Transport Layer Security |
| <b>UID</b> | Unique Identifier        |
| <b>USB</b> | Universal Serial Bus     |
| <b>VR</b>  | Value Representation     |


### **3.5 REFERENCES**

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <https://www.dicomstandard.org/>

## 4. NETWORKING

### 4.1 IMPLEMENTATION MODEL

#### 4.1.1 Application Data Flow



**Figure 4.1-1**  
**APPLICATION DATA FLOW DIAGRAM**



- The Verification SCU AE issues a C-ECHO to verify a DICOM connection to a remote AE. It is associated with the local real-world activity "Verify Connectivity". "Verify Connectivity" is performed via the Service Tool.
- The Storage SCU AE sends instances to a remote AE. It is associated with the local real-world activity "Send Instances". "Send Instances" is performed upon user request for specific instances selected.
- The MWM SCU AE receives worklist information from a remote AE. It is associated with the local real-world activity "Update Worklist". When the "Update Worklist" is performed the MWM SCU AE queries a remote AE for worklist items and provides the set of worklist items matching the query request. "Update Worklist" is performed manually or automatically.

## **4.1.2 Functional Definition of AEs**

### **4.1.2.1 Functional Definition of Verification SCU AE**

The Verification SCU AE issues a C-ECHO to verify a DICOM connection to a remote AE. It is performed via the Service Tool.


### **4.1.2.2 Functional Definition of Storage SCU AE**

The existence of a send-job queue entry with associated network destination will activate the Storage SCU AE. An association request is sent to the destination AE and upon successful negotiation of a Presentation Context the image transfer is started. If the image transfer fails, the Storage SCU AE will retry this send-job automatically.

### **4.1.2.3 Functional Definition of MWM SCU AE**

The MWM SCU AE attempts to download a worklist from a remote node. If the MWM SCU AE establishes an association to a remote AE, it will transfer patient's information and worklist items via the open association. The results will be displayed in a separate list. The patient's information will be used for the patient registration.

### 4.1.3 Sequencing of Real-World Activities



**Figure 4.1-2**  
**SEQUENCING CONSTRAINTS**

Under typical scheduled workflow conditions the sequencing constraints illustrated in Figure 4.1-2 apply:

1. Query Worklist
2. Receive Worklist of Modality Scheduled Procedure Steps (MSPS)
3. Select Workitem (MSPS) from Worklist
4. Store Acquired Instances

Other workflow situations (e.g. unscheduled procedure steps) will have other sequencing constraints. Some activities may be omitted according to situations.

## 4.2 AE SPECIFICATIONS

### 4.2.1 Verification SCU AE Specification

#### 4.2.1.1 SOP Classes

The Verification SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-1**  
**SOP CLASSES FOR THE VERIFICATION SCU AE**

| SOP Class Name | SOP Class UID     | SCU | SCP |
|----------------|-------------------|-----|-----|
| Verification   | 1.2.840.10008.1.1 | Yes | No  |

#### 4.2.1.2 Association Policies

##### 4.2.1.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-2**  
**DICOM APPLICATION CONTEXT FOR THE VERIFICATION SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

##### 4.2.1.2.2 Number of Associations

The Verification SCU AE initiates one association at a time.

**Table 4.2-3**  
**NUMBER OF ASSOCIATIONS INITIATED FOR THE VERIFICATION SCU AE**

|                                             |   |
|---------------------------------------------|---|
| Maximum number of simultaneous associations | 1 |
|---------------------------------------------|---|

##### 4.2.1.2.3 Asynchronous Nature

The Verification SCU AE does not support asynchronous communication (multiple outstanding transactions over a single association).

**Table 4.2-4**  
**ASYNCHRONOUS NATURE FOR THE VERIFICATION SCU AE**

|                                                         |   |
|---------------------------------------------------------|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---------------------------------------------------------|---|

##### 4.2.1.2.4 Implementation Identifying Information

The implementation information for the Verification SCU AE is:

**Table 4.2-5**  
**DICOM Implementation Class and Version FOR THE VERIFICATION SCU AE**

|                             |                                 |
|-----------------------------|---------------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.6.39.1000.1 |
| Implementation Version Name | CM_UL_DCM_V1.0                  |

4.2.1.3 Association Initiation Policy

4.2.1.3.1 Activity – Verify Connectivity

4.2.1.3.1.1 Description and Sequencing of Activities

The Verification SCU AE attempts to initiate a new association in order to issue a verification request (C-ECHO).




Figure 4.2-1  
SEQUENCING OF ACTIVITY – VERIFY CONNECTIVITY

A possible sequence of interactions between the Verification SCU AE and an Image Server (e.g. a storage or archive device supporting the Verification SOP Classes as an SCP) is illustrated in the figure above:

- 1. The Verification SCU AE opens an association with the Image Server.
- 2. The Verification SCU AE issues a verification request (C-ECHO) and the Image Server replies with a C-ECHO response (status success).
- 3. The Verification SCU AE closes the association with the Image Server.

4.2.1.3.1.2 Proposed Presentation Contexts

The Verification SCU AE will propose the Presentation Contexts shown in the following table:

Table 4.2-6  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY VERIFY CONNECTIVITY

| Presentation Context Table |                   |                           |                     |      |           |
|----------------------------|-------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                   | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID               | Name List                 | UID List            |      |           |
| Verification               | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                            |                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |

#### 4.2.1.3.1.3 SOP Specific Conformance for Verification SOP Class

The Verification SCU AE provides standard conformance to the Verification Service Class as an SCU.

The behavior of Verification SCU AE when encountering status codes in a C-ECHO response is summarized in the table below:

**Table 4.2-7**  
**VERIFICATION RESPONSE STATUS HANDLING BEHAVIOR**

| <b>Service Status</b> | <b>Further Meaning</b> | <b>Status Code</b> | <b>Behavior</b>                                                                    |
|-----------------------|------------------------|--------------------|------------------------------------------------------------------------------------|
| Success               | Success                | 0000               | The Verification SCU AE judges the remote AE is present and active on the network. |

The behavior of Verification SCU AE during communication failure is summarized in the table below:

**Table 4.2-8**  
**VERIFICATION COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>                                                   |
|--------------------------------------------------|-------------------------------------------------------------------|
| Timeout                                          | The association is aborted. The loading box disappears on the UI. |
| Association aborted by the SCP or network layers | The association is aborted. The loading box disappears on the UI. |

## 4.2.2 Storage SCU AE Specification

### 4.2.2.1 SOP Classes

The Storage SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-9**  
**SOP CLASSES FOR THE STORAGE SCU AE**

| SOP Class Name                       | SOP Class UID               | SCU | SCP |
|--------------------------------------|-----------------------------|-----|-----|
| Ultrasound Image Storage             | 1.2.840.10008.5.1.4.1.1.6.1 | Yes | No  |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 |     |     |

### 4.2.2.2 Association Policies

#### 4.2.2.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-10**  
**DICOM APPLICATION CONTEXT FOR THE STORAGE SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.2.2.2 Number of Associations

The Storage SCU AE can initiate up to an association at a time for each destination to which a transfer request is being processed in the active job queue list. Up to ten jobs, that instances will be sent to the different remote hosts, will be active at a time, the other remains pending until the active job is completed or failed.

**Table 4.2-11**  
**NUMBER OF ASSOCIATIONS INITIATED FOR THE STORAGE SCU AE**

|                                             |   |
|---------------------------------------------|---|
| Maximum number of simultaneous associations | 1 |
|---------------------------------------------|---|

#### 4.2.2.2.3 Asynchronous Nature

The Storage SCU AE does not support asynchronous communication (multiple outstanding transactions over a single association).

**Table 4.2-12**  
**ASYNCHRONOUS NATURE FOR THE STORAGE SCU AE**

|                                                         |   |
|---------------------------------------------------------|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---------------------------------------------------------|---|

#### 4.2.2.2.4 Implementation Identifying Information

The implementation information for the Storage SCU AE is:

**Table 4.2-13**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR THE STORAGE SCU AE**


|                             |                                 |
|-----------------------------|---------------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.6.39.1000.1 |
| Implementation Version Name | CM_UL_DCM_V1.0                  |

### 4.2.2.3 Association Initiation Policy

#### 4.2.2.3.1 Activity – Send Instances

##### 4.2.2.3.1.1 Description and Sequencing of Activities

The Storage SCU AE attempts to initiate a new association in order to issue a storage request (C-STORE). If the job contains multiple instances then multiple C-STORE requests will be issued over the same association. If the instance transfer fails, the Storage SCU AE will retry this send-job automatically.



**Figure 4.2-2**  
**SEQUENCING OF ACTIVITY – SEND INSTANCES**

A possible sequence of interactions between the Storage SCU AE and an Image Server (e.g. a storage or archive device supporting the Storage SOP Classes as an SCP) is illustrated in the Figure above:

1. The Storage SCU AE opens an association with the Image Server.
2. Acquired instances are transmitted to the Image Server using a storage request (C-STORE) and the Image Server replies with a C-STORE response (status success).
3. The Storage SCU AE closes the association with the Image Server.



**4.2.2.3.1.2 Proposed Presentation Contexts**

The Storage SCU AE will propose the Presentation Contexts in the following table that shows one Presentation Context Item per row:

**Table 4.2-14  
PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND INSTANCES**

| <b>Presentation Context Table</b>    |                             |                           |                        |             |                  |
|--------------------------------------|-----------------------------|---------------------------|------------------------|-------------|------------------|
| <b>Abstract Syntax</b>               |                             | <b>Transfer Syntax</b>    |                        | <b>Role</b> | <b>Ext. Neg.</b> |
| <b>Name</b>                          | <b>UID</b>                  | <b>Name List</b>          | <b>UID List</b>        |             |                  |
| Ultrasound Image Storage             | 1.2.840.10008.5.1.4.1.1.6.1 | Implicit VR Little Endian | 1.2.840.10008.1.2      | SCU         | None             |
|                                      |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1    |             |                  |
|                                      |                             | JPEG Lossy* <sup>1</sup>  | 1.2.840.10008.1.2.4.50 |             |                  |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | Implicit VR Little Endian | 1.2.840.10008.1.2      |             |                  |
|                                      |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1    |             |                  |
|                                      |                             | JPEG Lossy* <sup>1</sup>  | 1.2.840.10008.1.2.4.50 |             |                  |

\*1 JPEG Baseline (Process 1)

#### 4.2.2.3.1.3 SOP Specific Conformance for Storage SOP Classes

The Storage SCU AE provides standard conformance to the Storage Service Class as an SCU.

The behavior of Storage SCU AE when encountering status codes in a C-STORE response is summarized in the table below:

**Table 4.2-15**  
**STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning | Status Code           | Behavior                                                                                                                                     |
|----------------|-----------------|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Success        | Success         | 0000                  | The SCP has successfully stored the SOP Instance. If all SOP Instances in a send job have status success then the job is marked as complete. |
| *              | *               | Any other status code | The association is aborted and the send job is marked as failed.                                                                             |

The behavior of Storage SCU AE during communication failure is summarized in the table below:

**Table 4.2-16**  
**STORAGE COMMUNICATION FAILURE BEHAVIOR**

| Exception                                        | Behavior                                                                                                |
|--------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| Timeout                                          | The association is aborted and the send job is marked as failed. An error message will appear on the UI |
| Association aborted by the SCP or network layers | The send job is marked as failed. An error message will appear on the UI                                |

If the instance transfer fails, the Storage SCU AE will retry this send-job automatically (see Section 4.4.2).

The contents of Storage SOP Instances created by the Storage SCU AE conform to the IOD definitions and are described in section 8.1.

## 4.2.3 MWM SCU AE Specification

### 4.2.3.1 SOP Classes

The MWM SCU AE provides Standard Conformance to the following SOP Classes:

**Table 4.2-17**  
**SOP CLASSES FOR THE MWM SCU AE**

| SOP Class Name                             | SOP Class UID          | SCU | SCP |
|--------------------------------------------|------------------------|-----|-----|
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Yes | No  |

### 4.2.3.2 Association Policies

#### 4.2.3.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 4.2-18**  
**DICOM APPLICATION CONTEXT FOR THE MWM SCU AE**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 4.2.3.2.2 Number of Associations

The MWM SCU AE initiates one association at a time for a worklist request.

**Table 4.2-19**  
**NUMBER OF ASSOCIATIONS INITIATED FOR THE MWM SCU AE**

|                                             |   |
|---------------------------------------------|---|
| Maximum number of simultaneous associations | 1 |
|---------------------------------------------|---|

#### 4.2.3.2.3 Asynchronous Nature

The MWM SCU AE does not support asynchronous communication (multiple outstanding transactions over a single association).

**Table 4.2-20**  
**ASYNCHRONOUS NATURE FOR THE MWM SCU AE**

|                                                         |   |
|---------------------------------------------------------|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---------------------------------------------------------|---|

#### 4.2.3.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 4.2-21**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR THE MWM SCU AE**

|                             |                                 |
|-----------------------------|---------------------------------|
| Implementation Class UID    | 1.2.392.200036.9116.6.39.1000.1 |
| Implementation Version Name | CM_UL_DCM_V1.0                  |

### 4.2.3.3 Association Initiation Policy


#### 4.2.3.3.1 Activity – Update Worklist

##### 4.2.3.3.1.1 Description and Sequencing of Activities

The request for an "Update Worklist" is initiated by user interaction, i.e. pressing the buttons "Get Worklist" or automatically at the time of patient registration.

Upon initiation of the request, the MWM SCU AE will build an Identifier for the C-FIND request, will initiate an association to send the request and will wait for worklist responses. After retrieval of all responses, the MWM SCU AE will access the local database to add or update patient demographic data. The results will be displayed in a separate list.

The MWM SCU AE will initiate an association in order to issue a C-FIND request according to the Modality Worklist Information Model.



**Figure 4.2-3**  
**SEQUENCING OF ACTIVITY – UPDATE WORKLIST**

A possible sequence of interactions between the MWM SCU AE and a Department Scheduler (e.g. a device such as a RIS or HIS which supports the Modality Worklist SOP Class as an SCP) is illustrated in the Figure above:

1. The MWM SCU AE opens an association with the Department Scheduler
2. The MWM SCU AE sends a C-FIND request to the Department Scheduler containing the Worklist Query attributes.
3. The Department Scheduler returns a C-FIND response containing the requested attributes of the first matching worklist item.
4. The Department Scheduler returns another C-FIND response containing the requested attributes of the second matching worklist item.
5. The Department Scheduler returns another C-FIND response with status Success indicating that no further matching worklist items exist. This example assumes that only 2 worklist items match the Worklist Query.
6. The MWM SCU AE closes the association with the Department Scheduler.
7. The user selects a worklist item from the Worklist and prepares to acquire new instances.

#### 4.2.3.3.1.2 Proposed Presentation Contexts

The MWM SCU AE will propose Presentation Contexts shown in the following table:

**Table 4.2-22**  
**PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY UPDATE WORKLIST**

| Presentation Context Table                 |                        |                           |                     |      |          |
|--------------------------------------------|------------------------|---------------------------|---------------------|------|----------|
| Abstract Syntax                            |                        | Transfer Syntax           |                     | Role | Ext.Neg. |
| Name                                       | UID                    | Name List                 | UID List            |      |          |
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None     |
|                                            |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |          |

#### 4.2.3.3.1.3 SOP Specific Conformance for Modality Worklist SOP Class

The MWM SCU AE provides standard conformance to the Modality Worklist SOP Class as an SCU.

The behavior of the MWM SCU when encountering status codes in a Modality Worklist C-FIND response is summarized in the table below.

**Table 4.2-23**  
**MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR**

| Service Status | Further Meaning                                                                    | Status Code           | Behavior                                                                                                                                                                                                                                                          |
|----------------|------------------------------------------------------------------------------------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Success        | Matching is complete                                                               | 0000                  | The SCP has completed the matches. Worklist items are available for display or further processing.                                                                                                                                                                |
| Refused        | Out of Resources                                                                   | A700                  | The Association is aborted using A-ABORT and the worklist query is marked as failed.                                                                                                                                                                              |
| Failed         | Identifier does not match SOP Class                                                | A900                  |                                                                                                                                                                                                                                                                   |
| Failed         | Unable to Process                                                                  | Cxxx                  |                                                                                                                                                                                                                                                                   |
| Cancel         | Matching terminated due to Cancel request                                          | FE00                  | If the query was cancelled due to too many worklist items then the SCP has completed the matches. Worklist items are available for display or further processing. Otherwise, the Association is aborted using A-ABORT and the worklist query is marked as failed. |
| Pending        | Matches are continuing                                                             | FF00                  | The worklist item contained in the Identifier is collected for later display or further processing.                                                                                                                                                               |
| Pending        | Matches are continuing – Warning that one or more Optional Keys were not supported | FF01                  | The worklist item contained in the Identifier is collected for later display or further processing.                                                                                                                                                               |
| *              | *                                                                                  | Any other status code | The association is aborted using A-ABORT.                                                                                                                                                                                                                         |

The behavior of the MWM SCU AE during communication failure is summarized in the table below.

**Table 4.2-24**  
**MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR**

| <b>Exception</b>                                 | <b>Behavior</b>                                                                                                         |
|--------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
| Timeout                                          | The association is aborted using A-ABORT. The loading box disappears on the UI and existing worklists are not replaced. |
| Unsupported character sets                       |                                                                                                                         |
| Association aborted by the SCP or network layers | The association is aborted using A-ABORT. The loading box disappears on the UI and existing worklists are not replaced. |

Acquired instances will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

The table below provides a description of the MWM SCU AE Worklist Request Identifier and specifies the attributes that are copied into the instances. Unexpected attributes returned in a C-FIND response are ignored.

**Table 4.2-25  
WORKLIST REQUEST IDENTIFIER**

| Module Name<br>Attribute Name        | Tag         | VR | M   | R | D | IOD |
|--------------------------------------|-------------|----|-----|---|---|-----|
| <b>Scheduled Procedure Step</b>      |             |    |     |   |   |     |
| Scheduled Procedure Step Sequence    | (0040,0100) | SQ |     |   |   |     |
| >Modality                            | (0008,0060) | CS | S   | x | x | x   |
| >Scheduled Station AE Title          | (0040,0001) | AE |     | x | x |     |
| >Scheduled Procedure Step Start Date | (0040,0002) | DA | S,R | x | x |     |
| >Scheduled Procedure Step End Date   | (0040,0004) | DA |     | x |   |     |
| <b>Requested Procedure</b>           |             |    |     |   |   |     |
| Study Instance UID                   | (0020,000D) | UI |     | x |   |     |
| Requested Procedure ID               | (0040,1001) | SH |     | x |   |     |
| <b>Imaging Service Request</b>       |             |    |     |   |   |     |
| Accession Number                     | (0008,0050) | SH |     | x | x | x   |
| <b>Patient Identification</b>        |             |    |     |   |   |     |
| Patient's Name                       | (0010,0010) | PN | S   | x | x | x   |
| Patient ID                           | (0010,0020) | LO | S   | x | x | x   |
| <b>Patient Demographic</b>           |             |    |     |   |   |     |
| Patient's Birth Date                 | (0010,0030) | DA |     | x | x | x   |
| Patient's Sex                        | (0010,0040) | CS |     | x | x | x   |
| Patient's Age                        | (0010,1010) | AS |     | x | x | x   |
| Patient's Size                       | (0010,1020) | DS |     | x | x | x   |
| Patient's Weight                     | (0010,1030) | DS |     | x | x | x   |
| Patient's Address                    | (0010,1040) | LO |     | x |   |     |
| Patient's Telephone Numbers          | (0010,2154) | SH |     | x |   |     |

The above table should be read as follows:

- Module Name: The name of the associated module for supported worklist attributes.
- Attribute Name: Attributes supported to build the MWM SCU AE Worklist Request Identifier.
- Tag: DICOM tag for this attribute.
- VR: DICOM VR for this attribute.
- M: Matching keys for (automatic) Worklist Update.  
S: Single Value Matching  
R: Range Matching
- R: Return keys. An "x" will indicate that the MWM SCU AE will supply this attribute as Return Key with zero length for Universal Matching.
- D: Displayed keys. An "x" indicates that this worklist attribute is displayed to the user during a patient registration. For example, Patient Name will be displayed when registering the patient prior to an examination.
- IOD: An "x" indicates that this worklist attribute is included into all Object Instances created during performance of the related Procedure Step.

Notes: Specific Character Set (0008,0005) will be created if an extended or replacement character set is used in the matching keys.

Lists all schedules regardless of modality, but study always starts as US.

#### **4.2.3.4 Association Acceptance Policy**

The MWM SCU AE does not accept associations.



## **4.3 NETWORK INTERFACES**

### **4.3.1 Physical Network Interface**

The application is indifferent to the physical medium over which TCP/IP executes; which is dependent on the underlying operating system and hardware.

### **4.3.2 Additional Protocols**

When host names rather than IP addresses are used in the configuration properties to specify presentation addresses for remote AEs, the application is dependent on the name resolution mechanism of the underlying operating system.

### **4.3.3 IPv4 and IPv6 Support**

This product only supports IPv4 connections.

## 4.4 CONFIGURATION

### 4.4.1 AE Title/Presentation Address Mapping

#### 4.4.1.1 Local AE Titles

All local applications use the AE Titles and TCP/IP Ports configured via the service tool.

**Table 4.4-1**  
**AE TITLE CONFIGURATION TABLE**

| Application Entity | Default AE Title | Default TCP/IP Port |
|--------------------|------------------|---------------------|
| MWM SCU            | <User Specified> | Not Applicable      |
| Verification SCU   | <User Specified> |                     |
| Storage SCU        | <User Specified> |                     |

The default character repertoire excluding the highlighted characters can be used for the AE Titles:

**Table 4.4-2**  
**AE TITLE CHARACTER REPERTOIRE**

|      | 0  | 1 | 2 | 3 | 4  | 5 | 6 | 7 | 8 | 9 | A  | B   | C  | D  | E | F |
|------|----|---|---|---|----|---|---|---|---|---|----|-----|----|----|---|---|
| 0x00 |    |   |   |   |    |   |   |   |   |   | LF |     | FF | CR |   |   |
| 0x10 |    |   |   |   |    |   |   |   |   |   |    | ESC |    |    |   |   |
| 0x20 | SP | ! | " | # | \$ | % | & | ' | ( | ) | *  | +   | ,  | -  | . | / |
| 0x30 | 0  | 1 | 2 | 3 | 4  | 5 | 6 | 7 | 8 | 9 | :  | ;   | <  | =  | > | ? |
| 0x40 | @  | A | B | C | D  | E | F | G | H | I | J  | K   | L  | M  | N | O |
| 0x50 | P  | Q | R | S | T  | U | V | W | X | Y | Z  | [   | ¥  | ]  | ^ | _ |
| 0x60 | `  | a | b | c | d  | e | f | g | h | i | j  | k   | l  | m  | n | o |
| 0x70 | p  | q | r | s | t  | u | v | w | x | y | z  | {   |    | }  | ~ |   |

#### 4.4.1.2 Remote AE Title/Presentation Address Mapping

The AE Title, host name, port number and secured connection of each remote applications are configured using the service tool. The character repertoire of the AE Titles is listed in Table 4.4-2.

#### 4.4.2 Parameters

A large number of parameters related to acquisition and general operation can be configured using the service tool. The table below only shows those configuration parameters relevant to DICOM communication. See the Product's Service Manual for details on general configuration capabilities.

**Table 4.4-3**  
**CONFIGURATION PARAMETERS TABLE**

| Parameter                                                                                                      | Configurable<br>(Yes/No)<br>[Range] | Default Value                                       |
|----------------------------------------------------------------------------------------------------------------|-------------------------------------|-----------------------------------------------------|
| <b>General Parameters</b>                                                                                      |                                     |                                                     |
| Maximum PDU send/receive size                                                                                  | No                                  | 16384 bytes                                         |
| Time-out waiting for an acceptance or rejection response to an association request (Application Level Timeout) | No                                  | 300 sec                                             |
| Time-out waiting for a response to an association release request (Application Level Timeout)                  | No                                  | 300 sec                                             |
| Time-out waiting for completion of a TCP/IP connect request (Low-level timeout)                                | Yes<br>[1-9999]                     | 300 sec                                             |
| Time-out awaiting a response to a DIMSE request (Low-Level Timeout)                                            | No                                  | 300 sec                                             |
| Time-out for waiting for data between TCP/IP-packets (Low Level Timeout)                                       | Yes<br>[1-9999]                     | 300 sec                                             |
| <b>Storage SCU Parameters</b>                                                                                  |                                     |                                                     |
| Maximum number of simultaneously initiated associations by the Storage SCU AE                                  | No                                  | 10                                                  |
| Number of times a failed send job may be retried                                                               | No                                  | Forever, until the job succeeds or user cancels it. |
| <b>Modality Worklist SCU Parameters</b>                                                                        |                                     |                                                     |
| Maximum number of simultaneously initiated associations by the MWM SCU AE                                      | No                                  | 1                                                   |
| Maximum number of worklist items                                                                               | No                                  | 50                                                  |
| Query worklist for specific Scheduled Station AE Title                                                         | No                                  |                                                     |
| Query worklist for specific Modality                                                                           | No                                  | US                                                  |

## **5. MEDIA INTERCHANGE**

This product does not support Media Storage

Although this product does not support profiles such as AUG-GEN-CD, STD-GEN-DVD-JPEG, and STD-GEN-USB-JPEG, this product has the function to export in DICOM File format (PS 3.10).

## 6. SUPPORT OF CHARACTER SETS

This product supports the following character sets:

- ISO-IR 100 (Latin alphabet No.1)      Supplementary set of ISO 8859
- ISO-IR 144 (Cyrillic)                      Supplementary set of ISO 8859
- ISO-IR 192                                      Unicode in UTF-8
- GB18030                                        GB18030

Notes: If the MWM SCU AE receives worklist items that contain characters from unsupported character sets, it may abort the association using A-ABORT.

## 7. SECURITY

This product does not support any specific security measures.

It is assumed that the product is used within a secured environment. It is assumed that a secured environment includes at a minimum:

- a. Firewall or router protections to ensure that only approved external hosts have network access to the product.
- b. Firewall or router protections to ensure that the product only has network access to approved external hosts and services.
- c. Any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g. such as a Virtual Private Network (VPN))

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

## 8. ANNEXES

### 8.1 IOD CONTENTS

#### 8.1.1 Created SOP Instances

Table 8.1-1 specifies the attributes of an Ultrasound Image transmitted by the Storage SCU AE.

Table 8.1-2 specifies the attributes of an Ultrasound Multi-frame Image transmitted by the Storage SCU AE.

The following tables use a number of abbreviations. The abbreviations used in the "Presence of ..." column are:

|             |                                                                              |
|-------------|------------------------------------------------------------------------------|
| VNAP        | Value Not Always Present (attribute sent zero length if no value is present) |
| ANAP        | Attribute Not Always Present                                                 |
| ALWAYS      | Always Present                                                               |
| EMPTY       | Attribute is sent without a value                                            |
| Not Present | All attributes in this module are not present                                |

The abbreviations used in the "Source" column:

|        |                                                        |
|--------|--------------------------------------------------------|
| MWL    | the attribute value source is Modality Worklist        |
| USER   | the attribute value source is from user input          |
| AUTO   | the attribute value is generated automatically         |
| CONFIG | the attribute value source is a configurable parameter |

**8.1.1.1 US Image IOD**

**Table 8.1-1**  
**IOD OF CREATED US IMAGE SOP INSTANCES**

| <b>IE</b>          | <b>Module</b>              | <b>Reference</b> | <b>Presence of Module</b> |
|--------------------|----------------------------|------------------|---------------------------|
| Patient            | Patient                    | Table 8.1-3      | ALWAYS                    |
|                    | Clinical Trial Subject     | --               | Not Present               |
| Study              | General Study              | Table 8.1-4      | ALWAYS                    |
|                    | Patient Study              | Table 8.1-5      | ALWAYS                    |
|                    | Clinical Trial Study       | --               | Not Present               |
| Series             | General Series             | Table 8.1-6      | ALWAYS                    |
|                    | Clinical Trial Series      | --               | Not Present               |
| Frame of Reference | Frame of Reference         | --               | Not Present               |
|                    | Synchronization            | --               | Not Present               |
| Equipment          | General Equipment          | Table 8.1-7      | ALWAYS                    |
| Image              | General Image              | Table 8.1-8      | ALWAYS                    |
|                    | Image Pixel                | Table 8.1-9      | ALWAYS                    |
|                    | Contrast/bolus             | --               | Not Present               |
|                    | Palette Color Lookup Table | --               | Not Present               |
|                    | US Region Calibration      | Table 8.1-10     | ALWAYS                    |
|                    | US Image                   | Table 8.1-12     | ALWAYS                    |
|                    | Overlay Plane              | --               | Not Present               |
|                    | VOI LUT                    | --               | Not Present               |
|                    | SOP Common                 | Table 8.1-11     | ALWAYS                    |
|                    | Private Application        | --               | Not Present               |



**8.1.1.2 US Multi-frame Image IOD**

**Table 8.1-2**  
**IOD OF CREATED US MULTI-FRAME IMAGE SOP INSTANCES**

| <b>IE</b>          | <b>Module</b>              | <b>Reference</b> | <b>Presence of Module</b> |
|--------------------|----------------------------|------------------|---------------------------|
| Patient            | Patient                    | Table 8.1-3      | ALWAYS                    |
|                    | Clinical Trial Subject     | --               | Not Present               |
| Study              | General Study              | Table 8.1-4      | ALWAYS                    |
|                    | Patient Study              | Table 8.1-5      | ALWAYS                    |
|                    | Clinical Trial Study       | --               | Not Present               |
| Series             | General Series             | Table 8.1-6      | ALWAYS                    |
|                    | Clinical Trial Series      | --               | Not Present               |
| Frame of Reference | Frame of Reference         | --               | Not Present               |
|                    | Synchronization            | --               | Not Present               |
| Equipment          | General Equipment          | Table 8.1-7      | ALWAYS                    |
| Image              | General Image              | Table 8.1-8      | ALWAYS                    |
|                    | Image Pixel                | Table 8.1-9      | ALWAYS                    |
|                    | Contrast/bolus             | --               | Not Present               |
|                    | Cine                       | Table 8.1-13     | ALWAYS                    |
|                    | Multi-frame                | Table 8.1-14     | ALWAYS                    |
|                    | Frame Pointers             | --               | Not Present               |
|                    | Palette Color Lookup Table | --               | Not Present               |
|                    | US Region Calibration      | Table 8.1-10     | ALWAYS                    |
|                    | US Image                   | Table 8.1-15     | ALWAYS                    |
|                    | VOI LUT                    | --               | Not Present               |
|                    | SOP Common                 | Table 8.1-11     | ALWAYS                    |
|                    | Private Application        | --               | Not Present               |

### 8.1.1.3 Common Modules

**Table 8.1-3**  
**PATIENT MODULE OF CREATED SOP INSTANCES**

| Attribute Name           | Tag         | VR | Value                                             | Presence of Value | Source       |
|--------------------------|-------------|----|---------------------------------------------------|-------------------|--------------|
| Patient's Name           | (0010,0010) | PN |                                                   | VNAP              | MWL/<br>USER |
| Patient ID               | (0010,0020) | LO |                                                   | ALWAYS            | MWL/<br>USER |
| Patient's Birth Date     | (0010,0030) | DA |                                                   | VNAP              | MWL/<br>USER |
| Patient's Sex            | (0010,0040) | CS |                                                   | VNAP              | MWL/<br>USER |
| Patient Identity Removed | (0012,0062) | CS | Yes or NO                                         | ALWAYS            | USER         |
| De-identification Method | (0012,0063) | CS | DCM:113100/113105/113107/<br>113108/113109/113111 | ANAP              | USER         |

**Table 8.1-4**  
**GENERAL STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name             | Tag         | VR | Value | Presence of Value | Source |
|----------------------------|-------------|----|-------|-------------------|--------|
| Study Instance UID         | (0020,000D) | UI |       | ALWAYS            | AUTO   |
| Study Date                 | (0008,0020) | DA |       | ALWAYS            | AUTO   |
| Study Time                 | (0008,0030) | TM |       | ALWAYS            | AUTO   |
| Referring Physician's Name | (0008,0090) | PN |       | EMPTY             | USER   |
| Study ID                   | (0020,0010) | SH |       | ALWAYS            | AUTO   |
| Accession Number           | (0008,0050) | SH |       | VNAP              | USER   |
| Study Description          | (0008,1030) | LO |       | ALWAYS            | AUTO   |

**Table 8.1-5**  
**PATIENT STUDY MODULE OF CREATED SOP INSTANCES**

| Attribute Name   | Tag         | VR | Value | Presence of Value | Source       |
|------------------|-------------|----|-------|-------------------|--------------|
| Patient's Age    | (0010,1010) | AS |       | VNAP              | AUTO         |
| Patient's Size   | (0010,1020) | DS |       | VNAP              | MWL/<br>USER |
| Patient's Weight | (0010,1030) | DS |       | VNAP              | MWL/<br>USER |

**Table 8.1-6**  
**GENERAL SERIES MODULE OF CREATED SOP INSTANCES**

| Attribute Name      | Tag         | VR | Value | Presence of Value | Source |
|---------------------|-------------|----|-------|-------------------|--------|
| Modality            | (0008,0060) | CS | US    | ALWAYS            | AUTO   |
| Series Instance UID | (0020,000E) | UI |       | ALWAYS            | AUTO   |
| Series Number       | (0020,0011) | IS |       | ALWAYS            | AUTO   |
| Series Date         | (0008,0021) | DA |       | ALWAYS            | AUTO   |

|               |             |    |                                                                                                            |        |      |
|---------------|-------------|----|------------------------------------------------------------------------------------------------------------|--------|------|
| Series Time   | (0008,0031) | TM |                                                                                                            | ALWAYS | AUTO |
| Protocol Name | (0018,1030) | LO | Use the BodyMark category set in the preset.<br>= CARDIAC, GYN, GENERAL, ABD, OB, SMALL PARTS, URO, VESSEL | ALWAYS | AUTO |

**Table 8.1-7  
GENERAL EQUIPMENT MODULE OF CREATED SOP INSTANCES**

| Attribute Name            | Tag         | VR | Value                             | Presence of Value | Source |
|---------------------------|-------------|----|-----------------------------------|-------------------|--------|
| Manufacturer              | (0008,0070) | LO | CANON_MEC                         | ALWAYS            | AUTO   |
| Institution Name          | (0008,0080) | LO |                                   | ANAP              | AUTO   |
| Manufacturer's Model Name | (0008,1090) | LO | CUS-AAR00                         | ALWAYS            | AUTO   |
| Device Serial Number      | (0018,1000) | LO |                                   | ALWAYS            | AUTO   |
| Software Version          | (0018,1020) | LO | V1.0 SPxxxx*<br>(*4 digit number) | ALWAYS            | AUTO   |

**Table 8.1-8  
GENERAL IMAGE MODULE OF CREATED SOP INSTANCES**

| Attribute Name      | Tag         | VR | Value                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Presence of Value | Source |
|---------------------|-------------|----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|--------|
| Instance Number     | (0020,0013) | IS |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | ALWAYS            | AUTO   |
| Patient Orientation | (0020,0020) | CS |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | EMPTY             | AUTO   |
| Content Date        | (0008,0023) | DA |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | ALWAYS            | AUTO   |
| Content Time        | (0008,0033) | TM |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | ALWAYS            | AUTO   |
| Image Type          | (0008,0008) | CS | Value 1: Pixel Data Characteristics<br>"ORIGINAL" or "DERIVED"<br>Value 2: Patient Examination Characteristics<br>"PRIMARY" or "SECONDARY"<br>Value 3: use the BodyMark category set in the preset.<br>= CARDIAC, GYN, GENERAL, ABD, OB, SMALL PARTS, URO, VESSEL<br>Value 4: The following four hexadecimal digits, depending on the image mode. Do not set it for Secondary Capture images.<br>0001 2D Imaging<br>0002 M-Mode<br>0008 PW Doppler<br>0010 Color Doppler<br>0100 Color Power | ALWAYS            | AUTO   |

|                      |             |    |           |        |      |
|----------------------|-------------|----|-----------|--------|------|
| Acquisition Date     | (0008,0022) | DA |           | ALWAYS | AUTO |
| Acquisition Time     | (0008,0032) | TM |           | ALWAYS | AUTO |
| Image Comments       | (0020,4000) | LT |           | ANAP   | USER |
| Burned In Annotation | (0028,0301) | CS | YES or NO | ALWAYS | AUTO |

**Table 8.1-9**  
**IMAGE PIXEL MODULE OF CREATED SOP INSTANCES**

| Attribute Name             | Tag         | VR       | Value               | Presence of Value | Source |
|----------------------------|-------------|----------|---------------------|-------------------|--------|
| Samples per Pixel          | (0028,0002) | US       | 1 or 3              | ALWAYS            | AUTO   |
| Photometric Interpretation | (0028,0004) | CS       | RGB or YBR_FULL_422 | ALWAYS            | CONFIG |
| Planar Configuration       | (0028,0006) | US       | 0 or 1              | ANAP              | AUTO   |
| Rows                       | (0028,0010) | US       | [Rows], [Columns]   | ALWAYS            | AUTO   |
| Columns                    | (0028,0011) | US       | 1600, 1600          | ALWAYS            | AUTO   |
| Bits Allocated             | (0028,0100) | US       | 8                   | ALWAYS            | AUTO   |
| Bits Stored                | (0028,0101) | US       | 8                   | ALWAYS            | AUTO   |
| High Bit                   | (0028,0102) | US       | 7                   | ALWAYS            | AUTO   |
| Pixel Representation       | (0028,0103) | US       | 0                   | ALWAYS            | AUTO   |
| Pixel Data                 | (7FE0,0010) | OB or OW |                     | EMPTY             | AUTO   |

**Table 8.1-10**  
**US REGION CALIBRATION MODULE OF CREATED SOP INSTANCES**

| Attribute Name                 | Tag         | VR | Value | Presence of Value | Source |
|--------------------------------|-------------|----|-------|-------------------|--------|
| Sequence of Ultrasound Regions | (0018,6011) | SQ |       | ANAP              | AUTO   |
| >Region Spatial Format         | (0018,6012) | US |       | ALWAYS            | AUTO   |
| >Region Data Type              | (0018,6014) | US |       | ALWAYS            | AUTO   |
| >Region Flags                  | (0018,6016) | UL |       | ALWAYS            | AUTO   |
| >Region Location Min x0        | (0018,6018) | UL |       | ALWAYS            | AUTO   |
| >Region Location Min y0        | (0018,601A) | UL |       | ALWAYS            | AUTO   |
| >Region Location Max x1        | (0018,601C) | UL |       | ALWAYS            | AUTO   |
| >Region Location Max y1        | (0018,601E) | UL |       | ALWAYS            | AUTO   |
| >Physical Units X Direction    | (0018,6024) | US |       | ALWAYS            | AUTO   |
| >Physical Units Y Direction    | (0018,6026) | US |       | ALWAYS            | AUTO   |
| >Physical Delta X              | (0018,602C) | FD |       | ALWAYS            | AUTO   |
| >Physical Delta Y              | (0018,602E) | FD |       | ALWAYS            | AUTO   |

**Table 8.1-11**  
**SOP COMMON MODULE OF CREATED SOP INSTANCES**

| Attribute Name         | Tag         | VR | Value         | Presence of Value | Source |
|------------------------|-------------|----|---------------|-------------------|--------|
| Specific Character Set | (0008,0005) | CS | See Section 6 | ALWAYS            | AUTO   |
| Instance Creation Date | (0008,0012) | DA |               | ALWAYS            | AUTO   |
| Instance Creation Time | (0008,0013) | TM |               | ALWAYS            | AUTO   |

|                      |             |    |                                                                                                        |        |      |
|----------------------|-------------|----|--------------------------------------------------------------------------------------------------------|--------|------|
| Instance Creator UID | (0008,0014) | UI | 1.2.392.200036.9116.6.39.xx<br>xxxxxx* (*8 digit number)                                               | ALWAYS | AUTO |
| SOP Class UID        | (0008,0016) | UI | 1.2.840.10008.5.1.4.1.1.6.1<br>for US Image<br>1.2.840.10008.5.1.4.1.1.3.1<br>for US Multi-frame Image | ALWAYS | AUTO |
| SOP Instance UID     | (0008,0018) | UI |                                                                                                        | ALWAYS | AUTO |

#### 8.1.1.4 US Image Modules

**Table 8.1-12**  
**US IMAGE MODULE OF CREATED US IMAGE SOP INSTANCES**

| Attribute Name             | Tag         | VR             | Value                           | Presence of Value | Source |
|----------------------------|-------------|----------------|---------------------------------|-------------------|--------|
| Focus Depth                | (0018,5012) | DS             |                                 | ALWAYS            | AUTO   |
| Mechanical Index           | (0018,5022) | DS             |                                 | ALWAYS            | AUTO   |
| Transducer Type            | (0018,6031) | CS             |                                 | ALWAYS            | AUTO   |
| Samples Per Pixel          | (0028,0002) | US             | 1 or 3                          | ALWAYS            | AUTO   |
| Photometric Interpretation | (0028,0004) | CS             | RGB or YBR_FULL_422             | ALWAYS            | CONFIG |
| Planar Configuration       | (0028,0006) | US             | 0 or 1                          | ALWAYS            | AUTO   |
| Rows                       | (0028,0010) | US             | [Rows], [Columns]<br>1600, 1600 | ALWAYS            | AUTO   |
| Columns                    | (0028,0011) | US             |                                 | ALWAYS            | AUTO   |
| Bits Allocated             | (0028,0100) | US             | 8                               | ALWAYS            | AUTO   |
| Bits Stored                | (0028,0101) | US             | 8                               | ALWAYS            | AUTO   |
| High Bit                   | (0028,0102) | US             | 7                               | ALWAYS            | AUTO   |
| Pixel Representation       | (0028,0103) | US             | 0                               | ALWAYS            | AUTO   |
| Pixel Data                 | (7FE0,0010) | OB<br>or<br>OW |                                 | EMPTY             | AUTO   |

#### 8.1.1.5 US Multi-frame Image Modules

**Table 8.1-13**  
**CINE MODULE OF CREATED US MULTI-FRAME IMAGE SOP INSTANCES**

| Attribute Name                 | Tag         | VR | Value | Presence of Value | Source |
|--------------------------------|-------------|----|-------|-------------------|--------|
| Start Trim                     | (0008,2142) | IS |       | ALWAYS            | AUTO   |
| Stop Trim                      | (0008,2143) | IS |       | ALWAYS            | AUTO   |
| Recommended Display Frame Rate | (0008,2144) | IS |       | ALWAYS            | USER   |
| Cine Rate                      | (0018,0040) | IS |       | ALWAYS            | USER   |
| Effective Duration             | (0018,0072) | DS |       | ALWAYS            | AUTO   |
| Frame Time                     | (0018,1063) | DS |       | ALWAYS            | AUTO   |
| Frame Delay                    | (0018,1066) | DS |       | ALWAYS            | AUTO   |
| Actual Frame Duration          | (0018,1242) | IS |       | ALWAYS            | AUTO   |

**Table 8.1-14**  
**MULTI-FRAME MODULE OF CREATED US MULTI-FRAME IMAGE SOP INSTANCES**

| Attribute Name          | Tag         | VR | Value       | Presence of Value | Source |
|-------------------------|-------------|----|-------------|-------------------|--------|
| Number of Frames        | (0028,0008) | IS |             | ALWAYS            | AUTO   |
| Frame Increment Pointer | (0028,0009) | AT | <0018,1063> | ALWAYS            | AUTO   |

**Table 8.1-15**  
**US IMAGE MODULE OF CREATED US MULTI-FRAME IMAGE SOP INSTANCES**

| Attribute Name             | Tag         | VR | Value                           | Presence of Value | Source |
|----------------------------|-------------|----|---------------------------------|-------------------|--------|
| Focus Depth                | (0018,5012) | DS |                                 | ALWAYS            | AUTO   |
| Mechanical Index           | (0018,5022) | DS |                                 | ALWAYS            | AUTO   |
| Transducer Type            | (0018,6031) | CS |                                 | ALWAYS            | AUTO   |
| Samples per Pixel          | (0028,0002) | US | 1 or 3                          | ALWAYS            | AUTO   |
| Photometric Interpretation | (0028,0004) | CS | RGB or YBR_FULL_422             | ALWAYS            | AUTO   |
| Planar Configuration       | (0028,0006) | US | 0 or 1                          | ALWAYS            | AUTO   |
| Rows                       | (0028,0010) | US | [Rows], [Columns]<br>1600, 1600 | ALWAYS            | AUTO   |
| Columns                    | (0028,0011) | US |                                 | ALWAYS            | AUTO   |
| Bits Allocated             | (0028,0100) | US | 8                               | ALWAYS            | AUTO   |
| Bits Stored                | (0028,0101) | US | 8                               | ALWAYS            | AUTO   |
| High Bit                   | (0028,0102) | US | 7                               | ALWAYS            | AUTO   |
| Pixel Representation       | (0028,0103) | US | 0                               | ALWAYS            | AUTO   |
| Pixel Data                 | (7FE0,0010) | OB |                                 | EMPTY             | AUTO   |

### 8.1.1.6 Other Modules

The tables below show the attributes that extend the standard IODs of US Image and US Multi-frame Image.

**Table 8.1-16**  
**PATIENT DEMOGRAPHIC MODULE OF CREATED SOP INSTANCES**

| Attribute Name              | Tag         | VR | Value | Presence of Value | Source |
|-----------------------------|-------------|----|-------|-------------------|--------|
| Patient's Address           | (0010,1040) | LO |       | VNAP              | AUTO   |
| Patient's Telephone Numbers | (0010,2154) | SH |       | VNAP              | AUTO   |

### 8.1.2 Usage of Attributes from received IOD's

Not applicable.

### 8.1.3 Attribute Mapping

The tables below show the relationships between attributes received via Modality Worklist, stored in acquired images.

The cell content conventions should be read as follows:

Copy:

The value will be copied from a corresponding source attribute of another DICOM object, as defined by the table column.

Copy from: <DICOM attribute>: The source as specified in the referenced DICOM attribute will be used instead of using the DICOM attribute of the same row as the source.

Equal (internally generated): The value will be internally generated which may be used in more than one DICOM object.

**Table 8.1-17**

**SCHEDULED CASE - ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST AND IMAGE**

| Attribute Name   | Tag         | Modality Worklist | Image IOD |
|------------------|-------------|-------------------|-----------|
| Accession Number | (0008,0050) | Source            | Copy      |

#### **8.1.4 Coerced/Modified Fields**

Not applicable.



## **8.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES**

Not applicable.

## **8.3 CODED TERMINOLOGY AND TEMPLATES**

Not applicable.

## **8.4 GRAYSCALE IMAGE CONSISTENCY**

Not applicable.

## **8.5 STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES**

Not applicable.

## **8.6 PRIVATE TRANSFER SYNTAXES**

Not applicable.

## **8.7 STANDARD EXTENDED AND PRIVATE TEMPLATES**

Not applicable.

## **8.8 DICOM Security Profile Details**

Not applicable.