

**DICOM CONFORMANCE STATEMENT
STORAGE SCU & Q/R SCP & PRINT SCU
FOR
TOSHIBA SUPERCONDUCTING MRI SYSTEMS
FLEXART-series/VISART-series/
EXCELART-series/SECOND CONSOLE
(MIIMR0001EAA)**

TOSHIBA CORPORATION

© TOSHIBA CORPORATION 2000
ALL RIGHTS RESERVED

IMPORTANT !

- (1) No part of this manual may be copied or reprinted, in whole or in part, without written permission.
- (2) The contents of this manual are subject to change without prior notice and without our legal obligation.

Table of Contents

| | |
|------------------------------------------------------------------------------------|----------|
| 1 Introduction | 1 |
| 1.1 References | 1 |
| 1.2 Definitions | 1 |
| 1.3 Acronyms, Abbreviations and Symbols | 2 |
| 2 Implementation Model | 3 |
| 2.1 Application Data Flow Diagram | 3 |
| 2.2 Functional Definitions of AE's | 4 |
| 2.2.1 Export AE | 4 |
| 2.2.2 Import AE | 4 |
| 2.3 Sequencing of Real World Activities | 5 |
| 2.3.1 Features | 5 |
| 2.3.1.1 Manual send of image or study | 5 |
| 2.3.1.2 Manual Request of Print | 5 |
| 2.3.1.3 Response of Query/Retrieve | 5 |
| 2.3.2 Operation | 6 |
| 2.3.2.1 Manual send of image or study | 6 |
| 2.3.2.2 Manual Request of Print | 6 |
| 3 AE Specifications | 7 |
| 3.1 Export Specification | 7 |
| 3.1.1 Export Association Establishment Policies | 8 |
| 3.1.1.1 Export General | 8 |
| 3.1.1.2 Export Number of Associations | 8 |
| 3.1.1.3 Export Asynchronous Nature | 8 |
| 3.1.1.4 Export Implementation Identifying Information | 8 |
| 3.1.2 Export Association Initiation by Real-World Activity | 9 |
| 3.1.2.1 Export Real-World Activity - Storage | 9 |
| 3.1.2.1.1 Export Associated Real-World Activity - Storage | 9 |
| 3.1.2.1.2 Export Proposed Presentation Contexts - Storage | 9 |
| 3.1.2.2 Export Real-World Activity - Print | 10 |
| 3.1.2.2.1 Export Associated Real-World Activity - Print | 10 |
| 3.1.2.2.2 Export Proposed Presentation Contexts - Print | 10 |
| 3.1.3 Export Association Acceptance Policy | 10 |
| 3.2 Import Specification | 11 |
| 3.2.1 Import Association Establishment Policies | 11 |
| 3.2.1.1 Import General | 11 |
| 3.2.1.2 Import Number of Associations | 11 |
| 3.2.1.3 Import Asynchronous Nature | 11 |
| 3.2.1.4 Import Implementation Identifying Information | 11 |
| 3.2.2 Import Association Initiation by Real-World Activity | 12 |
| 3.2.3 Import Association Acceptance Policy | 12 |
| 3.2.3.1 Import Real-World Activity - Query/Retrieve (Find) | 12 |
| 3.2.3.1.1 Import Associated Real-World Activity - Query/Retrieve (Find) | 12 |
| 3.2.3.1.2 Import Presentation Context Table - Query/Retrieve (Find) | 12 |
| 3.2.3.1.3 Import Presentation Context Acceptance Criterion - Query/Retrieve (Find) | 13 |
| 3.2.3.1.4 Import Transfer Syntax Selection Policies - Query/Retrieve (Find) | 13 |
| 3.2.3.2 Import Real-World Activity - Query/Retrieve (Move) | 14 |
| 3.2.3.2.1 Import Associated Real-World Activity - Query/Retrieve (Move) | 14 |
| 3.2.3.2.2 Import Presentation Context Table - Query/Retrieve (Move) | 14 |
| 3.2.3.2.3 Import Presentation Context Acceptance Criterion - Query/Retrieve (Move) | 15 |

| | |
|-----------------------------------------------------------------------------|-----------|
| 3.2.3.2.4 Import Transfer Syntax Selection Policies - Query/Retrieve (Move) | 15 |
| 4 Communication Profiles | 16 |
| 4.1 Supported Communication Stacks | 16 |
| 4.2 OSI Stack | 16 |
| 4.3 TCP/IP Stack | 16 |
| 4.3.1 API | 16 |
| 4.3.2 Physical Media Support | 16 |
| 4.4 Point-to-Point Stack | 16 |
| 5 Extensions/Specializations/Privatizations | 17 |
| 6 Configuration | 18 |
| 6.1 AE Title/Presentation Address Mapping | 18 |
| 6.2 Configurable Parameters | 18 |
| 6.2.1 Time-out Value, Retry Count, Retry Interval | 18 |
| 6.3 Warning Status Criteria | 20 |
| 6.3.1 MR Image Storage | 20 |
| 6.3.1.1 C-STORE response | 20 |
| 6.3.2 Basic Grayscale Print Management | 20 |
| 6.3.2.1 Basic Film Session SOP Class | 20 |
| 6.3.2.1.1 N-CREATE response | 20 |
| 6.3.2.2 Basic Film Box SOP Class | 20 |
| 6.3.2.2.1 N-CREATE response | 20 |
| 6.3.2.2.2 N-ACTION response | 21 |
| 6.3.2.3 Basic Grayscale Image Box SOP Class | 21 |
| 6.3.2.3.1 N-SET response | 21 |
| 6.3.2.4 Printer SOP Class | 21 |
| 6.3.2.4.1 N-GET response | 21 |
| 6.4 Implementation Information and Maximum Reception PDU Size | 22 |
| 7 Support of Extended Character Sets | 23 |
| 8 MR Information Object Definition | 24 |
| 8.1 Entity Module Definitions | 24 |
| 8.1.1 MR IOD Modules | 24 |
| 8.2 Information Object Definitions | 24 |
| 8.2.1 Patient Module | 24 |
| 8.2.2 General Study Module | 25 |
| 8.2.3 Patient Study Module | 25 |
| 8.2.4 General Series Module | 26 |
| 8.2.5 Frame of Reference Module | 26 |
| 8.2.6 General Equipment Module | 27 |
| 8.2.7 General Image Module | 27 |
| 8.2.8 Image Plane Module | 28 |
| 8.2.9 Image Pixel Module | 28 |
| 8.2.10 Contrast/Bolus Module | 29 |
| 8.2.11 MR Image Module | 29 |
| 8.2.12 VOI LUT Module | 30 |
| 8.2.13 SOP Common Module | 30 |
| 9 DIMSE-Service and Attributes | 31 |
| 9.1 DIMSE-Services | 31 |

| | |
|----------------------------------------------------|-----------|
| 9.2 Basic Film Session SOP Class | 32 |
| 9.2.1 N-CREATE Attributes | 32 |
| 9.3 Basic Film Box SOP Class | 32 |
| 9.3.1 N-CREATE Attributes | 32 |
| 9.4 Basic Grayscale Image Box SOP Class | 33 |
| 9.4.1 N-SET Attributes | 33 |
| 9.5 Printer SOP Class | 34 |
| 9.5.1 N-EVENT-REPORT | 34 |
| 9.5.2 N-GET Attributes | 34 |
| | |
| 10 Search Keys | 35 |
| 10.1 Query/Retrieve SCP (C-FIND) | 35 |
| 10.1.1 Patient Root Q/R Information Model - FIND | 35 |
| 10.1.1.1 Patient Level | 35 |
| 10.1.1.2 Study Level | 35 |
| 10.1.1.3 Series Level | 35 |
| 10.1.1.4 Image Level | 35 |
| 10.1.2 Study Root Q/R Information Model - FIND | 36 |
| 10.1.2.1 Study Level | 36 |
| 10.1.2.2 Series Level | 36 |
| 10.1.2.3 Image Level | 36 |
| 10.1.3 Patient/Study Only Information Model - FIND | 36 |
| 10.1.3.1 Patient Level | 36 |
| 10.1.3.2 Study Level | 36 |

1 Introduction

This document is a DICOM Conformance Statement for Toshiba's Superconducting MRI Systems(FLEXART-series/VISART-series/EXCELART-series/SECOND CONSOLE). It is intended to provide the reader with the knowledge of how to integrate this product within a DICOM compliant hospital network. It details the DICOM Service Classes, Information Objects, and Communication Protocols which are supported by this product.

If the reader is unfamiliar with DICOM, it is recommended that they read the DICOM Specification(referenced below) prior to reading this conformance statement. Also note that this document is formatted according to the DICOM Specification, Part 2:Conformance.

1.1 References

- ACR-NEMA Digital Imaging and Communications in Medicine, DICOM V3.0.

1.2 Definitions

- **Association Establishment** - An Association Establishment is the first phase of communication between two DICOM Application Entities. The AEs use the Association Establishment to negotiate how data will be encoded and the type of data to be exchanged.
- **Called Application Entity Title** - The Called AE Title defines the intended receiver of an Association.
- **Calling Application Entity Title** - The Calling AE Title defines the requestor of an Association.
- **DICOM Message Service Element (DIMSE)** - A DIMSE defines the services and protocols utilized by an Application Entity to exchange messages.
- **Information Object Definition (IOD)** - An IOD is a data model which is an abstraction of real-world information. This data model defines the nature and attributes relevant to the class of real-world objects represented.
- **Service Class Provider (SCP)** - A Service Class Provider plays the "server" role to perform operations and invoke notifications during an Association. An example of a Storage Service Class Provider would be an image storage device. In this case, the image storage device is storing the image that was sent by a Service Class User.
- **Service Class User (SCU)** - A Service Class User plays the "client" role to invoke operations and perform notifications during an Association. An example of a Storage Service Class User would be an image acquisition device. In this case, the image acquisition device will create and send a DICOM image by requesting that a Service Class Provider store that image.
- **Service/Object Pair (SOP) Class** - A SOP Class is defined by the union of an Information Object Definition and a set of DIMSE Services. A DICOM Application Entity may support one or more SOP Classes. Each SOP Class is uniquely identified by a SOP Class UID.
- **SOP Instance** - A specific occurrence of a Information Object.
- **Transfer Syntax** - The Transfer Syntax is a set of encoding rules that allow DICOM Application Entities to negotiate the encoding techniques (e.g. data element structure, byte ordering, compression) they are able to support. The Transfer Syntax is negotiated during Association Negotiation.
- **Unique Identifier (UID)** - A Unique Identifier is a globally unique, ISO compliant, ASCII-numeric string. It guarantees uniqueness across multiple countries, sites, vendors and equipment.

1.3 Acronyms, Abbreviations and Symbols

- ACC American College of Cardiology
- ACR American College of Radiology
- ASCII American Standard Code for Information Interchange
- AE Application Entity
- ANSI American National Standards Institute
- CEN TC251 Comite Europeen de Normalisation - Technical Committee 251 - Medical Informatics

- DICOM Digital Imaging and Communications in Medicine
- DIMSE DICOM Message Service Element
- DIMSE-C DICOM Message Service Element - Composite
- DIMSE-N DICOM Message Service Element - Normalized
- HIS Hospital Information System
- HL7 Health Level 7
- IE Information Entity
- IOD Information Object Definition
- ISO International Standards Organization
- JIRA Japan Industries Association of Radiological Systems
- NEMA National Electrical Manufacturers Association
- OSI Open Systems Interconnection
- PDU Protocol Data Unit
- RIS Radiology Information System
- SCP Service Class Provider
- SCU Service Class User
- SOP Service-Object Pair
- TCP/IP Transmission Control Protocol/Internet Protocol
- UID Unique Identifier

2 Implementation Model

2.1 Application Data Flow Diagram




Figure 1

2.2 Functional Definitions of AE's

2.2.1 Export AE

Export AE is used to transmit images to a remote DICOM device. It therefore performs the following tasks:

- Builds DICOM MR Information Objects
- Establishes DICOM Association with remote DICOM device
- Performs storage of DICOM MR Information Objects to remote DICOM device

Export AE is used to transmit request for Print images to a remote DICOM device. It therefore performs the following tasks:

- Builds DICOM Basic Grayscale Print Objects
- Establishes DICOM Association with remote DICOM device
- Performs transmit of DICOM Basic Grayscale Print Objects to remote DICOM device

2.2.2 Import AE

Import AE is used to respond to requests of Query/Retrieve from a remote DICOM device.

2.3 Sequencing of Real World Activities

2.3.1 Features

2.3.1.1 Manual send of image or study

- Operator requests to send images after selecting the images to be transferred.
- When the transfer fails, Operator can manually attempt to resend the study at a later time.

2.3.1.2 Manual Request of Print

- The number of frames in the rows and columns on each film can be specified as desired, up to a total maximum of 48 frames per film.
- If an error occurs during printing, a request to retry printing is issued automatically.
- Print requests are placed on a queue, and are executed in the background.

2.3.1.3 Response of Query/Retrieve

- It returns the result of the search corresponding to the search request.
- Returns the requested images.

2.3.2 Operation

2.3.2.1 Manual send of image or study

The operation for manual image transferring is described below:

Step-1: Select the destination of image transfer.

Step-2: Select the image or study to be transferred.

Step-3: Request transfer.

2.3.2.2 Manual Request of Print

The operation for printing is described below:

Step-1: Enter the information for the film, and select each of the images to be printed.

Step-2: Execute the print request.

3 AE Specifications

3.1 Export Specification

Export AE provides Standard Conformance to the following DICOM SOP Classes as an SCU:

Table 1

| SOP Class Name | SOP Class UID |
|----------------------------------|---------------------------|
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 |
| Basic Grayscale Print Management | 1.2.840.10008.5.1.1.9 |

3.1.1 Export Association Establishment Policies

3.1.1.1 Export General

Export AE will utilize and understand the following Application Context Name:

Table 2

| | |
|--------------------------------|-----------------------|
| DICOM V3.0 Application Context | 1.2.840.10008.3.1.1.1 |
|--------------------------------|-----------------------|

Export AE supports a minimum PDU size of 16Kbytes and a maximum PDU size of 16Kbytes. The default value is set to 16Kbytes.

3.1.1.2 Export Number of Associations

Export AE can only establish one association at a time, independent of the number of destinations chosen.

3.1.1.3 Export Asynchronous Nature

Export AE allows a single outstanding operation on any association. Therefore, Export AE does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

3.1.1.4 Export Implementation Identifying Information

Export AE will specify the following Implementation Identifying Information:

- Implementation Class UID 1.2.392.200036.9116.4.1.12
- Implementation Version Name TM_MR_DCM_V1.2

3.1.2 Export Association Initiation by Real-World Activity

Export AE initiates an association when the following activity is chosen by the operator.

- "Manual send of image or study"
 - ◆ Storage - Create and store an MR image to a remote DICOM device
- "Manual Request of Print"
 - ◆ Print - Request print image to a remote DICOM device

Export AE also initiates an association when the following activity occurs.

- "Response of Query/Retrieve"
 - ◆ Storage - Create and store an MR image to a remote DICOM device

3.1.2.1 Export Real-World Activity - Storage

3.1.2.1.1 Export Associated Real-World Activity - Storage

Storage is executed by the MRI System after the operator's image transfer requests are queued.

3.1.2.1.2 Export Proposed Presentation Contexts - Storage

Export AE proposes the following Presentation Contexts shown below:

Table 3

| Presentation Context Table | | | | | |
|----------------------------|---------------------------|------------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.1.2.1.2.1 Export SOP Specific Conformance - MR Image Storage

- Export AE operation involves the following sequence of steps for each image transfer.
 - (1) Association establishment(requestor only)
 - (2) Data transfer(SCU only)
 - (3) Association release(requestor only)

Export AE judges that the transfer of one image succeeded when the result of (2) "Data transfer" is "Success" even if the result of (3) "Association release" is "Failure".

- MR Information Object Definition is described in chapter 8.

3.1.2.2 Export Real-World Activity - Print

3.1.2.2.1 Export Associated Real-World Activity - Print

Export AE performs DICOM printing to a destination device. If a communication or printing error occurs, the print operation is automatically retried several times.

3.1.2.2.2 Export Proposed Presentation Contexts - Print

Export AE proposes the following Presentation Contexts shown below:

Table 4

| Presentation Context Table | | | | | |
|----------------------------------|-----------------------|------------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Basic Grayscale Print Management | 1.2.840.10008.5.1.1.9 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.1.2.2.2.1 Export SOP Specific Conformance - Basic Grayscale Print Management

- Export AE operation involves the following sequence of steps for each request print image.
 - (1) Association establishment (requestor only)
 - (2) Print request (SCU only)
 - (3) Association release (requestor only)

Export AE judges that the request printing images succeeded when the result of (2) "Print request" is "Success" even if the result of (3) "Association release" is "Failure".

- DIMSE-Service and Attributes are described in chapter 9.

3.1.3 Export Association Acceptance Policy

Export AE does not accept any associations generated by remote applications.

3.2 Import Specification

Import AE provides Standard Conformance to the following DICOM SOP Classes as an SCP:

Table 5

| SOP Class Name | SOP Class UID |
|-------------------------------------------------------------|-----------------------------|
| Patient Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.1.1 |
| Study Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 |
| Patient/ Study Only Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.3.1 |
| Patient Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 |
| Study Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.2 |
| Patient/ Study Only Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.3.2 |

3.2.1 Import Association Establishment Policies

3.2.1.1 Import General

Import AE will utilize and understand the following Application Context Name:

Table 6

| | |
|--------------------------------|-----------------------|
| DICOM V3.0 Application Context | 1.2.840.10008.3.1.1.1 |
|--------------------------------|-----------------------|

Import AE supports a minimum PDU size of 16Kbytes and a maximum PDU size of 16Kbytes. The default value is set to 16Kbytes.

3.2.1.2 Import Number of Associations

Import AE can accept a maximum of five simultaneous associations.

3.2.1.3 Import Asynchronous Nature

Import AE allows a single outstanding operation on any association. Therefore, Import AE does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

3.2.1.4 Import Implementation Identifying Information

Import AE will specify the following Implementation Identifying Information:

- Implementation Class UID 1.2.392.200036.9116.4.1.12
- Implementation Version Name TM_MR_DCM_V1.2

3.2.2 Import Association Initiation by Real-World Activity

Import AE never initiates an association.

3.2.3 Import Association Acceptance Policy

Import AE accepts an association generated by remote applications.

- "Response of Query/Retrieve"
 - ◆ Query/Retrieve (Find) - Receive a search request from a remote DICOM device
 - ◆ Query/Retrieve (Move) - Receive a request for image transfer from a remote DICOM device

3.2.3.1 Import Real-World Activity - Query/Retrieve (Find)

3.2.3.1.1 Import Associated Real-World Activity - Query/Retrieve (Find)

When a request for a search is performed by a remote system, an association is accepted.

3.2.3.1.2 Import Presentation Context Table - Query/Retrieve (Find)

Import AE accepts the following Presentation Contexts shown below:

Table 7

| Presentation Context Table | | | | | |
|--------------------------------------------------|-----------------------------|---------------------------|-------------------|-------------|-----------------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Patient Root Q/R Information Model - FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| Study Root Q/R Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| Patient/ Study Only Q/R Information Model - FIND | 1.2.840.10008.5.1.4.1.2.3.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |

3.2.3.1.2.1 Import SOP Specific Conformance - Query/Retrieve (Find)

- The status is indicated by C-FIND-RSP as shown below:

Table 8

| Status | Meaning | Code |
|---------|---------------------------------------------------------------------------------------------------------------------------------------|-------|
| Cancel | Matching terminated due to Cancel request. | FE00H |
| Success | Matching is complete - No final identifier is supplied. | 0000H |
| Pending | Matches are continuing - Current Match is supplied and any Optional keys were supported in the same manner as Required keys. | FF00H |
| | Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this Identifier. | FF01H |

- If an error should occur, A-ABORT request is sent to a remote DICOM device.
- Search keys for Query/Retrieve SCP are described in chapter 10.1.

3.2.3.1.3 Import Presentation Context Acceptance Criterion - Query/Retrieve (Find)

Import AE accepts the Presentation Contexts listed in the Presentation Context Table.(Table 7)

3.2.3.1.4 Import Transfer Syntax Selection Policies - Query/Retrieve (Find)

Import AE supports only the Implicit VR Little Endian transfer syntax. It rejects any proposed Presentation Contexts which does not specify the default Implicit VR Little Endian transfer syntax.

3.2.3.2 Import Real-World Activity - Query/Retrieve (Move)**3.2.3.2.1 Import Associated Real-World Activity - Query/Retrieve (Move)**

When a request for the image transfer is made of a remote system, an association is accepted.

3.2.3.2.2 Import Presentation Context Table - Query/Retrieve (Move)

Import AE accepts the following Presentation Contexts shown below:

Table 9

| Presentation Context Table | | | | | |
|--------------------------------------------------|-----------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Patient Root Q/R Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| Study Root Q/R Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| Patient/ Study Only Q/R Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.3.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |

3.2.3.2.2.1 Import SOP Specific Conformance - Query/Retrieve (Move)

- The status is indicated by C-MOVE-RSP as shown below:

Table 10

| Status | Meaning | Code |
|---------|-----------------------------------------------------|-------|
| Refused | Out of Resources - Unable to perform Sub-operations | A702H |
| Success | Sub-operations complete - No Failures | 0000H |

- If an error should occur, A-ABORT request is sent to a remote DICOM device.

3.2.3.2.3 Import Presentation Context Acceptance Criterion - Query/Retrieve (Move)

Import AE accepts the Presentation Contexts listed in the Presentation Context Table.(Table 9)

3.2.3.2.4 Import Transfer Syntax Selection Policies - Query/Retrieve (Move)

Import AE supports only the Implicit VR Little Endian transfer syntax. It rejects any proposed Presentation Contexts which does not specify the default Implicit VR Little Endian transfer syntax.

4 Communication Profiles

4.1 Supported Communication Stacks

This product provides DICOM TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard.

4.2 OSI Stack

Not applicable to this product.

4.3 TCP/IP Stack

This product inherits its TCP/IP stack from the computer system upon which it executes.

4.3.1 API

Not applicable to this product.

4.3.2 Physical Media Support

This product is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the computer system upon which it executes.

4.4 Point-to-Point Stack

Not applicable to this product.

5 Extensions/Specializations/Privatizations

Not applicable to this product.

6 Configuration

For the MR Systems, the configuration can be set using the Online Setup interface.
 Note: Settings and changes are performed by Toshiba Service Personnel at the time of installation of the system.

6.1 AE Title/Presentation Address Mapping

Mapping from the AE titles to the presentation address are as follows:

- One port number and one AE title can be described for one host name.
- Each AE title is mapped to one port number.
- The MR Systems has following default values:

Local AE Title TM_MR_DCM_V1.0

"Response of Query/Retrieve"

Local Port No. 8500

6.2 Configurable Parameters

6.2.1 Time-out Value, Retry Count, Retry Interval

The time-out value, retry count , and retry interval in each status are shown below.




Figure 2

Table 11

| Status | Item | Time-out value | Retry count | Retry interval | Remarks |
|--------|-------------------------------------------------|-----------------------------------------|-------------|----------------|-------------------------------------------------------|
| S1 | Association establishment request waiting time | default:30 seconds range:1 to 10000 | Not set | Not set | Only one parameter can be set in Toshiba MRI Systems. |
| S2 | Association establishment response waiting time | default:30 seconds range:1 to 10000 | Not set | Not set | Only one parameter can be set in Toshiba MRI Systems. |
| S3 | Service request waiting time | default:30 seconds range:1 to 10000 | Not set | Not set | Only one parameter can be set in Toshiba MRI Systems. |
| S4 | Service response waiting time | default:300 seconds range:1 to 10000 | Not set | Not set | Only one parameter can be set in Toshiba MRI Systems. |
| S5 | Association release waiting time | default:30 seconds range:1 to 10000 | Not set | Not set | Only one parameter can be set in Toshiba MRI Systems. |

6.3 Warning Status Criteria

6.3.1 MR Image Storage

6.3.1.1 C-STORE response

If SUCCESS is set, this product judges that C-STORE request succeeded.

If FAIL is set, this product judges that C-STORE request failed.

Table 12

| Warning response | Default | Parameter setting range |
|-----------------------------------|---------|-------------------------|
| Coercion of Data Elements | FAIL | Not change |
| Data Set does not match SOP Class | FAIL | Not change |
| Elements discarded | FAIL | Not change |

6.3.2 Basic Grayscale Print Management

6.3.2.1 Basic Film Session SOP Class

6.3.2.1.1 N-CREATE response

If SUCCESS is set, this product judges that N-CREATE request succeeded.

If FAIL is set, this product judges that N-CREATE request failed.

Table 13

| Warning response | Default | Parameter setting range |
|---------------------------------|---------|-------------------------|
| Memory allocation not supported | FAIL | Not change |

6.3.2.2 Basic Film Box SOP Class

6.3.2.2.1 N-CREATE response

If SUCCESS is set, this product judges that N-CREATE request succeeded.

If FAIL is set, this product judges that N-CREATE request failed.

Table 14

| Warning response | Default | Parameter setting range |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|---------|-------------------------|
| Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. | FAIL | Not change |

6.3.2.2.2 N-ACTION response

If SUCCESS is set, this product judges that N-ACTION request succeeded.

If FAIL is set, this product judges that N-ACTION request failed.

Table 15

| Warning response | Default | Parameter setting range |
|--------------------------------------------------------------------------------------|---------|-------------------------|
| Film Box SOP Instance hierarchy does not contain Image Box SOP Instances(empty page) | FAIL | Not change |
| Image size is larger than image box size, the image has been demagnified. | FAIL | Not change |

6.3.2.3 Basic Grayscale Image Box SOP Class**6.3.2.3.1 N-SET response**

If SUCCESS is set, this product judges that N-SET request succeeded.

If FAIL is set, this product judges that N-SET request failed.

Table 16

| Warning response | Default | Parameter setting range |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|---------|-------------------------|
| Image size larger than image box size, the image has been demagnified. | FAIL | Not change |
| Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. | FAIL | Not change |

6.3.2.4 Printer SOP Class**6.3.2.4.1 N-GET response**

If SUCCESS is set, this product judges that N-GET request succeeded.

If FAIL is set, this product judges that N-GET request failed.

Table 17

| Warning response | Default | Parameter setting range |
|----------------------|---------|-------------------------|
| Attribute list error | FAIL | Not change |

6.4 Implementation Information and Maximum Reception PDU Size

The default values for the MRI System are used for the Implementation Class UID, the Implementation Version Name, and the Maximum length received. They cannot be changed.

Table 18

| Parameter | Default |
|--------------------------------------|----------------------------|
| Implementation Class UID | 1.2.392.200036.9116.4.1.12 |
| Implementation Version Name | TM_MR_DCM_V1.2 |
| Maximum length received (unit: byte) | 0x4000 |

7 Support of Extended Character Sets

This product supports the following character sets:

- ISO-IR 6 (default) ISO 646
- ISO-IR 87(Japanese) JIS X 0208(Kanji)

8 MR Information Object Definition

8.1 Entity Module Definitions

The information modules for the MR scanner are defined below.

8.1.1 MR IOD Modules

Table 19

| Information Entity | Module | Reference | Usage ^{*1} |
|--------------------|---------------------------|-----------|---------------------|
| Patient | Patient Module | 8.2.1 | M |
| Study | General Study Module | 8.2.2 | M |
| Study | Patient Study Module | 8.2.3 | U |
| Series | General Series Module | 8.2.4 | M |
| Frame of Reference | Frame of Reference Module | 8.2.5 | M |
| Equipment | General Equipment Module | 8.2.6 | M |
| Image | General Image Module | 8.2.7 | M |
| Image | Image Plane Module | 8.2.8 | M |
| Image | Image Pixel Module | 8.2.9 | M |
| Image | Contrast/bolus Module | 8.2.10 | C |
| Image | MR Image Module | 8.2.11 | M |
| Image | VOI LUT Module | 8.2.12 | U |
| Image | SOP Common Module | 8.2.13 | M |

*1:M=Mandatory, C=Conditional, U=User option

8.2 Information Object Definitions

8.2.1 Patient Module

Table 20

| Attribute Name | Tag | Type | Attribute Description |
|----------------------|-------------|------|--------------------------------|
| Patient's Name | (0010,0010) | 2 | Length=0 when no entry is made |
| Patient ID | (0010,0020) | 2 | Length=0 when no entry is made |
| Patient's Birth Date | (0010,0030) | 2 | Length=0 when no entry is made |
| Patient's Sex | (0010,0040) | 2 | Length=0 when no entry is made |
| Patient Comment | (0010,4000) | 3 | Not set when no entry is made |

8.2.2 General Study Module

Table 21

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------------------|-------------|-------------|--------------------------------|
| Study Instance UID | (0020,000D) | 1 | Always set |
| Study Date | (0008,0020) | 2 | Always set |
| Study Time | (0008,0030) | 2 | Always set |
| Referring Physician's Name | (0008,0090) | 2 | Length=0 when no entry is made |
| Study ID | (0020,0010) | 2 | Always set |
| Accession Number | (0008,0050) | 2 | Length=0 when no entry is made |
| Name of Physician(s) Reading Study | (0008,1060) | 3 | Not set when no entry is made |

8.2.3 Patient Study Module

Table 22

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|-------------------------------|
| Patient's Age | (0010,1010) | 3 | Not set when no entry is made |
| Patient's Size | (0010,1020) | 3 | Not set when no entry is made |
| Patient's Weight | (0010,1030) | 3 | Not set when no entry is made |

8.2.4 General Series Module

Table 23

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|-----------------------------------|
| Modality | (0008,0060) | 1 | Always set |
| Series Instance UID | (0020,000E) | 1 | Always set |
| Series Number | (0020,0011) | 2 | Always set |
| Laterality | (0020,0060) | 2C | Length=0 when no entry is made |
| Series Date | (0008,0021) | 3 | Not set when no data is available |
| Series Time | (0008,0031) | 3 | Not set when no data is available |
| Series Description | (0008,103E) | 3 | Always set |
| Operator's Name | (0008,1070) | 3 | Not set when no entry is made |
| Body Part Examined | (0018,0015) | 3 | Always set |
| Patient Position | (0018,5100) | 2C | Always set |

8.2.5 Frame of Reference Module

Table 24

| Attribute Name | Tag | Type | Attribute Description |
|------------------------------|-------------|-------------|------------------------------------|
| Position Reference Indicator | (0020,1040) | 2 | Length=0 when no data is available |
| Frame of Reference UID | (0020,0052) | 1 | Always set |

8.2.6 General Equipment Module

Table 25

| Attribute Name | Tag | Type | Attribute Description |
|---------------------------|-------------|-------------|-----------------------------------|
| Manufacturer | (0008,0070) | 2 | Always set |
| Institution Name | (0008,0080) | 3 | Not set when no data is available |
| Station Name | (0008,1010) | 3 | Not set when no data is available |
| Manufacturer's Model Name | (0008,1090) | 3 | Not set when no data is available |
| Device Serial Number | (0018,1000) | 3 | Not set when no data is available |
| Software Version(s) | (0018,1020) | 3 | Not set when no data is available |

8.2.7 General Image Module

Table 26

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|-----------------------------------|
| Image Type | (0008,0008) | 1 | Always set |
| Image Number | (0020,0013) | 2 | Always set |
| Patient Orientation | (0020,0020) | 2C | Always set |
| Acquisition Number | (0020,0012) | 3 | Always set |
| Acquisition Date | (0008,0022) | 3 | Not set when no data is available |
| Acquisition Time | (0008,0032) | 3 | Not set when no data is available |

8.2.8 Image Plane Module

Table 27

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------------|-------------|-------------|------------------------------|
| Pixel Spacing | (0028,0030) | 1 | Always set |
| Image Orientation (Patient) | (0020,0037) | 1 | Always set |
| Image Position (Patient) | (0020,0032) | 1 | Always set |
| Slice Thickness | (0018,0050) | 2 | Always set |
| Slice Location | (0020,1041) | 3 | Always set |

8.2.9 Image Pixel Module

Table 28

| Attribute Name | Tag | Type | Attribute Description |
|----------------------------|-------------|-------------|-----------------------------------|
| Samples per Pixel | (0028,0002) | 1 | Always set(1) |
| Photometric Interpretation | (0028,0004) | 1 | Always set("MONOCHROME2") |
| Rows | (0028,0010) | 1 | Always set |
| Columns | (0028,0011) | 1 | Always set |
| Bits Allocated | (0028,0100) | 1 | Always set(16) |
| Bits Stored | (0028,0101) | 1 | Always set(16) |
| High Bit | (0028,0102) | 1 | Always set(15) |
| Pixel Representation | (0028,0103) | 1 | Always set(1) |
| Pixel Data | (7FE0,0010) | 1 | Always set |
| Smallest Image Pixel Value | (0028,0106) | 3 | Not set when no data is available |
| Largest Image Pixel Value | (0028,0107) | 3 | Not set when no data is available |

8.2.10 Contrast/Bolus Module

Table 29

| Attribute Name | Tag | Type | Attribute Description |
|----------------------|-------------|------|--------------------------------|
| Contrast/Bolus Agent | (0018,0010) | 2 | Length=0 when no entry is made |

8.2.11 MR Image Module

Table 30

| Attribute Name | Tag | Type | Attribute Description |
|------------------------------|-------------|------|-----------------------------------|
| Image Type | (0008,0008) | 1 | Always set |
| Samples per Pixel | (0028,0002) | 1 | Always set(1) |
| Photometric Interpretation | (0028,0004) | 1 | Always set("MONOCHROME2") |
| Bits Allocated | (0028,0100) | 1 | Always set(16) |
| Scanning Sequence | (0018,0020) | 1 | Always set |
| Sequence Variant | (0018,0021) | 1 | Always set("NONE") |
| Scan Options | (0018,0022) | 2 | Length=0 when no entry is made |
| MR Acquisition Type | (0018,0023) | 2 | Length=0 when no entry is made |
| Sequence Name | (0018,0024) | 3 | Always set |
| Repetition Time | (0018,0080) | 2C | Always set |
| Echo Time | (0018,0081) | 2 | Always set |
| Echo Train Length | (0018,0091) | 2 | Length=0 when no entry is made |
| Inversion Time | (0018,0082) | 2C | Not set when no entry is made |
| Number of Averages | (0018,0083) | 3 | Always set |
| Imaging Frequency | (0018,0084) | 3 | Not set when no entry is made |
| Imaged Nucleus | (0018,0085) | 3 | Always set |
| Echo Number(s) | (0018,0086) | 3 | Always set |
| Trigger Time | (0018,1060) | 2C | Not set when no entry is made |
| Nominal Interval | (0018,1062) | 3 | Not set when no data is available |
| Acquisition Matrix | (0018,1310) | 3 | Always set |
| Flip Angle | (0018,1314) | 3 | Always set |
| Temporal Position Identifier | (0020,0100) | 3 | Not set when no data is available |

8.2.12 VOI LUT Module

Table 31

| Attribute Name | Tag | Type | Attribute Description |
|-----------------------|-------------|-------------|------------------------------|
| Window Center | (0028,1050) | 3 | Always set |
| Window Width | (0028,1051) | 1C | Always set |

8.2.13 SOP Common Module

Table 32

| Attribute Name | Tag | Type | Attribute Description |
|------------------------|-------------|-------------|-----------------------------------------|
| Specific Character Set | (0008,0005) | 1C | Not set when the system is English mode |
| SOP Class UID | (0008,0016) | 1 | Always set |
| SOP Instance UID | (0008,0018) | 1 | Always set |

9 DIMSE-Service and Attributes

9.1 DIMSE-Services

Table 33

| SOP Class | DIMSE Service Element | Usage SCU ^{*1} | Usage |
|------------------------------|-----------------------|-------------------------|----------|
| Basic Film Session SOP Class | N-CREATE | M | used |
| | N-SET | U | not used |
| | N-DELETE | U | used |
| | N-ACTION | U | not used |
| Basic Film Box SOP Class | N-CREATE | M | used |
| | N-SET | U | not used |
| | N-DELETE | U | used |
| | N-ACTION | M | used |
| Image Box SOP Class | N-SET | M | used |
| Printer SOP Class | N-EVENT-REPORT | M | used |
| | N-GET | U | used |

*1 : M = Mandatory, U = User option

9.2 Basic Film Session SOP Class

9.2.1 N-CREATE Attributes

Table 34

| Attribute Name | Tag | Usage | Attribute Description |
|--------------------|-------------|-------|-----------------------------------|
| Number of Copies | (2000,0010) | U | Always set |
| Print Priority | (2000,0020) | U | Always set("LOW") |
| Medium Type | (2000,0030) | U | Always set("CLEAR FILM") |
| Film Destination | (2000,0040) | U | Not set when no data is available |
| Film Session Label | (2000,0050) | U | Always set("TOSHIBA_MRI") |
| Memory Allocation | (2000,0060) | U | Not set when no data is available |

9.3 Basic Film Box SOP Class

9.3.1 N-CREATE Attributes

Table 35

| Attribute Name | Tag | Usage | Attribute Description |
|-------------------------------------|-------------|-------|-----------------------|
| Image Display Format | (2010,0010) | M | Always set |
| Film Orientation | (2010,0040) | U | Always set |
| Film Size ID | (2010,0050) | U | Always set |
| Magnification Type | (2010,0060) | U | Always set |
| Border Density | (2010,0100) | U | Always set |
| Empty Image Density | (2010,0110) | U | Always set |
| Trim | (2010,0140) | U | Always set |
| Referenced Film Session Sequence | (2010,0500) | M | Always set |
| >Referenced SOP Class UID | (0008,1150) | M | Always set |
| >Referenced SOP Instance UID | (0008,1155) | M | Always set |

9.4 Basic Grayscale Image Box SOP Class

9.4.1 N-SET Attributes

Table 36

| Attribute Name | Tag | Usage | Attribute Description |
|--------------------------------|-------------|-------|---------------------------|
| Image Position | (2020,0010) | M | Always set |
| Polarity | (2020,0020) | U | Always set |
| Magnification Type | (2010,0060) | U | Always set |
| Basic Grayscale Image Sequence | (2020,0110) | M | Always set |
| >Samples Per Pixel | (0028,0002) | M | Always set(1) |
| >Photometric Interpretation | (0028,0004) | M | Always set("MONOCHROME2") |
| >Rows | (0028,0010) | M | Always set |
| >Columns | (0028,0011) | M | Always set |
| >Pixel Aspect Ratio | (0028,0034) | MC | Always set |
| >Bits Allocated | (0028,0100) | M | Always set(16 or 8) |
| >Bits Stored | (0028,0101) | M | Always set |
| >High Bit | (0028,0102) | M | Always set |
| >Pixel Representation | (0028,0103) | M | Always set(0) |
| >Pixel Data | (7FE0,0010) | M | Always set |

9.5 Printer SOP Class

9.5.1 N-EVENT-REPORT

Table 37

| Event Type Name | Event Type ID | Attribute | Tag | Usage SCU/SCP |
|-----------------|---------------|----------------------------|-------------|---------------|
| NORMAL | 1 | | | |
| WARNING | 2 | Printer Name | (2110,0030) | U/U |
| | | Printer Status Information | (2110,0020) | U/M |
| FAILURE | 3 | Printer Name | (2110,0030) | U/U |
| | | Printer Status Information | (2110,0020) | U/M |

9.5.2 N-GET Attributes

Table 38

| Attribute Name | Tag | Usage SCU/SCP |
|---------------------------|-------------|---------------|
| Printer Status | (2110,0010) | U/M |
| Printer Status Info | (2110,0020) | U/M |
| Printer Name | (2110,0030) | U/U |
| Manufacturer | (0008,0070) | U/U |
| Manufacturer's Model Name | (0008,1090) | U/U |
| Device Serial Number | (0018,1000) | U/U |
| Software Version | (0018,1020) | U/U |
| Date of Last Calibration | (0018,1200) | U/U |
| Time of Last Calibration | (0018,1201) | U/U |

10 Search Keys

10.1 Query/Retrieve SCP (C-FIND)

The search keys used for the Query/Retrieve SCP(C-FIND) are shown.

10.1.1 Patient Root Q/R Information Model - FIND

10.1.1.1 Patient Level

Table 39

| Attribute Name | Tag | Type |
|----------------|-------------|------|
| Patient's Name | (0010,0010) | R |
| Patient ID | (0010,0020) | U |

10.1.1.2 Study Level

Table 40

| Attribute Name | Tag | Type |
|--------------------|-------------|------|
| Study Date | (0008,0020) | R |
| Study Time | (0008,0030) | R |
| Accession Number | (0008,0050) | R |
| Study ID | (0020,0010) | R |
| Study Instance UID | (0020,000D) | U |

10.1.1.3 Series Level

Table 41

| Attribute Name | Tag | Type |
|---------------------|-------------|------|
| Modality | (0008,0060) | R |
| Series Number | (0020,0011) | R |
| Series Instance UID | (0020,000E) | U |

10.1.1.4 Image Level

Table 42

| Attribute Name | Tag | Type |
|------------------|-------------|------|
| Image Number | (0020,0013) | R |
| SOP Instance UID | (0008,0018) | U |

10.1.2 Study Root Q/R Information Model - FIND**10.1.2.1 Study Level**

Table 43

| Attribute Name | Tag | Type |
|-----------------------|-------------|-------------|
| Patient's Name | (0010,0010) | R |
| Patient ID | (0010,0020) | R |
| Study Date | (0008,0020) | R |
| Study Time | (0008,0030) | R |
| Accession Number | (0008,0050) | R |
| Study ID | (0020,0010) | R |
| Study Instance UID | (0020,000D) | U |

10.1.2.2 Series Level

Attributes for the Series Level of the Study Root Query/Retrieve Information Model are the same as the Attributes for the Series Level of the Patient Root Query/Retrieve Information Model described in Section 10.1.1.3.

10.1.2.3 Image Level

Attributes for the Image Level of the Study Root Query/Retrieve Information Model are the same as the Attributes for the Image Level of the Patient Root Query/Retrieve Information Model described in Section 10.1.1.4.

10.1.3 Patient/Study Only Information Model - FIND**10.1.3.1 Patient Level**

Attributes for the Patient Level of the Patient/Study Only Query/Retrieve Information Model are the same as the Attributes for the Patient Level of the Patient Root Query/Retrieve Information Model described in Section 10.1.1.1.

10.1.3.2 Study Level

Attributes for the Study Level of the Patient/Study Only Query/Retrieve Information Model are the same as the Attributes for the Study Level of the Patient Root Query/Retrieve Information Model described in Section 10.1.1.2