

**DICOM
CONFORMANCE STATEMENT
FOR
MODEL TWS-2000 SERIES
MODEL TWS-3000 SERIES
MODEL TWS-5000 SERIES
MODEL STWS-005 SERIES
MODEL TIS-3000
MODEL TIS-4000
(MIIMS0002EAG)**

TOSHIBA MEDICAL SYSTEMS CORPORATION

© TOSHIBA MEDICAL SYSTEMS CORPORATION 1998-2008

ALL RIGHTS RESERVED

Trademarks

This document may include trademarks or registered trademarks of other companies.

IMPORTANT!

- (1) No part of this document may be copied or reprinted, in whole or in part, without written permission.
- (2) The contents of this document are subject to change without prior notice and without our legal obligation.

Table of Contents

	Page
1 Introduction	1
1.1 References	1
1.2 Definitions	1
1.3 Acronyms, Abbreviations, and Symbols	2
2 Implementation Model	3
2.1 Application Data Flow Diagram	3
2.2 Functional Definitions of AE's	4
2.2.1 Export AE	4
2.2.2 Import AE	4
2.3 Sequencing of Real World Activities	4
2.3.1 Features	4
2.3.2 Operation	5
3 AE Specifications	6
3.1 Export AE Specification	6
3.1.1 Export AE Association Establishment Policies	7
3.1.2 Export AE Association Initiation by Real-World Activity	8
3.1.3 Export AE Association Acceptance Policy	14
3.2 Import AE Specification	14
3.2.1 Import AE Association Establishment Policies	15
3.2.2 Import AE Association Initiation by Real-World Activity	15
3.2.3 Import AE Association Acceptance Policy	15
4 Communication Profiles	20
4.1 Supported Communication Stacks	20
4.2 OSI Stack	20
4.3 TCP/IP Stack	20
4.3.1 API	20
4.3.2 Physical Media Support	20
4.4 Point-to-Point Stack	20
5 Extensions/Specializations/Privateizations	21
5.1 Standard Extended/Specialized/Private SOP - Storage SCU	21
6 Configuration	22
6.1 AE Title/Presentation Address Mapping	22
6.2 Configurable Parameters	22
6.2.1 Time-out Value, Retry Count, Retry Interval	22
6.2.2 Warning Status Criteria	24
6.3 Implementation Information and Maximum Reception PDU Size	29
6.4 Default Transfer Syntax	30
6.4.1 Import AE	30

7	Support of Extended Character Sets	31
7.1	Query/Retrieve SCU	31
7.2	Storage SCU/SCP	31
7.3	Print SCU	31
8	Information Object Definition - Storage SCU	32
8.1	Entity Module Definitions	32
8.1.1	CR IMAGE IOD Modules	32
8.1.2	CT IMAGE IOD Modules	33
8.1.3	MR IMAGE IOD Modules	34
8.1.4	NM IMAGE IOD Modules	35
8.1.5	SC IMAGE IOD Modules	36
8.1.6	US IMAGE IOD Modules	37
8.1.7	US MULTI-FRAME IMAGE IOD Modules	38
8.1.8	XA IMAGE IOD Modules	39
8.1.9	XRF IMAGE IOD Modules	40
8.1.10	DX IMAGE IOD Modules	41
8.1.11	MG IMAGE IOD Modules	42
8.1.12	PET IMAGE IOD Modules	43
8.2	Information Object Definitions	44
8.2.1	Patient Module	44
8.2.2	General Study Module	44
8.2.3	Patient Study Module	45
8.2.4	General Series Module	45
8.2.5	CR Series Module	46
8.2.6	NM/PET Patient Orientation Module	46
8.2.7	Frame of Reference Module	47
8.2.8	US Frame of Reference Module	47
8.2.9	General Equipment Module	48
8.2.10	SC Equipment Module	48
8.2.11	General Image Module	49
8.2.12	Image Plane Module	50
8.2.13	Image Pixel Module	51
8.2.14	NM Image Pixel Module	52
8.2.15	Contrast/Bolus Module	53
8.2.16	Palette Color Lookup Table Module	54
8.2.17	US Region Calibration Module	55
8.2.18	Cine Module	57
8.2.19	Multi-Frame Module	57
8.2.20	NM Multi-Frame Module	58
8.2.21	Frame Pointers Module	59
8.2.22	Mask Module	59
8.2.23	Display Shutter Module	60
8.2.24	Device Module	60
8.2.25	Therapy Module	61
8.2.26	CR Image Module	62
8.2.27	CT Image Module	63
8.2.28	MR Image Module	64
8.2.29	NM Image Module	66
8.2.30	US Image Module	68
8.2.31	NM Isotope Module	71
8.2.32	NM Detector Module	73
8.2.33	NM TOMO Acquisition Module	74
8.2.34	NM Multi-Gated Acquisition Module	75
8.2.35	NM Phase Module	75
8.2.36	NM Reconstruction Module	76

8.2.37	SC Image Module	76
8.2.38	X-Ray Image Module	77
8.2.39	X-Ray Acquisition Module	78
8.2.40	X-Ray Collimator Module	79
8.2.41	X-Ray Table Module	79
8.2.42	XRF Positioner Module	80
8.2.43	XRF Tomography Acquisition Module	80
8.2.44	XA Positioner Module	80
8.2.45	Overlay Plane Module	81
8.2.46	Multi-frame Overlay Module	82
8.2.47	Curve Module	82
8.2.48	Modality LUT Module	83
8.2.49	VOI LUT Module	83
8.2.50	SOP Common Module	84
8.2.51	Curve Identification Module	85
8.2.52	Audio Module	86
8.2.53	Common Attribute Set for Code Sequence	86
8.2.54	Specimen Identification Module	87
8.2.55	DX Series Module	87
8.2.56	DX Anatomy Imaged Module	88
8.2.57	DX Image Module	89
8.2.58	DX Detector Module	90
8.2.59	DX Positioning Module	91
8.2.60	X-Ray Acquisition Dose Module	92
8.2.61	X-Ray Generation Module	93
8.2.62	X-Ray Filtration Module	93
8.2.63	X-Ray Grid Module	94
8.2.64	X-Ray Tomography Acquisition Module	94
8.2.65	Image Histogram Module	94
8.2.66	Acquisition Context Module	95
8.2.67	Mammography Series Module	95
8.2.68	Mammography Image Module	96
8.2.69	PET Series Module	97
8.2.70	PET Isotope Module	99
8.2.71	PET Multi-gated Acquisition Module	100
8.2.72	NM/PET Patient Orientation Module	100
8.2.73	PET Image Module	101
8.3	Private Data Elements	103
9	Information Object Definition - Storage SCP	106
9.1	Entity Module Definitions	106
10	Search Keys	107
10.1	Query/Retrieve SCU (C-Find)	107
10.1.1	Patient Root Q/R Information Model - Find	107
10.1.2	Study Root Q/R Information Model - Find	108
10.1.3	Patient/Study Only Q/R Information Model - Find	109
11	DIMSE-Service and Attributes	110
11.1	DIMSE-Services	110
11.2	Basic Film Session SOP Class	111
11.2.1	N-CREATE Attributes	111
11.3	Basic Film Box SOP Class	112
11.3.1	N-CREATE Attributes	112
11.4	Basic Grayscale Image Box SOP Class	113

11.4.1 N-SET Attributes	113
11.5 Basic Color Image Box SOP Class	114
11.5.1 N-SET Attributes	114
11.6 Printer SOP Class	115
11.6.1 N-GET Attributes	115

1 Introduction

This document is a DICOM Conformance Statement for Toshiba's TWS-2000 Series, TWS-3000 Series and TWS-5000 Series, STWS-005 Series, TIS-3000, TIS-4000 (hereinafter referred to as the TWS). It is intended to provide the reader with the knowledge of how to integrate this product within a DICOM compliant hospital network. It details the DICOM Service Classes, Information Objects, and Communication Protocols which are supported by this product.

If readers are unfamiliar with DICOM, it is recommended that they read the DICOM Specification (referenced below) prior to reading this conformance statement. Also note that this document is formatted according to the DICOM Specification, Part 2: Conformance.

1.1 References

- ACR-NEMA Digital Imaging and Communications in Medicine, DICOM V3.0.

1.2 Definitions

- **Association Establishment** - An Association Establishment is the first phase of communication between two DICOM Application Entities (AEs). The AEs use the Association Establishment to negotiate how data will be encoded and the type of data to be exchanged.
- **Called Application Entity Title** - The Called AE Title defines the intended receiver of an Association.
- **Calling Application Entity Title** - The Calling AE Title defines the requestor of an Association.
- **DICOM Message Service Element (DIMSE)** - A DIMSE defines the services and protocols utilized by an Application Entity to exchange messages.
- **Information Object Definition (IOD)** - An IOD is a data model which is an abstraction of real-world information. This data model defines the nature and attributes relevant to the class of real-world objects represented.
- **Service Class Provider (SCP)** - A Service Class Provider plays the "server" role to perform operations and invoke notifications during an Association. An example of a Storage Service Class Provider would be an image storage device. In this case, the image storage device stores the image that was sent by a Service Class User.
- **Service Class User (SCU)** - A Service Class User plays the "client" role to invoke operations and perform notifications during an Association. An example of a Storage Service Class User would be an image acquisition device. In this case, the image acquisition device creates and sends a DICOM image by requesting that a Service Class Provider store that image.
- **Service/Object Pair (SOP) Class** - An SOP Class is defined by the union of an Information Object Definition and a set of DIMSE Services.
A DICOM Application Entity may support one or more SOP Classes. Each SOP Class is uniquely identified by an SOP Class UID.
- **SOP Instance** - A specific occurrence of an Information Object.
- **Transfer Syntax** - The Transfer Syntax is a set of encoding rules that allow DICOM Application Entities to negotiate the encoding techniques (e.g. data element structure, byte ordering, compression) that they are able to support. The Transfer Syntax is negotiated during Association Negotiation.

- **Unique Identifier (UID)** - A Unique Identifier is a globally unique, ISO compliant, ASCII-numeric string. It guarantees uniqueness across multiple countries, sites, vendors, and equipment.

1.3 Acronyms, Abbreviations, and Symbols

- ACC American College of Cardiology
- ACR American College of Radiology
- ASCII American Standard Code for Information Interchange
- AE Application Entity
- ANSI American National Standards Institute
- CEN TC251 Comite Europeen de Normalisation - Technical Committee 251 - Medical Informatics
- DICOM Digital Imaging and Communications in Medicine
- DIMSE DICOM Message Service Element
- DIMSE-C DICOM Message Service Element - Composite
- DIMSE-N DICOM Message Service Element - Normalized
- HIS Hospital Information System
- HL7 Health Level 7
- IE Information Entity
- IOD Information Object Definition
- ISO International Standards Organization
- JIRA Japan Industries Association of Radiological Systems
- JIS Japanese Industrial Standards
- NEMA National Electrical Manufacturers Association
- MEDIS-DC the MEDical Information System Development Center
- OSI Open Systems Interconnection
- PDU Protocol Data Unit
- RIS Radiology Information System
- SCP Service Class Provider
- SCU Service Class User
- SOP Service-Object Pair
- TCP/IP Transmission Control Protocol/Internet Protocol
- UID Unique Identifier

2 Implementation Model

2.1 Application Data Flow Diagram




Figure 1

2.2 Functional Definitions of AE's

2.2.1 Export AE

Export AE is used to transmit Query and Retrieve requests to a remote DICOM device. It therefore performs the following tasks:

- Establishes DICOM Association with remote DICOM device
- Performs request of Query/Retrieve to remote DICOM device

Export AE is used to transmit images to a remote DICOM device. It therefore performs the following tasks:

- Establishes DICOM Association with remote DICOM device
- Performs storage of DICOM CR, CT, MR, NM, SC, US, US Multi-frame, XA, XRF, DX, MG, and PET images Information Objects to remote DICOM device

Export AE is used to transmit request for print images to a remote DICOM device. It therefore performs the following tasks:

- Builds DICOM Basic Grayscale/Color Print Objects
- Establishes DICOM Association with remote DICOM device
- Performs transmit of DICOM Basic Grayscale/Color Print Objects to remote DICOM device

2.2.2 Import AE

Import AE is used to respond to requests to verify that the TWS is present and active on the network and to receive CR, CT, MR, NM, SC, US, US Multi-frame, XA, XRF, DX, MG, and PET images from remote DICOM devices.

2.3 Sequencing of Real World Activities

2.3.1 Features

2.3.1.1 Request of Query/Retrieve

- Operator can use the search key defined as chapter 10, and can get a search list.
- The search list is displayed in units of patient, study, series, and image levels. Image transfer requests are also handled in units of patient, study, series, and image levels.
- If an error occurs when image data is being retrieved, an error message is displayed.

2.3.1.2 Manual Send of Image, Study or Patient Unit

- Operator requests to send data after selecting them from the Image, Study or Patient List.
- When the image transfer fails, display the Send Failed Image List. And operator can manually attempt to resend image at a later time.

2.3.1.3 Archiving of the Received Image to the Local File System

- The TWS receives CR, CT, MR, NM, SC, US, US Multi-frame, XA, XRF, DX, MG, and PET images from remote DICOM devices.
- The TWS archives the received images to the local file system.

2.3.1.4 Manual Print of Image, Series or Study Unit

- Operator requests to print images displayed on monitor(s).
- When the image transfer fails, display the Send Failed Image List. And operator can manually attempt to resend image at a later time.

2.3.2 Operation

2.3.2.1 Request of Query/Retrieve

The operation for a search request and an image transfer request are described below:

- Step-1: Enter the information, search key defined as chapter 10, for the items for which the operator wishes to search.
- Step-2: Execute the search request.
- Step-3: Select images from the results of the search.
- Step-4: Execute the image transfer request.

2.3.2.2 Manual Send of Image, Study or Patient Unit

The operation for sending images is described below:

- Step-1: Select the patient, study or the image to be sent.
- Step-2: Select the destination of image sending.
- Step-3: Request transfer.

2.3.2.3 Archiving of the Received Image to the Local File System

There is no specific operation for receiving and archiving images.

2.3.2.4 Manual Print of Image, Series or Study Unit

- Step-1: Select the images for Image, Series or Study Unit
- Step-2: Set up the parameters for printing images.
- Step-3: Select the destination of image printing.
- Step-4: Execute the print image request.

3 AE Specifications

3.1 Export AE Specification

Export AE provides Standard Conformance to the following DICOM SOP Classes as an SCU:

Table 1

SOP Class Name	SOP Class UID
CR Image Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Storage	1.2.840.10008.5.1.4.1.1.2
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
NM Image Storage	1.2.840.10008.5.1.4.1.1.20
SC Image Storage	1.2.840.10008.5.1.4.1.1.7
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1
US Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1
XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2
DX Image Storage	1.2.840.10008.5.1.4.1.1.1.1
MG Image Storage	1.2.840.10008.5.1.4.1.1.1.2
PET Image Storage	1.2.840.10008.5.1.4.1.1.128
Patient Root Query/Retrieve Information Model-Find	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve Information Model-Move	1.2.840.10008.5.1.4.1.2.1.2
Study Root Query/Retrieve Information Model-Find	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve Information Model-Move	1.2.840.10008.5.1.4.1.2.2.2
Patient/Study Only Query/Retrieve Information Model-Find	1.2.840.10008.5.1.4.1.2.3.1
Patient/Study Only Query/Retrieve Information Model-Move	1.2.840.10008.5.1.4.1.2.3.2
Basic Grayscale Print Management	1.2.840.10008.5.1.1.9
Basic Color Print Management	1.2.840.10008.5.1.1.18

The SOP Classes listed in Table2 indicate the SOP Classes regulated by the Basic Grayscale Print Management Meta SOP Class.

Table 2

SOP Class Name	SOP Class UID
Basic Film Session SOP Class	1.2.840.10008.5.1.1.1
Basic Film Box SOP Class	1.2.840.10008.5.1.1.2
Basic Grayscale Image Box SOP Class	1.2.840.10008.5.1.1.4
Printer SOP Class	1.2.840.10008.5.1.1.16

The SOP Classes listed in Table3 indicate the SOP Classes regulated by the Color Grayscale Print Management Meta SOP Class.

Table 3

SOP Class Name	SOP Class UID
Basic Film Session SOP Class	1.2.840.10008.5.1.1.1
Basic Film Box SOP Class	1.2.840.10008.5.1.1.2
Basic Color Image Box SOP Class	1.2.840.10008.5.1.1.4.1
Printer SOP Class	1.2.840.10008.5.1.1.16

3.1.1 Export AE Association Establishment Policies

3.1.1.1 Export AE General

Export AE will utilize and understand the following Application Context Name:

Table 4

DICOM V3.0 Application Context	1.2.840.10008.3.1.1.1
--------------------------------	-----------------------

Export AE supports a minimum PDU size of 16 Kbytes and a maximum PDU size of 64 Kbytes. The default value is set to 64 Kbytes.

3.1.1.2 Export AE Number of Associations

Export AE can establish two associations for FIND and MOVE, and one association for Storage and Print at once regardless of the number of the selected destinations.

3.1.1.3 Export AE Asynchronous Nature

Export AE allows a single outstanding operation on any association. Therefore, Export AE does not support asynchronous operation window negotiation other than the default as specified by the specification.

3.1.1.4 Export AE Implementation Identifying Information

Export AE will specify the following Implementation Identifying Information:

Table 5

Applicable Models	Implementation Class UID	Implementation Version Name
TWS-2100,TWS-2300P, TWS-2300L,TWS-2400A, TWS-2500P, TWS-2900, TWS-3100L, TWS-3300P, TWS-5100L, TWS-5200P, TWS-5300P, STWS-005, STWS-005Y, TIS-3000, TIS-4000	1.2.392.200036.9116.7.4.10	TM_OT_PTWS_1.0

3.1.2 Export AE Association Initiation by Real-World Activity

Export AE initiates an association when the following activity is chosen by the operator:

- "Manual Send of Image,Study or Patient Unit"
 - Storage - Store a CR, CT, MR, NM, SC, US, US Multi-frame,XA, XRF, DX, MG, or PET image to a remote DICOM device
- "Request of Query/Retrieve"
 - Find -Get an image list from a remote DICOM device
 - Move -Send an image transfer request to a remote DICOM device
- "Manual Print of Image,Series or Study Unit"
 - Print -Request print image to a remote DICOM device.

3.1.2.1 Export AE Real-World Activity – Storage

3.1.2.1.1 Export AE Associated Real-World Activity - Storage

The Storage is executed by the TWS after the operator requests the image transfer.

3.1.2.1.2 Export AE Proposed Presentation Contexts - Storage

Export AE proposes the following Presentation Contexts shown below:

Table 6

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
CR Image Storage	1.2.840.10008.5.1.4.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
NM Image Storage	1.2.840.10008.5.1.4.1.1.20	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
SC Image Storage	1.2.840.10008.5.1.4.1.1.7	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
US Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
DX Image Storage	1.2.840.10008.5.1.4.1.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
MG Image Storage	1.2.840.10008.5.1.4.1.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
PET Image Storage	1.2.840.10008.5.1.4.1.1.128	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

3.1.2.1.2.1 Export AE SOP Specific Conformance - Storage

Export AE operation involves the following sequence of steps for each image transfer.

- (1) Association establishment (requestor only)
- (2) Data transfer (SCU only)
- (3) Association release (requestor only)

Export AE judges that the transfer of one image succeeded when the result of (2) "Data transfer" is "Success" even if the result of (3) "Association release" is "Failure".

3.1.2.2 Export AE Real-World Activity – Query/Retrieve (Find)

3.1.2.2.1 Export AE Associated Real-World Activity - Query/Retrieve (Find)

The operator specifies a search and sends the request to a remote DICOM device. When the results of the search have been completely received, the association is released.

3.1.2.2.2 Export AE Proposed Presentation Contexts - Query/Retrieve (Find)

Export AE proposes the following presentation contexts shown below:

Table 7

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Q/R Info. Model-Find	1.2.840.10008.5.1.4.1.2.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
Study Root Q/R Info. Model-Find	1.2.840.10008.5.1.4.1.2.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
Patient/Study Only Q/R Info. Model-Find	1.2.840.10008.5.1.4.1.2.3.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

3.1.2.2.2.1 Export AE SOP Specific Conformance - Query/Retrieve (Find)

Export AE operation involves the following sequence of steps for each search request:

- (1) Association establishment (requestor only)
- (2) Query request (SCU only)
- (3) Association release (requestor only)

Export AE judges that query request succeeded when the result of (2) "Query request" is "Success", even if the result of (3) "Association release" is "Failure".

Search keys for the Query/Retrieve SCU are described in chapter 10.

3.1.2.3 Export AE Real-World Activity – Query/Retrieve (Move)

3.1.2.3.1 Export AE Associated Real-World Activity - Query/Retrieve (Move)

The operator specifies a patient, study, series, and image to retrieve and sends the request to a remote DICOM device. The association is released when all the image data corresponding to the request list has been received from the remote DICOM device, or when the operator makes cancel or abort request.

3.1.2.3.2 Export AE Proposed Presentation Contexts - Query/Retrieve (Move)

Export AE proposes the following presentation contexts shown below:

Table 8

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Q/R Info. Model-Move	1.2.840.10008.5.1.4.1.2.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
Study Root Q/R Info. Model-Move	1.2.840.10008.5.1.4.1.2.2.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
Patient/Study Only Q/R Info. Model-Move	1.2.840.10008.5.1.4.1.2.3.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

3.1.2.3.2.1 Export AE SOP Specific Conformance - Query/Retrieve (Move)

Export AE operation involves the following sequence of steps for each retrieve request:

- (1) Association establishment (requestor only)
- (2) Retrieve (image transfer) request (SCU only)
- (3) Association release (requestor only)

Export AE judges that Retrieve request succeeded when the result of (2) "Retrieve request" is Success, even if the result of (3) "Association release" is Failure.

3.1.2.4 Export Real-World Activity - Print

3.1.2.4.1 Export Associated Real-World Activity - Print

Export AE performs request print images manually to destination device.

3.1.2.4.2 Export Proposed Presentation Contexts - Print

Export AE proposes the following Presentation Contexts shown below:

Table 9

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Basic Grayscale Print Management	1.2.840.10008.5.1.1.9	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
Basic Color Print Management	1.2.840.10008.5.1.1.18	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

3.1.2.4.2.1 Export AE SOP Specific Conformance - Print

Export AE operation involves the following sequence of steps for each image transfer.

- (1) Association establishment(requestor only)
- (2) Data transfer(SCU only)
- (3) Association release(requestor only)

Export AE judges that the transfer of one image succeeded when the result of (2)"Data transfer" is "Success" even if the result of (3)"Association release" is "Failure".

DIMSE-Service and Attributes are described in chapter 11.

3.1.3 Export AE Association Acceptance Policy

Export AE does not accept any associations generated by remote applications.

3.2 Import AE Specification

Import AE provides Standard Conformance to the following DICOM SOP Classes as an SCP:

Table 10

SOP Class Name	SOP Class UID
Verification	1.2.840.10008.1.1
CR Image Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Storage	1.2.840.10008.5.1.4.1.1.2
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
NM Image Storage	1.2.840.10008.5.1.4.1.1.20
SC Image Storage	1.2.840.10008.5.1.4.1.1.7
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1
US Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1
XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2
DX Image Storage	1.2.840.10008.5.1.4.1.1.1.1
MG Image Storage	1.2.840.10008.5.1.4.1.1.1.2
PET Image Storage	1.2.840.10008.5.1.4.1.1.128

3.2.1 Import AE Association Establishment Policies

3.2.1.1 Import AE General

Import AE will utilize and understand the following Application Context Name:

Table 11

DICOM V3.0 Application Context	1.2.840.10008.3.1.1.1
--------------------------------	-----------------------

Import AE supports a minimum PDU size of 16 Kbytes and a maximum PDU size of 64 Kbytes. The default value is set to 64 Kbytes.

3.2.1.2 Import AE Number of Associations

Import AE supports up to three associations at a time.

3.2.1.3 Import AE Asynchronous Nature

Import AE allows a single outstanding operation on any association. Therefore, Import AE does not support asynchronous operation window negotiation other than the default as specified by the specification.

3.2.1.4 Import AE Implementation Identifying Information

Import AE will specify the following Implementation Identifying Information:

Table 12

Applicable Models	Implementation Class UID	Implementation Version Name
TWS-2100, TWS-2300P, TWS-2300L, TWS-2400A, TWS-2500P, TWS-2900, TWS-3100L, TWS-3300P, TWS-5100L, TWS-5200P, TWS-5300P, STWS-005, STWS-005Y, TIS-3000, TIS-4000	1.2.392.200036.9116.7.4.10	TM_OT_PTWS_1.0

3.2.2 Import AE Association Initiation by Real-World Activity

Import AE never initiates an association.

3.2.3 Import AE Association Acceptance Policy

When Import AE receives an association request, it will allow the following activities to be performed during that association:

- Verification - Allow a remote DICOM device to verify that the TWS is active on the DICOM network
- Storage - Allow a remote DICOM device to send a CR, CT, MR, NM, SC, US, US Multi-frame, XA, XRF, DX, MG or PET image to the TWS

3.2.3.1 Import AE Real-World Activity - Verification

3.2.3.1.1 Import AE Associated Real-World Activity - Verification

The TWS responds to the Verification made by a remote Verification SCU.

3.2.3.1.2 Import AE Presentation Context Table - Verification

Import AE accepts the following Presentation Contexts shown below:

Table 13

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None

3.2.3.1.2.1 Import AE SOP Specific Conformance - Verification

Import AE responds with the following status codes in response to a C-ECHO request.

Table 14

Service Status	Further Meaning	Protocol Codes	Description
Success	Success	0x0000	Operation performed properly

3.2.3.1.3 Import AE Presentation Context Acceptance Criterion - Verification

Import AE accepts the Presentation Contexts listed in the Presentation Context Table (Table 13).

3.2.3.1.4 Import AE Transfer Syntax Selection Policies - Verification

Import AE supports only the Implicit VR Little Endian transfer syntax. It rejects any proposed Presentation Context which does not specify the default Implicit VR Little Endian transfer syntax.

3.2.3.2 Import AE Real-World Activity - Storage

3.2.3.2.1 Import AE Associated Real-World Activity - Storage

The TWS receives image data sent by a remote Storage SCU, archives it to the local file system, and responds to the remote Storage SCU.

3.2.3.2.2 Import AE Presentation Context Table - Storage

Import AE accepts the following presentation contexts shown below:

Table 15

Presentation Context Table					
Abstract Syntax		Transfer Syntax			Extended
Name	UID	Name List	UID List	Role	Negotiation
CR Image Storage	1.2.840.10008.5.1.4.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
NM Image Storage	1.2.840.10008.5.1.4.1.1.20	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
SC Image Storage	1.2.840.10008.5.1.4.1.1.7	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
US Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Extended	
Name	UID	Name List	UID List	Role	Negotiation
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
DX Image Storage	1.2.840.10008.5.1.4.1.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
MG Image Storage	1.2.840.10008.5.1.4.1.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
PET Image Storage	1.2.840.10008.5.1.4.1.1.128	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None

3.2.3.2.2.1 Import AE SOP Specific Conformance - Storage

Import AE responds with the following status codes in response to a C-STORE request.

Import AE achieves Level 2 Conformance as described in Part 4 of the DICOM V3.0 Standard document.

Table 16

Service Status	Further Meaning	Protocol Codes	Description
Success	Success	0x0000	Operation performed properly.
Error	Data Set does not match SOP Class	0xA900	SOP Class UID does not match.
	Cannot understand	0xC000	Invalid data set or unsupported extended character sets. (See section 7 'Support of Extended Character Sets'.)
Refused	Out of Resources	0xA700	Local resource is insufficient.

When the service status is "Refused", check the free space in the local file system of the remote Storage SCU before resending.

- If the free space is insufficient, create enough space.
- If there is enough free space, wait till all large applications that are running are terminated.

3.2.3.3 Import AE Presentation Context Acceptance Criterion - Storage

Import AE accepts the Presentation Contexts listed in the Presentation Context Table (Table 15).

3.2.3.4 Import AE Transfer Syntax Selection Policies - Storage

Import AE accepts the Transfer Syntax listed in the Presentation Context Table (Table 15).

The selection priority of acceptable Transfer Syntax is the Default Transfer Syntax. See subsection 6.4.2.

4 Communication Profiles

4.1 Supported Communication Stacks

This product provides DICOM TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard.

4.2 OSI Stack

Not applicable to this product.

4.3 TCP/IP Stack

This product inherits its TCP/IP stack from the computer system upon which it executes.

4.3.1 API

Not applicable to this product.

4.3.2 Physical Media Support

This product is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the computer system upon which it executes.

4.4 Point-to-Point Stack

Not applicable to this product.

5 Extensions/Specializations/Privatizations

5.1 Standard Extended/Specialized/Private SOP - Storage SCU

Private Data Elements used in this product are listed in section 8.3

6 Configuration

For the TWS, the configuration can be set using the Network setup tool.

Note: Settings are performed by your Toshiba Service Personnel at the time of installation of the TWS.

6.1 AE Title/Presentation Address Mapping

The mapping from the AE titles to the presentation addresses is as follows:

- One port number and one AE title can be described for one host name.
- Each AE title is mapped to one port number.
- The TWS has following default values:
 - Local Port No. 5000
 - Local AE Title TM_OT_PTWS_1.0

6.2 Configurable Parameters

6.2.1 Time-out Value, Retry Count, Retry Interval

The time-out value, retry count, and retry interval in each status are shown below.




Figure 2

Table 17

Status	Item	Time-out value	Retry count	Retry interval	Remarks
S1	Association establishment request waiting time	Default: 60 seconds Range: 1 to 32767 seconds	Not set	Not set	Only one parameter can be set in the TWS.
S2	Association establishment response waiting time	Default: 60 seconds Range: 1 to 32767 seconds	Not set	Not set	Only one parameter can be set in the TWS.
S3	Service request waiting time	Default: 60 seconds Range: 1 to 32767 seconds	Not set	Not set	Only one parameter can be set in the TWS.
S4	Service response waiting time	Default: C-STORE=120seconds C-FIND=120 seconds C-MOVE=900 seconds Range: 1 to 32767 seconds	Not set	Not set	Can be set for each provided service.
S5	Association release waiting time	Default: 60 seconds Range: 1 to 32767 seconds	Not set	Not set	Only one parameter can be set in the TWS.

6.2.2 Warning Status Criteria

The warning status criteria can be set for each station and each service, for Export AE.

6.2.2.1 CR Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 18

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.2 CT Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 19

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.3 MR Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 20

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.4 NM Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 21

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.5 SC Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 22

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.6 US Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 23

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.7 US Multi-frame Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 24

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.8 XA Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 25

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.9 XRF Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 26

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.10 DX Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 27

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL

Elements discarded	FAIL	SUCCESS or FAIL
--------------------	------	-----------------

6.2.2.11 MG Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 28

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.12 PET Image Storage

If SUCCESS is set, the TWS judges that the image transfer succeeded.

If FAIL is set, the TWS judges that the image transfer failed.

Table 29

Warning response	Default	Parameter setting range
Coercion of Data Elements	FAIL	SUCCESS or FAIL
Data Set does not match SOP Class	FAIL	SUCCESS or FAIL
Elements discarded	FAIL	SUCCESS or FAIL

6.2.2.13 Query/Retrieve(Move)

If SUCCESS is set, the TWS judges that C-MOVE request succeeded.

If FAIL is set, the TWS judges that the C-MOVE request failed.

Table 30

Warning response	Default	Parameter setting range
Sub-operations Complete-One or more Failures	FAIL	SUCCESS or FAIL

6.2.2.14 Basic Grayscale Print Management

If SUCCESS is set, the TWS judges that the request printing images succeeded.

If FAIL is set, the TWS judges that the request printing images failed.

6.2.2.14.1 Basic Film Session SOP Class

6.2.2.14.1.1 N-CREATE response

Table 31

Warning response	Default	Parameter setting range
Memory allocation not supported.	FAIL	SUCCESS or FAIL

6.2.2.14.2 Basic Film Box SOP Class

6.2.2.14.2.1 N-CREATE response

Table 32

Warning response	Default	Parameter setting range
Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead.	FAIL	SUCCESS or FAIL

6.2.2.14.2.2 N-ACTION response

Table 33

Warning response	Default	Parameter setting range
Film Box SOP Instance hierarchy does not contain Image Box SOP Instances(empty page)	FAIL	SUCCESS or FAIL
Image size is larger than image box size, the image has been demagnified.	FAIL	SUCCESS or FAIL

6.2.2.14.3 Basic Grayscale Image Box SOP Class

6.2.2.14.3.1 N-SET response

Table 34

Warning response	Default	Parameter setting range
Image size larger than image box size, the image has been demagnified.	FAIL	SUCCESS or FAIL
Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead.	FAIL	SUCCESS or FAIL

6.2.2.15 Basic Color Print Management

If SUCCESS is set, the TWS judges that the request printing images succeeded.

If FAIL is set, the TWS judges that the request printing images failed.

6.2.2.15.1 Basic Film Session SOP Class

The warning status criteria are same as Basic Grayscale Print Management.

6.2.2.15.2 Basic Film Box SOP Class

The warning status criteria are same as Basic Grayscale Print Management.

6.2.2.15.3 Basic Color Image Box SOP Class

6.2.2.15.3.1 N-SET response

Table 35

Warning response item	Default	Parameter setting range
Image size larger than image box size, the image has been demagnified.	FAIL	SUCCESS or FAIL

6.3 Implementation Information and Maximum Reception PDU Size

The default values for the TWS are used for the Implementation Class UID, the Implementation Version name, and the Maximum length received. They cannot be changed.

Table 36

Parameter	Default
Implementation Class UID	1.2.392.200036.9116.7.4.10
Implementation Version Name	TM_OT_PTWS_1.0
Maximum length received	0x10000 (64 Kbytes)

6.4 Default Transfer Syntax

6.4.1 Import AE

Selection priority of acceptable Transfer Syntax is following Default Transfer Syntax:

Default = " Implicit VR Little Endian "

7 Support of Extended Character Sets

7.1 Query/Retrieve SCU

This product supports the following character sets:

- ISO-IR 6 (default) ISO 646

If Export AE receives result of a search that contains characters from an unsupported character set, a character may not be correctly displayed on a list.

7.2 Storage SCU/SCP

This product supports the following character sets:

- ISO-IR 6 (default) ISO 646
- ISO-IR 13 (Japanese) JIS X 0201 (Katakana)
- ISO-IR 14 (Japanese) JIS X 0201 (Romaji)
- ISO-IR 87 (Japanese) JIS X 0208 (Kanji)

If Import AE receives image data that contains characters from an unsupported character set, Import AE will respond with "Cannot understand" to the C-STORE request. (See subsection 3.2.3.2.2.1.)

7.3 Print SCU

This product supports the following character sets:

- ISO-IR 6 (default) ISO 646

8 Information Object Definition - Storage SCU

8.1 Entity Module Definitions

The information modules for TWS are defined below.

8.1.1 CR IMAGE IOD Modules

Table 37

Information Entity	Module	Reference	Usage¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
	CR Series Module	8.2.5	M
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	CR Image Module	8.2.26	M
	Overlay Plane Module	8.2.45	U
	Curve Module	8.2.47	U
	Modality LUT Module	8.2.48	U
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M

¹ M=Mandatory, C=Conditional, U=User option

8.1.2 CT IMAGE IOD Modules

Table 38

Information Entity	Module	Reference	Usage¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Frame of Reference	Frame of Reference Module	8.2.7	M
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Plane Module	8.2.12	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	CT Image Module	8.2.27	M
	Overlay Plane Module	8.2.45	U
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M

¹ M=Mandatory, C=Conditional, U=User option

8.1.3 MR IMAGE IOD Modules

Table 39

Information Entity	Module	Reference	Usage¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Frame of Reference	Frame of Reference Module	8.2.7	M
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Plane Module	8.2.12	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	MR Image Module	8.2.28	M
	Overlay Plane Module	8.2.45	U
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M

¹ M=Mandatory, C=Conditional, U=User option

8.1.4 NM IMAGE IOD Modules

Table 40

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
	NM/PET Patient Orientation Module	8.2.6	M
Frame of Reference	Frame of Reference Module	8.2.7	U
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	NM Image Pixel Module	8.2.14	M
	Multi-frame Module	8.2.19	M
	NM Multi-frame Module	8.2.20	M
	NM Image Module	8.2.29	M
	NM Isotope Module	8.2.31	M
	NM Detector Module	8.2.32	M
	NM TOMO Acquisition Module	8.2.33	C
	NM Multi-gated Acquisition Module	8.2.34	C
	NM Phase Module	8.2.35	C
	NM Reconstruction Module	8.2.36	C
	Overlay Plane Module	8.2.45	U
	Multi-frame Overlay Module	8.2.46	U
	Curve Module	8.2.47	U
VOI LUT Module	8.2.49	U	
SOP Common Module	8.2.50	M	

¹ M=Mandatory, C=Conditional, U=User option

8.1.5 SC IMAGE IOD Modules

Table 41

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Equipment	General Equipment Module	8.2.9	U
	SC Equipment Module	8.2.10	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	SC Image Module	8.2.37	M
	Overlay Plane Module	8.2.45	U
	Modality LUT Module	8.2.48	U
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M

¹ M=Mandatory, C=Conditional, U=User option

8.1.6 US IMAGE IOD Modules

Table 42

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Frame of Reference	Frame of Reference Module	8.2.7	U
	US Frame of Reference Module	8.2.8	C
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	Palette Color Lookup Table Module	8.2.16	C
	US Region Calibration Module	8.2.17	U
	US Image	8.2.30	M
	Overlay Plane Module	8.2.45	U
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M
Curve	Curve Module	8.2.47	M
	SOP Common Module	8.2.50	M
	Curve Identification Module	8.2.51	M
	Audio Module	8.2.52	U

¹ M=Mandatory, C=Conditional, U=User option

8.1.7 US MULTI-FRAME IMAGE IOD Modules

Table 43

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Frame of Reference	Frame of Reference Module	8.2.7	U
	US Frame of Reference Module	8.2.8	C
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	Cine Module	8.2.18	M
	Multi-frame Module	8.2.19	M
	Palette Color Lookup Table Module	8.2.16	C
	US Region Calibration Module	8.2.17	U
	US Image Module	8.2.30	M
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M
Curve	Curve Module	8.2.47	M
	SOP Common Module	8.2.50	M
	Curve Identification Module	8.2.51	M
	Audio Module	8.2.52	U

¹ M=Mandatory, C=Conditional, U=User option

8.1.8 XA IMAGE IOD Modules

Table 44

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	Cine Module	8.2.18	C
	Multi-frame Module	8.2.19	C
	Frame Pointers Module	8.2.21	U
	Mask Module	8.2.22	C
	Display Shutter Module	8.2.23	U
	Device Module	8.2.24	U
	Therapy Module	8.2.25	U
	X-Ray Image Module	8.2.38	M
	X-Ray Acquisition Module	8.2.39	M
	X-Ray Collimator Module	8.2.40	U
	X-Ray Table Module	8.2.41	C
	XA Positioner Module	8.2.44	M
	Overlay Plane Module	8.2.45	U
	Multi-frame Overlay Module	8.2.46	C
	Curve Module	8.2.47	U
Modality LUT Module	8.2.48	C/U	
VOI LUT Module	8.2.49	U	
SOP Common Module	8.2.50	M	

¹ M=Mandatory, C=Conditional, U=User option

8.1.9 XRF IMAGE IOD Modules

Table 45

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	C
	Cine Module	8.2.18	C
	Multi-frame Module	8.2.19	C
	Frame Pointers Module	8.2.21	U
	Mask Module	8.2.22	C
	X-Ray Image Module	8.2.38	M
	X-Ray Acquisition Module	8.2.39	M
	X-Ray Collimator Module	8.2.40	U
	Display Shutter Module	8.2.23	U
	Therapy Module	8.2.25	U
	Device Module	8.2.24	U
	X-Ray Table Module	8.2.41	U
	XRF Positioner Module	8.2.42	U
	XRF Tomo Acquisition Module	8.2.43	C
	Overlay Plane Module	8.2.45	U
	Multi-frame Overlay Module	8.2.46	C
	Curve Module	8.2.47	U
	Modality LUT Module	8.2.48	C/U
VOI LUT Module	8.2.49	U	
SOP Common Module	8.2.50	M	

¹ M=Mandatory, C=Conditional, U=User option

8.1.10 DX IMAGE IOD Modules

Table 46

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
	Specimen Identification Module	8.2.54	U
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
	DX Series Module	8.2.55	M
Frame of Reference	Frame of Reference Module	8.2.7	U
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	U
	Display Shutter Module	8.2.23	U
	Device Module	8.2.24	U
	Therapy Module	8.2.25	U
	DX Anatomy Imaged Module	8.2.56	M
	DX Image Module	8.2.57	M
	DX Detector Module	8.2.58	M
	X-Ray Collimator Module	8.2.40	U
	DX Positioning Module	8.2.59	U
	X-Ray Tomography Acquisition Module	8.2.43	U
	X-Ray Acquisition Dose Module	8.2.60	U
	X-Ray Generation Module	8.2.61	U
	X-Ray Filtration Module	8.2.62	U
	X-Ray Grid Module	8.2.63	U
	Overlay Plane Module	8.2.45	C
	Curve Module	8.2.47	U
	VOI LUT Module	8.2.49	C
	Image Histogram Module	8.2.65	U
Acquisition Context Module	8.2.66	M	
SOP Common Module	8.2.50	M	

¹ M=Mandatory, C=Conditional, U=User option

8.1.11 MG IMAGE IOD Modules

Table 47

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
	Specimen Identification Module	8.2.54	U
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
	DX Series Module	8.2.55	M
	Mammography Series Module	8.2.67	M
Frame of Reference	Frame of Reference Module	8.2.7	C
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Pixel Module	8.2.13	M
	Contrast/Bolus Module	8.2.15	U
	Display Shutter Module	8.2.23	U
	Device Module	8.2.24	U
	Therapy Module	8.2.25	U
	DX Anatomy Imaged Module	8.2.56	M
	DX Image Module	8.2.57	M
	DX Detector Module	8.2.58	M
	X-Ray Collimator Module	8.2.40	U
	DX Positioning Module	8.2.59	U
	X-Ray Tomography Acquisition Module	8.2.43	U
	X-Ray Acquisition Dose Module	8.2.60	U
	X-Ray Generation Module	8.2.61	U
	X-Ray Filtration Module	8.2.62	U
	X-Ray Grid Module	8.2.63	U
	Mammography Image Module	8.2.68	M
	Overlay Plane Module	8.2.45	C
	Curve Module	8.2.47	U
	VOI LUT Module	8.2.49	C
Image Histogram Module	8.2.65	U	
Acquisition Context Module	8.2.66	M	
SOP Common Module	8.2.50	M	

¹ M=Mandatory, C=Conditional, U=User option

8.1.12 PET IMAGE IOD Modules

Table 48

Information Entity	Module	Reference	Usage ¹
Patient	Patient Module	8.2.1	M
Study	General Study Module	8.2.2	M
	Patient Study Module	8.2.3	U
Series	General Series Module	8.2.4	M
	PET Series Module	8.2.69	M
	PET Isotope Module	8.2.70	M
	PET Multi-gated Acquisition Module	8.2.71	C
	NM/PET Patient Orientation Module	8.2.72	M
Frame of Reference	Frame of Reference Module	8.2.7	M
Equipment	General Equipment Module	8.2.9	M
Image	General Image Module	8.2.11	M
	Image Plane Module	8.2.12	M
	Image Pixel Module	8.2.13	M
	PET Image Module	8.2.73	M
	Overlay Plane Module	8.2.45	U
	VOI LUT Module	8.2.49	U
	SOP Common Module	8.2.50	M

¹ M=Mandatory, C=Conditional, U=User option

8.2 Information Object Definitions

8.2.1 Patient Module

Table 49

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	2	Length=0 when no entry is made
Patient ID	(0010,0020)	2	Length=0 when no entry is made
Patient's Birth Date	(0010,0030)	2	Length=0 when no entry is made
Patient's Sex	(0010,0040)	2	Length=0 when no entry is made
Referenced Patient Sequence	(0008,1120)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Patient's Birth Time	(0010,0032)	3	Not set when no entry is made
Other Patient IDs	(0010,1000)	3	Not set when no entry is made
Other Patient Names	(0010,1001)	3	Not set when no entry is made
Ethnic Group	(0010,2160)	3	Not set when no entry is made
Patient Comments	(0010,4000)	3	Not set when no entry is made

8.2.2 General Study Module

Table 50

Attribute Name	Tag	Type	Attribute Description
Study Instance UID	(0020,000D)	1	Always set
Study Date	(0008,0020)	2	Length=0 when no entry is made
Study Time	(0008,0030)	2	Length=0 when no entry is made
Referring Physician's Name	(0008,0090)	2	Length=0 when no entry is made
Study ID	(0020,0010)	2	Length=0 when no entry is made
Accession Number	(0008,0050)	2	Length=0 when no entry is made
Study Description	(0008,1030)	3	Not set when no entry is made
Physician(s) of Record	(0008,1048)	3	Not set when no entry is made
Name of Physician(s) Reading Study	(0008,1060)	3	Not set when no entry is made
Referenced Study Sequence	(0008,1110)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set

8.2.3 Patient Study Module

Table 51

Attribute Name	Tag	Type	Attribute Description
Admitting Diagnoses Description	(0008,1080)	3	Not set when no entry is made
Patient's Age	(0010,1010)	3	Not set when no entry is made
Patient's Size	(0010,1020)	3	Not set when no entry is made
Patient's Weight	(0010,1030)	3	Not set when no entry is made
Occupation	(0010,2180)	3	Not set when no entry is made
Additional Patient's History	(0010,21B0)	3	Not set when no entry is made

8.2.4 General Series Module

Table 52

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	Always set
Series Instance UID	(0020,000E)	1	Always set
Series Number	(0020,0011)	2	Length=0 when no entry is made
Laterality	(0020,0060)	2C	If the setting conditions are met, Length=0 when no entry is made
Series Date	(0008,0021)	3	Not set when no entry is made
Series Time	(0008,0031)	3	Not set when no entry is made
Performing Physicians' Name	(0008,1050)	3	Not set when no entry is made
Protocol Name	(0018,1030)	3	Not set when no entry is made
Series Description	(0008,103E)	3	Not set when no entry is made
Operators' Name	(0008,1070)	3	Not set when no entry is made
Referenced Study Component Sequence	(0008,1111)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Body Part Examined	(0018,0015)	3	Not set when no entry is made
Patient Position	(0018,5100)	2C	If the setting conditions are met, Length=0 when no entry is made
Smallest Pixel Value in Series	(0028,0108)	3	Not set when no entry is made
Largest Pixel Value in Series	(0028,0109)	3	Not set when no entry is made

8.2.5 CR Series Module

Table 53

Attribute Name	Tag	Type	Attribute Description
Body Part Examined	(0018,0015)	2	Length=0 when no entry is made
View Position	(0018,5101)	2	Length=0 when no entry is made
Filter Type	(0018,1160)	3	Not set when no entry is made
Collimator/grid Name	(0018,1180)	3	Not set when no entry is made
Focal Spot(s)	(0018,1190)	3	Not set when no entry is made
Plate Type	(0018,1260)	3	Not set when no entry is made
Phosphor Type	(0018,1261)	3	Not set when no entry is made

8.2.6 NM/PET Patient Orientation Module

Table 54

Attribute Name	Tag	Type	Attribute Description
Patient Orientation Code Sequence	(0054,0410)	2	Length=0 when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Patient Orientation Modifier Code Sequence	(0054,0412)	2C	If the setting conditions are met, Length=0 when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
Patient Gantry Relationship Code Sequence	(0054,0414)	2	Length=0 when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made

8.2.7 Frame of Reference Module

Table 55

Attribute Name	Tag	Type	Attribute Description
Frame of Reference UID	(0020,0052)	1	Always set
Position Reference Indicator	(0020,1040)	2	Length=0 when no entry is made

8.2.8 US Frame of Reference Module

Table 56

Attribute Name	Tag	Type	Attribute Description
Region Location Min x_0	(0018,6018)	1	Always set
Region Location Min y_0	(0018,601A)	1	Always set
Region Location Max x_1	(0018,601C)	1	Always set
Region Location Max y_1	(0018,601E)	1	Always set
Physical Units X Direction	(0018,6024)	1	Always set
Physical Units Y Direction	(0018,6026)	1	Always set
Physical Delta X	(0018,602C)	1	Always set
Physical Delta Y	(0018,602E)	1	Always set
Reference Pixel x_0	(0018,6020)	3	Not set when no entry is made
Reference Pixel y_0	(0018,6022)	3	Not set when no entry is made
Ref. Pixel Physical Value X	(0018,6028)	3	Not set when no entry is made
Ref. Pixel Physical Value Y	(0018,602A)	3	Not set when no entry is made

8.2.9 General Equipment Module

Table 57

Attribute Name	Tag	Type	Attribute Description
Manufacturer	(0008,0070)	2	Length=0 when no entry is made
Institution Name	(0008,0080)	3	Not set when no entry is made
Institution Address	(0008,0081)	3	Not set when no entry is made
Station Name	(0008,1010)	3	Not set when no entry is made
Institutional Department Name	(0008,1040)	3	Not set when no entry is made
Manufacturer's Model Name	(0008,1090)	3	Not set when no entry is made
Device Serial Number	(0018,1000)	3	Not set when no entry is made
Software Versions	(0018,1020)	3	Not set when no entry is made
Spatial Resolution	(0018,1050)	3	Not set when no entry is made
Date of Last Calibration	(0018,1200)	3	Not set when no entry is made
Time of Last Calibration	(0018,1201)	3	Not set when no entry is made
Pixel Padding Value	(0028,0120)	3	Not set when no entry is made

8.2.10 SC Equipment Module

Table 58

Attribute Name	Tag	Type	Attribute Description
Conversion Type	(0008,0064)	1	Always set
Modality	(0008,0060)	3	Not set when no entry is made
Secondary Capture Device ID	(0018,1010)	3	Not set when no entry is made
Secondary Capture Device Manufacturer	(0018,1016)	3	Not set when no entry is made
Secondary Capture Device Manufacturer's Model Name	(0018,1018)	3	Not set when no entry is made
Secondary Capture Device Software Version(s)	(0018,1019)	3	Not set when no entry is made
Video Image Format Acquired	(0018,1022)	3	Not set when no entry is made
Digital Image Format Acquired	(0018,1023)	3	Not set when no entry is made

8.2.11 General Image Module

Table 59

Attribute Name	Tag	Type	Attribute Description
Image Number	(0020,0013)	2	Length=0 when no entry is made
Patient Orientation	(0020,0020)	2C	If the setting conditions are met, Length=0 when no entry is made
Image Date	(0008,0023)	2C	If the setting conditions are met, Length=0 when no entry is made
Image Time	(0008,0033)	2C	If the setting conditions are met, Length=0 when no entry is made
Image Type	(0008,0008)	3	Always set
Acquisition Number	(0020,0012)	3	Not set when no entry is made
Acquisition Date	(0008,0022)	3	Not set when no entry is made
Acquisition Time	(0008,0032)	3	Not set when no entry is made
Referenced Image Sequence	(0008,1140)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Derivation Description	(0008,2111)	3	Not set when no entry is made
Source Image Sequence	(0008,2112)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Images in Acquisition	(0020,1002)	3	Not set when no entry is made
Image Comments	(0020,4000)	3	Not set when no entry is made
Lossy Image Compression	(0028,2110)	3	Not set when no entry is made

8.2.12 Image Plane Module**Table 60**

Attribute Name	Tag	Type	Attribute Description
Pixel Spacing	(0028,0030)	1	Always set
Image Orientation (Patient)	(0020,0037)	1	Always set
Image Position (Patient)	(0020,0032)	1	Always set
Slice Thickness	(0018,0050)	2	Length=0 when no entry is made
Slice Location	(0020,1041)	3	Not set when no entry is made

8.2.13 Image Pixel Module

Table 61

Attribute Name	Tag	Type	Attribute Description
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Rows	(0028,0010)	1	Always set
Columns	(0028,0011)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Pixel Representation	(0028,0103)	1	Always set
Pixel Data	(7FE0,0010)	1	Always set
Planar Configuration	(0028,0006)	1C	Always set when the setting conditions are met
Pixel Aspect Ratio	(0028,0034)	1C	Always set when the setting conditions are met
Smallest Image Pixel Value	(0028,0106)	3	Not set when no entry is made
Largest Image Pixel Value	(0028,0107)	3	Not set when no entry is made
Red Palette Color Lookup Table Descriptor	(0028,1101)	1C	Always set when the setting conditions are met
Green Palette Color Lookup Table Descriptor	(0028,1102)	1C	Always set when the setting conditions are met
Blue Palette Color Lookup Table Descriptor	(0028,1103)	1C	Always set when the setting conditions are met
Red Palette Color Lookup Table Data	(0028,1201)	1C	Always set when the setting conditions are met
Green Palette Color Lookup Table Data	(0028,1202)	1C	Always set when the setting conditions are met
Blue Palette Color Lookup Table Data	(0028,1203)	1C	Always set when the setting conditions are met

8.2.14 NM Image Pixel Module**Table 62**

Attribute Name	Tag	Type	Attribute Description
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Pixel Spacing	(0028,0030)	2	Length=0 when no entry is made

8.2.15 Contrast/Bolus Module

Table 63

Attribute Name	Tag	Type	Attribute Description
Contrast/Bolus Agent	(0018,0010)	2	Length=0 when no entry is made
Contrast/Bolus Agent Sequence	(0018,0012)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
Contrast/Bolus Route	(0018,1040)	3	Not set when no entry is made
Contrast/Bolus Administration Route Sequence	(0018,0014)	3	Not set when no entry is made
>Code Value	(0008,0100)	1C	Always set when the setting conditions are met
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Additional Drug Sequence	(0018,002A)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
Contrast/Bolus Volume	(0018,1041)	3	Not set when no entry is made
Contrast/Bolus Start Time	(0018,1042)	3	Not set when no entry is made
Contrast/Bolus Stop Time	(0018,1043)	3	Not set when no entry is made
Contrast/Bolus Total Dose	(0018,1044)	3	Not set when no entry is made
Contrast Flow Rate(s)	(0018,1046)	3	Not set when no entry is made
Contrast Flow Duration(s)	(0018,1047)	3	Not set when no entry is made
Contrast/Bolus Ingredient	(0018,1048)	3	Not set when no entry is made
Contrast/Bolus Ingredient Concentration	(0018,1049)	3	Not set when no entry is made

8.2.16 Palette Color Lookup Table Module

Table 64

Attribute Name	Tag	Type	Attribute Description
Red Palette Color Lookup Table Descriptor	(0028,1101)	1C	Always set when the setting conditions are met
Green Palette Color Lookup Table Descriptor	(0028,1102)	1C	Always set when the setting conditions are met
Blue Palette Color Lookup Table Descriptor	(0028,1103)	1C	Always set when the setting conditions are met
Palette Color Lookup Table UID	(0028,1199)	3	Not set when no entry is made
Red Palette Color Lookup Table Data	(0028,1201)	1C	Always set when the setting conditions are met
Green Palette Color Lookup Table Data	(0028,1202)	1C	Always set when the setting conditions are met
Blue Palette Color Lookup Table Data	(0028,1203)	1C	Always set when the setting conditions are met
Segmented Red Palette Color Lookup Table Data	(0028,1221)	1C	Always set when the setting conditions are met
Segmented Green Palette Color Lookup Table Data	(0028,1222)	1C	Always set when the setting conditions are met
Segmented Blue Palette Color Lookup Table Data	(0028,1223)	1C	Always set when the setting conditions are met

8.2.17 US Region Calibration Module

Table 65

Attribute Name	Tag	Type	Attribute Description
Sequence of Ultrasound Regions	(0018,6011)	1	Always set
>Region Location Min x_0	(0018,6018)	1	Always set
>Region Location Min y_0	(0018,601A)	1	Always set
>Region Location Max x_1	(0018,601C)	1	Always set
>Region Location Max y_1	(0018,601E)	1	Always set
>Physical Units X Direction	(0018,6024)	1	Always set
>Physical Units Y Direction	(0018,6026)	1	Always set
>Physical Delta X	(0018,602C)	1	Always set
>Physical Delta Y	(0018,602E)	1	Always set
>Reference Pixel x_0	(0018,6020)	3	Not set when no entry is made
>Reference Pixel y_0	(0018,6022)	3	Not set when no entry is made
>Ref. Pixel Physical Value X	(0018,6028)	3	Not set when no entry is made
>Ref. Pixel Physical Value Y	(0018,602A)	3	Not set when no entry is made
>Region Spatial Format	(0018,6012)	1	Always set
>Region Data Type	(0018,6014)	1	Always set
>Region Flags	(0018,6016)	1	Always set
>Pixel Component Organization	(0018,6044)	1C	Always set when the setting conditions are met
>Pixel Component Mask	(0018,6046)	1C	Always set when the setting conditions are met
>Pixel Component Range Start	(0018,6048)	1C	Always set when the setting conditions are met
>Pixel Component Range Stop	(0018,604A)	1C	Always set when the setting conditions are met
>Pixel Component Physical Units	(0018,604C)	1C	Always set when the setting conditions are met
>Pixel Component Data Type	(0018,604E)	1C	Always set when the setting conditions are met
>Number of Table Break Points	(0018,6050)	1C	Always set when the setting conditions are met
>Table of X Break Points	(0018,6052)	1C	Always set when the setting conditions are met
>Table of Y Break Points	(0018,6054)	1C	Always set when the setting conditions are met
>Number of Table Entries	(0018,6056)	1C	Always set when the setting conditions are met

Attribute Name	Tag	Type	Attribute Description
>Table of Pixel Values	(0018,6058)	1C	Always set when the setting conditions are met
>Table of Parameter Values	(0018,605A)	1C	Always set when the setting conditions are met
>Transducer Frequency	(0018,6030)	3	Not set when no entry is made
>Pulse Repetition Frequency	(0018,6032)	3	Not set when no entry is made
>Doppler Correction Angle	(0018,6034)	3	Not set when no entry is made
>Steering Angle	(0018,6036)	3	Not set when no entry is made
>Doppler Sample Volume X Position	(0018,6038)	3	Not set when no entry is made
>Doppler Sample Volume Y Position	(0018,603A)	3	Not set when no entry is made
>TM-Line Position x_0	(0018,603C)	3	Not set when no entry is made
>TM-Line Position y_0	(0018,603E)	3	Not set when no entry is made
>TM-Line Position x_1	(0018,6040)	3	Not set when no entry is made
>TM-Line Position y_1	(0018,6042)	3	Not set when no entry is made

8.2.18 Cine Module

Table 66

Attribute Name	Tag	Type	Attribute Description
Preferred Playback Sequencing	(0018,1244)	3	Not set when no entry is made
Frame Time	(0018,1063)	1C	Always set when the setting conditions are met
Frame Time Vector	(0018,1065)	1C	Always set when the setting conditions are met
Start Trim	(0008,2142)	3	Not set when no entry is made
Stop Trim	(0008,2143)	3	Not set when no entry is made
Recommended Display Frame Rate	(0008,2144)	3	Not set when no entry is made
Cine Rate	(0018,0040)	3	Not set when no entry is made
Frame Delay	(0018,1066)	3	Not set when no entry is made
Effective Series Duration	(0018,0072)	3	Not set when no entry is made
Actual Frame Duration	(0018,1242)	3	Not set when no entry is made

8.2.19 Multi-Frame Module

Table 67

Attribute Name	Tag	Type	Attribute Description
Number of Frames	(0028,0008)	1	Always set
Frame Increment Pointer	(0028,0009)	1	Always set

8.2.20 NM Multi-Frame Module

Table 68

Attribute Name	Tag	Type	Attribute Description
Frame Increment Pointer	(0028,0009)	1	Always set
Energy Window Vector	(0054,0010)	1C	Always set when the setting conditions are met
Number of Energy Windows	(0054,0011)	1	Always set
Detector Vector	(0054,0020)	1C	Always set when the setting conditions are met
Number of Detectors	(0054,0021)	1	Always set
Phase Vector	(0054,0030)	1C	Always set when the setting conditions are met
Number of Phases	(0054,0031)	1C	Always set when the setting conditions are met
Rotation Vector	(0054,0050)	1C	Always set when the setting conditions are met
Number of Rotations	(0054,0051)	1C	Always set when the setting conditions are met
R-R Interval Vector	(0054,0060)	1C	Always set when the setting conditions are met
Number of R-R Intervals	(0054,0061)	1C	Always set when the setting conditions are met
Time Slot Vector	(0054,0070)	1C	Always set when the setting conditions are met
Number of Time Slots	(0054,0071)	1C	Always set when the setting conditions are met
Slice Vector	(0054,0080)	1C	Always set when the setting conditions are met
Number of Slices	(0054,0081)	1C	Always set when the setting conditions are met
Angular View Vector	(0054,0090)	1C	Always set when the setting conditions are met
Time Slice Vector	(0054,0100)	1C	Always set when the setting conditions are met

8.2.21 Frame Pointers Module

Table 69

Attribute Name	Tag	Type	Attribute Description
Representative Frame Number	(0028,6010)	3	Not set when no entry is made
Frame Numbers of Interest(FOI)	(0028,6020)	3	Not set when no entry is made
Frame(s) of Interest Description	(0028,6022)	3	Not set when no entry is made

8.2.22 Mask Module

Table 70

Attribute Name	Tag	Type	Attribute Description
Mask Subtraction Sequence	(0028,6100)	1	Always set
> Mask Operation	(0028,6101)	1	Always set
> Applicable Frame Range	(0028,6102)	3	Not set when no entry is made
> Mask Frame Numbers	(0028,6110)	1C	Always set when the setting conditions are met
> Contrast Frame Averaging	(0028,6112)	3	Not set when no entry is made
> Mask Sub-pixel Shift	(0028,6114)	3	Not set when no entry is made
> TID Offset	(0028,6120)	2C	If the setting conditions are met, Length=0 when no entry is made
> Mask Operation Explanation	(0028,6190)	3	Not set when no entry is made
Recommended Viewing Mode	(0028,1090)	2	Length=0 when no entry is made

8.2.23 Display Shutter Module

Table 71

Attribute Name	Tag	Type	Attribute Description
Shutter Shape	(0018,1600)	1	Always set
Shutter Left Vertical Edge	(0018,1602)	1C	Always set when the setting conditions are met
Shutter Right Vertical Edge	(0018,1604)	1C	Always set when the setting conditions are met
Shutter Upper Horizontal Edge	(0018,1606)	1C	Always set when the setting conditions are met
Shutter Lower Horizontal Edge	(0018,1608)	1C	Always set when the setting conditions are met
Center of Circular Shutter	(0018,1610)	1C	Always set when the setting conditions are met
Radius of Circular Shutter	(0018,1612)	1C	Always set when the setting conditions are met
Vertices of the Polygonal Shutter	(0018,1620)	1C	Always set when the setting conditions are met

8.2.24 Device Module

Table 72

Attribute Name	Tag	Type	Attribute Description
Device Sequence	(0050,0010)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Device Length	(0050,0014)	3	Not set when no entry is made
> Device Diameter	(0050,0016)	3	Not set when no entry is made
> Device Diameter Units	(0050,0017)	2C	If the setting conditions are met, Length=0 when no entry is made
> Device Volume	(0050,0018)	3	Not set when no entry is made
> Inter-Marker Distance	(0050,0019)	3	Not set when no entry is made
> Device Description	(0050,0020)	3	Not set when no entry is made

8.2.25 Therapy Module

Table 73

Attribute Name	Tag	Type	Attribute Description
Interventional Therapy Sequence	(0018,0036)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Interventional Status	(0018,0038)	2	Length=0 when no entry is made
> Intervention Drug Code Sequence	(0018,0029)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Intervention Drug Start Time	(0018,0035)	3	Not set when no entry is made
> Intervention Drug Stop Time	(0018,0027)	3	Not set when no entry is made
> Administration Route Code Sequence	(0054,0302)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Therapy Description	(0018,0039)	3	Not set when no entry is made

8.2.26 CR Image Module**Table 74**

Attribute Name	Tag	Type	Attribute Description
KVP	(0018,0060)	3	Not set when no entry is made
Plate ID	(0018,1004)	3	Not set when no entry is made
Distance Source to Detector	(0018,1110)	3	Not set when no entry is made
Distance Source to Patient	(0018,1111)	3	Not set when no entry is made
Exposure Time	(0018,1150)	3	Not set when no entry is made
X-ray Tube Current	(0018,1151)	3	Not set when no entry is made
Exposure	(0018,1152)	3	Not set when no entry is made
Generator Power	(0018,1170)	3	Not set when no entry is made
Acquisition Device Processing Description	(0018,1400)	3	Not set when no entry is made
Acquisition Device Processing Code	(0018,1401)	3	Not set when no entry is made
Cassette Orientation	(0018,1402)	3	Not set when no entry is made
Cassette Size	(0018,1403)	3	Not set when no entry is made
Exposures on Plate	(0018,1404)	3	Not set when no entry is made
Relative X-ray Exposure	(0018,1405)	3	Not set when no entry is made
Sensitivity	(0018,6000)	3	Not set when no entry is made

8.2.27 CT Image Module

Table 75

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	Always set
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Rescale Intercept	(0028,1052)	1	Always set
Rescale Slope	(0028,1053)	1	Always set
KVP	(0018,0060)	2	Length=0 when no entry is made
Acquisition Number	(0020,0012)	2	Length=0 when no entry is made
Scan Options	(0018,0022)	3	Not set when no entry is made
Data Collection Diameter	(0018,0090)	3	Not set when no entry is made
Reconstruction Diameter	(0018,1100)	3	Not set when no entry is made
Distance Source to Detector	(0018,1110)	3	Not set when no entry is made
Distance Source to Patient	(0018,1111)	3	Not set when no entry is made
Gantry/Detector Tilt	(0018,1120)	3	Not set when no entry is made
Table Height	(0018,1130)	3	Not set when no entry is made
Rotation Direction	(0018,1140)	3	Not set when no entry is made
Exposure Time	(0018,1150)	3	Not set when no entry is made
X-ray Tube Current	(0018,1151)	3	Not set when no entry is made
Exposure	(0018,1152)	3	Not set when no entry is made
Filter Type	(0018,1160)	3	Not set when no entry is made
Generator Power	(0018,1170)	3	Not set when no entry is made
Focal Spot	(0018,1190)	3	Not set when no entry is made
Convolution Kernel	(0018,1210)	3	Not set when no entry is made

8.2.28 MR Image Module

Table 76

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	Always set
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Scanning Sequence	(0018,0020)	1	Always set
Sequence Variant	(0018,0021)	1	Always set
Scan Options	(0018,0022)	2	Length=0 when no entry is made
MR Acquisition Type	(0018,0023)	2	Length=0 when no entry is made
Repetition Time	(0018,0080)	2C	If the setting conditions are met, Length=0 when no entry is made
Echo Time	(0018,0081)	2	Length=0 when no entry is made
Echo Train Length	(0018,0091)	2	Length=0 when no entry is made
Inversion Time	(0018,0082)	2C	If the setting conditions are met, Length=0 when no entry is made
Trigger Time	(0018,1060)	2C	If the setting conditions are met, Length=0 when no entry is made
Sequence Name	(0018,0024)	3	Not set when no entry is made
Angio Flag	(0018,0025)	3	Not set when no entry is made
Number of Averages	(0018,0083)	3	Not set when no entry is made
Imaging Frequency	(0018,0084)	3	Not set when no entry is made
Imaged Nucleus	(0018,0085)	3	Not set when no entry is made
Echo Number	(0018,0086)	3	Not set when no entry is made
Magnetic Field Strength	(0018,0087)	3	Not set when no entry is made
Spacing Between Slices	(0018,0088)	3	Not set when no entry is made
Number of Phase Encoding Steps	(0018,0089)	3	Not set when no entry is made
Percent Sampling	(0018,0093)	3	Not set when no entry is made
Percent Phase Field of View	(0018,0094)	3	Not set when no entry is made
Pixel Bandwidth	(0018,0095)	3	Not set when no entry is made
Nominal Interval	(0018,1062)	3	Not set when no entry is made
Beat Rejection Flag	(0018,1080)	3	Not set when no entry is made
Low R-R Value	(0018,1081)	3	Not set when no entry is made
High R-R Value	(0018,1082)	3	Not set when no entry is made
Intervals Acquired	(0018,1083)	3	Not set when no entry is made
Intervals Rejected	(0018,1084)	3	Not set when no entry is made
PVC Rejection	(0018,1085)	3	Not set when no entry is made
Skip Beats	(0018,1086)	3	Not set when no entry is made

Attribute Name	Tag	Type	Attribute Description
Heart Rate	(0018,1088)	3	Not set when no entry is made
Cardiac Number of Images	(0018,1090)	3	Not set when no entry is made
Trigger Window	(0018,1094)	3	Not set when no entry is made
Reconstruction Diameter	(0018,1100)	3	Not set when no entry is made
Receiving Coil	(0018,1250)	3	Not set when no entry is made
Transmitting Coil	(0018,1251)	3	Not set when no entry is made
Acquisition Matrix	(0018,1310)	3	Not set when no entry is made
Phase Encoding Direction	(0018,1312)	3	Not set when no entry is made
Flip Angle	(0018,1314)	3	Not set when no entry is made
SAR	(0018,1316)	3	Not set when no entry is made
Variable Flip Angle Flag	(0018,1315)	3	Not set when no entry is made
dB/dt	(0018,1318)	3	Not set when no entry is made
Temporal Position Identifier	(0020,0100)	3	Not set when no entry is made
Number of Temporal Positions	(0020,0105)	3	Not set when no entry is made
Temporal Resolution	(0020,0110)	3	Not set when no entry is made

8.2.29 NM Image Module

Table 77

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	Always set
Image ID	(0054,0400)	3	Not set when no entry is made
Lossy Image Compression	(0028,2110)	1C	Always set when the setting conditions are met
Counts Accumulated	(0018,0070)	2	Length=0 when no entry is made
Acquisition Termination Condition	(0018,0071)	3	Not set when no entry is made
Table Height	(0018,1130)	3	Not set when no entry is made
Table Traverse	(0018,1131)	3	Not set when no entry is made
Actual Frame Duration	(0018,1242)	1C	Always set when the setting conditions are met
Count Rate	(0018,1243)	3	Not set when no entry is made
Processing Function	(0018,5020)	3	Not set when no entry is made
Corrected Image	(0028,0051)	3	Not set when no entry is made
Whole Body Technique	(0018,1301)	3	Not set when no entry is made
Scan Velocity	(0018,1300)	2C	If the setting conditions are met, Length=0 when no entry is made
Scan Length	(0018,1302)	2C	If the setting conditions are met, Length=0 when no entry is made
Referenced Overlay Sequence	(0008,1130)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Referenced Curve Sequence	(0008,1145)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Trigger Source or Type	(0018,1061)	3	Not set when no entry is made
Anatomic Region Sequence	(0008,2218)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Anatomic Region Modifier Sequence	(0008,2220)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
Primary Anatomic Structure Sequence	(0008,2228)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set

Attribute Name	Tag	Type	Attribute Description
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Primary Anatomic Structure Modifier Sequence	(0008,2230)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made

8.2.30 US Image Module

Table 78

Attribute Name	Tag	Type	Attribute Description
Samples Per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Planar Configuration	(0028,0006)	1C	Always set
Pixel Representation	(0028,0103)	1	Always set
Frame Increment Pointer	(0028,0009)	1C	Always set
Image Type	(0008,0008)	2	Length=0 when no entry is made
Lossy Image Compression	(0028,2110)	1C	Always set
Number of Stages	(0008,2124)	2C	If the setting conditions are met, Length=0 when no entry is made
Number of Views in Stage	(0008,212A)	2C	If the setting conditions are met, Length=0 when no entry is made
Ultrasound Color Data Present	(0028,0014)	3	Not set when no entry is made
Referenced Overlay Sequence	(0008,1130)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Referenced Curve Sequence	(0008,1145)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Stage Name	(0008,2120)	3	Not set when no entry is made
Stage Number	(0008,2122)	3	Not set when no entry is made
View Number	(0008,2128)	3	Not set when no entry is made
Number of Event Timers	(0008,2129)	3	Not set when no entry is made
Event Elapsed Time(s)	(0008,2130)	3	Not set when no entry is made
Event Timer Name(s)	(0008,2132)	3	Not set when no entry is made
Anatomic Region Sequence	(0008,2218)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
>Anatomic Region Modifier Sequence	(0008,2220)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made

Attribute Name	Tag	Type	Attribute Description
Primary Anatomic Structure Sequence	(0008,2228)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
>Primary Anatomic Structure Modifier Sequence	(0008,2230)	3	Not set when no entry is made
>>Code Value	(0008,0100)	1C	Always set when the sequence is set
>>Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>>Code Meaning	(0008,0104)	3	Not set when no entry is made
Transducer Position Sequence	(0008,2240)	3	Not set when no entry is made
>Code Value	(0008,0100)	1C	Always set when the sequence is set
>Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>Code Meaning	(0008,0104)	3	Not set when no entry is made
> Transducer Position Modifier Sequence	(0008,2242)	3	Not set when no entry is made
>>Code Value	(0008,0100)	1C	Always set when the sequence is set
>>Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>>Code Meaning	(0008,0104)	3	Not set when no entry is made
Transducer Orientation Sequence	(0008,2244)	3	Not set when no entry is made
>Code Value	(0008,0100)	1C	Always set when the sequence is set
>Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>Code Meaning	(0008,0104)	3	Not set when no entry is made
> Transducer Orientation Modifier Sequence	(0008,2246)	3	Not set when no entry is made
>>Code Value	(0008,0100)	1C	Always set when the sequence is set
>>Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>>Code Meaning	(0008,0104)	3	Not set when no entry is made
Trigger Time	(0018,1060)	3	Not set when no entry is made
Nominal Interval	(0018,1062)	3	Not set when no entry is made
Beat Rejection Flag	(0018,1080)	3	Not set when no entry is made
Low R-R Value	(0018,1081)	3	Not set when no entry is made
High R-R Value	(0018,1082)	3	Not set when no entry is made
Heart Rate	(0018,1088)	3	Not set when no entry is made
Output Power	(0018,5000)	3	Not set when no entry is made
Transducer Data	(0018,5010)	3	Not set when no entry is made
Transducer Type	(0018,6031)	3	Not set when no entry is made
Focus Depth	(0018,5012)	3	Not set when no entry is made
Preprocessing Function	(0018,5020)	3	Not set when no entry is made

Attribute Name	Tag	Type	Attribute Description
Mechanical Index	(0018,5022)	3	Not set when no entry is made
Bone Thermal Index	(0018,5024)	3	Not set when no entry is made
Cranial Thermal Index	(0018,5026)	3	Not set when no entry is made
Soft Tissue Thermal Index	(0018,5027)	3	Not set when no entry is made
Soft Tissue-focus Thermal Index	(0018,5028)	3	Not set when no entry is made
Soft Tissue-surface Thermal Index	(0018,5029)	3	Not set when no entry is made
Depth of Scan Field	(0018,5050)	3	Not set when no entry is made
Image Transformation Matrix	(0018,5210)	3	Not set when no entry is made
Image Translation Vector	(0018,5212)	3	Not set when no entry is made
Overlay Subtype	(60xx,0045)	3	Not set when no entry is made

8.2.31 NM Isotope Module

Table 79

Attribute Name	Tag	Type	Attribute Description
Energy Window Information Sequence	(0054,0012)	2	Length=0 when no entry is made
> Energy Window Name	(0054,0018)	3	Not set when no entry is made
> Energy Window Range Sequence	(0054,0013)	3	Not set when no entry is made
>> Energy Window Lower Limit	(0054,0014)	3	Not set when no entry is made
>> Energy Window Upper Limit	(0054,0015)	3	Not set when no entry is made
Radiopharmaceutical Information Sequence	(0054,0016)	2	Length=0 when no entry is made
> Radionuclide Code Sequence	(0054,0300)	2C	If the setting conditions are met, Length=0 when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Radiopharmaceutical Route	(0018,1070)	3	Not set when no entry is made
> Administration Route Code Sequence	(0054,0302)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Radiopharmaceutical Volume	(0018,1071)	3	Not set when no entry is made
> Radiopharmaceutical Start Time	(0018,1072)	3	Not set when no entry is made
> Radiopharmaceutical Stop Time	(0018,1073)	3	Not set when no entry is made
> Radionuclide Total Dose	(0018,1074)	3	Not set when no entry is made
> Calibration Data Sequence	(0054,0306)	3	Not set when no entry is made
>> Energy Window Number	(0054,0308)	1C	Always set when the sequence is set
>> Syringe Counts	(0018,1045)	3	Not set when no entry is made
>> Residual Syringe Counts	(0054,0017)	3	Not set when no entry is made
> Radiopharmaceutical	(0018,0031)	3	Not set when no entry is made
> Radiopharmaceutical Code Sequence	(0054,0304)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
Intervention Drug Information Sequence	(0018,0026)	3	Not set when no entry is made
> Intervention Drug Name	(0018,0034)	3	Not set when no entry is made

Attribute Name	Tag	Type	Attribute Description
> Intervention Drug Code Sequence	(0018,0029)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Administration Route Code Sequence	(0054,0302)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Intervention Drug Start Time	(0018,0035)	3	Not set when no entry is made
> Intervention Drug Stop Time	(0018,0027)	3	Not set when no entry is made
> Intervention Drug Dose	(0018,0028)	3	Not set when no entry is made

8.2.32 NM Detector Module

Table 80

Attribute Name	Tag	Type	Attribute Description
Detector Information Sequence	(0054,0022)	2	Length=0 when no entry is made
> Collimator/grid Name	(0018,1180)	3	Not set when no entry is made
> Collimator Type	(0018,1181)	2C	If the setting conditions are met, Length=0 when no entry is made
> Field of View Shape	(0018,1147)	3	Not set when no entry is made
> Field of View Dimension(s)	(0018,1149)	3	Not set when no entry is made
> Focal Distance	(0018,1182)	2C	If the setting conditions are met, Length=0 when no entry is made
> X Focus Center	(0018,1183)	3	Not set when no entry is made
> Y Focus Center	(0018,1184)	3	Not set when no entry is made
> Zoom Center	(0028,0032)	3	Not set when no entry is made
> Zoom Factor	(0028,0031)	3	Not set when no entry is made
> Center of Rotation Offset	(0018,1145)	3	Not set when no entry is made
> Gantry/Detector Tilt	(0018,1120)	3	Not set when no entry is made
> Distance Source to Detector	(0018,1110)	2C	If the setting conditions are met, Length=0 when no entry is made
> Start Angle	(0054,0200)	3	Not set when no entry is made
> Radial Position	(0018,1142)	3	Not set when no entry is made
> Image Orientation (Patient)	(0020,0037)	2C	If the setting conditions are met, Length=0 when no entry is made
> Image Position (Patient)	(0020,0032)	2C	If the setting conditions are met, Length=0 when no entry is made
> View Code Sequence	(0054,0220)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
>> View Angulation Modifier Code Sequence	(0054,0222)	2C	If the setting conditions are met, Length=0 when no entry is made
>>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>>> Code Meaning	(0008,0104)	3	Not set when no entry is made

8.2.33 NM TOMO Acquisition Module

Table 81

Attribute Name	Tag	Type	Attribute Description
Rotation Information Sequence	(0054,0052)	2	Length=0 when no entry is made
> Start Angle	(0054,0200)	1C	Always set when the sequence is set
> Angular Step	(0018,1144)	1C	Always set when the sequence is set
> Rotation Direction	(0018,1140)	1C	Always set when the sequence is set
> Scan Arc	(0018,1143)	1C	Always set when the sequence is set
> Actual Frame Duration	(0018,1242)	1C	Always set when the sequence is set
> Radial Position	(0018,1142)	3	Not set when no entry is made
> Distance Source to Detector	(0018,1110)	2C	If the setting conditions are met, Length=0 when no entry is made
> Number of Frames in Rotation	(0054,0053)	1C	Always set when the sequence is set
> Table Traverse	(0018,1131)	3	Not set when no entry is made
> Table Height	(0018,1130)	3	Not set when no entry is made
Type of Detector Motion	(0054,0202)	3	Not set when no entry is made

8.2.34 NM Multi-Gated Acquisition Module

Table 82

Attribute Name	Tag	Type	Attribute Description
Beat Rejection Flag	(0018,1080)	3	Not set when no entry is made
PVC Rejection	(0018,1085)	3	Not set when no entry is made
Skip Beats	(0018,1086)	3	Not set when no entry is made
Heart Rate	(0018,1088)	3	Not set when no entry is made
Gated Information Sequence	(0054,0062)	2C	If the setting conditions are met, Length=0 when no entry is made
> Trigger Time	(0018,1060)	3	Not set when no entry is made
> Framing Type	(0018,1064)	3	Not set when no entry is made
> Data Information Sequence	(0054,0063)	2C	If the setting conditions are met, Length=0 when no entry is made
>> Frame Time	(0018,1063)	1C	Always set when the sequence is set
>> Nominal Interval	(0018,1062)	3	Not set when no entry is made
>> Low R-R Value	(0018,1081)	3	Not set when no entry is made
>> High R-R Value	(0018,1082)	3	Not set when no entry is made
>> Intervals Acquired	(0018,1083)	3	Not set when no entry is made
>> Intervals Rejected	(0018,1084)	3	Not set when no entry is made
>> Time Slot Information Sequence	(0054,0072)	2C	If the setting conditions are met, Length=0 when no entry is made
>>> Time Slot Time	(0054,0073)	3	Not set when no entry is made

8.2.35 NM Phase Module

Table 83

Attribute Name	Tag	Type	Attribute Description
Phase Information Sequence	(0054,0032)	2C	If the setting conditions are met, Length=0 when no entry is made
> Phase Delay	(0054,0036)	1C	Always set when the sequence is set
> Actual Frame Duration	(0018,1242)	1C	Always set when the sequence is set
> Pause between Frames	(0054,0038)	1C	Always set when the sequence is set
> Number of Frames in Phase	(0054,0033)	1C	Always set when the sequence is set
> Trigger Vector	(0054,0210)	3	Not set when no entry is made
> Number of Triggers in Phase	(0054,0211)	1C	Always set when the sequence is set

8.2.36 NM Reconstruction Module**Table 84**

Attribute Name	Tag	Type	Attribute Description
Spacing Between Slices	(0018,0088)	2	Length=0 when no entry is made
Reconstruction Diameter	(0018,1100)	3	Not set when no entry is made
Convolution Kernel	(0018,1210)	3	Not set when no entry is made
Slice Thickness	(0018,0050)	2	Length=0 when no entry is made
Slice Location	(0020,1041)	3	Not set when no entry is made

8.2.37 SC Image Module**Table 85**

Attribute Name	Tag	Type	Attribute Description
Date of Secondary Capture	(0018,1012)	3	Not set when no entry is made
Time of Secondary Capture	(0018,1014)	3	Not set when no entry is made

8.2.38 X-Ray Image Module

Table 86

Attribute Name	Tag	Type	Attribute Description
Frame Increment Pointer	(0028,0009)	1C	Always set when the setting conditions are met
Lossy Image Compression	(0008,2110)	1C	Always set when the setting conditions are met
Image Type	(0008,0008)	1	Always set
Pixel Intensity Relationship	(0028,1040)	1	Always set
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Pixel Representation	(0028,0103)	1	Always set
Scan Options	(0018,0022)	3	Not set when no entry is made
Anatomic Region Sequence	(0008,2218)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Anatomic Region Modifier Sequence	(0008,2220)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
Primary Anatomic Structure Sequence	(0008,2228)	3	Not set when no entry is made
> Code Value	(0008,0100)	1C	Always set when the sequence is set
> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
> Code Meaning	(0008,0104)	3	Not set when no entry is made
> Primary Anatomic Structure Modifier Sequence	(0008,2230)	3	Not set when no entry is made
>> Code Value	(0008,0100)	1C	Always set when the sequence is set
>> Coding Scheme Designator	(0008,0102)	1C	Always set when the sequence is set
>> Code Meaning	(0008,0104)	3	Not set when no entry is made
R Wave Pointer	(0028,6040)	3	Not set when no entry is made
Reference Image Sequence	(0008,1140)	1C	Always set when the setting conditions are met
> Reference SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Reference SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Derivation Description	(0008,2111)	3	Not set when no entry is made

Attribute Name	Tag	Type	Attribute Description
Acquisition Device Processing Description	(0018,1400)	3	Not set when no entry is made
Calibration Image	(0050,0004)	3	Not set when no entry is made

8.2.39 X-Ray Acquisition Module

Table 87

Attribute Name	Tag	Type	Attribute Description
KVP	(0018,0060)	2	Length=0 when no entry is made
Radiation Setting	(0018,1155)	1	Always set
X-ray Tube Current	(0018,1151)	2C	If the setting conditions are met, Length=0 when no entry is made
Exposure Time	(0018,1150)	2C	If the setting conditions are met, Length=0 when no entry is made
Exposure	(0018,1152)	2C	If the setting conditions are met, Length=0 when no entry is made
Grid	(0018,1166)	3	Not set when no entry is made
Average Pulse Width	(0018,1154)	3	Not set when no entry is made
Radiation Mode	(0018,115A)	3	Not set when no entry is made
Type of Filters	(0018,1161)	3	Not set when no entry is made
Intensifier Size	(0018,1162)	3	Not set when no entry is made
Field of View Shape	(0018,1147)	3	Not set when no entry is made
Field of View Dimension(s)	(0018,1149)	3	Not set when no entry is made
Imager Pixel Spacing	(0018,1164)	3	Not set when no entry is made
Focal Spot	(0018,1190)	3	Not set when no entry is made
Image Area Dose Product	(0018,115E)	3	Not set when no entry is made

8.2.40 X-Ray Collimator Module

Table 88

Attribute Name	Tag	Type	Attribute Description
Collimator Shape	(0018,1700)	1	Always set
Collimator Left Vertical Edge	(0018,1702)	1C	Always set when the setting conditions are met
Collimator Right Vertical Edge	(0018,1704)	1C	Always set when the setting conditions are met
Collimator Upper Horizontal Edge	(0018,1706)	1C	Always set when the setting conditions are met
Collimator Lower Horizontal Edge	(0018,1708)	1C	Always set when the setting conditions are met
Center of Circular Collimator	(0018,1710)	1C	Always set when the setting conditions are met
Radius of Circular Collimator	(0018,1712)	1C	Always set when the setting conditions are met
Vertices of the Polygonal Collimator	(0018,1720)	1C	Always set when the setting conditions are met

8.2.41 X-Ray Table Module

Table 89

Attribute Name	Tag	Type	Attribute Description
Table Motion	(0018,1134)	2	Length=0 when no entry is made
Table Vertical Increment	(0018,1135)	2C	If the setting conditions are met, Length=0 when no entry is made
Table Longitudinal Increment	(0018,1137)	2C	If the setting conditions are met, Length=0 when no entry is made
Table Lateral Increment	(0018,1136)	2C	If the setting conditions are met, Length=0 when no entry is made
Table Angle	(0018,1138)	3	Not set when no entry is made

8.2.42 XRF Positioner Module

Table 90

Attribute Name	Tag	Type	Attribute Description
Distance Source to Detector	(0018,1110)	3	Not set when no entry is made
Distance Source to Patient	(0018,1111)	3	Not set when no entry is made
Estimated Radiographic Magnification Factor	(0018,1114)	3	Not set when no entry is made
Column Angulation	(0018,1450)	3	Not set when no entry is made

8.2.43 XRF Tomography Acquisition Module

Table 91

Attribute Name	Tag	Type	Attribute Description
Tomo Layer Height	(0018,1460)	1	Always set
Tomo Angle	(0018,1470)	3	Not set when no entry is made
Tomo Time	(0018,1480)	3	Not set when no entry is made

8.2.44 XA Positioner Module

Table 92

Attribute Name	Tag	Type	Attribute Description
Distance Source to Patient	(0018,1111)	3	Not set when no entry is made
Distance Source to Detector	(0018,1110)	3	Not set when no entry is made
Estimated Radiographic Magnification Factor	(0018,1114)	3	Not set when no entry is made
Positioner Motion	(0018,1500)	2C	If the setting conditions are met, Length=0 when no entry is made
Positioner Primary Angle	(0018,1510)	2	Length=0 when no entry is made
Positioner Secondary Angle	(0018,1511)	2	Length=0 when no entry is made
Positioner Primary Angle Increment	(0018,1520)	2C	If the setting conditions are met, Length=0 when no entry is made
Positioner Secondary Angle Increment	(0018,1521)	2C	If the setting conditions are met, Length=0 when no entry is made
Detector Primary Angle	(0018,1530)	3	Not set when no entry is made
Detector Secondary Angle	(0018,1531)	3	Not set when no entry is made

8.2.45 Overlay Plane Module

Table 93

Attribute Name	Tag	Type	Attribute Description
Overlay Rows	(60xx,0010)	1	Always set
Overlay Columns	(60xx,0011)	1	Always set
Overlay Type	(60xx,0040)	1	Always set
Overlay Origin	(60xx,0050)	1	Always set
Overlay Bits Allocated	(60xx,0100)	1	Always set
Overlay Bit Position	(60xx,0102)	1	Always set
Overlay Data	(60xx,3000)	1C	Always set when the setting conditions are met
Overlay Description	(60xx,0022)	3	Not set when no entry is made
Overlay Subtype	(60xx,0045)	3	Not set when no entry is made
Overlay Label	(60xx,1500)	3	Not set when no entry is made
ROI Area	(60xx,1301)	3	Not set when no entry is made
ROI Mean	(60xx,1302)	3	Not set when no entry is made
ROI Standard Deviation	(60xx,1303)	3	Not set when no entry is made
Overlay Descriptor - Gray (retired)	(60xx,1100)	3	Set when the received image contains this entry
Overlay Descriptor - Red (retired)	(60xx,1101)	3	Set when the received image contains this entry
Overlay Descriptor - Green (retired)	(60xx,1102)	3	Set when the received image contains this entry
Overlay Descriptor - Blue (retired)	(60xx,1103)	3	Set when the received image contains this entry
Overlays - Gray (retired)	(60xx,1200)	3	Set when the received image contains this entry
Overlays - Red (retired)	(60xx,1201)	3	Set when the received image contains this entry
Overlays - Green (retired)	(60xx,1202)	3	Set when the received image contains this entry
Overlays - Blue (retired)	(60xx,1203)	3	Set when the received image contains this entry

8.2.46 Multi-frame Overlay Module

Table 94

Attribute Name	Tag	Type	Attribute Description
Number of Frames in Overlay	(60xx,0015)	1	Always set
Image Frame Origin	(60xx,0051)	3	Not set when no entry is made

8.2.47 Curve Module

Table 95

Attribute Name	Tag	Type	Attribute Description
Curve Dimensions	(50xx,0005)	1	Always set
Number of Points	(50xx,0010)	1	Always set
Type of Data	(50xx,0020)	1	Always set
Data Value Representation	(50xx,0103)	1	Always set
Curve Data	(50xx,3000)	1	Always set
Curve Description	(50xx,0022)	3	Not set when no entry is made
Axis Units	(50xx,0030)	3	Not set when no entry is made
Axis Labels	(50xx,0040)	3	Not set when no entry is made
Minimum Coordinate Value	(50xx,0104)	3	Not set when no entry is made
Maximum Coordinate Value	(50xx,0105)	3	Not set when no entry is made
Curve Range	(50xx,0106)	3	Not set when no entry is made
Curve Data Descriptor	(50xx,0110)	1C	Always set when the setting conditions are met
Coordinate Start Value	(50xx,0112)	1C	Always set when the setting conditions are met
Coordinate Step Value	(50xx,0114)	1C	Always set when the setting conditions are met
Curve Label	(50xx,2500)	3	Not set when no entry is made
Referenced Overlay Sequence	(50xx,2600)	3	Not set when no entry is made
> Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
> Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
> Referenced Overlay Group	(50xx,2610)	1C	Always set when the sequence is set

8.2.48 Modality LUT Module

Table 96

Attribute Name	Tag	Type	Attribute Description
Modality LUT Sequence	(0028,3000)	3	Not set when no entry is made
> LUT Descriptor	(0028,3002)	1C	Always set when the sequence is set
> LUT Explanation	(0028,3003)	3	Not set when no entry is made
> Modality LUT Type	(0028,3004)	1C	Always set when the sequence is set
> LUT Data	(0028,3006)	1C	Always set when the sequence is set
Rescale Intercept	(0028,1052)	1C	Always set when the setting conditions are met
Rescale Slope	(0028,1053)	1C	Always set when the setting conditions are met
Rescale Type	(0028,1054)	1C	Always set when the setting conditions are met

8.2.49 VOI LUT Module

Table 97

Attribute Name	Tag	Type	Attribute Description
VOI LUT Sequence	(0028,3010)	3	Not set when no entry is made
>LUT Descriptor	(0028,3002)	1C	Always set when the sequence is set
>LUT Explanation	(0028,3003)	3	Not set when no entry is made
>LUT Data	(0028,3006)	1C	Always set when the sequence is set
Window Center	(0028,1050)	3	Not set when no entry is made
Window Width	(0028,1051)	1C	Always set when the setting conditions are met
Window Center & Width Explanation	(0028,1055)	3	Not set when no entry is made

8.2.50 SOP Common Module**Table 98**

Attribute Name	Tag	Type	Attribute Description
SOP Class UID	(0008,0016)	1	Always set
SOP Instance UID	(0008,0018)	1	Always set
Specific Character Set	(0008,0005)	1C	Always set when the setting conditions are met
Instance Creation Date	(0008,0012)	3	Not set when no entry is made
Instance Creation Time	(0008,0013)	3	Not set when no entry is made
Instance Creator UID	(0008,0014)	3	Not set when no entry is made

8.2.51 Curve Identification Module

Table 99

Attribute Name	Tag	Type	Attribute Description
Curve Number	(0020,0024)	2	Length=0 when no entry is made
Curve Date	(0008,0025)	3	Not set when no entry is made
Curve Time	(0008,0035)	3	Not set when no entry is made
Referenced Image Sequence	(0008,1140)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the setting conditions are met
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the setting conditions are met
Referenced Overlay Sequence	(0008,1130)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the setting conditions are met
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the setting conditions are met
Referenced Curve Sequence	(0008,1145)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the setting conditions are met
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the setting conditions are met

8.2.52 Audio Module

Table 100

Attribute Name	Tag	Type	Attribute Description
Audio Type	(50xx,2000)	1	Always set
Audio Sample Format	(50xx,2002)	1	Always set
Number of Channels	(50xx,2004)	1	Always set
Number of Samples	(50xx,2006)	1	Always set
Sample Rate	(50xx,2008)	1	Always set
Total Time	(50xx,200A)	1	Always set
Audio Sample Data	(50xx,200C)	1	Always set
Referenced Image Sequence	(0008,1140)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the setting conditions are met
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the setting conditions are met
Audio Comments	(50xx,200E)	3	Not set when no entry is made

8.2.53 Common Attribute Set for Code Sequence (Invoked as “Code Sequence Macro”)

Table 101

Attribute Name	Tag	Type	Attribute Description
Code Value	(0008,0100)	1C	Always set when the setting conditions are met
Coding Scheme Designator	(0008,0102)	1C	Always set when the setting conditions are met
Coding Scheme Version	(0008,0103)	1C	Always set when the setting conditions are met
Code Meaning	(0008,0104)	1C	Always set when the setting conditions are met
Context Identifier	(0008,010F)	3	Not set when no entry is made
Mapping Resource	(0008,0105)	1C	Always set when the setting conditions are met
Context Group Version	(0008,0106)	1C	Always set when the setting conditions are met
Code Set Extension Flag	(0008,010B)	3	Not set when no entry is made
Context Group Local Version	(0008,0107)	1C	Always set when the setting conditions are met
Private Coding Scheme Creator UID	(0008,010C)	3	Not set when no entry is made
Code Set Extension Creator UID	(0008,010D)	1C	Always set when the setting conditions are met

8.2.54 Specimen Identification Module

Table 102

Attribute Name	Tag	Type	Attribute Description
Specimen Accession Number	(0040,050A)	1	Always set
Specimen Sequence	(0040,0550)	2	Not set when no entry is made
>Specimen Identifier	(0040,0551)	2C	Always set when the setting conditions are met
>Specimen Type Code Sequence	(0040,059A)	2C	Always set when the setting conditions are met
>>Include 'Code Sequence Macro' Table 101			
>Slide Identifier	(0040,06FA)	2C	Always set when the setting conditions are met

8.2.55 DX Series Module

Table 103

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	Always set
Referenced Study Component Sequence	(0008,1111)	1C	Always set when the setting conditions are met
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Presentation Intent Type	(0008,0068)	1	Always set

8.2.56 DX Anatomy Imaged Module

Table 104

Attribute Name	Tag	Type	Attribute Description
Image Laterality	(0020,0062)	1	Always set
Anatomic Region Sequence	(0008,2218)	2	Length=0 when no entry is made
>Include 'Code Sequence Macro' Table 101			
>Anatomic Region Modifier Sequence	(0008,2220)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
Primary Anatomic Structure Sequence	(0008,2228)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
>Primary Anatomic Structure Modifier Sequence	(0008,2230)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			

8.2.57 DX Image Module

Table 105

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	Always set
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Pixel Representation	(0028,0103)	1	Always set
Pixel Intensity Relationship	(0028,1040)	1	Always set
Pixel Intensity Relationship Sign	(0028,1041)	1	Always set
Rescale Intercept	(0028,1052)	1	Always set
Rescale Slope	(0028,1053)	1	Always set
Rescale Type	(0028,1054)	1	Always set
Presentation LUT Shape	(2050,0020)	1	Always set
Lossy Image Compression	(0028,2110)	1	Always set
Lossy Image Compression Ratio	(0028,2112)	1C	Always set when the setting conditions are met
Derivation Description	(0008,2111)	3	Not set when no entry is made
Acquisition Device Processing Description	(0018,1400)	3	Not set when no entry is made
Acquisition Device Processing Code	(0018,1401)	3	Not set when no entry is made
Patient Orientation	(0020,0020)	1	Always set
Calibration Image	(0050,0004)	3	Not set when no entry is made
Burned In Annotation	(0028,0301)	1	Always set
VOI LUT Sequence	(0028,3010)	1C	Always set when the setting conditions are met
>LUT Descriptor	(0028,3002)	1C	Always set when the sequence is set
>LUT Explanation	(0028,3003)	3	Not set when no entry is made
>LUT Data	(0028,3006)	1C	Always set when the sequence is set
Window Center	(0028,1050)	1C	Always set when the setting conditions are met
Window Width	(0028,1051)	1C	Always set when the setting conditions are met
Window Center & Width Explanation	(0028,1055)	3	Not set when no entry is made

8.2.58 DX Detector Module

Table 106

Attribute Name	Tag	Type	Attribute Description
Detector Type	(0018,7004)	2	Length=0 when no entry is made
Detector Configuration	(0018,7005)	3	Not set when no entry is made
Detector Description	(0018,7006)	3	Not set when no entry is made
Detector Mode	(0018,7008)	3	Not set when no entry is made
Detector ID	(0018,700A)	3	Not set when no entry is made
Date of Last Detector Calibration	(0018,700C)	3	Not set when no entry is made
Time of Last Detector Calibration	(0018,700E)	3	Not set when no entry is made
Exposures on Detector Since Last Calibration	(0018,7010)	3	Not set when no entry is made
Exposures on Detector Since Manufactured	(0018,7011)	3	Not set when no entry is made
Detector Time Since Last Exposure	(0018,7012)	3	Not set when no entry is made
Detector Active Time	(0018,7014)	3	Not set when no entry is made
Detector Activation Offset From Exposure	(0018,7016)	3	Not set when no entry is made
Detector Binning	(0018,701A)	3	Not set when no entry is made
Detector Conditions Nominal Flag	(0018,7000)	3	Not set when no entry is made
Detector Temperature	(0018,7001)	3	Not set when no entry is made
Sensitivity	(0018,6000)	3	Not set when no entry is made
Field of View Shape	(0018,1147)	3	Not set when no entry is made
Field of View Dimension(s)	(0018,1149)	3	Not set when no entry is made
Field of View Origin	(0018,7030)	1C	Always set when the setting conditions are met
Field of View Rotation	(0018,7032)	1C	Always set when the setting conditions are met
Field of View Horizontal Flip	(0018,7034)	1C	Always set when the setting conditions are met
Imager Pixel Spacing	(0018,1164)	1	Always set
Detector Element Physical Size	(0018,7020)	3	Not set when no entry is made
Detector Element Spacing	(0018,7022)	3	Not set when no entry is made
Detector Active Shape	(0018,7024)	3	Not set when no entry is made
Detector Active Dimension(s)	(0018,7026)	3	Not set when no entry is made
Detector Active Origin	(0018,7028)	3	Not set when no entry is made

8.2.59 DX Positioning Module

Table 107

Attribute Name	Tag	Type	Attribute Description
Projection Eponymous Name Code Sequence	(0018,5104)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
Patient Position	(0018,5100)	3	Not set when no entry is made
View Position	(0018,5101)	3	Not set when no entry is made
View Code Sequence	(0054,0220)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
>View Modifier Code Sequence	(0054,0222)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
Patient Orientation Code Sequence	(0054,0410)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
> Patient Orientation Modifier Code Sequence	(0054,0412)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
Patient Gantry Relationship Code Sequence	(0054,0414)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
Distance Source to Patient	(0018,1111)	3	Not set when no entry is made
Distance Source to Detector	(0018,1110)	3	Not set when no entry is made
Estimated Radiographic Magnification Factor	(0018,1114)	3	Not set when no entry is made
Positioner Type	(0018,1508)	2	Length=0 when no entry is made
Positioner Primary Angle	(0018,1510)	3	Not set when no entry is made
Positioner Secondary Angle	(0018,1511)	3	Not set when no entry is made
Detector Primary Angle	(0018,1530)	3	Not set when no entry is made
Detector Secondary Angle	(0018,1531)	3	Not set when no entry is made
Column Angulation	(0018,1450)	3	Not set when no entry is made
Table Type	(0018,113A)	3	Not set when no entry is made
Table Angle	(0018,1138)	3	Not set when no entry is made
Body Part Thickness	(0018,11A0)	3	Not set when no entry is made
Compression Force	(0018,11A2)	3	Not set when no entry is made

8.2.60 X-Ray Acquisition Dose Module

Table 108

Attribute Name	Tag	Type	Attribute Description
KVP	(0018,0060)	3	Not set when no entry is made
X-Ray Tube Current	(0018,1151)	3	Not set when no entry is made
X-Ray Tube Current in μ A	(0018,8151)	3	Not set when no entry is made
Exposure Time	(0018,1150)	3	Not set when no entry is made
Exposure Time in μ S	(0018,8150)	3	Not set when no entry is made
Exposure	(0018,1152)	3	Not set when no entry is made
Exposure in μ As	(0018,1153)	3	Not set when no entry is made
Distance Source to Detector	(0018,1110)	3	Not set when no entry is made
Distance Source to Patient	(0018,1111)	3	Not set when no entry is made
Image Area Dose Product	(0018,115E)	3	Not set when no entry is made
Body Part Thickness	(0018,11A0)	3	Not set when no entry is made
Relative X-Ray Exposure	(0018,1405)	3	Not set when no entry is made
Entrance Dose	(0040,0302)	3	Not set when no entry is made
Entrance Dose in mGy	(0040,8302)	3	Not set when no entry is made
Exposed Area	(0040,0303)	3	Not set when no entry is made
Distance Source to Entrance	(0040,0306)	3	Not set when no entry is made
Comments on Radiation Dose	(0040,0310)	3	Not set when no entry is made
X-Ray Output	(0040,0312)	3	Not set when no entry is made
Half Value Layer	(0040,0314)	3	Not set when no entry is made
Organ Dose	(0040,0316)	3	Not set when no entry is made
Organ Exposed	(0040,0318)	3	Not set when no entry is made
Anode Target Material	(0018,1191)	3	Not set when no entry is made
Filter Material	(0018,7050)	3	Not set when no entry is made
Filter Thickness Minimum	(0018,7052)	3	Not set when no entry is made
Filter Thickness Maximum	(0018,7054)	3	Not set when no entry is made
Rectification Type	(0018,1156)	3	Not set when no entry is made

8.2.61 X-Ray Generation Module

Table 109

Attribute Name	Tag	Type	Attribute Description
KVP	(0018,0060)	3	Not set when no entry is made
X-Ray Tube Current	(0018,1151)	3	Not set when no entry is made
X-Ray Tube Current in μA	(0018,8151)	3	Not set when no entry is made
Exposure Time	(0018,1150)	3	Not set when no entry is made
Exposure Time in μS	(0018,8150)	3	Not set when no entry is made
Exposure	(0018,1152)	3	Not set when no entry is made
Exposure in μAs	(0018,1153)	3	Not set when no entry is made
Exposure Control Mode	(0018,7060)	3	Not set when no entry is made
Exposure Control Mode Description	(0018,7062)	3	Not set when no entry is made
Exposure Status	(0018,7064)	3	Not set when no entry is made
Phototimer Setting	(0018,7065)	3	Not set when no entry is made
Focal Spot	(0018,1190)	3	Not set when no entry is made
Anode Target Material	(0018,1191)	3	Not set when no entry is made
Rectification Type	(0018,1156)	3	Not set when no entry is made

8.2.62 X-Ray Filtration Module

Table 110

Attribute Name	Tag	Type	Attribute Description
Filter Type	(0018,1160)	3	Not set when no entry is made
Filter Material	(0018,7050)	3	Not set when no entry is made
Filter Thickness Minimum	(0018,7052)	3	Not set when no entry is made
Filter Thickness Maximum	(0018,7054)	3	Not set when no entry is made

8.2.63 X-Ray Grid Module

Table 111

Attribute Name	Tag	Type	Attribute Description
Grid	(0018,1166)	3	Not set when no entry is made
Grid Absorbing Material	(0018,7040)	3	Not set when no entry is made
Grid Spacing Material	(0018,7041)	3	Not set when no entry is made
Grid Thickness	(0018,7042)	3	Not set when no entry is made
Grid Pitch	(0018,7044)	3	Not set when no entry is made
Grid Aspect Ratio	(0018,7046)	3	Not set when no entry is made
Grid Period	(0018,7048)	3	Not set when no entry is made
Grid Focal Distance	(0018,704C)	3	Not set when no entry is made

8.2.64 X-Ray Tomography Acquisition Module

Table 112

Attribute Name	Tag	Type	Attribute Description
Tomography Layer Height	(0018,1460)	1	Always set
Tomography Angle	(0018,1470)	3	Not set when no entry is made
Tomography Time	(0018,1480)	3	Not set when no entry is made
Tomography Type	(0018,1490)	3	Not set when no entry is made
Tomography Class	(0018,1491)	3	Not set when no entry is made
Number of Tomosynthesis Source Images	(0018,1495)	3	Not set when no entry is made

8.2.65 Image Histogram Module

Table 113

Attribute Name	Tag	Type	Attribute Description
Image Histogram Module	(0060,3000)	1	Always set
>Histogram Number of Bins	(0060,3002)	1C	Always set when the sequence is set
>Histogram First Bin Value	(0060,3004)	1C	Always set when the sequence is set
>Histogram Last Bin Value	(0060,3006)	1C	Always set when the sequence is set
>Histogram Bin Width	(0060,3008)	1C	Always set when the sequence is set
>Histogram Explanation	(0060,3010)	3	Not set when no entry is made
>Histogram Data	(0060,3020)	1C	Always set when the sequence is set

8.2.66 Acquisition Context Module**Table 114**

Attribute Name	Tag	Type	Attribute Description
Acquisition Context Sequence	(0040,0555)	2	Length=0 when no entry is made
>Concept-name Code Sequence	(0040,A043)	1C	Always set when the sequence is set
>>Include 'Code Sequence Macro' Table 101			
>Referenced Frame Numbers	(0040,A136)	1C	Always set when the sequence is set
>Numeric Value	(0040,A30A)	1C	Always set when the sequence is set
>Measurement Units Code Sequence	(0040,08EA)	1C	Always set when the sequence is set
>>Include 'Code Sequence Macro' Table 101			
>Date	(0040,A121)	1C	Always set when the sequence is set
>Time	(0040,A122)	1C	Always set when the sequence is set
>Person Name	(0040,A123)	1C	Always set when the sequence is set
>Text Value	(0040,A160)	1C	Always set when the sequence is set
>Concept Code Sequence	(0040,A168)	1C	Always set when the sequence is set
>>Include 'Code Sequence Macro' Table 101			
Acquisition Context Description	(0040,0556)	3	Not set when no entry is made

8.2.67 Mammography Series Module**Table 115**

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	1	Always set

8.2.68 Mammography Image Module**Table 116**

Attribute Name	Tag	Type	Attribute Description
Positioner Type	(0018,1508)	1	Always set
Positioner Primary Angle	(0018,1510)	3	Not set when no entry is made
Positioner Secondary Angle	(0018,1511)	3	Not set when no entry is made
Image Laterality	(0020,0062)	1	Always set
Organ Exposed	(0040,0318)	1	Always set
Implant Present	(0028,1300)	3	Not set when no entry is made
Partial View	(0028,1350)	3	Not set when no entry is made
Partial View Description	(0028,1351)	3	Not set when no entry is made
Anatomic Region Sequence	(0008,2218)	1	Always set
>Include 'Code Sequence Macro' Table 101			
View Code Sequence	(0054,0220)	1	Always set
>Include 'Code Sequence Macro' Table 101			
>View Modifier Code Sequence	(0054,0222)	2	Length=0 when no entry is made
>>Include 'Code Sequence Macro' Table 101			

8.2.69 PET Series Module

Table 117

Attribute Name	Tag	Type	Attribute Description
Series Date	(0008,0021)	1	Always set
Series Time	(0008,0031)	1	Always set
Units	(0054,1001)	1	Always set
Counts Source	(0054,1002)	1	Always set
Series Type	(0054,1000)	1	Always set
Reprojection Method	(0054,1004)	2C	If the setting conditions are met, Length=0 when no entry is made
Number of R-R Intervals	(0054,0061)	1C	Always set when the setting conditions are met
Number of Time Slots	(0054,0071)	1C	Always set when the setting conditions are met
Number of Time Slices	(0054,0101)	1C	Always set when the setting conditions are met
Number of Slices	(0054,0081)	1	Always set
Corrected Image	(0028,0051)	2	Length=0 when no entry is made
Randoms Correction Method	(0054,1100)	3	Not set when no entry is made
Attenuation Correction Method	(0054,1101)	3	Not set when no entry is made
Scatter Correction Method	(0054,1105)	3	Not set when no entry is made
Decay Correction	(0054,1102)	1	Always set
Reconstruction Diameter	(0018,1100)	3	Not set when no entry is made
Convolution Kernel	(0018,1210)	3	Not set when no entry is made
Reconstruction Method	(0054,1103)	3	Not set when no entry is made
Detector Lines of Response Used	(0054,1104)	3	Not set when no entry is made
Acquisition Start Condition	(0018,0073)	3	Not set when no entry is made
Acquisition Start Condition Data	(0018,0074)	3	Not set when no entry is made
Acquisition Termination Condition	(0018,0071)	3	Not set when no entry is made
Acquisition Termination Condition Data	(0018,0075)	3	Not set when no entry is made
Field of View Shape	(0018,1147)	3	Not set when no entry is made
Field of View Dimensions	(0018,1149)	3	Not set when no entry is made
Gantry/Detector Tilt	(0018,1120)	3	Not set when no entry is made
Gantry/Detector Slew	(0018,1121)	3	Not set when no entry is made
Type of Detector Motion	(0054,0202)	3	Not set when no entry is made
Collimator Type	(0018,1181)	2	Length=0 when no entry is made

Collimator/Grid Name	(0018,1180)	3	Not set when no entry is made
Axial Acceptance	(0054,1200)	3	Not set when no entry is made
Axial Mash	(0054,1201)	3	Not set when no entry is made
Transverse Mash	(0054,1202)	3	Not set when no entry is made
Detector Element Size	(0054,1203)	3	Not set when no entry is made
Coincidence Window Width	(0054,1210)	3	Not set when no entry is made
Energy Window Range Sequence	(0054,0013)	3	Not set when no entry is made
>Energy Window Lower Limit	(0054,0014)	3	Not set when no entry is made
>Energy Window Upper Limit	(0054,0015)	3	Not set when no entry is made
Secondary Counts Type	(0054,1220)	3	Not set when no entry is made

8.2.70 PET Isotope Module**Table 118**

Attribute Name	Tag	Type	Attribute Description
Radiopharmaceutical Information Sequence	(0054,0016)	2	Length=0 when no entry is made
>Radionuclide Code Sequence	(0054,0300)	2	Length=0 when no entry is made
>>Include 'Code Sequence Macro' Table 101			
>Radiopharmaceutical Route	(0018,1070)	3	Not set when no entry is made
>Administration Route Code Sequence	(0054,0302)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
>Radiopharmaceutical Volume	(0018,1071)	3	Not set when no entry is made
>Radiopharmaceutical Start Time	(0018,1072)	3	Not set when no entry is made
>Radiopharmaceutical Stop Time	(0018,1073)	3	Not set when no entry is made
>Radionuclide Total Dose	(0018,1074)	3	Not set when no entry is made
>Radionuclide Half Life	(0018,1075)	3	Not set when no entry is made
>Radionuclide Positron Fraction	(0018,1076)	3	Not set when no entry is made
>Radiopharmaceutical Specific Activity	(0018,1077)	3	Not set when no entry is made
>Radiopharmaceutical	(0018,0031)	3	Not set when no entry is made
>Radiopharmaceutical Code Sequence	(0054,0304)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
Intervention Drug Information Sequence	(0018,0026)	3	Not set when no entry is made
>Intervention Drug Name	(0018,0034)	3	Not set when no entry is made
>Intervention Drug Code Sequence	(0018,0029)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
>Intervention Drug Start Time	(0018,0035)	3	Not set when no entry is made
>Intervention Drug Stop Time	(0018,0027)	3	Not set when no entry is made
>Intervention Drug Dose	(0018,0028)	3	Not set when no entry is made

8.2.71 PET Multi-gated Acquisition Module**Table 119**

Attribute Name	Tag	Type	Attribute Description
Beat Rejection Flag	(0018,1080)	2	Length=0 when no entry is made
Trigger Source or Type	(0018,1061)	3	Not set when no entry is made
PVC Rejection	(0018,1085)	3	Not set when no entry is made
Skip Beats	(0018,1086)	3	Not set when no entry is made
Heart Rate	(0018,1088)	3	Not set when no entry is made
Framing Type	(0018,1064)	3	Not set when no entry is made

8.2.72 NM/PET Patient Orientation Module**Table 120**

Attribute Name	Tag	Type	Attribute Description
Patient Orientation Code Sequence	(0054,0410)	2	Length=0 when no entry is made
>Include 'Code Sequence Macro' Table 101			
> Patient Orientation Modifier Code Sequence	(0054,0412)	2C	If the setting conditions are met, Length=0 when no entry is made
>>Include 'Code Sequence Macro' Table 101			
Patient Gantry Relationship Code Sequence	(0054,0414)	2	Length=0 when no entry is made
>Include 'Code Sequence Macro' Table 101			

8.2.73 PET Image Module

Table 121

Attribute Name	Tag	Type	Attribute Description
Image Type	(0008,0008)	1	Always set
Samples per Pixel	(0028,0002)	1	Always set
Photometric Interpretation	(0028,0004)	1	Always set
Bits Allocated	(0028,0100)	1	Always set
Bits Stored	(0028,0101)	1	Always set
High Bit	(0028,0102)	1	Always set
Rescale Intercept	(0028,1052)	1	Always set
Rescale Slope	(0028,1053)	1	Always set
Frame Reference Time	(0054,1300)	1	Always set
Trigger Time	(0018,1060)	1C	Always set when the setting conditions are met
Frame Time	(0018,1063)	1C	Always set when the setting conditions are met
Low R-R Value	(0018,1081)	1C	Always set when the setting conditions are met
High R-R Value	(0018,1082)	1C	Always set when the setting conditions are met
Lossy Image Compression	(0028,2110)	1C	Always set when the setting conditions are met
Image Index	(0054,1330)	1	Always set
Acquisition Date	(0008,0022)	2	Length=0 when no entry is made
Acquisition Time	(0008,0032)	2	Length=0 when no entry is made
Actual Frame Duration	(0018,1242)	2	Length=0 when no entry is made
Nominal Interval	(0018,1062)	3	Not set when no entry is made
Intervals Acquired	(0018,1083)	3	Not set when no entry is made
Intervals Rejected	(0018,1084)	3	Not set when no entry is made
Primary (Prompts) Counts Accumulated	(0054,1310)	3	Not set when no entry is made
Secondary Counts Accumulated	(0054,1311)	3	Not set when no entry is made
Slice Sensitivity Factor	(0054,1320)	3	Not set when no entry is made
Decay Factor	(0054,1321)	1C	Always set when the setting conditions are met
Dose Calibration Factor	(0054,1322)	3	Not set when no entry is made
Scatter Fraction Factor	(0054,1323)	3	Not set when no entry is made
Dead Time Factor	(0054,1324)	3	Not set when no entry is made
Referenced Overlay Sequence	(0008,1130)	3	Not set when no entry is made

>Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Referenced Curve Sequence	(0008,1145)	3	Not set when no entry is made
>Referenced SOP Class UID	(0008,1150)	1C	Always set when the sequence is set
>Referenced SOP Instance UID	(0008,1155)	1C	Always set when the sequence is set
Anatomic Region Sequence	(0008,2218)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
>Anatomic Region Modifier Sequence	(0008,2220)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			
Primary Anatomic Structure Sequence	(0008,2228)	3	Not set when no entry is made
>Include 'Code Sequence Macro' Table 101			
>Primary Anatomic Structure Modifier Sequence	(0008,2230)	3	Not set when no entry is made
>>Include 'Code Sequence Macro' Table 101			

8.3 Private Data Elements

Private data elements for the TWS are defined below.

Table 122

Attribute Name	Tag	Type	VR	VM	Attribute Description
Private Creator	(7019,00XX)	1C* ¹	LO	1	"TOSHIBA_MEC_OT3"
Annotation Version No.	(7019,XX01)	3	DS	1	Not set when no entry is made
Annotation Sequence	(7019,XX02)	3	SQ	1	Not set when no entry is made
> Private Creator	(7019,00XX)	1C* ¹	LO	1	"TOSHIBA_MEC_OT3"
> Annotation Object Number	(7019,XX04)	3	SS	1	Not set when no entry is made
> Annotation Type	(7019,XX05)	3	SS	1	Not set when no entry is made
> Annotation Font Size	(7019,XX06)	3	SS	1	Not set when no entry is made
> Annotation Font Color	(7019,XX07)	3	UL	1	Not set when no entry is made
> Display Position X of Result of Measure	(7019,XX08)	3	SS	1	Not set when no entry is made
> Display Position Y of Result of Measure	(7019,XX09)	3	SS	1	Not set when no entry is made
> Line Width	(7019,XX0A)	3	SS	1	Not set when no entry is made
> Line Color	(7019,XX0B)	3	UL	1	Not set when no entry is made
> Measurement Type	(7019,XX0C)	3	SS	1	Not set when no entry is made
> Character String Font Size	(7019,XX0D)	3	SS	1	Not set when no entry is made
> Character String Font Color	(7019,XX0E)	3	UL	1	Not set when no entry is made
> Line Position Sequence	(7019,XX10)	3	SQ	1	Not set when no entry is made
>> Private Creator	(7019,00XX)	1C* ¹	LO	1	"TOSHIBA_MEC_OT3"
>> Line Position x0	(7019,XX11)	3	SS	1	Not set when no entry is made
>> Line Position y0	(7019,XX12)	3	SS	1	Not set when no entry is made
>> Line Position x1	(7019,XX13)	3	SS	1	Not set when no entry is made
>> Line Position y1	(7019,XX14)	3	SS	1	Not set when no entry is made
>> Click Position X	(7019,XX15)	3	SS	1	Not set when no entry is made
>> Click Position Y	(7019,XX16)	3	SS	1	Not set when no entry is made
> ROI Center Position X	(7019,XX20)	3	SS	1	Not set when no entry is made
> ROI Center Position Y	(7019,XX21)	3	SS	1	Not set when no entry is made
> Length X from Center of ROI	(7019,XX22)	3	SS	1	Not set when no entry is made
> Length Y from Center of ROI	(7019,XX23)	3	SS	1	Not set when no entry is made
> Rotation Angle of ROI	(7019,XX24)	3	FL	1	Not set when no entry is made
> Vertex Number of Polygonal ROI	(7019,XX2A)	3	SS	1	Not set when no entry is made

Attribute Name	Tag	Type	VR	VM	Attribute Description
> Vertex Position of Polygonal ROI Sequence	(7019,XX2B)	3	SQ	1	Not set when no entry is made
>> Private Creator	(7019,00XX)	1C* ¹	LO	1	"TOSHIBA_MEC_OT3"
>> Vertex Position X	(7019,XX2C)	3	SS	1	Not set when no entry is made
>> Vertex Position Y	(7019,XX2D)	3	SS	1	Not set when no entry is made
> Edit Box Position X	(7019,XX30)	3	SS	1	Not set when no entry is made
> Edit Box Position Y	(7019,XX31)	3	SS	1	Not set when no entry is made
> Edit Box Character Set	(7019,XX32)	3	CS	1	Not set when no entry is made
> Edit Box String Length	(7019,XX33)	3	SS	1	Not set when no entry is made
> Edit Box String	(7019,XX34)	3	ST	1	Not set when no entry is made
> Start X Of Arrow	(7019,XX61)	3	SS	1	Not set when no entry is made
> Start Y Of Arrow	(7019,XX62)	3	SS	1	Not set when no entry is made
> End Y Of Arrow	(7019,XX63)	3	SS	1	Not set when no entry is made
> End X Of Arrow	(7019,XX64)	3	SS	1	Not set when no entry is made
> End X Of Arrow With Character	(7019,XX65)	3	SS	1	Not set when no entry is made
> End Y Of Arrow With Character	(7019,XX66)	3	SS	1	Not set when no entry is made
> Character Code Of Arrow With Character	(7019,XX67)	3	CS	1	Not set when no entry is made
> Character Length Of Arrow With Character	(7019,XX68)	3	SS	1	Not set when no entry is made
> Character String Of Arrow With Character	(7019,XX69)	3	ST	1	Not set when no entry is made
Image Type for MEDIS-DC media format	(7019,XX40)	3	CS	1	Not set when no entry is made
Information Change Sequence	(7019,XX45)	3	SQ	1	Not set when no entry is made
> Operators' Name	(0008,1070)	3	PN	1-n	Not set when no entry is made
> Private Creator	(7019,00XX)	1C* ¹	LO	1	"TOSHIBA_MEC_OT3"
> Information Change Date	(7019,XX46)	3	DA	1	Not set when no entry is made
> Information Change Time	(7019,XX47)	3	TM	1	Not set when no entry is made
> Changed Data Element	(XXXX,XXXX)	3	*	1	Not set when no entry is made
Pixel Size	(7019,XX50)	3	DS	2	Not set when no entry is made
Original Image SOP Instance UID	(7019,XX51)	3	UI	1	Not set when no entry is made
Original Image Series Instance UID	(7019,XX52)	3	UI	1	Not set when no entry is made
Copy Image Creation Date	(7019,XX53)	3	DA	1	Not set when no entry is made
Copy Image Creation Time	(7019,XX54)	3	TM	1	Not set when no entry is made

*1:Always set when the TWS creates private data.

9 Information Object Definition - Storage SCP

9.1 Entity Module Definitions

The acceptable information objects are the same as those defined at section 8"Information Object Definition - Storage SCU", it is recommended that the remote Storage SCU set the following tags:

Table 123

Attribute Name	Module	Tag	Type	Reasons
Patient's Name	Patient	(0010,0010)	2	To archive in MEDIS-DC media format* ¹
Study Date	General Study	(0008,0020)	2	To archive in MEDIS-DC media format* ¹
Study Time	General Study	(0008,0030)	2	To archive in MEDIS-DC media format* ¹
Institution Name	General Equipment	(0008,0080)	3	To archive in MEDIS-DC media format* ¹
Patient ID	Patient	(0010,0020)	2	To archive in MEDIS-DC media format* ¹
Patient's Birth Date	Patient	(0010,0030)	2	To archive in MEDIS-DC media format* ¹
Patient's Sex	Patient	(0010,0040)	2	To archive in MEDIS-DC media format* ¹

*1: MEDIS-DC media format is provided by MEDIS-DC.(MDS A 0008-1995)

10 Search Keys

10.1 Query/Retrieve SCU (C-Find)

The search keys used for the Query/Retrieve SCU (C-Find) are shown.

10.1.1 Patient Root Q/R Information Model - Find

10.1.1.1 Patient Level SCU Request

Table 124

Attribute Name	Tag	Type	Attribute Description
Patient's Name	(0010,0010)	R	Specified search condition
Patient ID	(0010,0020)	U	Specified search condition

Type: U=Unique Key, R=Required Key

10.1.1.2 Study Level SCU Request

Table 125

Attribute Name	Tag	Type	Attribute Description
Study Date	(0008,0020)	R	Specified search condition
Study Time	(0008,0030)	R	Specified search condition
Accession Number	(0008,0050)	R	Always set Length=0
Study ID	(0020,0010)	R	Specified search condition
Study Instance UID	(0020,000d)	U	Specified search condition

Type: U=Unique Key, R=Required Key

10.1.1.3 Series Level SCU Request

Table 126

Attribute Name	Tag	Type	Attribute Description
Modality	(0008,0060)	R	Specified search condition
Series Number	(0020,0011)	R	Specified search condition
Series Instance UID	(0020,000e)	U	Specified search condition

Type: U=Unique Key, R=Required Key

10.1.1.4 Image Level SCU Request

Table 127

Attribute Name	Tag	Type	Attribute Description
Image Number	(0020,0013)	R	Specified search condition
SOP Instance UID	(0008,0018)	U	Specified search condition

Type: U=Unique Key, R=Required Key

10.1.2 Study Root Q/R Information Model - Find

10.1.2.1 Study Level SCU Request

Table 128

Attribute Name	Tag	Type	Attribute Description
Study Date	(0008,0020)	R	Specified search condition
Study Time	(0008,0030)	R	Specified search condition
Accession Number	(0008,0050)	R	Always set Length=0
Patient's Name	(0010,0010)	R	Specified search condition
Patient ID	(0010,0020)	R	Specified search condition
Study ID	(0020,0010)	R	Specified search condition
Study Instance UID	(0020,000d)	U	Specified search condition

Type: U=Unique Key, R=Required Key

10.1.2.2 Series Level SCU Request

Attributes for the Series Level of the Study Root Query/Retrieve Information Model are the same as the attributes for the Series Level of the Patient Root Query/Retrieve Information Model described in 10.1.1.3

10.1.2.3 Image Level SCU Request

Attributes for the Image Level of the Study Root Query/Retrieve Information Model are the same as the attributes for the Image Level of the Patient Root Query/Retrieve Information Model described in 10.1.1.4

10.1.3 Patient/Study Only Q/R Information Model - Find

10.1.3.1 Patient Level SCU Request

Attributes for the Patient Level of the Patient/Study Only Query/Retrieve Information Model are the same as the attributes for the Patient Level of the Patient Root Query/Retrieve Information Model described in 10.1.1.1

10.1.3.2 Study Level SCU Request

Attributes for the Study Level of the Patient/Study Only Query/Retrieve Information Model are the same as the attributes for the Study Level of the Patient Root Query/Retrieve Information Model described in 10.1.1.2

11 DIMSE-Service and Attributes

The attribute listed in the following tables represent a small set of the possible attributes which could be supported by each services.

11.1 DIMSE-Services

Table 129

SOP Class	DIMSE-Service	Reference	Usage SCU/SCP *1	Usage
Basic Film Session SOP Class	N-CREATE	11.2.1	M/M	Used
Basic Film Box SOP Class	N-CREATE	11.3.1	M/M	Used
	N-ACTION	-	M/M	Used
Basic Grayscale Image Box SOP Class	N-SET	11.4.1	M/M	Used
Basic Color Image Box SOP Class	N-SET	11.5.1	M/M	Used
Printer SOP Class	N-EVENT-REPORT	-	M/M	Used
	N-GET	11.6.1	U/M	Used

*1 : M=Mandatory, U=User Option

11.2 Basic Film Session SOP Class

11.2.1 N-CREATE Attributes

Table 130

Attribute name	Tag	Usage SCU/SCP	Attributes
Number of Copies	(2000,0010)	U/M	Always set. (1)
Print Priority	(2000,0020)	U/M	Always set. (default:"MED")
Medium Type	(2000,0030)	U/M	Always set. (default: "BLUE FILM")
Film Destination	(2000,0040)	U/M	16bytes maximum, or no set.
Memory Allocation	(2000,0060)	U/U	12bytes maximum, or no set.

11.3 Basic Film Box SOP Class

11.3.1 N-CREATE Attributes

Table 131

Attribute name	Tag	Usage SCU/SCP	Attributes
Image Display Format	(2010,0010)	M/M	Always set. ("STANDARD\C,R")
Film Orientation	(2010,0040)	U/M	Always set.
Film Size ID	(2010,0050)	U/M	Always set.
Magnification Type	(2010,0060)	U/M	Always set.
Smoothing Type	(2010,0080)	U/U	If it defines, always set.
Border Density	(2010,0100)	U/U	Always set.
Empty Image Density	(2010,0110)	U/U	Always set.
Min Density	(2010,0120)	U/U	0 to 65535, or no set.
Max Density	(2010,0130)	U/M	0 to 65535, or no set.
Trim	(2010,0140)	U/U	Always set
Configuration Information	(2010,0150)	U/M	1024bytes maximum, or no set.
Referenced Film Session Sequence	(2010, 0500)	M/M	Always set.
>Referenced SOP Class UID	(0008,1150)	M/M	Always set.
>Referenced SOP Instance UID	(0008,1155)	M/M	Always set.

11.4 Basic Grayscale Image Box SOP Class

11.4.1 N-SET Attributes

Table 132

Attribute name	Tag	Usage SCU/SCP	Attributes
Image position	(2020, 0010)	M/M	Always set.
Polarity	(2020, 0020)	U/M	Always set. ("NORMAL")
Magnification Type	(2010, 0060)	U/U	16bytes maximum, or no set.
Smoothing Type	(2010, 0080)	U/U	16bytes maximum, or no set.
Basic Grayscale Image Sequence	(2020, 0110)	M/M	Always set.
>Samples Per Pixel	(0028, 0002)	M/M	Always set. (1)
>Photometric Interpretation	(0028, 0004)	M/M	Always set. ("MONOCHROME2")
>Rows	(0028, 0010)	M/M	Always set.
>Columns	(0028, 0011)	M/M	Always set.
>Pixel Aspect Ratio	(0028, 0034)	MC/M	Always set. (1\1)
>Bits Allocated	(0028, 0100)	M/M	Always set. (8)
>Bits Stored	(0028, 0101)	M/M	Always set. (8)
>High Bit	(0028, 0102)	M/M	Always set. (7)
>Pixel Representation	(0028, 0103)	M/M	Always set. (0)
>Pixel Data	(7FE0, 0010)	M/M	Always set.

11.5 Basic Color Image Box SOP Class

11.5.1 N-SET Attributes

Table 133

Attribute name	Tag	Usage SCU/SCP	Attributes
Image position	(2020, 0010)	M/M	Always set.
Polarity	(2020, 0020)	U/M	Always set. ("NORMAL")
Magnification Type	(2010, 0060)	U/U	16bytes maximum, or no set.
Smoothing Type	(2010, 0080)	U/U	16bytes maximum, or no set.
Basic Color Image Sequence	(2020, 0111)	M/M	Always set.
>Samples Per Pixel	(0028, 0002)	M/M	Always set. (3)
>Photometric Interpretation	(0028, 0004)	M/M	Always set. ("RGB")
>Planar Configuration	(0028,0006)	M/M	Always set (1)
>Rows	(0028, 0010)	M/M	Always set.
>Columns	(0028, 0011)	M/M	Always set.
>Pixel Aspect Ratio	(0028, 0034)	MC/M	Always set. (1\1)
>Bits Allocated	(0028, 0100)	M/M	Always set. (8)
>Bits Stored	(0028, 0101)	M/M	Always set. (8)
>High Bit	(0028, 0102)	M/M	Always set. (7)
>Pixel Representation	(0028, 0103)	M/M	Always set. (0)
>Pixel Data	(7FE0, 0010)	M/M	Always set.

11.6 Printer SOP Class

11.6.1 N-GET Attributes

Table 134

Attribute name	Tag	Usage SCU/SCP	Attributes
Printer Status	(2110,0010)	U/M	
Printer Status Info	(2110,0020)	U/M	
Printer Name	(2110,0030)	U/U	
Manufacturer	(0008,0070)	U/U	
Manufacturer Model Name	(0008,1090)	U/U	