

**DICOM  
CONFORMANCE STATEMENT  
FOR  
DIAGNOSTIC ULTRASOUND SYSTEM  
APLIO MODEL SSA-770A**

**TOSHIBA CORPORATION**

**IMPORTANT !**

- (1) No part of this manual may be copied or reprinted, in whole or in part, without written permission.
- (2) The contents of this manual are subject to change without prior notice and without our legal obligation.

## Table Of Contents

|         |   |   |
|---------|---|---|
| 1       | INTRODUCTION.....   | 1 |
| 1.1     | REFERENCE.....  | 1 |
| 1.2     | DEFINITIONS.....  | 1 |
| 1.3     | ACRONYMS, ABBREVIATIONS AND SYMBOLS.....                                    | 3 |
| 2       | IMPLEMENTATION MODEL.....   | 4 |
| 2.1     | MOD MEDIUM.....   | 4 |
| 2.1.1   | APPLICATION DATA FLOW DIAGRAM.....  | 4 |
| 2.1.2   | FUNCTIONAL DEFINITIONS OF AE <sub>s</sub> .....                             | 5 |
| 2.1.3   | SEQUENCING OF REAL WORLD ACTIVITIES.....                                    | 5 |
| 2.1.3.1 | FEATURES.....   | 5 |
| 2.1.3.2 | OPERATIONS.....   | 6 |
| 2.1.4   | FILE META INFORMATION FOR INPLEMENTATION CLASS UID AND<br>VERSION NAME..... | 6 |
| 3       | AE SPECIFICATIONS.....  | 7 |
| 3.1     | MEDIUM AE SPECIFICAT ION.....   | 7 |
| 3.1.1   | FILE META INFORMATION FOR THE APPLICATION ENTITY.....                       | 7 |
| 3.1.2   | REAL WORLD ACTIVITY - REMOVABLE MEDIA.....                                  | 7 |
| 4       | AUGMENTATED AND PRIVATE APPLICATION PROFILES.....                           | 8 |
| 4.1     | AUGMENTED APPLICATION PROFILES.....   | 8 |
| 4.2     | SOP CLASS AUGMENTATIONS.....  | 8 |
| 4.3     | DIRECTORY AUGMENTATIONS.....  | 8 |
| 4.4     | OTHER AUGMENTATIONS.....  | 8 |
| 4.5     | PRIVATE APPLICATION PROFILES.....   | 9 |

5 EXTENSIONS/SPECIALIZATIONS/PRIVATIZATIONS ..... 9

6 CONFIGURATION ..... 9

7 SUPPORT OF EXTENDED CHARACTER SETS ..... 9

8 MOD STORAGE INFORMATION OBJECT DEFINITION ..... 10

8.1 ENTRY MODULE DEFINITIONS ..... 10

8.1.1 DICOM FILE META INFORMATION ..... 10

8.1.1.1 DICOM FILE META INFORMATION OF DIRECTORY IOD ..... 10

8.1.1.2 DICOM FILE META INFORMATION IMAGE IOD ..... 11

8.1.2 BASIC DIRECTORY INFORMATION OBJECT ..... 12

8.1.2.1 FILE-SET IDENTIFICATION MODULE ..... 12

8.1.2.2 DIRECTORY INFORMATION MODULE ..... 12

8.1.3 DEFINITION OF SPECIFIC DIRECTORY RECORDS ..... 13

8.1.3.1 PATIENT KEYS ..... 13

8.1.3.2 STUDY KEYS ..... 13

8.1.3.3 SERIES KEYS ..... 14

8.1.3.4 IMAGE KEYS ..... 15

## 1 INTRODUCTION

This document is a DICOM Conformance Statement for Toshiba's Diagnostic Ultrasound System APLIO MODEL SSA-770A. It is intended to provide the reader with the knowledge of how to integrate this product within a DICOM compliant hospital archive. It details the DICOM Service Classes, Information Objects, and Communication Protocols which are supported by this product as follows:

- MOD Medium Storage Service Class (FSC/FSR/FSU)

If the reader is unfamiliar with DICOM, it is recommended that they read the DICOM Specification prior to reading this conformance statement. Also note that this document is formatted according to the DICOM Specification, Part 2: Conformance.

### 1.1 REFERENCE

ACR-NEMA Digital Imaging and Communications in Medicine, DICOM V3.0.

### 1.2 DEFINITIONS

- **Application Profile:** A Media Storage Application Profile defines a selection of choices at the various layers of the DICOM Media Storage Model which are applicable to a specific need or context in which the media interchange is intended to be performed.
- **DICOM File Service:** The DICOM File Service specifies a minimum abstract view of files to be provided by the Media Format Layer. Constraining access to the content of files by the Application Entities through such a DICOM File Service boundary ensures Media Format and Physical Media independence.
- **DICOM File:** A DICOM File is a File with a content formatted according to the requirements of this Part of the DICOM Standard. In particular such files shall contain, the File Meta Information and a properly formatted Data Set.
- **DICOMDIR File:** A unique and mandatory DICOM File within a File-set which contains the Media Storage Directory SOP Class. This File is given a single component File ID, DICOMDIR.
- **File:** A File is an ordered string of zero or more bytes, where the first byte is at the beginning of the file and the last byte at the end of the File. Files are identified by a unique File ID and may be written, read and/or deleted.
- **File ID:** Files are identified by a File ID which is unique within the context of the File-set they belong to. A set of ordered File ID Components (up to a maximum of eight) forms a File ID.

- **File ID Component:** A string of one to eight characters of a defined character set.
- **File Meta Information:** The File Meta Information includes identifying information on the encapsulated Data Set. It is a mandatory header at the beginning of every DICOM File.
- **File-set:** A File-set is a collection of DICOM Files (and possibly non-DICOM Files) that share a common naming space within which File IDs are unique.
- **File-set Creator:** An Application Entity that creates the DICOMDIR File (see section 8.6) and zero or more DICOM Files.
- **File-set Reader:** An Application Entity that accesses one or more files in a File-set.
- **File-set Updater:** An Application Entity that accesses Files, creates additional Files, or deletes existing Files in a File-set. A File-set Updater makes the appropriate alterations to the DICOMDIR file reflecting the additions or deletions.
- **DICOM File Format:** The DICOM File Format provides a means to encapsulate in a File the Data Set representing a SOP Instance related to a DICOM Information Object.
- **Media Format:** Data structures and associated policies which organizes the bit streams defined by the Physical Media format into data file structures and associated file directories.
- **Media Storage Model:** The DICOM Media Storage Model pertains to the data structures used at different layers to achieve interoperability through media interchange.
- **Physical Media:** A piece of material with recording capabilities for streams of bits. Characteristics of a Physical Media include form factor, mechanical characteristics, recording properties and rules for recording and organizing bit streams in accessible structures
- **Information Object Definition (IOD)** - An IOD is a data model which is an abstraction of real-world information. This data model defines the nature and attributes relevant to the class of real-world objects represented.
- **Service/Object Pair (SOP) Class** - A SOP Class is defined by the union of an Information Object Definition and a set of DIMSE Services. A DICOM Application Entity may support one or more SOP Classes. Each SOP Class is uniquely identified by a SOP Class UID.
- **SOP Instance** - A specific occurrence of a Information Object.
- **Transfer Syntax** - The Transfer Syntax is a set of encoding rules that allow DICOM Application Entities to negotiate the encoding techniques (e.g. data element structure, byte ordering, compression) they are able to support. The Transfer Syntax is negotiated during Association Negotiation.
- **Unique Identifier (UID)** - A Unique Identifier is a globally unique, ISO compliant, ASCII-numeric string. It guarantees uniqueness across multiple countries, sites, vendors and equipment.

### 1.3 ACRONYMS, ABBREVIATIONS AND SYMBOLS

The following acronyms and abbreviations are used in this document.

- ACC American College of Cardiology
- ACR American College of Radiology
- ASCII American Standard Code for Information Interchange
- AE Application Entity
- ANSI American National Standards Institute
- CEN TC251 Comite Europeen de Normalisation - Technical Committee  
251 - Medical Informatics
- DICOM Digital Imaging and COmmunications in Medicine
- DIMSE DICOM Message Service Element
- DIMSE-C DICOM Message Service Element-Composite
- DIMSE-N DICOM Message Service Element-Normalized
- FSC File-Set Creator
- FSR File-Set Reader
- FSU File-Set Updater
- HIS Hospital Information System
- HL7 Health Level 7
- IE Information Entity
- IOD Information Object Definition
- ISO International Standard Organization
- JIRA Japan Industries Association of Radiological Systems
- NEMA National Electrical Manufacturers Association
- RIS Radiology Information System
- SOP Service Object Pair
- UID Unique Identifier

## 2 IMPLEMENTATION MODEL

### 2.1 MOD MEDIUM

The DICOM archive application serves as an interface with the MOD offline medium device. It incorporates the offline media directory into the browser and copies SOP Instances to a medium or retrieves SOP Instances from a medium to local storage. The DICOM Archive application supports Standard MOD media. The FSU role updates new SOP Instances only to media with preexisting File-sets conforming to the Application Profiles supported. The contents of DICOMDIR will be temporarily stored in the Archive Database.

#### 2.1.1 APPLICATION DATA FLOW DIAGRAM




Figure 1



### **2.1.2 FUNCTIONAL DEFINITIONS OF AE<sub>s</sub>**

Medium AE is used to create/read/modify studies/images to/from an offline DICOM MOD Medium. It therefore performs the following tasks:

- Builds DICOM SC, US Information Objects.
- CREATES a DICOMDIR file that represents the contents of the DICOM SC, US Information Objects to be recorded.
- RECORDS DICOM SC, US Information Objects from local storage to a MOD medium.
- READS the DICOMDIR file that represents the contents of the data as recorded.
- DISPLAYS the ordered list of studies/images, identifying information, or icon images.
- READS the selected studies/images from a MOD medium and displays them on the monitor of the DICOM Reader.
- READS the File-set of the MOD medium and writes it to the local storage of DICOM Reader.
- MODIFIES the DICOMDIR file and the studies/images.

### **2.1.3 SEQUENCING OF REAL WORLD ACTIVITIES**

The DICOM Archive application will not perform updates before the directory information of DICOMDIR is read completely.

#### **2.1.3.1 FEATURES**

- Operator requests storage of study/images to a MOD medium after selecting studies from the Study List.
- Operator requests immediate and automatic storage of an image after capture.
- Operator requests retrieval of study/images from a MOD medium to the local disk.
- Storage requests are placed in a queue and are executed in the background.

### 2.1.3.2 OPERATIONS

The operations for manual image transfer/modification are described below:

Operation-1

Step-1: Display the image to be transferred/modified.

Step-2: Request transfer/modification.

Operation-2

Step-1: Select the study to be transferred/modified.

Step-2: Request transfer/modification.

### 2.1.4 FILE META INFORMATION FOR IMPLEMENTATION CLASS UID AND VERSION NAME

Medium File Meta Information will specify the following Implementation Identifying Information

Table 1

| System            | Implementation Class UID                | Implementation Version Name | File Meta Information Version |
|-------------------|---|-----------------------------|-------------------------------|
| APLIO<br>SSA-770A | 1.2.392.200036.9116.7.8.10.<br>46.6.1.1 | TM_APLIO_1.0                | Version 1                     |

### 3 AE SPECIFICATIONS

#### 3.1 MEDIUM AE SPECIFICATION

Medium Activity : Medium AE provides Conformance to DICOM Interchange of the Media Storage Service Class. The Application Profiles and roles are listed in the following table.

Table 2

| Application Profiles Supported       | Real-World Activity            | Roles | SC Option   |
|--------------------------------------|--------------------------------|-------|-------------|
| AUG-US-ID-MF-MOD128 (128MB 90mm MOD) | Transfer of SC, US Examination | FSC   | Interchange |
| AUG-US-ID-MF-MOD230 (230MB 90mm MOD) | Modifier of SC, US Examination | FSU   | Interchange |
| AUG-US-ID-MF-MOD540 (540MB 90mm MOD) | Review of SC, US Examination   | FSR   | Interchange |
| AUG-US-ID-MF-MOD650 (650MB 90mm MOD) |                                |       |             |
| AUG-US-ID-MF-MOD12 (1.2GB 90mm MOD)  |                                |       |             |

##### 3.1.1 FILE META INFORMATION FOR THE APPLICATION ENTITY

The Source Application Entity Title is set by the user in the configuration file.

##### 3.1.2 REAL WORLD ACTIVITY - REMOVABLE MEDIA

A DICOM conformant Magneto-Optical Disk (MOD) is created when a non-conformant MOD is inserted into APLIO and one or more DICOM Exams are transferred to the MOD. When Exams are first transferred, their files are added to the MOD in DICOM Part 10 format and a valid DICOMDIR is created and saved to the MOD. APLIO can add images to an existing DICOM conformant MOD and update its DICOMDIR. APLIO can be a File-set Reader and a File-set Updater.

## 4 AUGMENTATED AND PRIVATE APPLICATION PROFILES

### 4.1 AUGMENTED APPLICATION PROFILES

Augmented Application Profiles cover from Standard Application Profiles in addition to dealing with SC IOD.

Table 3

| Augmented Application Profiles   | Standard Application Profiles |
|----------------------------------|-------------------------------|
| AUG-US-ID-MF-MOD128              | STD-US-ID-MF-MOD128           |
| AUG-US-ID-MF-MOD230              | STD-US-ID-MF-MOD230           |
| AUG-US-ID-MF-MOD540              | STD-US-ID-MF-MOD540           |
| AUG-US-ID-MF-MOD650 <sup>1</sup> | STD-US-ID-MF-MOD650           |
| AUG-US-ID-MF-MOD12 <sup>1</sup>  | STD-US-ID-MF-MOD12            |

<sup>1</sup> AUG-US-ID-MF-MOD650 and AUG-US-ID-MF-MOD12 are replaced from 130mm MOD to 90mm MOD

### 4.2 SOP CLASS AUGMENTATIONS

Aforementioned Application Profiles sport following SOP Class UID and Transfer Syntax.

Table 4

| Information Object Definition   | SOP Class UID             | Transfer Syntax                        | Transfer Syntax UID |
|---------------------------------|---------------------------|--|---------------------|
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Explicit VR Little Endian Uncompressed | 1.2.840.10008.1.2.1 |

### 4.3 DIRECTORY AUGMENTATIONS

Not applicable to this product

### 4.4 OTHER AUGMENTATIONS

Not applicable to this product

#### **4.5 PRIVATE APPLICATION PROFILES**

Not applicable to this product

#### **5 EXTENSIONS/SPECIALIZATIONS/PRIVATIZATIONS**

Not applicable to this product.

#### **6 CONFIGURATION**

Not applicable to this product.

#### **7 SUPPORT OF EXTENDED CHARACTER SETS**

ISO-IR 100 (Latin 1 alphabet No.1) Supplementary set of ISO8859

## 8 MOD STORAGE INFORMATION OBJECT DEFINITION

### 8.1 ENTRY MODULE DEFINITIONS

Table 5

| Object                                   | Module/Key                     | Reference         | Usage <sup>1</sup> |
|--|--------------------------------|-------------------|--------------------|
| DICOM File Meta Information              |                                | 8.1.1.1 / 8.1.1.2 | M                  |
| Basic Directory Information Object       | File-set Identification Module | 8.1.2.1           | M                  |
|  | Directory Information Module   | 8.1.2.2           | U                  |
| Definition of Specific Directory Records | PATIENT keys                   | 8.1.3.1           | M                  |
|  | STUDY keys                     | 8.1.3.2           | M                  |
|  | SERIES keys                    | 0                 | M                  |
|  | IMAGE keys                     | 8.1.3.4           | M                  |

<sup>1</sup> M=Mandatory, C=Conditional, U=User option

#### 8.1.1 DICOM FILE META INFORMATION

##### 8.1.1.1 DICOM FILE META INFORMATION OF DIRECTORY IOD

Table 6

| Attribute Name                 | Tag          | Type | Attribute Description               |
|--------------------------------|--------------|------|-------------------------------------|
| File Preamble                  | No Tag       | 1    | Always set                          |
| DICOM Prefix                   | No Tag       | 1    | Always set ("DICM")                 |
| Group Length                   | (0002, 0000) | 1    | Always set                          |
| File Meta Information Version  | (0002, 0001) | 1    | Always set (0001H)                  |
| Media Storage SOP Class UID    | (0002, 0002) | 1    | Always set ("1.2.840.10008.1.3.10") |
| Media Storage SOP Instance UID | (0002, 0003) | 1    | Always set                          |
| Transfer Syntax UID            | (0002, 0010) | 1    | Always set                          |
| Implementation Class UID       | (0002, 0012) | 1    | Always set                          |
| Implementation Version Name    | (0002,0013)  | 3    | Always set ("TM_APLIO_1.0")         |

### 8.1.1.2 DICOM FILE META INFORMATION IMAGE IOD

Table 7

| Attribute Name                 | Tag          | Type | Attribute Description       |
|--------------------------------|--------------|------|-----------------------------|
| File Preamble                  | No Tag       | 1    | Always set                  |
| DICOM Prefix                   | No Tag       | 1    | Always set ("DICM")         |
| Group Length                   | (0002, 0000) | 1    | Always set                  |
| File Meta Information Version  | (0002, 0001) | 1    | Always set (0001H)          |
| Media Storage SOP Class UID    | (0002, 0002) | 1    | Always set                  |
| Media Storage SOP Instance UID | (0002, 0003) | 1    | Always set                  |
| Transfer Syntax UID            | (0002, 0010) | 1    | Always set                  |
| Implementation Class UID       | (0002, 0012) | 1    | Always set                  |
| Implementation Version Name    | (0002, 0013) | 3    | Always set ("TM_APLIO_1.0") |

## 8.1.2 BASIC DIRECTORY INFORMATION OBJECT

### 8.1.2.1 FILE-SET IDENTIFICATION MODULE

Table 8

| Attribute Name | Tag          | Type | Attribute Description |
|----------------|--------------|------|-----------------------|
| File-set ID    | (0004, 1130) | 2    | Always set            |

### 8.1.2.2 DIRECTORY INFORMATION MODULE

Table 9

| Attribute Name  | Tag          | Type | Attribute Description   |
|---|--------------|------|---|
| Offset of the First Directory Record of the Root Directory Entity | (0004, 1200) | 1    | Always set  |
| Offset of the Last Directory Record of the Root Directory Entity  | (0004, 1202) | 1    | Always set  |
| File-set Consistency Flag   | (0004, 1212) | 1    | Always set  |
| Directory Record Sequence   | (0004, 1220) | 2    | Set if present in image object.   |
| >Offset of the Next Directory Record                              | (0004, 1400) | 1C   | Set if present in image object.   |
| >Record In-use Flag   | (0004, 1410) | 1C   | Always set (FFFFH)  |
| >Offset of Referenced Lower-Level Directory Entity                | (0004, 1420) | 1C   | Set if present in image object.   |
| >Directory Record Type  | (0004, 1430) | 1C   | Always set<br>"PATIENT",<br>"STUDY",<br>"SERIES",<br>"IMAGE",<br>"CURVE", or<br>"PRIVATE" |
| >Referenced File ID   | (0004, 1500) | 1C   | Contains the filename in media for the Directory Records of Type IMAGE and PRIVATE        |
| >Referenced SOP Class UID in File                                 | (0004, 1510) | 1C   | For the Directory Record Types IMAGE and PRIVATE  |
| >Referenced SOP Instance UID in File                              | (0004, 1511) | 1C   | For the Directory Record Types IMAGE and PRIVATE  |
| >Referenced Transfer Syntax UID in                                | (0004, 1512) | 1C   | For the Directory Record Types  |



|                         |           |  |                   |
|-------------------------|-----------|--|-------------------|
| File                    |           |  | IMAGE and PRIVATE |
| > Record Selection Keys | see below |  |                   |

### 8.1.3 DEFINITION OF SPECIFIC DIRECTORY RECORDS

#### 8.1.3.1 PATIENT KEYS

Table 10

| Attribute Name         | Tag          | Type | Attribute Description           |
|------------------------|--------------|------|---------------------------------|
| Specific Character Set | (0008,0005)  | 1C   | Always set ("ISO_IR 100")       |
| Patient's Name         | (0010, 0010) | 2    | Set if present in image object. |
| Patient ID             | (0010, 0020) | 1    | Always set                      |
| Patient's Birth Date   | (0010,0030)  | 3    | Set if present in image object. |
| Patient's Sex          | (0010,0040)  | 3    | Set if present in image object. |
| Patient's Age          | (0010,1010)  | 3    | Set if present in image object. |
| Patient's Size         | (0010,1020)  | 3    | Set if present in image object. |
| Patient's Weight       | (0010,1030)  | 3    | Set if present in image object. |
| Patient's Address      | (0010,1040)  | 3    | Set if present in image object. |
| Patient Comment        | (0010,4000)  | 3    | Set if present in image object. |

#### 8.1.3.2 STUDY KEYS

Table 11

| Attribute Name                  | Tag          | Type | Attribute Description                                |
|---------------------------------|--------------|------|--|
| Specific Character Set          | (0008,0005)  | 1C   | Always set ("ISO_IR 100")                            |
| Study Date                      | (0008, 0020) | 1    | Always set   |
| Study Time                      | (0008, 0030) | 1    | Always set   |
| Accession Number                | (0008, 0050) | 2    | Set if present in image object.                      |
| Referring Physician's Name      | (0008,0090)  | 3    | Set if present in image object.                      |
| Study Description               | (0008, 1030) | 2    | Set if present in image object.                      |
| Admitting Diagnosis Description | (0008,1080)  | 3    | Set if present in image object.                      |
| Study Instance UID              | (0020, 000D) | 1C   | Always set   |
| Study ID                        | (0020, 0010) | 1    | Generates automatically, if not present. Value = "-" |

### 8.1.3.3 SERIES KEYS

Table 12

| Attribute Name              | Tag          | Type | Attribute Description           |
|-----------------------------|--------------|------|---------------------------------|
| Specific Character Set      | (0008,0005)  | 1C   | Always set ("ISO_IR 100")       |
| Series Date                 | (0008,0021)  | 3    | Set if present in image object. |
| Modality                    | (0008, 0060) | 1    | Always set                      |
| Institution Name            | (00080080)   | 3    | Set if present in image object. |
| Institution Address         | (00080081)   | 3    | Set if present in image object. |
| Performing Physician's Name | (0008,1050)  | 3    | Set if present in image object. |
| Operator's Name             | (0008,1070)  | 3    | Set if present in image object. |
| Series Instance UID         | (0020, 000E) | 1    | Always set                      |
| Series Number               | (0020, 0011) | 1    | Always set                      |

## 8.1.3.4 IMAGE KEYS

Table 13

| Attribute Name              | Tag          | Type | Attribute Description           |
|-----------------------------|--------------|------|---------------------------------|
| Specific Character Set      | (0008,0005)  | 1C   | Always set ("ISO_IR 100")       |
| Image Type                  | (0008, 0008) | 3    | Always set                      |
| SOP Class UID               | (0008,0016)  | 3    | Always set                      |
| SOP Instance UID            | (0008,0018)  | 3    | Always set                      |
| Image Date                  | (0008,0023)  | 3    | Set if present in image object. |
| Image Time                  | (0008,0033)  | 3    | Set if present in image object. |
| Instance Number             | (0020, 0013) | 1    | Always set                      |
| Rows                        | (0028, 0010) | 3    | Always set                      |
| Columns                     | (0028, 0011) | 3    | Always set                      |
| Calibration Image           | (0050,0004)  | 3    | Always set                      |
| Pixel Spacing               | (0028, 0030) | 3    | Set if present in image object. |
| Icon Image Sequence         | (0088,0200)  | 3    | Always set                      |
| >Samples per Pixel          | (0028,0002)  |      | 1                               |
| >Photometric Interpretation | (0028,0004)  |      | "MONOCHROME2"                   |
| >Rows                       | (0028,0010)  |      | 64                              |
| >Columns                    | (0028,0011)  |      | 64                              |
| >Bits Allocated             | (0028,0100)  |      | 8                               |
| >Bits Stored                | (0028,0101)  |      | 8                               |
| >High Bit                   | (0028,0102)  |      | 7                               |
| >Pixel Representation       | (0028,0103)  |      | 0 (Unsigned)                    |
| >Pixel Data                 | (7FE0,0010)  |      | Icon Image pixel data           |