

**PRODUCT SECURITY PATCH LIST
FOR
ULTRASOUND IMAGE ANALYSIS PROGRAM
(2G984-0051EN)**

CANON MEDICAL SYSTEMS CORPORATION

Trademarks

Windows® is a registered trademark of Microsoft Corporation in the United States and other countries.

This document may also include trademarks or registered trademarks of other companies.

Note that the trademark symbol "™" and the registered trademark symbol "®" may or may not be used in this document.

IMPORTANT!

1. No part of this document may be copied or reprinted, in whole or in part, without written permission.
2. The contents of this document are subject to change without prior notice and without legal obligation.
3. The contents of this document are correct to the best of our knowledge. Please inform us of any ambiguous or erroneous descriptions, missing information, etc.

REVISION RECORD					
REV.	DATE (YYYY-MM)	REASON /AUTHOR	PAGE CHANGED	SER. No.	DOC. PRODUCT.
INI.	2025-05	Ookouchi	-----		NASU

Security Patch Package name

NA

Applicable Products for the Security Patch

CUW-U001S	V1.0 SP0003* - SP0004*
	V2.0 SP0000* - SP0005*

Steps to Apply the Security Patch

Please follow the instruction of the administrator of the facility or organization. For details, please contact the customer service for each system.

Verified Security Patch List

KB	MSRC	explanation(en-us)	URL/File Name	url
KB2565063_x64	MS11-025	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package (KB2565063)	vc_redist_x64.exe	https://catalog.s.download.windowssupdate.com/msdownload/update/software/secu/2011/07/vcredist_x64_15d032d669078aa6f0f7fd1cbf4115a070bd034d.exe
VS2013_x64	n/a	Microsoft Visual C++ 2013 Update 5 Redistributable Package	vc_redist_x64.exe	https://download.microsoft.com/download/c/c/2/cc2df5f8-4454-44b4-802d-5ea68d086676/vcredist_x64.exe
VS2019_x86	n/a	Microsoft Visual C++ 2015-2019 Redistributable (x86)	vc_redist.x86.exe	https://download.visualstudio.microsoft.com/download/pr/48431a06-59c5-4b63-a102-20b66a521863/CAA38FD474164A38AB47AC1755C8CCCA5CCFACFA9A874F62609E6439924E87EC/VC_redist.x86.exe
VS2019_x64	n/a	Microsoft Visual C++ 2015-2019 Redistributable (x64)	vc_redist.x64.exe	https://download.visualstudio.microsoft.com/download/pr/48431a06-59c5-4b63-a102-20b66a521863/4B5890EB1AEFD F8DFA3234B5032147EB90F050C5758A80901B201AE969780107/VC_redist.x64.exe
KB5055523	n/a	2025-04 Cumulative Update for Windows 11 Version 24H2 for x64-based Systems (KB5055523)	windows11.0-kb5055523-x64.msu	https://catalog.sfdelivery.microsoft.com/filestreamingservice/files/beae415c-d56f-477c-9a3a-3aa4336890f6/public/windows11.0-kb5055523-x64_b1df8c7b11308991a9c45ae3fba6caa0e2996157.msu

KB	MSRC	explanation(en-us)	URL/File Name	url
KB5054979	n/a	2025-04 Cumulative Update for .NET Framework 3.5 and 4.8.1 for Windows 11, version 24H2 for x64 (KB5054979)	windows11.0-kb5054979-x64-ndp481.msu	https://catalog.sf.dl.delivery.mp.microsoft.com/filestreamingservice/files/02cd9c83-8312-424d-9a06-b042095804a8/public/windows11.0-kb5054979-x64-ndp481_8e2f730bc747de0f90aaee95d4862e4f88751c07.msu
Edge_x64_135.0.3179.54	n/a	Microsoft Edge-Stable Channel Version 135 Update for x64 based Editions (Build 135.0.3179.54)	microsoftedgeenterprisex64_bbbbee981ef38bda2759e4ec8855121348882ac6e.cab	https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/updt/2025/04/microsoftedgeenterprisex64_bbbbee981ef38bda2759e4ec8855121348882ac6e.cab
KB890830_v5.132	n/a	Windows Malicious Software Removal Tool x64 - v5.132 (KB890830)	windows-kb890830-x64-v5.132.exe	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/uprl/2025/02/windows-kb890830-x64-v5.132_99e2f6f80f1b15eaff192caf421e2c4c139796d.exe



LEGAL MANUFACTURER

CANON MEDICAL SYSTEMS CORPORATION

1385, SHIMOISHIGAMI, OTAWARA-SHI, TOCHIGI 324-8550, JAPAN