

**PRODUCT SECURITY PATCH LIST
FOR
MRI SCANNER
(2G984-0140EN)**

CANON MEDICAL SYSTEMS CORPORATION

Trademarks

Windows® is a registered trademark of Microsoft Corporation in the United States and other countries.

This document may also include trademarks or registered trademarks of other companies.

Note that the trademark symbol "™" and the registered trademark symbol "®" may or may not be used in this document.

IMPORTANT!

1. No part of this document may be copied or reprinted, in whole or in part, without written permission.
2. The contents of this document are subject to change without prior notice and without legal obligation.
3. The contents of this document are correct to the best of our knowledge. Please inform us of any ambiguous or erroneous descriptions, missing information, etc.

REVISION RECORD

REV.	DATE (YYYY-MM)	REASON /AUTHOR	PAGE CHANGED	SER. No.	DOC. PRODUCT.
INI.	2026-03	Shishido	-----		NASU

Security Patch Package name

SPZM45-20085_202512MR_hotfix_win10_LTSC2019_x64_enu

Applicable Products for the Security Patch

Vantage Galan 3T	MRT-3020	V6.0 SP2015*
		V7.0 SP0070*
		V8.0 SP0040* - V8.0 SP0044*
		V8.0 SP0050*
		V9.0 SP0020* - V9.0 SP0022*
Vantage Orian	MRT-1550	V9.0 SP0031* - V9.0 SP0055*
		V6.0 SP2015* - V6.0 SP2025*
		V7.0 SP0070* - V7.0 SP0075*
		V8.0 SP0040* - V8.0 SP0044*
		V8.0 SP0050*
Vantage Elan	MRT-2020	V9.0 SP0020* - V9.0 SP0022*
		V9.0 SP0030* - V9.0 SP0055*
		V6.0 SP1080* - V6.0 SP1085*
Vantage Titan 3T	MRT-3010	V8.1 SP0000* - V8.1 SP0032*
Vantage Titan	MRT-1510 MRT-1504	V8.1 SP0050* - V8.1 SP0055*
		V6.0 SP2032*
Vantage Titan	MRT-1510 MRT-1504	V6.0 SP1080* - V6.0 SP1085*

Steps to Apply the Security Patch

Please follow the instruction of applying patches defined by each of the product Security Patch Package. For details, please contact the customer service for each system.

Verified Security Patch List

KB	MSRC	explanation(en-us)	URL/File Name	url
KB2565063_x86	MS11-025	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package (KB2565063)	vc redistrib_x86.exe	https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/07/vcredist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe
KB2565063_x64	MS11-025	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package (KB2565063)	vc redistrib_x64.exe	https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/07/vcredist_x64_15d032d669078aa6f07fd1cbf4115a070bd034d.exe
VS2013_x86	n/a	Visual C++ Redistributable Packages for Visual Studio 2013	vc redistrib_x86.exe	https://download.visualstudio.microsoft.com/download/pr/10912115/88b8f4396e3edaa58c1350db34cdd751/vcredist_x86.exe

KB	MSRC	explanation(en-us)	URL/File Name	url
VS2013_x64	n/a	Visual C++ Redistributable Packages for Visual Studio 2013	vc redistrib_x64.exe	https://download.visualstudio.microsoft.com/download/pr/10912049/3831217601a8a928c5119d826bafd8f7/vcredist_x64.exe
VC2015-2022x86	n/a	Microsoft Visual C++ 2015-2022 Redistributable (x86) - 14.31.31103	VC_redist.x86.exe	https://download.visualstudio.microsoft.com/download/pr/8e32d7eb-5130-4dc8-9c3e-5891f375e112/B7AE307237F869E09F7413691A2CD1944357B5CEE28049C0A0D3430B47BB3EDC/VC_redist.x86.exe
VC2015-2022x64	n/a	Microsoft Visual C++ 2015-2022 Redistributable (x64) - 14.31.31103	VC_redist.x64.exe	https://download.visualstudio.microsoft.com/download/pr/d22ecb93-6eab-4ce1-89f3-97a816c55f04/37ED59A66699C0E5A7EBEEF7352D7C1C2ED5EDE7212950A1B0A8EE289AF4A95B/VC_redist.x64.exe
VC2015-2019x86	n/a	Microsoft Visual C++ 2015-2019 Redistributable (x86)	VC_redist.x86.exe	https://download.visualstudio.microsoft.com/download/pr/8ecb9800-52fd-432d-83ee-d6e037e96cc2/50A3E92ADE4C2D8F310A2812D46322459104039B9DEADB7FDD483B5C697C0C8/VC_redist.x86.exe
VC2015-2019x64	n/a	Microsoft Visual C++ 2015-2019 Redistributable (x64)	VC_redist.x64.exe	https://download.visualstudio.microsoft.com/download/pr/89a3b9df-4a09-492e-8474-8f92c115c51d/B1A32C71A6B7D5978904FB223763263EA5A7EB23B2C44A0D60E90D234AD99178/VC_redist.x64.exe
KB5005112	n/a	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems (KB5005112)	windows10.0-kb5005112-x64_81d09dc6978520e1a6d44b3b15567667f83eba2c.msu	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2021/08/windows10.0-kb5005112-x64_81d09dc6978520e1a6d44b3b15567667f83eba2c.msu
ASP.NETCore_6.0.27	n/a	2024-02 ASP.NET Core 6.0 Runtime (v6.0.27) - Windows x64 Installer	aspnetcore-runtime-6.0.27-win-x64.exe	https://download.visualstudio.microsoft.com/download/pr/856b04b7-f893-4fb1-9205-052413fde70f/09996e15acebe136113a3aa77b28fe5e/aspnetcore-runtime-6.0.27-win-x64.exe
DotNetRuntime_6.0.35	n/a	2024-10 .NET 6.0 Runtime (v6.0.35) - Windows x64 Installer	dotnet-runtime-6.0.35-win-x64.exe	https://download.visualstudio.microsoft.com/download/pr/c4f65621-b36b-46a9-8380-d5b660bef27e/0185fd72055dcdca86166b99add71686/dotnet-runtime-6.0.35-win-x64.exe
KB5021953_dotnet-runtime-3.1.32	n/a	2022-12 .NET Core 3.1.32 Security Update for x64 Server (KB5021953)	dotnet-runtime-3.1.32-win-x64.exe	https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2022/12/dotnet-runtime-3.1.32-win-x64_b319156006fcf5a0bc0cadfd536cbafa99b219b2.exe
KB5021953_aspnetcore-runtime-3.1.32	n/a	2022-12 .NET Core 3.1.32 Security Update for x64 Server (KB5021953)	aspnetcore-runtime-3.1.32-win-x64.exe	https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2022/12/aspnetcore-runtime-3.1.32-win-x64_d89589de054d33694272729bc5d713c02f2cc6af.exe

KB	MSRC	explanation(en-us)	URL/File Name	url
KB890830_v5.138	n/a	Windows Malicious Software Removal Tool x64 - v5.138 (KB890830)	windows-kb890830-x64-v5.138.exe	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/uprl/2025/12/windows-kb890830-x64-v5.138_dc0e0b97d146a92303d881bf9e2912299cbbbac5.exe
KB5066143	n/a	2025-10 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64 (KB5066143)	windows10.0-kb5066143-x64_461230ac1884e9691f54936b91ff04e46b80c2bc.msu	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2025/09/windows10.0-kb5066143-x64_461230ac1884e9691f54936b91ff04e46b80c2bc.msu
KB5074975	n/a	2025-12 Cumulative Update for Windows 10 Version 1809 for x64-based Systems (KB5074975)	windows10.0-kb5074975-x64_acfdfe5d5ea12c7ec55ca3cb48730c84748c6dbe.msu	https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/updt/2025/12/windows10.0-kb5074975-x64_acfdfe5d5ea12c7ec55ca3cb48730c84748c6dbe.msu
KB5019181	n/a	2023-03 Update for Windows 10 Version 1809 for x64-based Systems (KB5019181)	windows10.0-kb5019181-x64_01173052877dc5762f73069a958c108afd1125be.msu	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/updt/2023/02/windows10.0-kb5019181-x64_01173052877dc5762f73069a958c108afd1125be.msu
ASP.NETCore_8.0.21	n/a	2025-10 ASP.NET Core 8.0 Runtime (v8.0.21) - Windows x64 Installer	aspnetcore-runtime-8.0.21-win-x64.exe	https://builds.dotnet.microsoft.com/dotnet/aspnetcore/Runtime/8.0.21/aspnetcore-runtime-8.0.21-win-x64.exe
DotNetRuntime_8.0.24	n/a	2026-02 .NET 8.0 Runtime (v8.0.24) - Windows x64 Installer	dotnet-runtime-8.0.24-win-x64.exe	https://builds.dotnet.microsoft.com/dotnet/Runtime/8.0.24/dotnet-runtime-8.0.24-win-x64.exe



LEGAL MANUFACTURER

CANON MEDICAL SYSTEMS CORPORATION

1385, SHIMOISHIGAMI, OTAWARA-SHI, TOCHIGI 324-8550, JAPAN