

**PRODUCT SECURITY PATCH LIST  
FOR  
X-RAY CT SCANNER  
(2G984-0157EN)**

Canon Inc.

## Trademarks

Windows® is a registered trademark of Microsoft Corporation in the United States and other countries.

This document may also include trademarks or registered trademarks of other companies.

Note that the trademark symbol "™" and the registered trademark symbol "®" may or may not be used in this document.

### **IMPORTANT!**

1. No part of this document may be copied or reprinted, in whole or in part, without written permission.
2. The contents of this document are subject to change without prior notice and without legal obligation.
3. The contents of this document are correct to the best of our knowledge. Please inform us of any ambiguous or erroneous descriptions, missing information, etc.

<b>REVISION RECORD</b>						
REV.	DATE (YYYY-MM)	REASON	/AUTHOR	PAGE CHANGED	SER. No.	DOC. PRODUCT.
INI.	2026-05		Abe	-----		BU

Security Patch Package name

SPZX76-20103\_202603CT\_hotfix\_Server2016\_x64\_enu

Applicable Products for the Security Patch

Aquilion Prime SP V10.2, V10.4, V10.17, Base V16.2SP0001 is not applicable	TSX-303B
Aquilion Lightning V10.2, V10.4, V10.17, Base V16.2SP0001 is not applicable	TSX-036A
Aquilion ONE V10.0, V10.4, V10.11, V10.12, V10.16, Base V16.5SP0001 is not applicable	TSX-305A/ TSX-306A
Aquilion Exceed LB V10.9, V10.18	TSX-202A
Aquilion LB V6.6, V6.7	TSX-201A

Steps to Apply the Security Patch

Please follow the instruction of applying patches defined by each of the product Security Patch Package. For details, please contact the customer service for each system.

Verified Security Patch List

KB	MSRC	explanation(en-us)	URL/File Name	url
KB2538242_x86	MS11-025	Microsoft Visual C++ 2005 Service Pack 1 Redistributable Package MFC Security Update	vc redistrib_x86.exe	<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/06/vcredist_x86_b8fab0bb7f62a24ddf77b19cd9a1451abd7b847.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/06/vcredist_x86_b8fab0bb7f62a24ddf77b19cd9a1451abd7b847.exe</a>
KB2538242_x64	MS11-025	Microsoft Visual C++ 2005 Service Pack 1 Redistributable Package MFC Security Update	vc redistrib_x64.exe	<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/06/vcredist_x64_ee916012783024dac67fc606457377932c826f05.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/06/vcredist_x64_ee916012783024dac67fc606457377932c826f05.exe</a>
KB2538243_x86	MS11-025	Security Update for Microsoft Visual C++ 2008 Service Pack 1 Redistributable Package (KB2538243)	vc redistrib_x86.exe	<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/05/vcredist_x86_470640aa4bb7db8e69196b5edb0010933569e98d.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/05/vcredist_x86_470640aa4bb7db8e69196b5edb0010933569e98d.exe</a>
KB2538243_x64	MS11-025	Security Update for Microsoft Visual C++ 2008 Service Pack 1 Redistributable Package (KB2538243)	vc redistrib_x64.exe	<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/05/vcredist_x64_a7c83077b8a28d409e36316d2d7321fa0ccdb7e8.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/05/vcredist_x64_a7c83077b8a28d409e36316d2d7321fa0ccdb7e8.exe</a>
KB2565063_x86	MS11-025	Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package MFC Security Update	vc redistrib_x86.exe	<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/07/vcredist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secur/2011/07/vcredist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe</a>

KB	MSRC	explanation(en-us)	URL/File Name	url
KB2565063_x64	MS11-025	Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package MFC Security Update	vc_redist_x64.exe	<a href="https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcredist_x64_15d032d669078aa6f0f7fd1cbf4115a070bd034d.exe">https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcredist_x64_15d032d669078aa6f0f7fd1cbf4115a070bd034d.exe</a>
VS2017_x86	n/a	Microsoft Visual C++ Visual Studio 2017	vc_redist.x86.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/12319034/ccd261eb0e095411af3b306273231b68/VC_redist.x86.exe">https://download.visualstudio.microsoft.com/download/pr/12319034/ccd261eb0e095411af3b306273231b68/VC_redist.x86.exe</a>
VS2017_x64	n/a	Microsoft Visual C++ Visual Studio 2017	vc_redist.x64.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/11100230/15ccb3f02745c7b206ad10373cbca89b/VC_redist.x64.exe">https://download.visualstudio.microsoft.com/download/pr/11100230/15ccb3f02745c7b206ad10373cbca89b/VC_redist.x64.exe</a>
KB5075902	n/a	2026-02 Servicing Stack Update for Windows Server 2016 for x64-based Systems (KB5075902)	windows10.0-kb5075902-x64_a840f48f81e79f45131128f9f823633a59f57879.msu	<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2026/02/windows10.0-kb5075902-x64_a840f48f81e79f45131128f9f823633a59f57879.msu">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2026/02/windows10.0-kb5075902-x64_a840f48f81e79f45131128f9f823633a59f57879.msu</a>
KB5078938	n/a	2026-03 Cumulative Update for Windows Server 2016 for x64-based Systems (KB5078938)	windows10.0-kb5078938-x64_4f2644d71dfb3459b6777dfb1a3a4e16cf3877f5.msu	<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2026/03/windows10.0-kb5078938-x64_4f2644d71dfb3459b6777dfb1a3a4e16cf3877f5.msu">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2026/03/windows10.0-kb5078938-x64_4f2644d71dfb3459b6777dfb1a3a4e16cf3877f5.msu</a>
KB890830_v5.139	n/a	Windows Malicious Software Removal Tool x64 - v5.139 (KB890830)	windows-kb890830-x64-v5.139.exe	<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/uprl/2026/02/windows-kb890830-x64-v5.139_c6a3e646b705898985bbd3b88ab2f672b84a8784.exe">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/uprl/2026/02/windows-kb890830-x64-v5.139_c6a3e646b705898985bbd3b88ab2f672b84a8784.exe</a>
KB5019182	n/a	2023-04 Update for Windows Server 2016 for x64-based Systems (KB5019182)	windows10.0-kb5019182-x64_9dda0f01f337746913602d2c7269f5aa9cea9235.msu	<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/updt/2023/04/windows10.0-kb5019182-x64_9dda0f01f337746913602d2c7269f5aa9cea9235.msu">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/updt/2023/04/windows10.0-kb5019182-x64_9dda0f01f337746913602d2c7269f5aa9cea9235.msu</a>



LEGAL MANUFACTURER

**Canon Inc.**

---

30-2, SHIMOMARUKO 3-CHOME, OHTA-KU, TOKYO 146-8501, JAPAN