

**PRODUCT SECURITY PATCH LIST  
FOR  
MRI SCANNER  
(2G984-0166EN)**

**CANON INC.**

## Trademarks

Windows® is a registered trademark of Microsoft Corporation in the United States and other countries.

This document may also include trademarks or registered trademarks of other companies.

Note that the trademark symbol "™" and the registered trademark symbol "®" may or may not be used in this document.

### **IMPORTANT!**

1. No part of this document may be copied or reprinted, in whole or in part, without written permission.
2. The contents of this document are subject to change without prior notice and without legal obligation.
3. The contents of this document are correct to the best of our knowledge. Please inform us of any ambiguous or erroneous descriptions, missing information, etc.

<b>REVISION RECORD</b>						
REV.	DATE (YYYY-MM)	REASON	/AUTHOR	PAGE CHANGED	SER. No.	DOC. PRODUCT.
INI.	2026-06		Kikuchi	-----		BU

Security Patch Package name

SPZM45-20088\_202605MR\_hotfix\_win10\_LTSC2019\_x64\_enu

Applicable Products for the Security Patch

Vantage Galan 3T	MRT-3020	V6.0 SP2015*
		V7.0 SP0070*
		V8.0 SP0040* - V8.0 SP0055*
		V9.0 SP0020* - V9.0 SP0022*
		V9.0 SP0031* - V9.0 SP0055*
Vantage Orian	MRT-1550	V6.0 SP2015* - V6.0 SP2025*
		V7.0 SP0070* - V7.0 SP0075*
		V8.0 SP0040* - V8.0 SP0050*
		V9.0 SP0020* - V9.0 SP0022*
		V9.0 SP0030* - V9.0 SP0055*
Vantage Elan	MRT-2020	V6.0 SP1080* - V6.0 SP1085*
		V8.1 SP0000* - V8.1 SP0032*
		V8.1 SP0050* - V8.1 SP0055*
Vantage Titan 3T	MRT-3010	V6.0 SP2032*
Vantage Titan	MRT-1510	V6.0 SP1080* - V6.0 SP1085*
	MRT-1504	

Steps to Apply the Security Patch

Please follow the instruction of applying patches defined by each of the product Security Patch Package. For details, please contact the customer service for each system.

Verified Security Patch List

KB	MSRC	explanation(en-us)	URL/File Name	url
KB2565063_x86	MS11-025	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package (KB2565063)	vc redistrib_x86.exe	<a href="https://catalog.s.download.windowupdate.com/msdownload/update/software/secu/2011/07/vcredist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe">https://catalog.s.download.windowupdate.com/msdownload/update/software/secu/2011/07/vcredist_x86_28c54491be70c38c97849c3d8cfbfd0d3c515cb.exe</a>
KB2565063_x64	MS11-025	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package (KB2565063)	vc redistrib_x64.exe	<a href="https://catalog.s.download.windowupdate.com/msdownload/update/software/secu/2011/07/vcredist_x64_15d032d669078aa6f0f7fd1cbf4115a070bd034d.exe">https://catalog.s.download.windowupdate.com/msdownload/update/software/secu/2011/07/vcredist_x64_15d032d669078aa6f0f7fd1cbf4115a070bd034d.exe</a>
VS2013_x86	n/a	Visual C++ Redistributable Packages for Visual Studio 2013	vc redistrib_x86.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/10912115/88b8f4396e3edaa58c1350db34cdd751/vcredist_x86.exe">https://download.visualstudio.microsoft.com/download/pr/10912115/88b8f4396e3edaa58c1350db34cdd751/vcredist_x86.exe</a>

KB	MSRC	explanation(en-us)	URL/File Name	url
VS2013_x64	n/a	Visual C++ Redistributable Packages for Visual Studio 2013	vc_redist_x64.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/10912049/3831217601a8a928c5119d826bafd8f7/vc_redist_x64.exe">https://download.visualstudio.microsoft.com/download/pr/10912049/3831217601a8a928c5119d826bafd8f7/vc_redist_x64.exe</a>
VC2015-2022x86	n/a	Microsoft Visual C++ 2015-2022 Redistributable (x86) - 14.31.31103	VC_redist.x86.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/8e32d7eb-5130-4dc8-9c3e-5891f375e112/B7AE307237F869E09F7413691A2CD1944357B5CEE28049C0A0D3430B47BB3EDC/VC_redist.x86.exe">https://download.visualstudio.microsoft.com/download/pr/8e32d7eb-5130-4dc8-9c3e-5891f375e112/B7AE307237F869E09F7413691A2CD1944357B5CEE28049C0A0D3430B47BB3EDC/VC_redist.x86.exe</a>
VC2015-2022x64	n/a	Microsoft Visual C++ 2015-2022 Redistributable (x64) - 14.31.31103	VC_redist.x64.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/d22ecb93-6eab-4ce1-89f3-97a816c55f04/37ED59A66699C0E5A7EBEEF7352D7C1C2ED5EDE7212950A1B0A8EE289AF4A95B/VC_redist.x64.exe">https://download.visualstudio.microsoft.com/download/pr/d22ecb93-6eab-4ce1-89f3-97a816c55f04/37ED59A66699C0E5A7EBEEF7352D7C1C2ED5EDE7212950A1B0A8EE289AF4A95B/VC_redist.x64.exe</a>
VC2015-2019x86	n/a	Microsoft Visual C++ 2015-2019 Redistributable (x86)	VC_redist.x86.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/8ecb9800-52fd-432d-83ee-d6e037e96cc2/50A3E92ADE4C2D8F310A2812D46322459104039B9DEADB7FDD483B5C697C0C8/VC_redist.x86.exe">https://download.visualstudio.microsoft.com/download/pr/8ecb9800-52fd-432d-83ee-d6e037e96cc2/50A3E92ADE4C2D8F310A2812D46322459104039B9DEADB7FDD483B5C697C0C8/VC_redist.x86.exe</a>
VC2015-2019x64	n/a	Microsoft Visual C++ 2015-2019 Redistributable (x64)	VC_redist.x64.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/89a3b9df-4a09-492e-8474-8f92c115c51d/B1A32C71A6B7D5978904FB223763263EA5A7EB23B2C44A0D60E90D234AD99178/VC_redist.x64.exe">https://download.visualstudio.microsoft.com/download/pr/89a3b9df-4a09-492e-8474-8f92c115c51d/B1A32C71A6B7D5978904FB223763263EA5A7EB23B2C44A0D60E90D234AD99178/VC_redist.x64.exe</a>
KB5005112	n/a	2021-08 Servicing Stack Update for Windows 10 Version 1809 for x64-based Systems (KB5005112)	windows10.0-kb5005112-x64_81d09dc6978520e1a6d44b3b15567667f83eba2c.msu	<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2021/08/windows10.0-kb5005112-x64_81d09dc6978520e1a6d44b3b15567667f83eba2c.msu">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2021/08/windows10.0-kb5005112-x64_81d09dc6978520e1a6d44b3b15567667f83eba2c.msu</a>
ASP.NETCore_6.0.27	n/a	2024-02 ASP.NET Core 6.0 Runtime (v6.0.27) - Windows x64 Installer	aspnetcore-runtime-6.0.27-win-x64.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/856b04b7-f893-4fb1-9205-052413fde70f/099996e15acebe136113a3aa77b28fe5e/aspnetcore-runtime-6.0.27-win-x64.exe">https://download.visualstudio.microsoft.com/download/pr/856b04b7-f893-4fb1-9205-052413fde70f/099996e15acebe136113a3aa77b28fe5e/aspnetcore-runtime-6.0.27-win-x64.exe</a>
DotNetRuntime_6.0.35	n/a	2024-10 .NET 6.0 Runtime (v6.0.35) - Windows x64 Installer	dotnet-runtime-6.0.35-win-x64.exe	<a href="https://download.visualstudio.microsoft.com/download/pr/c4f65621-b36b-46a9-8380-d5b660bef27e/0185fd72055dcdca86166b99add71686/dotnet-runtime-6.0.35-win-x64.exe">https://download.visualstudio.microsoft.com/download/pr/c4f65621-b36b-46a9-8380-d5b660bef27e/0185fd72055dcdca86166b99add71686/dotnet-runtime-6.0.35-win-x64.exe</a>
KB5021953_dotnet-runtime-3.1.32	n/a	2022-12 .NET Core 3.1.32 Security Update for x64 Server (KB5021953)	dotnet-runtime-3.1.32-win-x64.exe	<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2022/12/dotnet-runtime-3.1.32-win-x64_b319156006fcf5a0bc0cadfd536cbafa99b219b2.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2022/12/dotnet-runtime-3.1.32-win-x64_b319156006fcf5a0bc0cadfd536cbafa99b219b2.exe</a>
KB5021953_aspnetcore-runtime-3.1.32	n/a	2022-12 .NET Core 3.1.32 Security Update for x64 Server (KB5021953)	aspnetcore-runtime-3.1.32-win-x64.exe	<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2022/12/aspnetcore-runtime-3.1.32-win-x64_d89589de054d33694272729bc5d713c02f2cc6af.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2022/12/aspnetcore-runtime-3.1.32-win-x64_d89589de054d33694272729bc5d713c02f2cc6af.exe</a>

KB	MSRC	explanation(en-us)	URL/File Name	url
KB890830_v5.141	n/a	Windows Malicious Software Removal Tool x64 - v5.141 (KB890830)	windows-kb890830-x64-v5.141.exe	<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/upr1/2026/05/windows-kb890830-x64-v5.141_28e375ae266a48518b85561b124c6d21b82d84a6.exe">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/upr1/2026/05/windows-kb890830-x64-v5.141_28e375ae266a48518b85561b124c6d21b82d84a6.exe</a>
KB5087061	n/a	2026-05 Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10 Version 1809 for x64 (KB5087061)	windows10.0-kb5087061-x64_25b4ce6e32862a6288609b4471ff1cd8e07d8e3f.msu	<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2026/04/windows10.0-kb5087061-x64_25b4ce6e32862a6288609b4471ff1cd8e07d8e3f.msu">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/secu/2026/04/windows10.0-kb5087061-x64_25b4ce6e32862a6288609b4471ff1cd8e07d8e3f.msu</a>
KB5087538	n/a	2026-05 Cumulative Update for Windows 10 Version 1809 for x64-based Systems (KB5087538)	windows10.0-kb5087538-x64_86844c656b842c1c8bf11e6153e4edfb4e914853.msu	<a href="https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2026/05/windows10.0-kb5087538-x64_86844c656b842c1c8bf11e6153e4edfb4e914853.msu">https://catalog.s.download.windowsupdate.com/c/msdownload/update/software/secu/2026/05/windows10.0-kb5087538-x64_86844c656b842c1c8bf11e6153e4edfb4e914853.msu</a>
KB5019181	n/a	2023-03 Update for Windows 10 Version 1809 for x64-based Systems (KB5019181)	windows10.0-kb5019181-x64_01173052877dc5762f73069a958c108afd1125be.msu	<a href="https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/updt/2023/02/windows10.0-kb5019181-x64_01173052877dc5762f73069a958c108afd1125be.msu">https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/updt/2023/02/windows10.0-kb5019181-x64_01173052877dc5762f73069a958c108afd1125be.msu</a>
ASP.NETCore_8.0.25	n/a	2026-03 ASP.NET Core 8.0 Runtime (v8.0.25) - Windows x64 Installer	aspnetcore-runtime-8.0.25-win-x64.exe	<a href="https://builds.dotnet.microsoft.com/dotnet/aspnetcore/Runtime/8.0.25/aspnetcore-runtime-8.0.25-win-x64.exe">https://builds.dotnet.microsoft.com/dotnet/aspnetcore/Runtime/8.0.25/aspnetcore-runtime-8.0.25-win-x64.exe</a>
DotNetRuntime_8.0.26	n/a	2026-04 .NET 8.0 Runtime (v8.0.26) - Windows x64 Installer	dotnet-runtime-8.0.26-win-x64.exe	<a href="https://builds.dotnet.microsoft.com/dotnet/Runtime/8.0.26/dotnet-runtime-8.0.26-win-x64.exe">https://builds.dotnet.microsoft.com/dotnet/Runtime/8.0.26/dotnet-runtime-8.0.26-win-x64.exe</a>



LEGAL MANUFACTURER

**CANON INC.**

---

30-2, SHIMOMARUKO 3-CHOME, OHTA-KU, TOKYO 146-8501, JAPAN