

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	RapideyeCore™	V3.1 SP0000J	11 Oct. 2016
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Technical Framework (Rev 15.0)			
Access to Radiology Information	Image Manager/Archive	None	
Scheduled Workflow(SWF.b)	Image Manager/Archive	None	
Consistent Presentation of Images	Image Manager/Archive	None	
Mammography Image	Image Manager/Archive	None	
Evidence Documents	Image Manager/Archive	None	
Radiation Exposure Monitoring	Image Manager/Archive	None	
Key Image Note	Image Manager/Archive	None	
Import Reconciliation Workflow	Image Manager/Archive	None	
Patient Information Reconciliation	Image Manager/Archive	None	
IT Infrastructure IHE IT Infrastructure Technical Framework (Rev 13.0)			
Patient Administration Management	Patient Encounter Consumer	None	
<u>Internet address for vendor's IHE information:</u>			
Japan: https://www.medical.canon/Interoperability/IHE/JP			
Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

“ RapideyeCore“ is a trademark of Canon Medical Systems Corporation.