

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	Ultimax-i, ZEXIRA, KALARE, Plessart ZERO, Plessart VIVO, Plessart EX8 / HDR-08A	V1.31.01 or later	16 Jan.2017
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 15.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
<u>Internet address for vendor's IHE information:</u>			
Japan: https://www.medical.canon/Interoperability/IHE/JP			
Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

"Ultimax", "ZEXIRA", "KALARE", "Plessart" are trademarks of Canon Medical Systems Corporation.