

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	Aquilion™ Aquilion ONE™ Astelion™ Alexion™ TSX-101A,TSX-301A, TSX-301B,TSX-301C, TSX-302A, TSX-032A,TSX-033A, TSX-303A, TSX-034A,TSX-201A, TSX-035A, TSX-305A, TSX-036A, TSX-303B	V7.0 or later	May 2017
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 15.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Charge Posting	Acquisition Modality	PPS Exception Management Modality Group Case	
Presentation of Grouped Procedures	Acquisition Modality	No options defined	
Evidence Documents	Acquisition Modality	TID 3900 Cardiovascular Analysis Report Calcium Score Evidence	
Portable Data for Imaging	Portable Media Creator	None	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Cardiology IHE Cardiology Technical Framework (Rev 5.0)			
Evidence Documents	Acquisition Modality	TID 3900 Cardiovascular Analysis Report Calcium Score Evidence	
<u>Internet address for vendor's IHE information:</u>			
Japan: https://www.medical.canon/Interoperability/IHE/JP			
Global: https://www.medical.canon/Interoperability/IHE/EN			

Links to Standards Conformance Statements for the Implementation		
HL7	None	
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org

“Aquilion”, “ Aquilion ONE”, “Astelion” and “Alexion” are trademarks of Canon Medical Systems Corporation.