

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	Vantage Orian™ MRT-1550 Vantage Galan™ 3T MRT-3020	V4.5 and V5.0	14 Mar 2018
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 11.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Modality Group Case Japanese Character Set	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	DVD Media * No CD-R but DVD-R support * No Compression only	
IT Infrastructure IHE IT Infrastructure Technical Framework (Rev 11.0)			
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: Japan: https://www.medical.canon/Interoperability/IHE/JP Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net/	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org/	

"Vantage Orian" and "Vantage Galan" are trademarks of Canon Medical Systems Corporation.