

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Systems Corporation	Vantage Orian™ MRT-1550 Vantage Galan™ 3T MRT-3020	V4.5, V5.0 and V6.0	April 2019
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 11.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Modality Group Case Japanese Character Set	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	DVD Media * No CD-R but DVD-R support * No Compression only	
<b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 11.0)			
Consistent Time	Time Client	None	
<u>Internet address for vendor's IHE information:</u> Japan: <a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a> Global: <a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a> Global: <a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>	

"Vantage Orian" and "Vantage Galan" are trademarks of Canon Medical Systems Corporation.