

IHE Integration Statement			
Vendor	Product Name	Version	Date
CANON INC.	Viamo sv7 CUS-VSV7	V1.5	Jun. 2026
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 16.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Assisted Acquisition Protocol Setting Modality Group Case Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Cardiology IHE Cardiology Technical Framework (Rev 5.0)			
Echocardiography Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management	
IT Infrastructure IHE IT Infrastructure Technical Framework (Rev 14.0)			
Consistent Time	Time Client	None	
<u>Internet address for vendor’s IHE information:</u> Japan: https://www.medical.canon/Interoperability/IHE/JP Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

HL7® is the registered trademark of Health Level Seven International and the use of this trademark does not constitute an endorsement by HL7.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.