

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Systems Corporation	Alphenix™ series / DFP-8000	V8.0 or later	June 2018,
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 16)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
<b>Cardiology</b> IHE Cardiology Technical Framework (Rev 2.1)			
Cardiac Catheterization Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
<b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 3.0)			
Consistent Time	Time Client	None	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a>			
Global: <a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a> Global: <a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	

“ Alphenix“ is a trademark of Canon Medical Systems Corporation.