

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Systems Corporation	Aquilion Start™ TSX-037A	V9.0 or later	Mar, 2019
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 15.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Charge Posting	Acquisition Modality	PPS Exception Management Modality Group Case	
Presentation of Grouped Procedures	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a>			
Global: <a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a> Global: <a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	

"Aquilion Start" is a trademark of Canon Medical Systems Corporation.