

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	Vantage Elan MRT-2020 Vantage Titan MRT-1504,MRT-1510, MRT-2004, MRT-3010 Vantage Orian MRT-1550 Vantage Galan MRT-3020	V6.0 and V7.0	27 February 2023
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 11.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Modality Group Case Japanese Character Set	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	DVD Media	
IT Infrastructure IHE IT Infrastructure Technical Framework (Rev 11.0)			
Consistent Time	Time Client	None	
<u>Internet address for vendor's IHE information:</u> Japan: https://www.medical.canon/Interoperability/IHE/JP Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net/	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org/	

HL7® is the registered trademark of Health Level Seven International and the use of this trademark does not constitute an endorsement by HL7.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.