

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Systems Corporation	UltraExtend™ NX CUW-U001S	V2.0	1 December 2022
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 16.0)			
Scheduled Workflow	Image Display	None	
Patient Information Reconciliation	Image Display	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Evidence Documents	Image Display	No options defined TID 5000 OB-GYN Ultrasound Procedure Report	
<b>Cardiology</b> IHE Cardiology Technical Framework (Rev 5.0)			
Echocardiography Workflow	Image Display	None	
Evidence Documents	Image Display	Echo Evidence TID 5200 Echocardiography Procedure Report	
Stress Testing Workflow	Image Display	None	
<b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 17.0)			
Audit Trail and Node Authentication	Secure Application Secure Node	None	
Consistent Time	Time Client	None	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a>			
Global: <a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a> Global: <a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	

"UltraExtend" is a trademark of Canon Medical Systems Corporation.