

IHE Integration Statement			
Vendor	Product Name	Version	Date
CANON INC.	UltraExtend NX CUW-U001S	V2.0	Jun. 2026
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework (Rev 16.0)			
Scheduled Workflow	Image Display	None	
Patient Information Reconciliation	Image Display	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Evidence Documents	Image Display	No options defined TID 5000 OB-GYN Ultrasound Procedure Report	
Cardiology IHE Cardiology Technical Framework (Rev 5.0)			
Echocardiography Workflow	Image Display	None	
Evidence Documents	Image Display	Echo Evidence TID 5200 Echocardiography Procedure Report	
Stress Testing Workflow	Image Display	None	
IT Infrastructure IHE IT Infrastructure Technical Framework (Rev 17.0)			
Audit Trail and Node Authentication	Secure Application Secure Node	None	
Consistent Time	Time Client	None	
<u>Internet address for vendor’s IHE information:</u> Japan: https://www.medical.canon/Interoperability/IHE/JP Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

HL7® is the registered trademark of Health Level Seven International and the use of this trademark does not constitute an endorsement by HL7.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.