

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Systems Corporation	Aplio™ <i>flex</i> CUS-AFL00 Aplio™ <i>go</i> CUS-AGG00	V1.0	30 Jan. 2022
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 16.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Assisted Acquisition Protocol Setting Modality Group Case Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Evidence Documents	Acquisition Modality	No options defined TID 5000 OB-GYN Ultrasound Procedure Report	
<b>Cardiology</b> IHE Cardiology Technical Framework (Rev 5.0)			
Echocardiography Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management	
Evidence Documents	Acquisition Modality	Echo Evidence TID 5100 Vascular Ultrasound Report TID 5200 Echocardiography Procedure Report	
Stress Testing Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management	
<b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 17.0)			
Audit Trail and Node Authentication	Secure Application Secure Node	Radiology Audit Trail STX: BCP195 TLS Secure Transport Connection -TLS 1.2 Floor ATX: TLS... FQDN Validation of Server Certificate	
Consistent Time	Time Client	None	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a>			
Global: <a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a>			

<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	None	
<b>DICOM</b>	Japan: <a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a> Global: <a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>

“Aplio” is a trademark of Canon Medical Systems Corporation.