

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Canon Medical Systems Corporation	Aquilion™ TSX-307A	V1.2 or later	16 June, 2022
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b> IHE Radiology Technical Framework (Rev 15.0)			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Charge Posting	Acquisition Modality	PPS Exception Management Modality Group Case	
Presentation of Grouped Procedures	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
<b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework (Rev 17.0)			
Audit Trail and Node Authentication	Secure Application Secure Node	STX: TLS 1.2 Floor using BCP195 ATX: TLS Syslog ATX: UDP Syslog	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a>			
Global: <a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a> Global: <a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	

“Aquilion” is a trademark of Canon Medical Systems Corporation.