

| IHE Integration Statement                                          |                                   |                                                                                                                                                                            |              |
|--------------------------------------------------------------------|-----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| Vendor                                                             | Product Name                      | Version                                                                                                                                                                    | Date         |
| CANON INC.                                                         | Aplio <i>beyond</i> CUS-ABE00     | V1.0                                                                                                                                                                       | Jun.<br>2026 |
| Integration Profiles Implemented                                   | Actors Implemented                | Options Implemented                                                                                                                                                        |              |
| <b>Radiology</b> IHE Radiology Technical Framework                 |                                   |                                                                                                                                                                            |              |
| Scheduled Workflow                                                 | Acquisition Modality              | Patient Based Worklist Query<br>Broad Worklist Query<br>PPS Exception Management<br>Assisted Acquisition Protocol Setting<br>Modality Group Case<br>Japanese Character Set |              |
| Patient Information Reconciliation                                 | Acquisition Modality              | No options defined                                                                                                                                                         |              |
| Portable Data for Imaging                                          | Portable Media Creator            | None                                                                                                                                                                       |              |
| Evidence Documents                                                 | Acquisition Modality              | No options defined<br>TID 5000 OB-GYN<br>Ultrasound Procedure Report                                                                                                       |              |
| <b>Cardiology</b> IHE Cardiology Technical Framework               |                                   |                                                                                                                                                                            |              |
| Echocardiography Workflow                                          | Acquisition Modality              | Patient Based Worklist Query<br>Broad Worklist Query<br>PPS Exception Management<br>Stress Echo                                                                            |              |
| Evidence Documents                                                 | Acquisition Modality              | Echo Evidence<br>TID 5100 Vascular<br>Ultrasound Report<br>TID 5200 Echocardiography<br>Procedure Report<br>TID 3300 Stress Testing<br>Report                              |              |
| Stress Testing Workflow                                            | Acquisition Modality              | Patient Based Worklist Query<br>Broad Worklist Query<br>PPS Exception Management<br>Stress Echo                                                                            |              |
| <b>IT Infrastructure</b> IHE IT Infrastructure Technical Framework |                                   |                                                                                                                                                                            |              |
| Audit Trail and Node Authentication                                | Secure Application<br>Secure Node | Radiology Audit Trail<br>STX: BCP195 TLS Secure<br>Transport Connection -TLS<br>1.3 Floor<br>ATX: TLS...<br>FQDN Validation of Server<br>Certificate                       |              |
| Consistent Time                                                    | Time Client                       | None                                                                                                                                                                       |              |

|                                                                                                                              |                                                                                                                                                                                                                                                                     |                                                                      |
|------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------|
| <b><u>Internet address for vendor's IHE information:</u></b>                                                                 |                                                                                                                                                                                                                                                                     |                                                                      |
| Japan:<br><a href="https://www.medical.canon/Interoperability/IHE/JP">https://www.medical.canon/Interoperability/IHE/JP</a>  |                                                                                                                                                                                                                                                                     |                                                                      |
| Global:<br><a href="https://www.medical.canon/Interoperability/IHE/EN">https://www.medical.canon/Interoperability/IHE/EN</a> |                                                                                                                                                                                                                                                                     |                                                                      |
| <b>Links to Standards Conformance Statements for the Implementation</b>                                                      |                                                                                                                                                                                                                                                                     |                                                                      |
| <b>HL7</b>                                                                                                                   | None                                                                                                                                                                                                                                                                |                                                                      |
| <b>DICOM</b>                                                                                                                 | Japan:<br><a href="https://www.medical.canon/Interoperability/DICOM/JP">https://www.medical.canon/Interoperability/DICOM/JP</a><br>Global:<br><a href="https://www.medical.canon/Interoperability/DICOM/EN">https://www.medical.canon/Interoperability/DICOM/EN</a> |                                                                      |
| <b>Links to general information on IHE</b>                                                                                   |                                                                                                                                                                                                                                                                     |                                                                      |
| In North America:<br><a href="http://www.ihe.net">http://www.ihe.net</a>                                                     | In Europe:<br><a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>                                                                                                                                                                                     | In Japan:<br><a href="http://www.ihe-j.org">http://www.ihe-j.org</a> |

HL7® is the registered trademark of Health Level Seven International and the use of this trademark does not constitute an endorsement by HL7.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.