

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	DoseXross ERIS-5KDXA	V1.2 or later	01 Feb. 2024
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology IHE Radiology Technical Framework			
Radiation Exposure Monitoring (REM)	Dose Info Consumer	None	
IT Infrastructure IHE IT Infrastructure Technical Framework			
Consistent Time	Time Client	None	
<u>Internet address for vendor's IHE information:</u> Japan: https://www.medical.canon/Interoperability/IHE/JP Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: https://www.ihe.net	In Europe: https://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

HL7® is the registered trademark of Health Level Seven International and the use of this trademark does not constitute an endorsement by HL7.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.