

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	Ultimax-i DREX-UI80 (Model UDR-10A/U1)	V2.7	25.08.2025
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology	IHE Radiology Technical Framework		
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
IT Infrastructure	IHE IT Infrastructure Technical Framework		
Audit Trail and Node Authentication	Secure Node	None	
Consistent Time	Time Client	None	
<u>Internet address for vendor's IHE information:</u>			
Japan: https://www.medical.canon/Interoperability/IHE/JP			
Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: https://www.ihe.net	In Europe: https://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

HL7® is the registered trademark of Health Level Seven International and the use of this trademark does not constitute an endorsement by HL7.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its Standards publications relating to digital communications of medical information.