

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Systems Corporation	RADREX™-i series / DRAD-3000E	V3.60 or Later	6 Feb. 2012
This product implements all transactions required in the IHE Technical Framework (Rev 10.0) and IHE IT Infrastructure Technical Framework (Rev 7.0) to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting Billing and Material Management Broad Worklist Query Modality Group Case Patient Based Worklist Query PPS Exception Management	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
IT Infrastructure			
Consistent Time	Time Client	None	
<u>Internet address for vendor's IHE information:</u>			
Japan: https://www.medical.canon/Interoperability/IHE/JP			
Global: https://www.medical.canon/Interoperability/IHE/EN			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: https://www.medical.canon/Interoperability/DICOM/JP Global: https://www.medical.canon/Interoperability/DICOM/EN		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

“ RADREX“ is a trademark of Canon Medical Systems Corporation.