

IHE Integration Statement			
Vendor	Product Name	Version	Date
Toshiba Medical Systems Corporation	Aquilion™ 32 64 TSX-101A	V3.00	31 May.2006
This product implements all transactions required in the IHE Technical Framework (Rev 6.0) to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Charge Posting	Acquisition Modality	PPS Exception Management Modality Group Case	
Presentation of Grouped Procedures	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
<u>Internet address for vendor's IHE information:</u>			
Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe/ http://www.toshiba-medical.co.jp/tmd/english/ihe/ (International page)			
North America: http://medical.toshiba.com/services/iheintegrationstatements.htm			
Europe: http://www.toshiba-europe.com/medical/medical.asp?PageID=153			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom/ http://www.toshiba-medical.co.jp/tmd/english/dicom/ (International page) North America: http://medical.toshiba.com/services/dicomconformancstatements.htm Europe: http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532		
Links to general information on IHE			
In North America: http://www.ihe.net/	In Europe: http://www.ihe-europe.org/	In Japan: http://www.ihe-j.org/	

“ Aquilion“ is a trademark of Toshiba Medical Systems Corporation.