

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Toshiba Medical Systems Corporation	Aquilion™ ,Alexion™ , Aquilion ONE™ , Aquilion™ PRIME, Astelion™ TSX-101A, TSX-301A, TSX-301B, TSX-302A, TSX-032A, TSX-033A, TSX-034A,TSX-201A	V4.70 or Later	17 Oct.2013
This product implements all transactions required in the IHE Technical Framework (Rev 12.0), IHE Cardiology Technical Framework (Rev 5.0) to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b>			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Acquisition Modality	No options defined	
Charge Posting	Acquisition Modality	PPS Exception Management Modality Group Case	
Presentation of Grouped Procedures	Acquisition Modality	No options defined	
Evidence Documents	Acquisition Modality	No options defined System generates TID 3900 Cardiovascular Analysis Report Calcium Score Evidence	
Portable Data for Imaging	Portable Media Creator	None	
Radiation Exposure Monitoring	Acquisition Modality	No options defined	
<b>Cardiology</b>			
Evidence Documents	Acquisition Modality	No options defined System generates TID 3900 Cardiovascular Analysis Report Calcium Score Evidence	

<b>Internet address for vendor's IHE information:</b>		
Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/ihe">http://www.toshiba-medical.co.jp/tmd/products/ihe</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe">http://www.toshiba-medical.co.jp/tmd/english/ihe</a> (International page)		
North America: <a href="http://medical.toshiba.com/services/enterprise-integration/ihe-integration/index.php">http://medical.toshiba.com/services/enterprise-integration/ihe-integration/index.php</a>		
Europe: <a href="http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements">http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	None	
<b>DICOM</b>	Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/dicom">http://www.toshiba-medical.co.jp/tmd/products/dicom</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom">http://www.toshiba-medical.co.jp/tmd/english/dicom</a> (International page) North America: <a href="http://medical.toshiba.com/service/dicom/">http://medical.toshiba.com/service/dicom/</a> Europe: <a href="http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements">http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>

Aquilion™ and Alexion™, AquilionONE™ and Astelion™ is a trademark of Toshiba Medical Systems Corporation.