

| IHE Integration Statement | | | |
|---|--|---|-------------|
| Vendor | Product Name | Version | Date |
| Toshiba Medical Systems Corporation | Aquilion™ TSX-101A,TSX-201A Activion™ TSX-031A | V3.10 or later | 01 Jun.2011 |
| This product implements all transactions required in the IHE Technical Framework (Rev 10.0), IHE Cardiology Technical Framework (Rev 3.0) and IHE IT Infrastructure Technical Framework (Rev 7.0) to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Radiology | | | |
| Scheduled Workflow | Acquisition Modality | Patient Based Worklist Query Broad Worklist Query Assisted Acquisition Protocol Setting Modality Group Case PPS Exception Management Billing and Material Management Japanese Character Set | |
| Patient Information Reconciliation | Acquisition Modality | No options defined | |
| Consistent Presentation of Images | Acquisition Modality | No options defined | |
| Charge Posting | Acquisition Modality | PPS Exception Management Modality Group Case | |
| Presentation of Grouped Procedures | Acquisition Modality | No options defined | |
| Portable Data for Imaging | Portable Media Creator | No options defined | |
| <u>Internet address for vendor's IHE information:</u> | | | |
| Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe http://www.toshiba-medical.co.jp/tmd/english/ihe (International page) | | | |
| North America: http://medical.toshiba.com/service/ihe/index.php | | | |
| Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| HL7 | None | | |
| DICOM | Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom http://www.toshiba-medical.co.jp/tmd/english/dicom (International page) North America: http://medical.toshiba.com/service/dicom/index.php Europe: http://www.toshiba-medical.eu/en/Our-Product-Range/CT/Compliance-Statements | | |
| Links to general information on IHE | | | |
| In North America: http://www.ihe.net | In Europe: http://www.ihe-europe.net | In Japan: http://www.ihe-j.org | |

“Aquilion” and “Activion” is a trademark of Toshiba Medical Systems Corporation.