

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Toshiba Medical Systems Corporation	Vantage Titan™ 3T MRT-3010 Vantage Titan™ MRT-1504, MRT-2004 EXCELART Vantage™ Atlas MRT-1503,MRT-2003	V1.35 or later	4 Apr. 2015
This product implements all transactions required in the IHE Technical Framework (Rev 11.0), IHE Cardiology Technical Framework (Rev 4.0) and IHE IT Infrastructure Technical Framework (Rev 9.0) to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b>			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Modality Group Case Japanese Character Set	
Consistent Presentation of Images	Acquisition Modality	No options defined	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/ihe/">http://www.toshiba-medical.co.jp/tmd/products/ihe/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe/">http://www.toshiba-medical.co.jp/tmd/english/ihe/</a> (International page)			
North America: <a href="http://medical.toshiba.com/services/enterprise-integration/ihe-integration/index.php">http://medical.toshiba.com/services/enterprise-integration/ihe-integration/index.php</a>			
Europe: <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe/">http://www.toshiba-medical.co.jp/tmd/english/ihe/</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/dicom/">http://www.toshiba-medical.co.jp/tmd/products/dicom/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom/">http://www.toshiba-medical.co.jp/tmd/english/dicom/</a> (International page)		
	North America: <a href="http://medical.toshiba.com/services/enterprise-integration/dicom-conformance/index.php">http://medical.toshiba.com/services/enterprise-integration/dicom-conformance/index.php</a>		
	Europe: <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom/">http://www.toshiba-medical.co.jp/tmd/english/dicom/</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>	

"Vantage Titan" and "EXCELART Vantage" are trademark of Toshiba Medical Systems Corporation.