

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Toshiba Medical Systems Corporation	Aplio™ / SSA-700A/770A	V7.00	1 <sup>st</sup> .Jun.2006
This product implements all transactions required in the IHE Technical Framework (Rev 6.0), and IHE Cardiology Technical Framework (Rev 2.0) to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b>			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
<b>Cardiology</b>			
Evidence Documents	Acquisition Modality	Echo Evidence	
<b><u>Internet address for vendor's IHE information:</u></b>			
Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/ihe/">http://www.toshiba-medical.co.jp/tmd/products/ihe/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe/">http://www.toshiba-medical.co.jp/tmd/english/ihe/</a> (International page)			
North America: <a href="http://medical.toshiba.com/services/iheintegrationstatements.htm">http://medical.toshiba.com/services/iheintegrationstatements.htm</a>			
Europe: <a href="http://www.toshiba-europe.com/medical/medical.asp?PageID=153">http://www.toshiba-europe.com/medical/medical.asp?PageID=153</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/dicom/">http://www.toshiba-medical.co.jp/tmd/products/dicom/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom/">http://www.toshiba-medical.co.jp/tmd/english/dicom/</a> (International page)		
	North America: <a href="http://medical.toshiba.com/services/dicomconformancestatements.htm">http://medical.toshiba.com/services/dicomconformancestatements.htm</a>		
	Europe: <a href="http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532">http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe-europe.org/">http://www.ihe-europe.org/</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>	

"Aplio" is a trademark of Toshiba Medical Systems Corporation.