

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Toshiba Medical Systems Corporation	Aplio™ XG /SSA-790A	V2.00	8.Mar.2007
This product implements all transactions required in the IHE Technical Framework (Rev 7.0), IHE Cardiology Technical Framework (Rev 2.1), IHE Cardiology Technical Framework Supplement 2006-2007, and IHE IT Infrastructure Technical Framework (Rev 3.0) to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Radiology</b>			
Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Portable Data for Imaging	Portable Media Creator	None	
Evidence Documents	Acquisition Modality	Echo Evidence TID 5000 OB-GYN Ultrasound Procedure Report	
Teaching file and Clinical trial Export	Export Selector	Additional Teaching File Information Delay for Reason	
<b>Cardiology</b>			
Evidence Documents	Acquisition Modality	Echo Evidence TID 5100 Vascular Ultrasound Report TID 5200 Echocardiography Procedure Report	
Echocardiography Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo	
Stress Testing Workflow	Acquisition Modality	Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo	
<b>IT Infrastructure</b>			
Consistent Time	Time Client	None	

<b>Internet address for vendor's IHE information:</b>		
Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/ihe/">http://www.toshiba-medical.co.jp/tmd/products/ihe/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/ihe/">http://www.toshiba-medical.co.jp/tmd/english/ihe/</a> (International page)		
North America: <a href="http://medical.toshiba.com/services/iheintegrationstatements.htm">http://medical.toshiba.com/services/iheintegrationstatements.htm</a>		
Europe: <a href="http://www.toshiba-europe.com/medical/medical.asp?PageID=153">http://www.toshiba-europe.com/medical/medical.asp?PageID=153</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	None	
<b>DICOM</b>	Japan: <a href="http://www.toshiba-medical.co.jp/tmd/products/dicom/">http://www.toshiba-medical.co.jp/tmd/products/dicom/</a> <a href="http://www.toshiba-medical.co.jp/tmd/english/dicom/">http://www.toshiba-medical.co.jp/tmd/english/dicom/</a> (International page) North America: <a href="http://medical.toshiba.com/services/dicomconformancestatements.htm">http://medical.toshiba.com/services/dicomconformancestatements.htm</a> Europe: <a href="http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532">http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Europe: <a href="http://www.ihe-europe.org/">http://www.ihe-europe.org/</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>

“Aplio” is a trademark of Toshiba Medical Systems Corporation.