

| IHE Integration Statement | | | |
|--|---------------------------|--|-------------|
| Vendor | Product Name | Version | Date |
| Toshiba Medical Systems Corporation | Xario™ XG /SSA-680A | V1.00 | 15.Oct.2007 |
| This product implements all transactions required in the IHE Technical Framework (Rev 7.0), IHE Cardiology Technical Framework (Rev 2.1), IHE Cardiology Technical Framework Supplement 2006-2007, and IHE IT Infrastructure Technical Framework (Rev 3.0) to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Radiology | | | |
| Scheduled Workflow | Acquisition Modality | Patient Based Worklist Query Broad Worklist Query PPS Exception Management | |
| Patient Information Reconciliation | Acquisition Modality | No options defined | |
| Portable Data for Imaging | Portable Media Creator | None | |
| Evidence Documents | Acquisition Modality | Echo Evidence TID 5000 OB-GYN Ultrasound Procedure Report | |
| Teaching file and Clinical trial Export | Export Selector | Additional Teaching File Information Delay for Reason | |
| Cardiology | | | |
| Evidence Documents | Acquisition Modality | Echo Evidence TID 5100 Vascular Ultrasound Report TID 5200 Echocardiography Procedure Report | |
| Echocardiography Workflow | Acquisition Modality | Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo | |
| Stress Testing Workflow | Acquisition Modality | Patient Based Worklist Query Broad Worklist Query PPS Exception Management Stress Echo | |
| IT Infrastructure | | | |
| Consistent Time | Time Client | None | |

| | | |
|---|--|--|
| Internet address for vendor's IHE information: | | |
| Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe/ http://www.toshiba-medical.co.jp/tmd/english/ihe/ (International page) | | |
| North America: http://medical.toshiba.com/services/iheintegrationstatements.htm | | |
| Europe: http://www.toshiba-europe.com/medical/medical.asp?PageID=153 | | |
| Links to Standards Conformance Statements for the Implementation | | |
| HL7 | None | |
| DICOM | Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom/ http://www.toshiba-medical.co.jp/tmd/english/dicom/ (International page) North America: http://medical.toshiba.com/services/dicomconformancestatements.htm Europe: http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532 | |
| Links to general information on IHE | | |
| In North America: http://www.ihe.net/ | In Europe: http://www.ihe-europe.org/ | In Japan: http://www.ihe-j.org/ |

“Xario” is a trademark of Toshiba Medical Systems Corporation.