

| IHE Integration Statement | | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|-------------|
| Vendor | Product Name | Version | Date |
| Toshiba Medical Systems Corporation | DI Station™ ADR-2000A | V9.3.1 | 3 Feb. 2006 |
| This product implements all transactions required in the IHE Technical Framework (Rev 6.0) to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Radiology | | | |
| Scheduled Workflow | Acquisition Modality | Patient Based Worklist Query Broad Worklist Query Japanese Character Set | |
| Patient Information Reconciliation | Acquisition Modality | No options defined | |
| Portable Data for Imaging | Portable Media Creator | None | |
| Internet address for vendor's IHE information: | | | |
| Japan: http://www.toshiba-medical.co.jp/tmd/products/ihe/ http://www.toshiba-medical.co.jp/tmd/english/ihe/ (International page) | | | |
| North America: http://medical.toshiba.com/services/iheintegrationstatements.htm | | | |
| Europe: http://www.toshiba-europe.com/medical/medical.asp?PageID=153 | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| HL7 | None | | |
| DICOM | Japan: http://www.toshiba-medical.co.jp/tmd/products/dicom/ http://www.toshiba-medical.co.jp/tmd/english/dicom/ (International page) | | |
| | North America: http://medical.toshiba.com/services/dicomconformancstatements.htm | | |
| | Europe: http://www.toshiba-europe.com/medical/medicallr.asp?PageID=1532 | | |
| Links to general information on IHE | | | |
| In North America: http://www.ihe.net/ | In Europe: http://www.ihe-europe.org/ | In Japan: http://www.ihe-j.org/ | |

“DI Station“ is a trademark of Toshiba Medical Systems Corporation.