Canon

No. 2G984-0212EN

PRODUCT SECURITY PATCH LIST
FOR
ULTRASOUND IMAGE ANALYSIS PROGRAM

(2G984-0212EN**)**

CANON INC.

© CANON INC. 2026
ALL RIGHTS RESERVED

Trademarks

Windows® is a registered trademark of Microsoft Corporation in the United States and other countries.
This document may also include trademarks or registered trademarks of other companies.
Note that the trademark symbol "™" and the registered trademark symbol "®" may or may not be used in this document.

IMPORTANT!

1. No part of this document may be copied or reprinted, in whole or in part, without written permission.
2. The contents of this document are subject to change without prior notice and without legal obligation.
3. The contents of this document are correct to the best of our knowledge. Please inform us of any ambiguous or erroneous descriptions, missing information, etc.

REVISION RECORD

REV.	DATE (YYYY-MM)	REASON /AUTHOR	PAGE CHANGED	SER. No.	DOC. PRODUCT.
INI.	2026-10	Konno	-------		BU

Security Patch Package name

NA

Applicable Products for the Security Patch

CUW-U001S	V1.0 SP0003* - SP0006*
	V2.0 SP0000* - SP0005*

Steps to Apply the Security Patch

Please follow the instruction of the administrator of the facility or organization. For details, please contact the customer service for each system.

Verified Security Patch List

KB	MSRC	explanation(en-us)	URL/File Name	url
KB2565063_x64	MS11-025	Security Update for Microsoft Visual C++ 2010 Service Pack 1 Redistributable Package (KB2565063)	vcredist_x64.exe	https://catalog.s.download.windowsupdate.com/msdownload/update/software/secu/2011/07/vcredist_x64_15d032d669078aa6f0f7fd1cbf4115a070bd034d.exe
VS2013_x64	n/a	Microsoft Visual C++ 2013 Update 5 Redistributable Package	vcredist_x64.exe	https://download.microsoft.com/download/c/c/2/cc2df5f8-4454-44b4-802d-5ea68d086676/vcredist_x64.exe
VS2022_x86	n/a	Microsoft Visual C++ 2015-2022 Redistributable (x86)	VC_redist.x86.exe	https://download.visualstudio.microsoft.com/download/pr/9c69db26-cda4-472d-bdae-f2b87f4a0177/A32DD41EAAB0C5E1EAA78BE3C0BB73B48593DE8D97A7510B97DE3FD993538600/VC_redist.x86.exe
VS2022_x64	n/a	Microsoft Visual C++ 2015-2022 Redistributable (x64)	VC_redist.x64.exe	https://download.visualstudio.microsoft.com/download/pr/1754ea58-11a6-44ab-a262-696e194ce543/3642E3F95D50CC193E4B5A0B0FFBF7FE2C08801517758B4C8AEB7105A091208A/VC_redist.x64.exe
KB5129195	n/a	2026-09 Cumulative Update for Windows 11, version 24H2 for x64-based Systems (KB5129195) (26100.9457)	windows11.0-kb5129195-x64_ed361878ec2b56a7dfdb8f256565263a7fdc8eaa.msu	https://catalog.sf.dl.delivery.mp.microsoft.com/filestreamingservice/files/471c8f54-b9fc-491d-9944-234d5db009f0/public/windows11.0-kb5129195-x64_ed361878ec2b56a7dfdb8f256565263a7fdc8eaa.msu

KB	MSRC	explanation(en-us)	URL/File Name	url
KB5126052	n/a	2026-09 Cumulative Update for .NET Framework 3.5 and 4.8.1 for Windows 11, version 24H2 for x64 (KB5126052)	windows11.0-kb5126052-x64-ndp481_082cfd5816914c794b261ea68bf6c91419fc32a8.msu	https://catalog.sf.dl.delivery.mp.microsoft.com/filestreamingservice/files/aa6d5327-d621-4af8-8340-5a291e597f6b/public/windows11.0-kb5126052-x64-ndp481_082cfd5816914c794b261ea68bf6c91419fc32a8.msu
Edge_x64_154.0.4258.48	n/a	Microsoft Edge-Stable Channel Version 154 Update for x64 based Editions (Build 154.0.4258.48)	microsoftedgeenterprisex64.cab	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/updt/2026/09/microsoftedgeenterprisex64_f6b4d865aa6e24c4b27cc4d8f00a2e8cc506b21b.cab
KB890830_v5.145	n/a	Windows Malicious Software Removal Tool x64 - v5.145 (KB890830)	windows-kb890830-x64-v5.145.exe	https://catalog.s.download.windowsupdate.com/d/msdownload/update/software/uprl/2026/09/windows-kb890830-x64-v5.145_860b493566a29282606ed075513294c239b1a545.exe
MSTeams-x64.msix_26246.1604.5133.838	n/a	n/a	MSTeams-x64.msix	https://statics.teams.cdn.office.net/production-windows-x64/enterprise/webview2/lkg/MSTeams-x64.msix
teamsbootstrapper.exe	n/a	n/a	teamsbootstrapper.exe	https://statics.teams.cdn.office.net/production-teamsprovision/lkg/teamsbootstrapper.exe

LEGAL MANUFACTURER



CANON INC.

30-2, SHIMOMARUKO 3-CHOME, OHTA-KU, TOKYO 146-8501, JAPAN